

Softform Premier Active Cushion



Static Dynamic Antidecubitus Cushion



The unique design of the Invacare **Softform Premier Active Cushion** features an alternating air insert beneath a castellated foam insert. The design of the foam insert follows the same principle of the **Softform Premier Mattress**. The cushion retains its properties as an at risk static cushion, but should a patient require stepping up to a dynamic surface, then a discreet pump can be fitted to the air insert, transferring the cushion into an alternating surface delivering additional levels of pressure relief. When the **Active Cushion** is accompanied by a dynamic mattress, patients can be offered a 24hr support surface.

How it works:



To step up a patient, retract the pipes from the cushion.



Fit pipes to the Active pump. Turn pump on – no additional adjustments required.



Patient is now on an alternating surface.



Compatible with Active Mattress

The cushion uses the same pump as the **Premier Active mattress**. When both products are combined, patients are offered a 24 hr support system and continuity of care.



Sensory software within the pump can gauge patient weight

Pressure can be maintained at a level appropriate to individual patients. Fully alarmed in case of power failure.

An innovative static dynamic cushion








The air insert contained within the *Invacare Softform Premier Active Cushion* allows stepping up and down between static and dynamic alternating systems with minimum patient disruption. The air cells work on a two up, two down function beneath the insert, delivering a very low impact at the point of change, thereby delivering high levels of comfort in addition to improved pressure reduction.

Air insert within Softform Premier Foam

Allows "two up, two down" alternating cycle stepping up and down between static and dynamic alternating systems.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Width	Length	Height	Total product weight	Max. user weight	Max. cleaning temperature	Warranty
						
460 mm	510 mm	152 mm	1.5 kg	108 kg	80°	Cover: 12 months Foam: 12 months Pump: 24 months Air Insert: 12 months

Softform Premier Active Cushion


Manufactured to European Standards: Foam: BS 5852: Crib 5; Cover BS7175: Crib 5.

Colour

Fuchsia



QR Code



Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Benkenstrasse 260
4108 Witterswil
Switzerland

