français

Outil de mesure KSL

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
Measuring RW/RW messen/Mesurer RW

KSL Measuring Tool

The KSL measuring tool is a device designed to determine the tippiness of your wheelchair.

Measuring the Rearwheel Position

- Position the handle ① of the measuring tool ③ approx 2cm above the seat in the centre of the backrest tube (SHh + 2cm). It must be in-plane ② with the backrest (see *Top view*)
- Adjust the lines marked Pos. 1 – 5 to the centre of the wheel (see *Side view*).
- Adjust the template so the air bubble in the spirit level is centered between both aid lines.
- Read off the template which position the rear wheel axle has been adjusted to. Pos. 5 equals the minimum and Pos.1 the maximum tippiness.

If the rear wheel axle is between two positions we advise to round up to the next highest value, to receive a higher stability.

In addition the measuring tool can be used to measure the backrest angle of your wheelchair.

Measuring the Backrest Angle (RW)

- Remove the seat cushion and backrest upholstery.
- Lay the tool flush with the seat frame. The handle should be on the seat area, the angle lines (74° - 90°) pointing to the backrest (see *Measuring RW*).
- Determine which of the angle lines (74° - 90°) runs parallel to the backrest tube. This is the backrest angle (RW) of your wheelchair.

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