Yes, you can.®

Invacare® **Birdie™**

Invacare® **Birdie™ Plus**

- EN** *User's Manual*
- FR** *Manuel d'utilisation*
- ES** *Manual del usuario*
- DE** *Gebrauchsanweisung*
- NO** *Brukerveiledning*
- FI** *Käyttöönje*
- NL** *Gebruiksaanwijzing*
- DA** *Brugermanual*
- SV** *Bruksanvisning*
- IT** *Manuale d'uso*
- PT** *Manual de Utilização*

EN	<i>User's Manual</i>	4
FR	<i>Manuel d'utilisation</i>	20
ES	<i>Manual del usuario</i>	36
DE	<i>Gebrauchsanweisung</i>	52
NO	<i>Brukervejledning</i>	68
FI	<i>Kaäyttöönje</i>	84
NL	<i>Gebruiksaanwijzing</i>	100
DA	<i>Brugermanual</i>	116
SV	<i>Bruksanvisning</i>	132
IT	<i>Manuale d'uso</i>	148
PT	<i>Manual de Utilizaçao</i>	164

UK *User's Manual*

Stamp of the Distributor

Table of Contents

1. General	6
2. Main parts of the hoist.	8
3. Unpacking, assembly and disassembly of the hoist	9
4. Operating the hoist	10
5. Emergency lowering, emergency lifting and emergency stop .	13
6. Accessories	14
7. Maintenance and cleaning	14
8. Waste disposal.	15
9. Trouble-shooting	16
10. Technical specifications	17
11. Symbols	19

Birdie™, Birdie™ Plus and Birdie™ Compact

Thank you for choosing Invacare®'s **Birdie™** mobile hoist. Invacare®'s **Birdie™** are especially designed to lift the patient to and from wheelchairs, beds, toilets and floors. **Birdie™** hoists offer easy handling and increased comfort within the nursing sector. The design of **Birdie™** provides a longer reach, wider leg span and a wider lifting range.

The choice of either electrical or manual leg span combined with the large range of spreader bars and slings available from Invacare® makes it easy to adjust the hoist for optimum patient care.

Throughout the entire production process, our materials and components are quality controlled by the operators. A final test is made when the product is fully assembled. If the product does not correspond to the quality demands of Invacare®, it will not pass quality control.

In the event of a problem in connection with the delivered product, please contact your local Invacare® supplier.

Invacare® will take no responsibility if the product is used or assembled in any way other than stated in this user manual. Only accessories mentioned in this user manual may be used on the **Birdie™** mobile hoists. This product must only be operated by qualified staff who has received the necessary instruction and training.

Please read this user manual carefully before you are using the mobile hoist.

1. General

- **Birdie™, Birdie™ Plus and Birdie™ Compact** are CE-marked in accordance with directive 93/42/EEC concerning medical devices.
- **Birdie™, Birdie™ Plus and Birdie™ Compact** have been tested and approved according to EN/ISO 10535.
- Control unit and motors have been approved according to EN 60601.
- **Birdie™, Birdie™ Plus and Birdie™ Compact** have undergone a risk analysis according to EN/ISO 14971.
- Control unit, motors and hand control are IPX4 protected.
- Lifting capacity: **Birdie™** and **Birdie™ Plus** 170 kg, **Birdie™ Compact** 150 kg.
- The hoist has an expected lifetime of 8 years.
- The hoist will not operate when in charge mode.

Disconnect the charger from the hoist before moving or using the hoist. Make sure the mains cable is not squeezed or damaged when assembly or disassembly of the hoist or when moving the hoist.

If the functions of the hoist change, please see the section "Maintenance and cleaning". This product must be serviced and maintained by qualified personnel.

Please be aware of the following possible risks when using the hoists:

The hoist must be carefully supervised if the patient is a child.

Please note that it is more difficult to manoeuvre the hoist on carpeted or uneven/textured surfaces as opposed to wooden or flat surfaces; this is due to the starting forces required to initiate the movement that may create additional friction exerted on the castors.

It is recommended not to move the user with the jib at maximum height. Take great care when moving the hoist on carpeted, wet, slippery, rough or uneven surfaces. Never use on slopes.

Caution! *Never attempt to lift a patient outside the base area of the hoist. When lifting a patient from the floor, the patient must always be located between the two hoist legs, and never outside this area. To avoid injury, utmost care must be taken when raising and lowering a patient.*

*Invacare® accepts no liability for any use, change or assembly of the **Birdie™**, **Birdie™ Plus** and **Birdie™ Compact** hoists other than as stated in this user manual.*

There is a risk of squeezing and entrapment during operation and transport of the hoist.

Never exceed the maximum recommended lifting capacity.

It is important that the cables do not become entrapped or jammed in moving parts and castors during operation and transport.

Important! *Invacare® recommend that operation of any mobile lifting equipment is carried out by personnel with prior professional training in moving and handling techniques.*

Always choose the sling design and size according to the patient's weight, size and physical ability. Always ensure, that the sling is applied correctly. Never leave a patient in the sling unattended.


Please ensure that the hoist is not exposed to water, as this can cause damage.



Locking pin must be located into the base position when transporting the hoist. Failure to do so may cause injury during handling.

2. Main parts of the hoist

Birdie™ Plus

3. Unpacking, assembly and disassembly of the hoist

- The packaging is developed for optimal protection of the hoist.
- If you return the hoist, it must be packed as delivered in the original box.
- You can order replacement boxes from Invacare®.

Unpacking


1. Should the packaging be damaged upon receipt, each part of the hoist must be examined for visible defects. In the event of any damage, please contact Invacare®.
2. Carefully take out and identify all the parts in the box.
3. The box contains the following parts:
 - 1 hoist
 - 1 mains cable
 - 1 user's manual
 - 1 spreader bar
 - 1 manual leg spreader handle (optional)
 - 1 sling (optional)

To prevent self-discharging of the battery or accidental operation of the hoist, the emergency stop is activated during shipping. Please remember to deactivate the emergency stop and charge the battery prior to first use.

Assembling of the hoist

The hoists are easy to assemble without the use of tools.

1. Perform unpacking and assembly operation at floor level
2. Take out all the parts of the box (picture 1-4).
3. Apply castor brakes before assembling hoist.
4. Remove the locking pin and raise the mast of the hoist by pulling the handle bar (picture 5-6).
5. Insert the locking pin (picture 7 and 12-13 below). Ensure that the locking pin is correctly inserted as the picture 12.
6. Loosen the spreader bar by pulling it downwards (picture 8).
7. Remove the pipe pin. Assemble the motor and the jib by using the pipe pin (picture 9-10).
8. Correctly assembled hoist (picture 11).



Activate the emergency stop before assembly or disassembly to prevent entrapment/squeezing.
Remove optional leg spreader handle before disassembly.

There is a possible risk of squeezing of both limbs and wires during assembly and disassembly of the hoist. Take the utmost care when lifting components during assembly - some parts are heavy. Always remember to adopt the correct lifting position.

Ensure that the locking pin is correctly inserted.

Disassembling of the hoist

1. Lower the jib and narrow both legs completely, and activate the emergency stop button.
2. Apply castor brakes before disassembling hoist.
3. Remove the pipe pin and the motor piston from the jib, reinstall the pipe pin in the piston end, and lock the motor into the clips on the mast.
4. Attach the spreader bar into the welded fork on the mast.
5. Remove the locking pin from the base of the mast, release the safety latch, lower the mast, and relocate the locking pin into the mast near the suspension axle of the mast.

The hoist can now be located in the packaging box, pulled on the rear wheels, or parked in an upright position with the mast/jib assembly pointing upwards

4. Operating the hoist

Lifting capacity

Birdie™ and **Birdie™ Plus**: 170 kg.

Birdie™ Compact: 150 kg.

Operation of the manual leg spread

The manual leg span is operated by the 2 pedals on the leg section of the hoist. The pedals are operated from one side to the other to adjust the distance between the legs of the hoist.

The hand control has 2 buttons – one to raise and one to lower the jib.

1. Stand behind the hoist and take a firm hold of the handle bar using both hands. Using your foot, press the left pedal on the base to increase the distance between the legs.
2. Using your foot, press the right pedal to decrease the distance between the legs.


Operation of the electric legs

The distance between the legs is controlled by 2 electric motors controlled by the hand control.

The hand control has 4 buttons: 2 for raising and lowering the jib and 2 for operating the electric legs.

Using the hand control

1. Press "arrow up" on the hand control to raise the jib.
2. Press arrow down on the hand control to lower the jib.
3. Operation of electric legs.

Operating forces of buttons: Max. 5 N.

LED diodes on the control unit

Birdie™ and **Birdie™ Compact** - the upper yellow diode will blink during charging, and switch to continuous light when fully charged. The lower green diode will light continuously when the control unit is connected to the mains, and light up when any button on the remote control is pressed, or when the electric emergency lowering is activated.

Birdie™ Plus - the left green diode will light continuously when the control unit is connected to the mains. The yellow right diode will light continuously during charging, and go out when the battery is fully charged. There is no indication when hand control buttons are pressed, or when the emergency lowering/raising system is activated.

Brakes

1. Stand behind the hoist and take a good hold of the handle bar.
2. Press the pedals downwards by feet on the back castors to activate the brake.
3. Press the pedals upwards by feet to release the brake.

Moving the hoist

1. Stand behind the hoist and take a good hold of the handle bar.
2. Release the brakes of the back castors.
3. Now the hoist can be pushed or pulled to the desired location.

Note: To negotiate raised obstacles during movement of the unloaded hoist, it is recommended to pull the hoist backwards. This operation should only be performed when the hoist is in the unloaded condition.

Tilt pedal

The tilt pedal is used for tilting when the hoist must be moved from one room to another, to force a footstep. Place one foot on the tilt pedal, and pull the empty hoist rearwards to a suitable angle.

WARNING! Never use the tilt pedal with a patient in the hoist.

Mounting the spreader bar

The spreader bar is easy mounted and dismantled on the carabineer hook:


Open the carabineer hook with one finger, and mount or dismount the spreader bar.

Use only spreader bars made for the **Birdie™** mobile hoists.

Make sure the spreader bar is suitable for the patient and the actual lift or transfer required.

Carabineer correctly positioned on the lifting arm, ready for use.

Correctly attached spreader bar.

Easy opening of the carabineer by pushing the safety catch backwards with one finger and attachment of the spreader bar.

Fitting the sling

Slings designed for 2- or 4-point spreader bars can be used with this mobile hoist, please refer to Invacare®'s sling brochure for details.

Note: The size of the sling and the width of the spreader bar should correspond.

1. Choose the most appropriate sling for the patient. The carer must be trained in the use of hoists before performing any transfers.
2. Carefully inspect the hoist and the sling. If defective parts are found, the damaged parts must be replaced.
3. Only use slings suited for the hoist and the patient.
4. Place the patient in the sling as described in the sling manual. Pay special attention to the fact that the sling is correctly applied.
5. Lower the jib and mount the straps of the sling. If the patient is in a sitting position, the jib must be lowered to chest height and the hoist must be moved as close as possible to the patient.
6. Make sure that there are no obstructions when raising the jib. Ensure that brakes are not engaged before lifting or lowering.

Caution! Do not lift or lower a patient with the brakes applied. Always let the hoist find the correct centre of gravity. Make sure that the lifting area is between the legs of the hoist. **Never** operate the hoist when the patient is outside the base and leg area. Always be careful when lifting or lowering a patient. Never exceed the maximum recommended lifting capacity.

Charging the battery

The hoists are equipped with an internal charger. It is recommended to charge the batteries regularly to ensure optimal use of the hoist and prolong the life of the batteries. Furthermore, it is recommended to charge the batteries before first use.

The control unit is equipped with a sound signal, which will beep when operating with low battery capacity. It is recommended to charge the batteries as soon as the sound signal is heard. In addition, **Birdie™ Plus** is equipped with a display, indicating the approximately remaining battery capacity.

How to charge the batteries:

1. Connect the mains cable to the control unit and plug it in. It takes about 4 hours to charge the batteries. Charging must take place in a room with good air ventilation. The charger stops automatically when the batteries are fully charged.
2. Remember to disconnect the charger cable before using the hoist again.

The emergency stop must not be activated - otherwise it will be impossible to charge the batteries. While charging takes place the hoist cannot be used. Do not use or move the hoist without unplugging from the socket outlet. Do not attempt to use the hoist if the battery housing is damaged. Replace a damaged battery housing before further use.

5. Emergency lowering, emergency lifting and emergency stop

Electric emergency lowering Birdie™, Birdie™ Compact and Birdie™ Plus

If the hand control fails, the jib can be lowered by using the circular switch for emergency lowering. This is located at the front of the control unit. The jib will be lowered as long as the button is pressed. On Birdie™ Plus the button should be operated by means of a sharp object, such as a pencil.

Birdie™ and Birdie™ Compact


Emergency lowering

Emergency lifting

Birdie™ Plus

Jumbo Care Control Unit

Electric emergency lifting - Birdie™ Plus

If the hand control fails, the jib can be raised by using the button for emergency lifting. This is located at the front of the control unit, and should be operated with a sharp object, such as a pencil. The jib will be lifted as long as the button is pressed.

Manual emergency lowering - Birdie™ and Birdie™ Plus

In case of partial or total power failure, or if the battery runs down while using the hoist, Birdie™ and Birdie™ Plus is equipped with a manual emergency lowering system located at the bottom of the motor. The weight is preset for 75 kg. If the patient's weight is more or less than 75 kg, it might be necessary to adjust the manual emergency lowering.


Adjustment of the manual emergency lowering - Birdie™ and Birdie™ Plus

If the lowering speed has to be adjusted according to the patient's weight, then the screw in the red release handle must be adjusted. Loosen to increase the speed and tighten to decrease the speed.

Note: Manual emergency lowering is only possible when a patient is sitting in the hoist.

Emergency stop

In case of emergency: Press the red button for emergency stop. The emergency stop can be reset by turning the stop button

Note! The mechanical emergency lowering system will only operate when there is a patient in the hoist. The mechanical emergency lowering system has to be adjusted according to the patient's weight.

6. Accessories

It is compulsory to use original Invacare® spare parts which you can buy through any Invacare® dealer. For repair, please contact your local Invacare®dealer.

Sling : See the sling brochure

7. Maintenance and cleaning

At normal daily operation, a service check-up should take place every year, according to the maintenance chart. When performing annual or regular maintenance, all parts designed to carry load must be as a minimum tested with maximum load. All safety features must be checked according to EN ISO 10535: 2006 Annex B.

LOLER statement.

The UK Health and Safety Executive's Lifting Operations and Lifting Equipment Regulations 1998, require any equipment that is used in the workplace to lift a load be subject to safety inspection on a six monthly basis. Please refer to the HSE web site for guidance www.hse.gov.uk

Jumbo care functions.

In add for a safe maintenance, the JUMBO Care offers five (5) main functions:

- 1) Service intervals*
- 2) Service information on PC*
- 3) Actuator life time*
- 4) Resetting of service interval*
- 5) Actuator exchange*

(For more information, see technical Jumbo care user manual or ask your local Invacare dealer)

It is the responsibility of the person responsible for the equipment to ensure that LOLER regulations are adhered to.

The hoist must be stored at normal room temperature. If it is stored in a damp, cold or wet environment then the motor and other mounting parts may be prone to corrosion.

- The hoist have been designed so that it requires a minimum amount of maintenance.*
- It is recommended that the batteries are charged frequently, if possible every night, as it will extend the lifetime of the batteries and will ensure high performance.*
- It is recommended that the hoist is cleaned after use.*
- Please check that the spreader bar connections, control unit and motors are in good working order and not damaged in any way. If there is any damage to the hoist, refrain from using it, and contact your local Invacare®dealer before further use.*

Wear and damage of loading parts

It is necessary to maintain and check all parts that are exposed to static or dynamic strain, such as the sling, the spreader bar and the turning points, do not have fractures, are fragile, askew or damaged.

Frequent use of the hoist will cause the carabiner to ware under the contact points of the spreader bar (b) and the mast pin (a). Do not use the hoist if the thickness of the carabineer measures less than 6 mm at these points. (See maintenance chart page 15)

Damaged parts must be replaced before further use.

Cleaning

- To prevent cross-infection, the hoists must be cleaned after use.*
- The hoist must be wiped with a moist firmly wrung cloth with ordinary household disinfectants.*
- Only use officially approved disinfection detergents.*
- Dry the hoist carefully after cleaning.*
- Never use acids, alkaline or solvents for cleaning the hoist.*

Motors, control unit and mounting parts can be destroyed if the hoist is cleaned in any other way than stated above.

Maintenance chart

Only personnel having received the necessary instruction or training by Invacare® must perform service and maintenance on Birdie™, Birdie™ Plus and Birdie™ Compact.										
Serial number (located on the hoist base):										
Date:	//	//	//	//	//	//	//	//	//	//
Initials:										
Visual inspection of all hoist parts										
Inspection of the spreader bar and carabiner - welding and assemblies										
Inspection of carabineer* (the thickness "a" and "b" mustn't measure less than 6 mm)			"a" =	"a" =	"a" =	"a" =	"a" =	"a" =	"a" =	"a" =
			"b" =	"b" =	"b" =	"b" =	"b" =	"b" =	"b" =	"b" =
Inspection of the motors - undamaged										
Inspection of cables and plugs - undamaged										
Inspection of emergency stop, electrical emergency lowering and -raising (only Birdie™ Plus)										
Inspection of the castors - attachment, braking ability and free rolling										

Invacare® offers courses in service and maintenance of the hoist.

For all servicing and maintenance requirements, please contact your local Invacare®dealer.

8. Waste disposal

This product has been supplied from an environmentally aware manufacturer that complies with the Waste Electrical and Electronic Equipment (WEEE) Directive 2002/96/CE.

This product may contain substances that could be harmful to the environment if disposed of in places (landfills) that are not appropriate according to legislation.

The 'crossed out wheelie bin' symbol is placed on this product to encourage you to recycle wherever possible.

Please be environmentally responsible and recycle this product through your recycling facility at the end of it's life.

Pb

Hg

Cd

9. Trouble-shooting


Only personnel having received the necessary instruction or training by Invacare® must perform service and maintenance on **Birdie™**, **Birdie™ Plus** and **Birdie™ Compact**.

Symptom	Possible cause	Remedy
<i>The mast seems loose</i>	<i>The mast axis screw has become loose</i>	<i>Insert the locking pin into the mast axis screw, and tighten the screw</i>
<i>Castors are noisy</i>	<i>Dust or dirt in the castors</i>	<i>Clean the castors</i>
<i>Hinge between jib and mast is noisy during operation</i>	<i>Missing lubrication</i>	<i>Lubricate the hinge</i>
<i>The motor is not running</i>	<i>Hand control or motor is not connected</i>	<i>Insert the plugs fully</i>
	<i>No power on the battery</i>	<i>Charge or replace the battery</i>
	<i>The red emergency stop is activated</i>	<i>Turn the red button clockwise to deactivate</i>
<i>Motor noise, but no movement in the piston rod</i>	<i>Motor is damaged</i>	<i>Replace the motor</i>
<i>The piston rod only moves in and not out</i>	<i>Motor is damaged</i>	<i>Replace the motor</i>
<i>The control unit emits a beeping sound during lifting, and the jib motor stops (Birdie™ and Birdie™ Compact). The jib motor stops (Birdie™ Plus)</i>	<i>Max. load is exceeded</i>	<i>Reduce the load (and the hoist will function normally)</i>

Contact your dealer if the above does not solve your problems.

Lubricate with medically clean oil, e.g. Kemtura Kem Lub KEM-WO 50, order no. 813239.

10. Technical specifications

	General specifications		Manual and electric Birdie™		Birdie™ Plus	Manual Birdie™ Compact	
			Ø 75 mm	Ø 100 mm	Ø 100 mm + Ø 125 mm	Ø 75 mm	Ø 100 mm
	Castor						
	Maximum lifting capacity	170 kg		170 kg		150 kg	
	Lowest position (min. height of CSP**) 3*	r'	65,5 cm	67 cm	67 cm	73 cm	74,5 cm
	Lifting area (height range)		46,5 - 168,5 cm	47,5 - 170 cm	47,5 - 170 cm	51,5 - 158,5 cm	53 - 160 cm
	Maximum height of CSP**	k'	189 cm	190,5 cm	190,5 cm	180 cm	181,5 cm
	Leg length		115 cm		116,5 cm	100 cm	
	Minimum internal width	r'	58 cm		58 cm	46 cm	
	Total width (closed), external measure		65,5 cm	68 cm	68 cm	53 cm	55 cm
	Internal width at maximum reach	q'	91 cm		91 cm	78 cm	
	Total width (open) internal measure		104 cm		104 cm	89 cm	
	Total width (open) centre to centre of castor	p'	108,5 cm		108,5 cm	93 cm	
	Turning radius		140 cm		140 cm	107 cm	
	Weight, mast, incl. battery, excl. spreader bar		21 kg		21,5 kg	17,5 kg	
	Total weight incl. spreader bar		42 kg		44 kg	36 kg	
	Weight, leg section *		19 kg		21 kg	16,5 kg	
	Height to upper edge of legs		10 cm	11,5 cm	11,5 cm	10 cm	11,5 cm
	Min. free height		2 cm	3,5 cm	3,5 cm	2 cm	3,5 cm
	Min. space for patient (to motor) in top position		34 cm		34 cm	30 cm	
	Maximum reach at 60 cm	a'	66,5 cm		66,5 cm	54,5 cm	
	Maximum reach from base	b'	66,5 cm		66,5 cm	54,5 cm	
	Reach from base with legs spread to 70 cm	c'	24 cm		24 cm	41,5 cm	
	Hoisting range	m'	123,5 cm		123,5 cm	107 cm	
	Voltage output		24 V ⁼⁼⁼ max. 250 VA				
	Voltage supply		100 -240 V ===50/60 Hz				
	Max. current input		Max. 200 mA		Max. 400 mA		Max. 200 mA
	Operating temperature.		5° to 40° C				
	Air humidity		20% to 90% @ 30° C - not condensing				
	Sound pressure		45-50 dB(A)				
	Atmospheric pressure		700 to 1060 hPa				
	Protection class, control unit		IPX4				
	Protection class, hand control		IPX4				
	Protection class, motor		IPX4				
	Insulation class		II - type B				
	Working ability		40 full lifts without battery charge with batteries at 50% of full capacity				
	Intermittents		10%, max. 2 minutes/18 minutes				
	Battery capacity		2,9 Ah		2,9 Ah + display	2,9 Ah	
	Manual emergency lowering		Yes		Yes	No	
	Electric emergency lowering/ lifting		Yes/No		Yes/Yes	Yes/No	

All measurements are given including a 45 cm 2 point spreader bar, except dimensions marked with *.
Invacare® reserves the right to change the measurements without warning.

** CSP: Central Suspension Point

1: Highest position, 2: Maximum reach position, 3: Lowest position

11. Symbols

The patient is not separated from the ground and the chassis:

Direct current:

Alternating current

Double insulated:

X kg max. load (SWL) = (Patient + sling + spreader bar).

X = 150 kg (**Birdie™ Compact**)

X = 170 kg (**Birdie™ and Birdie Plus™**)

The product should be reused where possible:

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Invacare® A/S
Invacare® Aquatec GmbH
Invacare® Poirier SAS
Invacare® Ireland Ltd
Invacare® AB
Invacare® B.V.