microAIR®

This manual **MUST** be given to the user of the product.
BEFORE using this product, read this manual and save for future reference.

Invacare® MA50CE Series

MA50CE Alternating Pressure System User Manual	3	EN
MA50CE Wechseldrucksystem Gebrauchsanweisung	31	DE
Système à pression alternée MA50CE Manuel d'utilisation'	64	FR
Sistema de presión alterna MA50CE Manual del usuario	94	ES
Sistema a pressione alternata MA50CE Manuale d'uso	124	IT
Sistema de Alternância de Pressão MA50CE Manual de utilização	155	PT
MA50CE växeltryckssystem Bruksanvisning	185	SV
MA50CE vekseltrykksystem Brukerveiledning	216	NO
MA50CE-vekseltrykksystem Brugsanvisning	244	DA

Yes, you can.®

© 2012 Invacare Corporation. All rights reserved. Reproduction, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

User Manual

Invacare® MA50CE Series

This manual MUST be given to the user of the product.

BEFORE using this product, read this manual and save for future reference.

MA50CE Alternating Pressure System

1	GENERAL	4
	Symbols	4
2	OVERVIEW	5
	Label Location.....	5
	Typical Product Parameters	6
	Intended Use.....	9
3	GENERAL GUIDELINES	10
	General Guidelines	10
	EMC Information.....	14
4	OPERATION	15
	Installing the Mattress Overlay System (MA50CE).....	15
	Installing the Power Unit.....	16
	Connecting the Hose.....	16
	Connecting the Power Cord.....	18
	Using the Front Panel.....	19
	Powering Up the System.....	23
	Placing the Patient on the Overlay.....	23
	Transferring Patient From/To a Trolley	24
	Transferring Patient From/To a Wheelchair	24
	Preparing for CPR Procedure.....	25
	About Power Outage and Transportation.....	25
5	MAINTENANCE AND TROUBLESHOOTING	26
	Cleaning and Care	26
	Troubleshooting.....	30

I General

This product complies with Directive 93/42/EEC concerning medical devices.
The launch date of this product is stated in the CE declaration of conformity.

I.1 Symbols

Warnings

Signal words are used in this manual and apply to hazards or unsafe practices which could result in personal injury or property damage. See the information below for definitions of the signal words.

DANGER

Danger indicates an imminently hazardous situation which, if not avoided, will result in death or serious injury.

CAUTION

Caution indicates a potentially hazardous situation which, if not avoided, may result in property damage or minor injury or both.

WARNING

Warning indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

Gives useful tips, recommendations and information for efficient, trouble-free use.

IMPORTANT

Indicates a hazardous situation that could result in damage to property if it is not avoided.

2 Overview

2.1 Label Location

ETL LISTED
CONFORMS TO
IEC STD 60601-1
IEC STD 60601-1-2
EN STD 60601-1
ANSI/AAMI STD ES60601-1
CERTIFIED TO
CAN/CSA STD C22.2 NO. 60601-1

Fuse Label

Refer to Typical Product Parameters on page 6 for icon definitions


 CARELESSNESS CAUSES FIRE GENERAL CLEANING Wipe down using soap & water/mild detergent solution or remove cover for laundering at 60 deg C. Do not use bleach. FOR DECONTAMINATION REFER TO YOUR INFECTION CONTROL SPECIALIST. Invacare Operations UK Ltd Pencoed, Bridgend CF35 5AQ Tel: +44 1656 776200 Fax: +44 1656 776201 www.invacare.com Manufactured for: Invacare Operations UK Ltd Made in the USA	 CARELESSNESS CAUSES FIRE GENERAL CLEANING Wipe down using soap & water/mild detergent solution or remove cover for laundering at 60 deg C. Do not use bleach. FOR DECONTAMINATION REFER TO YOUR INFECTION CONTROL SPECIALIST. Invacare Operations UK Ltd Pencoed, Bridgend CF35 5AQ Tel: +44 1656 776200 Fax: +44 1656 776201 www.invacare.com Manufactured for: Invacare Operations UK Ltd Made in the USA
OR	
MA50CE Overlay Size: 200x88x13cm Date: _____ Serial No: _____ Q.C.: _____  Cover please BS 7176 Crib 5 Cigarette and EN957 - 2 Match Equivalent. HMI No : 69681 Minimum User Weight : 15kgs	MA50CE Overlay Size: 200x83x13cm Date: _____ Serial No: _____ Q.C.: _____  Cover please BS 7176 Crib 5 Cigarette and EN957 - 2 Match Equivalent. HMI No : 69680 Minimum User Weight : 15kgs

INVACARE CORPORATION	
Alternating Pressure microAir 50CE	
⚠ DANGER-EXPLOSION HAZARD: DO NOT use in the presence of flammable anesthetics. RISQUE D'EXPLOSION; NE PAS employer en présence d'anesthésiques inflammables. ⚠ WARNING: Equipment should be connected to a properly grounded receptacle. Risk of electrical shock. DO NOT remove back. Disconnect air hose before administering CPR.	 5 055050 401418 > MODEL / REF: MA50CE S/N: 12D KP500001
K-0 Auto Sensing Input 90/120 VAC, 60Hz, 1A 220/240 VAC, 50Hz, 0.5A IEC 60601-1-2:2007 3rd Ed.	MANUFACTURED FOR: Invacare Operations UK Ltd, Pencoed, Bridgend, CF35 5AQ UK. Tel: +44 1656 776200 Made In U.S.A.

2 OVERVIEW

Refer to [Using the Front Panel](#) on page 19, for button descriptions

2.2 Typical Product Parameters

Electrical Parameters

If the unit is unplugged at the power inlet it will be disconnected from the main power.

INPUT VOLTAGE AC:	230 V
INPUT FREQUENCY:	50 Hz
CURRENT:	1 A
MAXIMUM POWER CONSUMPTION:	5 ± 3 W
CIRCUIT PROTECTION:	Dual fused, 250 V, 1 A time delay fuses
FUSE TYPE:	Bussmann #GMD-I-R
BREAKING CAPACITY (BRK. CAP.):	@125 VAC is 10kA, @ 250 VAC is 35A

MODE OF OPERATION:	Continuous
COMPLIANCE:	Class I equipment, conforms to UL STD 60601-I ETL certified as complying with CAN/CSA STD C22.2 NO. 60601-I

Environmental Parameters

Disposal Requirements: This equipment should be disposed of at your local Recycling Center (Non-hazardous waste) when it has reached the end of its service life.

OPERATING CONDITIONS AMBIENT TEMPERATURE: RELATIVE HUMIDITY: ATMOSPHERIC PRESSURE:	10° - 40° C 30% - 75% Non-Condensing 70 - 106 kPa
STORAGE AND SHIPPING CONDITIONS AMBIENT TEMPERATURE: RELATIVE HUMIDITY: ATMOSPHERIC PRESSURE:	-40° to 70° C 10% - 100% Non-Condensing 50 - 106 kPa

Performance Parameters

WEIGHT CAPACITY PRESSURE ZONE: MAXIMUM FLOW: MAXIMUM FLOW PRESSURE: MAXIMUM FLOW TIMER: SUPPORT SURFACE INFLATION TIME:	180 Kg 2 8 ± 4 LPM 35 ± 5 mmHg 45 minutes 20 - 45 minutes
PATIENT COMFORT CONTROL PRESSURES SOFT PRESSURE: FIRM PRESSURE: CYCLE TIME:	8 ± 4 mmHg 32 ± 4 mmHg 5, 10, 15, 20 minutes

2 OVERVIEW

PATIENT CONTACT:	Control unit and mattress have Latex-Free components
MINIMUM CAPACITY:	15 kgs
NUMBER OF AIR CUSHIONS:	16

Mechanical Parameters

	CONTROL UNIT
DIMENSIONS (L X W X H):	23 cm x 13 cm x 13 cm
WEIGHT:	2.2 Kg
POWER CORD:	3m - 4.2 m Long, Hospital Grade
CONNECTION:	Two 6 mm Flow Couplings
PACKAGING:	1 piece/box
AIR FILTER:	None

	OVERLAY
DIMENSIONS (L X W X H):	200 cm 88 cm x 13 cm 200 cm 83 cm x 13 cm
WEIGHT:	8 lbs (3.5 kg)
AIR CELL DIMENSION (L X W X H):	83 cm x 12.5 cm x 12.5 cm 88 cm x 12.5 cm x 12.5 cm

	DO NOT pierce or cut		DO NOT bleach
	DO NOT dry clean.		User weight limit of 180 Kg (397 lbs)
	DO NOT put near flame.		Tumble dry low heat setting

	Wash at or below 95° C (203° F)		DO NOT iron
	Consult Instructions for Use		Line Dry

2.3 Intended Use

This alternating pressure overlay and pump are intended to be used in conjunction with an appropriately sized bed frame and mattress, as part of an overall pressure ulcer program of care.

This product is designed to be used as an overlay only with a mattress used as a base. Use without a foam mattress under the alternating pressure overlay is prohibited.

Ensure that articles between the support surface and patient are minimized.

Bed sheets must be loosely fitted, with creases smoothed out. Care must always be taken to ensure that the support surface in contact with the user is kept free from crumbs and other food debris, and that drip cables, stents, and other foreign objects do not become entrapped between the user and the pressure reducing surface of the mattress, as this may result in the development of pressure ulcers.

The introduction of certain third party products between the mattress surface and the user may reduce or impede the clinical effectiveness of this product.

'Third party products' may include, but are not limited to items including under blankets, plastic sheets and sheepskins, etc. Heated over blankets must only be used in consultation with a suitably qualified healthcare professional, as an increase in temperature can increase the risk of developing pressure ulcers.

3 GENERAL GUIDELINES

3 General Guidelines

3.1 General Guidelines

DANGER: RISK OF INJURY OR DEATH

Use **ONLY Invacare** rails, mattresses, bed extenders and other accessories with Invacare bed products.

Bed accessories designed by other manufacturers have **NOT** been tested by Invacare.

Use of **NON-Invacare** bed accessories may result in injury or death.

WARNING

DO NOT use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as owner's manuals, service manuals or instruction sheets supplied with this product or optional equipment. If you are unable to understand the warnings, cautions or instructions, contact a healthcare professional, dealer or technical personnel before attempting to use this equipment - otherwise, injury or damage may occur.

DO NOT modify this equipment without authorization of the manufacturer.

ACCESSORIES WARNINGS

Invacare products are specifically designed and manufactured for use in conjunction with Invacare accessories. Accessories designed by other manufacturers have not been tested by Invacare and are not recommended for use with Invacare products.

NOTICE

THE INFORMATION CONTAINED IN THIS DOCUMENT IS SUBJECT TO CHANGE WITHOUT NOTICE.

Check all parts for shipping damage. Test before using. In case of damage, DO NOT use. Contact Invacare/Carrier for further instruction.

Contraindications

WARNING

ALWAYS consult the patient's physician before using the MA50CE system.

Installation

WARNING

The MA50CE series is recommended to be installed on medical bed frames with bed side or assist rails. To properly install this overlay system onto a mattress, refer to [Installing the Mattress Overlay System \(MA50CE\)](#) on page 15. It is preferred that the rails be in the raised position whenever a patient is on the bed. Health care professionals assigned to each case should make the final determination whether side or assist rails are warranted after assessing patient risks of entrapment and falls in accordance with applicable patient restraint legislation or facility interpretation of such legislation.

Controls on the footboard may be obstructed by the power unit on a few bed frames. It may be necessary to relocate the power unit. Refer to [Installing the Power Unit](#) on page 16.

Check that air hoses and power cord are clear of moving bed components before placing a patient on the bed. Operate all bed frame motorized functions through their full range of motion to be certain that there is no pulling, interference or pinching.

3 GENERAL GUIDELINES

Fire Hazard

DANGER

Smoking

DO NOT SMOKE while using this device. This system uses room air for circulation through the overlay. A cigarette can burn a hole in the bed surface and cause damage to the overlay. Also, patient clothing, bed sheets, etc. may be combustible and cause a fire. Failure to observe this warning can result in severe fire, property damage and cause physical injury or death.

Smoking by visitors in the room will contaminate the system. Therefore, visitor smoking is NOT permitted.

Anesthesia Equipment

There is an explosion risk if used with flammable anesthetics.

Oxygen

There is a possible fire hazard when used with oxygen administering equipment other than nasal mask or half bed tent type. The oxygen tent should NOT extend below mattress support level.

Entrapment May Occur

WARNING

Patient entrapment with bed side rails may cause injury or death. Mattress **MUST** fit bed frame and side rails snugly to prevent patient entrapment. Follow the manufacturer's instructions. Monitor patient frequently. Read and understand the Owner's/Operator's Manual prior to using this equipment. Invacare product manuals are available at www.invacare.com or your dealer.

Proper patient assessment and monitoring, and proper maintenance and use of equipment is required to reduce the risk of entrapment. Variations in bed rail dimensions, and mattress thickness, size or density could increase the risk of entrapment. Refer to the owner's manuals for beds and rails for additional product and safety information.

Assist rails with dimensions different from the original equipment supplied or specified by the bed manufacturer may not be interchangeable and may result in entrapment or other injury.

After any adjustments, repair or service and before use, make sure all attaching hardware is tightened securely.

3 GENERAL GUIDELINES

Electrical

DANGER - Risk of Electrical Shock

DO NOT remove cover. Refer to qualified service personnel.

WARNING

Before performing any maintenance to the power unit, disconnect the power cord from the wall outlet.

DO NOT insert items into any openings of the control unit. Doing so may cause fire or electric shock by shorting the internal components.

The control unit **MUST** be kept away from all heat sources and radiators during operation.

Connect the equipment to properly grounded wall outlet using 3m - 4.2m hospital grade power cord provided with the product.

Grounding reliability depends upon a properly grounded receptacle.

DO NOT expose the electronic control unit to any liquid, otherwise, severe electrical hazard could result.

ALWAYS use fuses of the same ratings as specified in the [Typical Product Parameters](#) on page 6. Using fuses with higher ratings could result in injury and/or damage to the product

Handle and transport the control unit carefully. **DO NOT** drop otherwise damage to the unit may result.

DO NOT place any objects or items such as blankets on or cover the control unit.

3.2 EMC Information

There is no known risk of adverse effects on the MA50CE control unit/pump caused by other electromagnetic devices, present at the time of treatment. Similarly, there is no known risk of the MA50CE control unit/pump causing adverse effects on other electromagnetic devices.

4 Operation

CAUTION

The control unit and overlay on the MA50CE series are designed to be used as a system. **DO NOT** replace overlays or control units with other models or other brands. Otherwise, damage to the system may occur. Contact your supplier to get the correct replacement if needed.

4.1 Installing the Mattress Overlay System (MA50CE)

The overlay system comes with the following straps - head end elastic strap, foot end elastic strap and two buckle straps.

1. If necessary, lower the side rails to facilitate installation of the overlay.
2. Place the overlay directly over the existing foam mattress.
3. Loop the head end elastic strap around the head end of the existing foam mattress.
4. Loop the foot end elastic strap around the foot end of the existing foam mattress.
5. Secure the overlay to the existing foam mattress using the two buckle straps.

Ensure that the side rails can move freely. Corner Rings and cable management system help manage sheets and power cord to keep neat and safe.



FIGURE I Installing the Mattress Overlay System (MA50CE)

4 OPERATION

4.2 Installing the Power Unit

1. Swing out the bed hooks on the back of the control unit.
2. Hang the control unit on the footboard with the control panel facing the caregiver.

Control unit should be facing away from the patient

If the bed does not have a footboard, place the control unit on a flat surface, leaving room for the hose to hang down.

FIGURE 2 Installing the Power Unit

4.3 Connecting the Hose

CAUTION

Ensure that the hose connecting the control unit to the overlay is routed such that it cannot be stepped on, kinked, squeezed or otherwise damaged.

The power cord to the control unit should be positioned to avoid a trip hazard and/or damage to the cord. Invacare recommends placing the cord through the cable management system located along the side of the mattress and attaching it to an electrical outlet at the head of the bed.

i For this procedure, refer to FIGURE 3 on page 17.

1. Locate the hose at the foot end of the overlay.
2. Locate the control unit connectors on the right side of the control unit.
3. Squeeze and hold the tabs on the hose connectors.
4. Push the hose connectors onto the control connectors until an audible click is heard.

i The audible click indicates that the hose connectors are properly engaged with the control connectors.


FIGURE 3 Connecting the Hose

4 OPERATION

4.4 Connecting the Power Cord

WARNING

DO NOT alter plug to fit a non-conforming outlet. Instead, have an electrician install a properly grounded outlet. Failure to use the correct plug and outlet can result in a potential safety hazard.

CAUTION

Ensure that the power cord of the control unit is not pinched, or has any objects placed on it. Ensure the power cord has been secured in the cable management system located along the side of the overlay. Ensure it is not located where it can be stepped on or tripped over.

1. Examine the hospital grade power cord supplied with the control unit.
2. Perform one of the following:
 - If the plug shows signs of damage or the cord insulation is damaged/pinched - Call your supplier for a replacement hospital grade cord.
 - If the plug is not damaged - Plug the end of the supplied hospital grade power cord into the power outlet on the side of the control unit.
3. Plug the other end of the plug into a properly grounded outlet on the wall.

Make sure the control units power inlet connection is positioned to easily disconnect the power cord from the unit.

4.5 Using the Front Panel

DETAIL “A” -
FRONT PANEL


FIGURE 4 Using the Front Panel

Power Button

- I. To turn the control unit on or off, press and release the Standby/Power button ().

Once the unit is plugged in, an Amber LED on the control unit is lit indicating that the system is in Standby mode. Once the Standby/Power button is pressed and released, a green led illuminates indicating that the control unit is on.

If the power comes on after a power outage, the control unit will go through its system initialization routine for a few seconds and then resume the desired function.

4 OPERATION

Static/Alternating Pressure Button

1. Press the Alternating Pressure button () to select between Static Mode and one of four Alternating Pressure Times.

The alternating pressure times indicate the frequency of deflation or inflation of half of the air cushions (even or odd numbered). For example, 10 minutes is selected. In this example, the even numbered air cushions in the mattress will deflate, while the odd numbered cushions remain at constant pressure. When 10 minutes has elapsed, the even numbered air cushions will inflate while the odd numbered cushions remain at constant pressure. When 10 minutes has elapsed again, the odd numbered air cushions deflate while the even numbered cushions remain at constant pressure. The air cushions continue to inflate or deflate like this as long as an alternating pressure time has been selected. At all times, half of the air cushions remain at a constant pressure when an alternating pressure time has been selected.

In Static Mode, all air cushions in the overlay are maintained at constant pressure. The Static LED is illuminated Green.

Adjusting the Minimum Pressure Setting (A/P Low)

The minimum pressure setting during Alternating Pressure is set to 0 (2 mmHg) at the factory. This means that the overlay will deflate as much as possible during the deflation time when using 10, 15 or 20 minute Alternating Pressure times. The setting can be changed so the overlay will only deflate to 25%, 50%, or 75% of the selected pressure setting.

1. Ensure the Standby LED is lit.
2. Press and hold the Max Inflate and Fowler buttons.
3. Read the LCD Display. The LCD Display will flash:
 - 0 if the current setting is 0 (2 mmHg).
 - 2 if the current setting is 25% of the selected pressure.
 - 5 if the current setting is 50% of the selected pressure.
 - 7 if the current setting is 75% of the selected pressure.
4. Press the Static/Alternating Pressure button to toggle between the four settings.
5. Press the Power button to save the setting.

Firm/Soft Buttons

- I. Select comfort pressure settings by pressing Firm () or Soft () button.
 - Soft Button - Pressing this button reduces the pressure setting in the mattress.
 - Firm Button - Pressing this button increase the pressure setting in the mattress.

The patient comfort pressure ranges from Soft (level 0 = 6 ± 4 mmHg) to Firm (level 9 = 32 ± 6 mmHg). The Comfort Control LED displays the patient comfort pressure levels from 0 to 9 and provides a guide to the caregiver to set approximate comfort pressure level depending on the patient weight. If the patient's weight to height ratio is above average, increase the pressure setting by approximately 20%.

Fowler Button

Use when articulating the bed to a sitting position by raising the head.

- I. Press the Fowler button () to activate the patient fowler mode.

The Fowler LED illuminates when in this mode. When this mode is activated, the control unit increases the pressure in the overlay to prevent the patient from bottoming out.

Max Inflate Button

- I. Press the Max Inflate button () to select the Max Inflate mode and inflate the overlay rapidly to maximum firmness (pressurized to 35 ± 6 mmHg). A series of beeps sound every three minutes as a reminder that the Max Inflate mode is active. After 20 minutes, the Max Inflate mode deactivates and the control unit defaults to the previous setting.

It takes 20-45 minutes for the overlay to inflate fully (inflation time depends on size of overlay).

Max Inflate mode can be manually disengaged by pressing the Max Inflate button. This will deactivate the Max Inflate LED.

It is recommended that Max Inflate setting be used during patient ingress/egress, patient wound care, patient turning or patient cleaning.

4 OPERATION

2. Press the Max Inflate button again to manually disengage Max Inflate mode.

This will deactivate the Max Inflate LED.

To protect the patient from lying on a firm surface too long, the pressure will automatically return to normal range after 20 minutes.

Lock Button

1. Press the Lock button () to lock out all control unit functions, including the On button, to prevent any tampering with the settings. Press the button for approximately 3-5 seconds for the Lock LED to activate.

The Lock LED illuminates when in this mode.

Alarm Silence Button

1. Press the Alarm Silence button () to silence the alarm that sounds in the event of power failure or when the hose is disconnected from the control unit.

The Alarm LED illuminates when in this mode.

Power Fail LED

In the event of power outage, an alarm sounds and Power Fail LED flashes AMBER. The control unit has internal memory and retains the previous settings during the power outage.

During a power outage, the overlay retains the air as long as the overlay is connected to the control unit.

Low Pressure LED

In the event that the overlay hose disconnects, an alarm sounds and Low Pressure LED flashes AMBER. Once the low pressure problem is fixed, the control unit resumes operation in the previously set mode.

4.6 Powering Up the System

For this procedure, refer to FIGURE 4 on page 19.

1. Turn on the power to the system by pressing the Power button on the control unit.

Once the button is released, a green LED illuminates when the unit is on.

4.7 Placing the Patient on the Overlay

1. Press the Max Inflate button to turn on the control unit to maximum flow.

In this mode, the Max Inflate LED lights up.

2. After the overlay is fully inflated, place the patient on the overlay.
3. Ensure that the patient's feet are towards the foot end of the overlay (the end with the hose).
4. Center the patient on the bed from side-to-side and head-to-foot.

Special positioning may be required with contracted patients to provide comfortable positions.

5. After placing the patient, make certain no objects will fall under the patient, such as feeding tubes, IVs etc.
6. Wait five minutes for the overlay pressure to stabilize.
7. Once the overlay inflates to its normal size, set the comfort pressure to the desired comfort level.
8. Wait five minutes for the overlay pressure to stabilize.
9. If the patient feels that the bed is too soft/hard, press the Soft or Firm button to adjust the comfort settings.
10. Use a regular pillow to help support and stabilize the patient's head.

4 OPERATION

4.8 Transferring Patient From/To a Trolley

WARNING

ALWAYS engage the wheel locks of the bed and the wheel locks of the trolley before transferring the patient between the bed and the trolley.

1. Engage the wheel locks of the bed. Refer to the owner's manual provided with the bed.
2. Engage the wheel locks of the trolley. Refer to the owner's manual provided with the trolley.
3. Press the Max Inflate or Mode button to achieve maximum overlay pressure.
4. Raise or lower the bed to match the trolley height. Refer to the owner's manual provided with the bed.
5. When the overlay has reached maximum firmness, perform one of the following:
 - Bed to Trolley Transfer - Slide the patient onto the trolley.
 - Trolley to Bed Transfer - Slide the patient onto the bed.

4.9 Transferring Patient From/To a Wheelchair

WARNING

ALWAYS engage the wheel locks of the bed and the wheel locks of the wheelchair before transferring the patient between the bed and the wheelchair.

1. Engage the wheel locks of the bed. Refer to the owner's manual provided with the bed.
2. Engage the wheel locks of the wheelchair, if applicable. Refer to the owner's manual provided with the wheelchair.
3. Press the Max Inflate or Mode button to achieve maximum overlay pressure.
4. Raise or lower the bed to match the wheelchair height. Refer to the owner's manual provided with the bed.
5. When the overlay has reached maximum firmness, perform one of the following:
 - Bed to Wheelchair Transfer - Slide the patient onto the wheelchair.
 - Wheelchair to Bed Transfer - Slide the patient onto the bed.

4.10 Preparing for CPR Procedure

For this procedure, refer to FIGURE 5.

1. Press and hold the tabs on hose connector while pulling the hose from the control unit.
2. Twist the CPR connector at the foot end of the overlay and all air will be evacuated in the torso and chest region of the overlay within 18 seconds.


FIGURE 5 Preparing for CPR Procedure

4.11 About Power Outage and Transportation

If the hoses remain connected to the control unit, the overlay retains air during a power outage, during transportation, when the control unit is unplugged or when the control unit is turned Off. If the hose becomes disconnected or damaged, the overlay deflates.

5 Maintenance and Troubleshooting

5.1 Cleaning and Care

WARNING

Before cleaning or disassembling the MA50CE series, check the underside of the overlay folds for sharp objects such as scissors, needles, etc. These objects should be removed and discarded before proceeding with further cleaning or disassembly. Failure to follow these recommendations may result in personal injury or damage to the product.

Because of potential risk of infectious exposure, cleaning with the patient on the bed is not recommended.

All equipment should be inspected. Any item that is visibly soiled with the patient's blood or other body fluids should be properly cleaned or removed.

Laundry workers should treat all soiled bedding as if it were contaminated with pathogenic microorganisms.

Cleaning the Control Unit

- I. Perform the following steps to clean the control unit and hose fittings:
 - A. Wipe all controls, chassis and hose fittings with a quaternary disinfectant solution.
 - B. Using a nylon brush, gently clean all crevices as they can harbor microorganisms.
 - C. Air dry all treated surfaces.

Cleaning Overlay Components (air cushions, overlay base, etc.):

1. Using a clean sponge or paper towel, wipe down the overlay components with a dilute detergent solution of quaternary cleaner disinfectant or other germicidal detergent solution.
2. Wipe the overlay components with a clean dry cloth.

CAUTION

Laundry workers should always wash their hands before working with clean bedding.

DO NOT overload the machine.

DO NOT use chlorine bleach because it may damage the fabric coating.

High air temperatures will damage the fabrics and void the Invacare warranty.

Cleaning Covers

(Cleaning the covers removes contaminants such as dust and organic matter.)

1. Remove all covers for laundering.
2. Launder the covers in water (recommended temperature to be at or below 95° C.)

! IMPORTANT

Washing at higher temperatures may cause shrinkage.

Drying the Covers

1. Perform one of the following:
 - Hang overlay covers from a line or bar and drip dry in a clean indoor environment.
 - Tumble dry on a low setting.

! IMPORTANT

Tumble dry setting must not exceed 40 ° C.

DO NOT tumble for longer than 10 minutes.

Washing at higher temperatures may cause shrinkage.

Dry covers thoroughly before re-fitting to the overlay.

5 MAINTENANCE AND TROUBLESHOOTING

Disinfecting (Reducing the number of microorganisms)

The overlay may be disinfected by spraying or wiping with tested, approved disinfectants. The Robert Koch Institute has a list of currently approved disinfectants at <http://www.rki.de>.

Please contact your hygiene specialist in the event of contamination.

! **IMPORTANT**

Ensure that any residual detergent has been removed prior to disinfection.

Light Soilage

- Wipe down the cover using a 0.1% Chlorine Solution (1,000 ppm).
- Rinse the cover with clean water using a single use non-abrasive cloth.
- Dry the cover thoroughly.

Heavy Soilage

Where the mattress is badly soiled, Invacare recommends cleaning with a dilute cleaning solution at or below 95° C in the washing machine.

! **IMPORTANT**

1% Chlorine Solution used on a regular basis can diminish the life of the cover if not rinsed and dried properly.

DO NOT use granules.

- Clean up all spillages of bodily fluids (i.e. blood, urine, faeces sputum, wound exudate and all other bodily secretions) as soon as possible using a 1% Chlorine Solution (10,000 ppm).

Large spillages of blood should be absorbed and removed with paper towels first and then cleaned as above.

- Rinse the cover with clean water using a single use non-abrasive cloth and dry thoroughly.

! **IMPORTANT**

Polyurethane coated fabrics can absorb liquids for short periods causing a temporary change to the polyurethane characteristics. The mattress cover swells temporarily and is more vulnerable to physical damage for a period after it is completely surface dried, by which time it will revert to its previous state.

- Dry cover thoroughly.

Storing the System

1. Ensure that the control unit is off and disconnect the power cord from the wall outlet.
2. Clean the system. Refer to [Cleaning and Care](#) on page 26.
3. Disconnect the air hose connector(s) from the control unit and allow air to vent from the overlay.
4. Gently roll up the overlay with minimal handling and agitation.
5. Ensure the top cover surface is inside the roll beginning at the head section.
6. Store the unit, keeping the overlay with the control unit.

5 MAINTENANCE AND TROUBLESHOOTING

5.2 Troubleshooting

PROBLEM	CAUSE	SOLUTION
Overlay not inflating Not alternating properly	Overlay hose disconnected	Connect hose connectors and lock them in place
	Air hose kinked or split	Straighten hose or replace split hose
	Major leak in the air cushions or overlay pad	Replace leaking air cushions or overlay pad
	Kinked or split manifold	Unkink manifold or replace split manifold
	Has power and fuse is good, control unit does not come on	Check all connections. Unplug from power and wait 30 seconds before connection to power source. If problem persists, send control unit back to factory for repair
	Not alternating, solenoid malfunction	Send control unit back to factory for repair
	No air, pump malfunction	Send unit for repair
	CPR connector left open	Twist CPR connector to close
No power	Control unit off	Check power source and turn on unit
	Power cord disconnected	Connect power cord to the power source
	No power in the power source	Check power source has power and turn on unit.
	Power outage	Wait till the power source has power
	Blown fuse	Replace blown fuse with an equivalent fuse

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
INVACARE

CARELESSNESS CAUSES FIRE
GENERAL CLEANING
Wipe down using soap & water/mild detergent solution or remove cover for handwashing in HC deg C. Do not use bleach.
FOR DECONTAMINATION REFER TO YOUR INFECTION CONTROL SPECIALIST.

Invacare Operations UK Ltd,
Pencoed, Bridgend CF35 5AQ
Tel: +44 1656 776201
Fax: +44 1656 776201
www.invacare.com

Manufactured for:
Invacare Operations UK Ltd
Made in the USA.

MA50CE Overlay
Size: 200x85x13cm Date: _____
Serial No: _____
Q.C.: _____

Cover passes BS 7175 Crib 5
Complies Passes EN 507-1
Cigarette and EN597 - 2 Match Equivalent.

HMI No : 69681
Minimum User Weight : 15kgs

INVACARE

CARELESSNESS CAUSES FIRE
GENERAL CLEANING
Wipe down using soap & water/mild detergent solution or remove cover for handwashing in HC deg C. Do not use bleach.
FOR DECONTAMINATION REFER TO YOUR INFECTION CONTROL SPECIALIST.

Invacare Operations UK Ltd,
Pencoed, Bridgend CF35 5AQ
Tel: +44 1656 776201
Fax: +44 1656 776201
www.invacare.com

Manufactured for:
Invacare Operations UK Ltd
Made in the USA.

MA50CE Overlay
Size: 200x85x13cm Date: _____
Serial No: _____
Q.C.: _____

Cover passes BS 7175 Crib 5
Complies Passes EN 507-1
Cigarette and EN597 - 2 Match Equivalent.

HMI No : 69680
Minimum User Weight : 15kgs

INVACARE CORPORATION
Alternating Pressure **microAir 50CE**

DANGER-EXPLOSION HAZARD: DO NOT use in the presence of flammable anesthetics.
RISQUE D'EXPLOSION: NE PAS employer en presence d'anesthésiques inflammables.
WARNING: Equipment should be connected to a properly grounded receptacle. Risk of electrical shock. DO NOT remove back. Disconnect air hose before administering CPR.

5 055050 401418 >

MODEL / REF: MA50CE

19/04/2012 **S/N: 12D KP500001**

MANUFACTURED FOR: Invacare Operations UK Ltd, Pencoed, Bridgend, CF35 5AQ UK. Tel: +44 1656 776200

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 CARELESSNESS CAUSES FIRE GENERAL CLEANING Wipe down using soap & water/mild detergent solution or remove cover for laundering at 60 deg C. Do not use bleach. FOR DECONTAMINATION REFER TO YOUR INFECTION CONTROL SPECIALIST. Invacare Operations UK Ltd Pencoed, Bridgend CF35 5AQ Tel: +44 1656 776200 Fax: +44 1656 776201 www.invacare.com Manufactured for: Invacare Operations UK Ltd Made in the USA	 CARELESSNESS CAUSES FIRE GENERAL CLEANING Wipe down using soap & water/mild detergent solution or remove cover for laundering at 60 deg C. Do not use bleach. FOR DECONTAMINATION REFER TO YOUR INFECTION CONTROL SPECIALIST. Invacare Operations UK Ltd Pencoed, Bridgend CF35 5AQ Tel: +44 1656 776200 Fax: +44 1656 776201 www.invacare.com Manufactured for: Invacare Operations UK Ltd Made in the USA
MA50CE Overlay Size: 200x68x13cm Date: _____ Serial No: _____ Q.C.: _____	MA50CE Overlay Size: 200x68x13cm Date: _____ Serial No: _____ Q.C.: _____
 Cover passes BS 7176 Crib 1 Composite Passes EN 207 - 1 Cigarette and EN597 - 2 Match Equivalent. HMI No : 69681 Minimum User Weight : 15kgs	 Cover passes BS 7176 Crib 1 Composite Passes EN 207 - 1 Cigarette and EN597 - 2 Match Equivalent. HMI No : 69680 Minimum User Weight : 15kgs

 INVACARE CORPORATION Alternating Pressure microAir 50CE		 5 055050 401418 >
⚠ DANGER-EXPLOSION HAZARD: DO NOT use in the presence of flammable anesthetics. ⚠ RISQUE D'EXPLOSION: NE PAS employer en presence d'anesthésiques inflammables. ⚠ WARNING: Equipment should be connected to a properly grounded receptacle. Risk of electrical shock. DO NOT remove back. Disconnect air hose before administering CPR.		
P/N: 400375 Rev: 03  19/04/2012	S/N: 12D KP500001	MANUFACTURED FOR: Invacare Operations UK Ltd, Pencoed, Bridgend, CF35 5AQ UK. Tel: +44 1656 776200 Made In U.S.A.

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
INVACARE

CARELESSNESS CAUSES FIRE
GENERAL CLEANING
Wipe down using soap & water/mild detergent solution or remove cover for transferring at 95 deg C. Do not use bleach.
FOR DECONTAMINATION REFER TO YOUR INFECTION CONTROL SPECIALIST.

Invacare Operations UK Ltd
Pencoed, Bridgend CF35 5AQ
Tel: +44 1656 776200
Fax: +44 1656 776201
www.invacare.com

Manufactured for:
Invacare Operations UK Ltd
Made in the USA

MA50CE Overlay
Size: 200x83x13cm Date: _____
Serial No: _____
Q.C.: _____

Cover passes BS 7175 Crb 1
Cigarette and EN597 - 2 Match Equivalent.

HMI No : 69681
Minimum User Weight : 15kgs

INVACARE

CARELESSNESS CAUSES FIRE
GENERAL CLEANING
Wipe down using soap & water/mild detergent solution or remove cover for transferring at 95 deg C. Do not use bleach.
FOR DECONTAMINATION REFER TO YOUR INFECTION CONTROL SPECIALIST.

Invacare Operations UK Ltd
Pencoed, Bridgend CF35 5AQ
Tel: +44 1656 776200
Fax: +44 1656 776201
www.invacare.com

Manufactured for:
Invacare Operations UK Ltd
Made in the USA

MA50CE Overlay
Size: 200x83x13cm Date: _____
Serial No: _____
Q.C.: _____

Cover passes BS 7175 Crb 1
Cigarette and EN597 - 2 Match Equivalent.

HMI No : 69680
Minimum User Weight : 15kgs

INVACARE CORPORATION
Alternating Pressure **microAir 50CE**

⚠ DANGER-EXPLOSION HAZARD: DO NOT use in the presence of flammable anesthetics. RISQUE D'EXPLOSION: NE PAS employer en presence d'anesthésiques inflammables.

⚠ WARNING: Equipment should be connected to a properly grounded receptacle. Risk of electrical shock. DO NOT remove back. Disconnect air hose before administering CPR.

5 055050 401418 >

MODEL / REF: MA50CE

S/N: 19/04/2012 12D KP500001


MANUFACTURED FOR:
Invacare Operations UK Ltd, Pencoed, Bridgend, CF35 5AQ UK. Tel: +44 1656 776200

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
 CARELESSNESS CAUSES FIRE GENERAL CLEANING Wipe down using soap & water/mild detergent solution or remove cover for laundering at 60 deg C. Do not use bleach. FOR DECONTAMINATION REFER TO YOUR INFECTION CONTROL SPECIALIST. Invacare Operations UK Ltd, Pencoed, Bridgend CF35 5AQ Tel: +44 1656 776200 Fax: +44 1656 776201 www.invacare.com Manufactured for: Invacare Operations UK Ltd Made in the USA	 CARELESSNESS CAUSES FIRE GENERAL CLEANING Wipe down using soap & water/mild detergent solution or remove cover for laundering at 60 deg C. Do not use bleach. FOR DECONTAMINATION REFER TO YOUR INFECTION CONTROL SPECIALIST. Invacare Operations UK Ltd, Pencoed, Bridgend CF35 5AQ Tel: +44 1656 776200 Fax: +44 1656 776201 www.invacare.com Manufactured for: Invacare Operations UK Ltd Made in the USA
MA50CE Overlay Size: 200x85x13cm Date: _____ Serial No: _____ Q.C.: _____	MA50CE Overlay Size: 200x85x13cm Date: _____ Serial No: _____ Q.C.: _____
 Cover passes BS 7175 Crib 5 Complies Passes EN 527-1 Cigarette and EN597 - 2 Match Equivalent. HMI No : 69681 Minimum User Weight : 15kgs	 Cover passes BS 7175 Crib 5 Complies Passes EN 527-1 Cigarette and EN597 - 2 Match Equivalent. HMI No : 69680 Minimum User Weight : 15kgs

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2.1 Plassering av etikett


ETL US
Intertek
2003311

ETL LISTED
CONFORMS TO
IEC STD 60601-1
EN STD 60601-1-2
ANSI/AAMI STD ES60601-1
CERTIFIED TO
CAN/CSA STD C22.2 NO. 60601-1

Sikringsetikett

Se Vanlige produktparametere på side 219 for ikonforklaringer

<p>INVACARE</p> <p>CARELESSNESS CAUSES FIRE GENERAL CLEANING Wipe down using soap & water/mild detergent solution or remove cover for laundering at 90 deg C. Do not use bleach. FOR DECONTAMINATION REFER TO YOUR INFECTION CONTROL SPECIALIST.</p> <p>Invacare Operations UK Ltd Pencoed, Bridgend CF35 3AQ Tel: +44 1656 776200 Fax: +44 1656 776201 www.invacare.com</p> <p>Manufactured for: Invacare Operations UK Ltd Made in the USA</p> <p>MA50CE Overlay Size: 200x83x13cm Date: _____ Serial No.: _____ Q.C.: _____</p> <p>Icons: No fire, No open flames, No smoking, No sparks, No heat, No electrical equipment, No liquids, No sharp objects, No sharp objects, No sharp objects, No sharp objects, No sharp objects, CE</p> <p>HMI No : 69681 Minimum User Weight : 15kgs</p>	<p>INVACARE</p> <p>CARELESSNESS CAUSES FIRE GENERAL CLEANING Wipe down using soap & water/mild detergent solution or remove cover for laundering at 90 deg C. Do not use bleach. FOR DECONTAMINATION REFER TO YOUR INFECTION CONTROL SPECIALIST.</p> <p>Invacare Operations UK Ltd Pencoed, Bridgend CF35 3AQ Tel: +44 1656 776200 Fax: +44 1656 776201 www.invacare.com</p> <p>Manufactured for: Invacare Operations UK Ltd Made in the USA</p> <p>ELLER</p> <p>MA50CE Overlay Size: 200x83x13cm Date: _____ Serial No.: _____ Q.C.: _____</p> <p>Icons: No fire, No open flames, No smoking, No sparks, No heat, No electrical equipment, No liquids, No sharp objects, No sharp objects, No sharp objects, No sharp objects, No sharp objects, CE</p> <p>Cover passes BS 7175 Crib 5 Composite Purses EN 957 - 1 Cigarette and EN957 - 2 Match Equivalents.</p> <p>HMI No : 69680 Minimum User Weight : 15kgs</p>
---	---

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Invacare AG, Benkenstrasse 260, CH-4108 Witterswil
Invacare Limited, Pencoed Technology Park, Pencoed,