Invacare® Lynx SERVICE MANUAL







This manual contains information on:

Troubleshooting

Maintenance

Repair

Edition: 22.05.07

	Mobitec Mobilitätshilfen GmbH	≅	+43 - 6232 - 55 35 0	
\mathbf{A}	Herzog Odilostrasse 101	– Fax:	+43 - 6232 - 55 35 4	
	A-5310 Mondsee		office@mobitec-austria.com	
	Austria	@:	austria@invacare.com	
	Austria	www:	www.mobitec-austria.com	
		******	www.mobiled additia.com	
(B)	Invacare® n.v.	~	+32 - (0)50 - 83 10 10	
<u> </u>	Autobaan 22	Fax:	+32 - (0)50 - 83 10 11	
	B-8210 Loppem (Brugge)	@:	belgium@invacare.com	
	Belgium	www:	www.invacare.be	
(CH)	Mobitec Rehab AG	~	+41 - (0)61 - 48 77 08 0	
<u>UII</u>	Benkenstraße 260	Fax:	+41 - (0)61 - 48 77 08 1	
	CH-4108 Witterswil	@:	office@mobitec-rehab.ch	
	Switzerland	@:	switzerland@invacare.com	
		www:	www.mobitec-rehab.ch	
	Invacare Aquatec	~	0 75 62 / 7 00 - 251	
(n)	Alemannenstraße 10	Fax	08 00 / 6 73 81 72	
	88316 Isny	@:	info@invacare-aquatec.de	
	Deutschland	www:	www.invacare-aquatec.de	
(BIZ)	Invacare® A/S	☎ (Kundeservice):	+45 - (0)36 - 90 00 00	
UK	Sdr. Ringvej 39	Fax (Kundeservice):	+45 - (0)36 - 90 00 01	
	DK-2605 Brøndby	@:	denmark@invacare.com	
	Danmark	www:	www.invacare.dk	
	Invacare® SA	☎ :	+34 - (0)972 - 49 32 00	
	c/ Areny, s/n	Fax:	+34 - (0)972 - 49 32 20	
	Poligon Industrial de Celrà	@:	contactsp@invacare.com	
	17460 Celrà (Girona)	WWW:	www.invacare.es	
	ESPAÑA			
	Invacare® Poirier SAS	* :	+33 - (0)247 - 62 64 66	
	Route de St Roch	Fax:	+33 - (0)247 - 42 12 24	
	F-37230 Fondettes	@:	contactfr@invacare.com	
	France	www:	www.invacare.fr	
(CP)	Invacare® Ltd		+44 - (0)1656 - 664 321	
(GB)	South Road		e): +44 - (0)1656 - 667 532	
	Bridgend Industrial Estate	@:	uk@invacare.com	
	Mid Glamorgan - CF31-3PY	@:	eire@invacare.com	
	United Kingdom	WWW:	www.invacare.co.uk	
	Onited Kingdoni	V V V V V V .	WWW.IIIWacaic.co.uk	

	Invacare Mecc San s.r.l.		+39 - 0445 - 38 00 59	
\cup	Via Dei Pini, 62	Fax:	+39 - 0445 - 38 00 34	
	I - 36016 Thiene (VI)	@:	italia@invacare.com	
	ITALIA	WWW:	www.invacare.it	
(IE)	Invacare Ireland	*	+353 - 18 10 70 84	
	Unit 5 Seatown Business Campus	Fax:	+353 - 18 10 70 85	
	Seatown Rd, Swords	@:	eire@invacare.com	
	Dublin			
	Ireland			
	Invacare® AS	☎ (Kundeservice):	+47 - (0)22 57 95 10	
(N)	Grensesvingen 9	Fax (Kundeservice):	+47 - (0)22 57 95 01	
	Postboks 6230	@:	norway@invacare.com	
	N-0603 Oslo	WWW:	www.invacare.no	
	Norge			
(NL)	Invacare® B.V.	2 :	+31 - (0)318 - 69 57 57	
ME	Celsiusstraat 46	Fax:	+31 - (0)318 - 69 57 58	
	NL-6716 BZ Ede	@:	csede@invacare.com	
	The Netherlands	www:	www.invacare.nl	
	Invacare® PORTUGAL Lda	2 :	+351-225105946	
	Rua Senhora de Campanhã 105	Fax:	+351-225105739	
	P-4369-001 Porto	@:	portugal@invacare.com	
	PORTUGAL	WWW:	www.invacare.pt	



Återförsäljare:

Invacare® AB



Fagerstagatan 9

S-163 91 Spånga

Sverige

MÖLNDAL

WWW:

@:

@:

Invacare® Deutschland GmbH

Kleiststraße 49

Tillverkare:

D-32457 Porta Westfalica

Deutschland

+46 - (0) 31 - 86 36 00 Fax: +46 - (0) 31 - 86 36 06

+46 - (0) 8 761 70 90

+46 - (0) 8 761 81 08

sweden@invacare.com

finland@invacare.com

www.invacare.se

@: ginvacare@invacare.com

LANDSKRONA

☎ (Kundtjänst):

Fax (Kundtjänst):

2 +46 - (0) 418 - 285 40 Fax: +46 - (0) 418 - 180 89

@: linvacare@invacare.com

OSKARSHAMN

+46 - (0) 491 - 101 40

Fax: +46 - (0) 491 - 101 80

@: oinvacare@invacare.com

4

Table of Contents

Cha	pter		Page
1	INTRODUCTION	6	
1.1	General information	6	
1.2	Notes on transport	6	
1.3	Important symbols in this manual	6	
2	SAFETY AND ASSEMBLY INSTRUCTIONS	8	
2.1	Before any inspection or repair work	8	
2.2	General safety information and notes on assembly / disassembly	8	
3	TIGHTENING TORQUES	10	
4	SERVICE PLAN (1X ANNUALLY)	11	
5	OPERATIONAL FAULTS	12	
5.1 5.1	Operational faults .1 Drive fault diagnosis	12 12	
5.2	Error Codes and Diagnostic Codes	13	
6	REPAIR WORK	15	
6.1	General warning information on installation work	15	
6.2	Replacing the drive motor	15	
6.3	Replacing the electronics	18	
6.4	Updating software	20	
6.5	Replacing the potentiometer	21	
6.6	Replacing operating console	23	
6.7	Replacing front axle	25	

1 Introduction

1.1 General information

- All maintenance and overhaul work must be carried out in accordance with these repair instructions.
- Please observe all safety instructions.
- Information about operation or about general maintenance and care work should be taken from the electric vehicle Operating Manual.
- You can find information about ordering spare parts in the spare parts catalogue.
- Use only genuine Invacare® spare parts. Using parts from any other source will void the warranty!
- We reserve the right to make any alterations on the grounds of technical improvements.
- The electric vehicle may only be maintained and overhauled by qualified personnel.
- The minimum requirement for service technicians is relevant training, such as in the cycle or orthopaedic mechanics fields, or suitably long-term job experience.
 - Experience and knowledge of electrical measuring devices (Multimeter) is also a requirement.
 - Special Invacare training sessions are recommended.
- Alterations to the electric vehicle which occur as a result of incorrectly or improperly executed maintenance or overhaul work lead to the exclusion of all liability on the part of INVACARE.
- If you have any problems or questions please contact INVACARE SERVICE.

1.2 Notes on transport

- If the electric vehicle has to be shipped back to the manufacturer for major repairs, you should always use the original packaging for transport.
- You should also include as accurate a fault description as possible.

1.3 Important symbols in this manual



WARNING!

This symbol warns you of danger!

Always follow the instructions to avoid injury to the user or damage to the product!



EXPLOSION HAZARD!

This symbol warns you of an explosion hazard, which can be caused by excessive tyre pressure in a pneumatic tyre!

Always follow the instructions to avoid injury to the user or damage to the product!



BURN HAZARD!

This symbol warns you of burns due, for example, to leaking battery acid!

• Always follow the instructions to avoid injury to the user or damage to the product!



NOTE:

This symbol identifies general information which is intended to simplify working with your product and which refers to special functions.



Requirements:

 This symbol identifies a list of various tools, components and items which you will need in order to carry out certain work.

2 Safety and assembly instructions

These safety instructions are intended to prevent accidents during work and it is imperative that they are observed.

2.1 Before any inspection or repair work

- Read and observe this repair manual and the associated operating manual!
- Observe the minimum requirements for carrying out the work (see chapter entitled "General information")!

2.2 General safety information and notes on assembly / disassembly



Danger of injury by crushing!

- Please note the heavy weight of some components. This applies especially to removal of drive units and batteries!
- Prop up the lifted electric vehicle with appropriate supports before starting the disassembly or assembly!



Danger of fire and burns due to electrical short-circuit!

- The electric vehicle must be switched off before removal of voltage-carrying components! To do this, disconnect the batteries!
- When making measurements on voltage-carrying components, avoid short-circuiting the contacts. Danger of fire and combustion!



Danger of injury and damage to the vehicle can result from incorrect or incomplete maintenance!

- Only ever use tools which are undamaged in good condition!
- Some moving parts have Teflon bushings! These parts must never be lubricated with grease!
- Never use standard nuts instead of self-locking nuts!
- Always use correctly dimensioned washers or spacers!
- Cable binders which have been cut off during disassembly should be replaced with new ones during reassembly!
- After completing maintenance work and before operating the electric vehicle, make sure all fixations are correctly secured! Check all parts for correct interlocking1
- Only operate the electric vehicle with correct tyre pressure (see Technical Specifications)!
- Check electrical components for correct functioning, incorrect polarity of cables can result in damage to the electronics!
- As a last check, always carry out a test-drive!



Notes

Mark all current settings for the electric vehicle (seat, armrests, backrest etc.), and the cable connecting plugs associated, before any removals. This makes reassembly easier. All plugs are fitted with mechanical safety devices which prevent release of the connecting plugs during operation. To release the connecting plugs the safety devices must be pressed in. When reassembling, ensure that these safety devices are correctly engaged.



WARNING: Any alteration to the drive programme can influence vehicle handling and the tipping stability of the electric vehicle!

- Alterations to the drive programme may only be carried out by trained Invacare® dealers!
- Invacare® supplies all electric vehicles from the factory with a standard drive programme. Invacare® can only assume a warranty for the safe vehicle handling of the electric vehicle – in particular tipping stability - for this standard drive programme!

3 Tightening torques

The tightening torques stated in the following table are dependent on the thread diameters for the nuts and bolts for which no special values are determined. All values apply to dry and grease-free threads.

Thread	M4	M5	M6	M8	M10	M12	M14	M16
Tightening torque in Nm ±10%	3 Nm	6 Nm	10 Nm	25 Nm	49 Nm	80 Nm	120 Nm	180 Nm

Caution: All other nuts or plastic connectors not noted here must be tightened FINGERTIGHT!

4 Service plan (1x annually)

Component	Check	Remedy	Notes	✓
Seat:	Welded seams, fixings and upholstery	Tighten screws, replace parts if damaged		
Frames (chassis) / battery mounting	 Check fixings, welded seams and battery mounting Check battery fixing straps 	Tighten screws, replace parts if damaged		
Wheel suspension and wheels	Check drive wheels for tight fit	Tighten hub nuts, replace if necessary		
	Check steering wheels for tight fit, float and side play	Adjust / replace	See "Replacing front axle" on page 25	
	Pneumatic tyre	Repair or replace if damaged	See operating manual	
Drive units, coupling mechanism	 Check functions in drive and push modes Check clutch mechanism 	 Replace motor if necessary Tighten screws / nuts, adjust or replace if necessary 	See "Replacing the drive motor" on page 15	
Lighting	Check function Check cable	Replace bulbs or cable if necessary		
Batteries	Check batteries for damage	Replace batteries if necessary	See operating manual	
	Check battery voltage	charge batteries	See operating manual	
	Check contacts and terminals	Clean contacts and terminals		
Drive electronics	 Status display flashing Fixing Cable, connecting plug Drive lever function Power supply 	 Evaluate blinking code Tighten, replace Replace Replace drive lever Replace cable, connecting plug or console 		
Drive program	Check drive electronics program version Newer version available?	Update software	See "Updating software" on page 20	

5 Operational Faults

5.1 Operational faults

Please proceed as follows if you have any problems:

- First assess the possible cause of the problem using the following table.
- Check the drive console status display. Evaluate the blink error code.
- Carry out the necessary checks and repairs as recommended in the following table.

5.1.1 Drive fault diagnosis

PROBLEM	OTHER SYMPTOMS	POSSIBLE CAUSE	SOLUTION	Documentation
Scooter will not start	Status display on operating console does not illuminate	Batteries possibly defective	Replace batteries	See operating manual
		Batteries possibly over- discharged	Pre-charge batteries	See operating manual
		Power supply to operating console possibly interrupted	Check master fuse	See operating instructions for main fuse position
			Check cable between modules for loose connections or damage	See "Replacing the electronics" on page 18
		Operating console possibly defective	Replace operating console	See "Replacing operating console" on page 23
	Operating unit status display blinking	Various causes	Assess error code	See "Error Codes and Diagnostic Codes" on page 13

PROBLEM	OTHER SYMPTOMS	POSSIBLE CAUSE	SOLUTION	Documentation
Scooter judders in drive mode	None	Drive motor(s) possibly defective	Replace motor(s)	See "Replacing the drive motor" on page 15
		Drive lever potentiometer possibly defective	Replacing the potentiometer	See "Replacing the potentiometer" on page 21
Batteries not being charged	None	Batteries possibly defective	Replace batteries	See operating manual
	LEDs blinking on charging unit	Charging device possibly defective	Replace charging unit	See charging unit operating manual
Scooter runs too slowly	None	Operating console possibly defective	Replace operating console	See "Replacing operating console" on page 23
		Batteries possibly defective	Replace batteries	See operating manual

5.2 Error Codes and Diagnostic Codes

Blink code	Fault	Consequence for the Scooter	Comments
1	Battery must be charged	Continues to drive	The batteries are discharged. Charge the battery as soon as possible.
2	Battery voltage too low	Stops driving	The batteries are depleted. Charge batteries.
			If you switch the Scooter off for a few minutes, the battery can often recuperate to such a stage that a short journey is still possible. You should only do this in emergency, however, because this causes the batteries to become excessively discharged!
3	Battery voltage too high	Stops driving	The battery voltage is too high. If the battery charger is connected, disconnect it from the Scooter.
			The electronic system charges the batteries when running downhill and when braking. This fault is caused when the battery voltage becomes too high during this process. Switch the Scooter off and on again. Switch on the lights to use up some power.

Blink code	Fault	Consequence for the Scooter	Comments
4	Power time exceeded	Stops driving	The Scooter has used up too much power for a long period, probably because the motor has been overloaded or the Scooter has attempted to work against insurmountable resistance. Switch off the Scooter, wait a few minutes and then switch it on again.
			The electronic system has determined a motor short-circuit. Check the wiring harness for short-circuits and check the motor
5	Brake failure	Stops driving	Ensure that the disengaging lever is in the engaged position.
			 There is a defect in the braking coil or in the cabling. Check the wheel lock and cabling for open or short-circuited circuitry.
6	No neutral position when switching Scooter on.	Stops driving	Drive lever is not in neutral when the keyswitch was turned. Put the drive lever in neutral, turn the power off and then turn on again.
			 It may be necessary to recalibrate the drive lever.
7	Fault in speed potentiometer	Stops driving	The drive lever electronics could be faulty or incorrectly connected.
8	Motor voltage error	Stops driving	The motor or its cabling is defective
9	Miscellaneous internal fault	Stops driving	The electronic system may have an internal error. Check all cables and connections.
10	Push/freewheel mode error	Movement braked	The Scooter has exceeded the permissible maximum speed during pushing or freewheeling. Switch the electronics system off and on again.

6 Repair Work

6.1 General warning information on installation work



CAUTION: Risk of damage to the vehicle! Collisions can be caused if shim rings are removed from the drive wheels during installation work!

 Shim rings are frequently placed between the drive shaft and the wheel hub to compensate tolerances. Collisions can be caused if these shim rings are removed and not re-installed! Please install all shim rings in exactly the same positions they were in before dismantling.

6.2 Replacing the drive motor



Note:

When disassembling, note the position of small parts such as screws and washers. Put small parts down so that they can be reassembled in the right sequence.



Note:

When connecting the motor and magnetic brake cable, ensure that the plug is in the right position. The motor rotation direction will be reversed if the connection is wrong.




Requirements:


- Open spanner, 13 mm.
- Phillips screwdriver
- Rubber hammer

From the operating manual:


- Loosen rear wheel fixing bolts.
- Remove the seat
- · Remove the battery pack.
- · Remove the drive unit.
- Further work:
- Support the drive unit with a wooden wedge.




- Loosen the fixing bolt (1) from the handle and remove it.
- Remove handle.
- Loosen the fixing bolt (2) on the rear panelling and remove it.
- Remove the rear panelling.




 Unplug the drive motor plug M+/ M- (1) and the magnetic brake plug B+/ B- (2) from the socket strip.




- Fold the drive forwards onto the locking mechanics.
- Remove wheels (1).



- Loosen motor fixing bolts (1) on both sides and remove.
- Remove the motor from the vehicle frame.



- Check rubber bearings (1) on both sides of the motor suspension. Replace damaged rubber bearings.
- Reassembly of the motor takes place in reverse order
- Complete reassembly of vehicle.
- Check all functions (trial run).



6.3 Replacing the electronics



CAUTION! Burn and scald hazard if power cable is short-circuited!

 Interrupt the power supply by removing the batteries before you start work on the electronic system.



Note:

The electronic system is supplied with a standard drive program. If you have carried out customer-specific modifications to the drive program, you will have to make these changes again after installing the new electronics module.



CAUTION: Any modifications made to the drive program can affect the drive behaviour and tipping stability!

- Changes to the drive program may only be carried out by trained Invacare® specialist dealers!
- Invacare can only provide a guarantee for safe driving especially tipping stability for unmodified standard drive programs.




Requirements:


- Phillips screwdriver
- To adapt the drive program: programming software or hand programming device and ACS electronics system installation manual, available from Invacare[®].
- From the operating manual:
- Remove the seat
- Remove the battery pack.

• Further work:


- Unlock the expanding rivet (1) on the rubber mat fixing by pressing the central clamping prong.
- Remove the expanding rivet and the rubber mat.




- Remove the electronics cover screws (1).
- Fold up the electronics cover (2) and remove it.




- Remove all plugs (1) from electronic system.
- Remove screws (2) on both sides of the electronics system.
- Replace the electronics system.
- Fit the electronics cover panel.



• Press the expanding rivet clamping prong (1) back in (2).



- Locate the rubber mat.
- Insert expanding rivet (1) and secure by pressing in clamping prong.
- Modify the drive program.
- Complete reassembly of vehicle.
- Check all vehicle functions (trial run).



6.4 Updating software

Invacare is continuously carrying out further development and improvements to drive programs for electric vehicles. For this reason, you should always check whether the drive program version number is up-to-date when carrying out any repairs or regular maintenance. If a newer version is available, the drive program must be updated. The procedure for updating the drive program is described in the wizard software operating manual.



Note:

The electronic system is supplied with a standard drive program. If you have carried out customer-specific modifications to the drive program, you will have to make these changes again after installing the new electronics module.



CAUTION: Any modifications made to the drive program can affect the drive behaviour and tipping stability!

- Changes to the drive program may only be carried out by trained Invacare® specialist dealers!
- Invacare can only provide a guarantee for safe driving especially tipping stability for unmodified standard drive programs.



Requirements:

- Dynamic[®] Wizard software
- Operating manual for Wizard software
- Further requirements, such as a minimum system configuration for the PC used for programming, required programming cables etc. can be taken from the Wizard software operating manual.

6.5 Replacing the potentiometer



CAUTION! Burn and scald hazard if power cable is short-circuited!

• Interrupt the power supply by removing the batteries before you start work on the electronic system.



CAUTION! Risk of accidents!

• The new potentiometer must be set to the middle position before inserting. Potentiometers which are not in the middle position can result in dangerous driving situations!



Requirements:

- Phillips screwdriver
- Allen key open-ended spanner 2 mm
- Digital multimeter
- From the operating manual:
- Remove the seat
- Remove the battery pack.
- Further work:
- Remove screws (1) below operating console.




• Pull the operating console upwards out of the front fixing (1).




- Remove drive lever (1).
- Separate the potentiometer connecting plug (2).
- Remove potentiometer (3).
- Set the potentiometer to the middle position (see below).
- Reassembly takes place in the reverse order.
- Complete reassembly of vehicle.
- Check all operating console displays.

Setting the potentiometer to the middle position:

- Loosen the bolts (1) on the drive lever mounting.
- Twist the potentiometer shaft (2) until the same resistance (approx. 5kOhm) is measured between the central connection cable (3) and both outputs (4).
- Retighten the bolts on the drive lever mounting.
- Check settings and repeat the adjustment if necessary.





6.6 Replacing operating console



CAUTION! Burn and scald hazard if power cable is short-circuited!

• Interrupt the power supply by removing the batteries before you start work on the electronic system.



Requirements:

- Phillips screwdriver
- From the operating manual:
- Remove the seat
- Remove the battery pack.
- Further work:
- Remove screws (1) below operating console.



 Pull the operating console upwards out of the front fixing.



- Remove plugs (3) from operating console.
- Replace operating console.
- Reassembly takes place in the reverse order.
- Complete reassembly of vehicle.
- Check all operating console displays.



6.7 Replacing front axle



Note:

When disassembling, note the position of small parts such as screws and washers. Put small parts down so that they can be reassembled in the right sequence.



Requirements:

- Open spanner, 9 mm.
- Open spanner, 10 mm.
- Open spanner, 12 mm.
- Open spanner, 13 mm.

From the operating manual:


- Remove the seat
- · Remove batteries.
- · Remove the drive unit.
- Remove front wheels.

Further work:

 Fold the steering column inwards and place the chassis on its side.



- Loosen the steering link nuts (1) and remove them (note information below).
- Remove steering link from deflector.
- Loosen the axle suspension nut (2) and remove it.
 In doing so, secure the axle suspension bolts (3) against turning using an openended spanner.
- Pull the bolts (3) out and remove the axle casing.



 When loosening the nuts (1), secure the steering link head against turning using an open-ended spanner.



- Reassembly of the axle takes place in reverse order
- Complete reassembly of vehicle.