Bedienungsanleitung
Operating Instructions
Mode d'emploi
Instrucciones de uso
Istruzioni per l'uso
Gebruiksaanwijzing
Bruksanvisning

*AQUATEC MINOR**
*AQUATEC CLASSIC**

AQUATEC MAJOR

AQUATEC®

(D)

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Important References:

For proper operation of your **AQUATEC** bath lift, please ensure the following conditions have been met before operating:

- The equipment should only be used if the distance from the tub spout to the lift at its highest position is at least 5 inches (12 cm).
- The lift must not be connected to pressureless reservoirs or flow heaters. Also do not connect to equipment with over-pressure regulators. There is a danger of the equipment bursting. **Always have a specialist connect the lift.**
- Water pressure required is 3–6 bar (21 psi minimum), dependent upon the weight being lifted (see pressure load charts).
- **During and after bathing, always shut off tub spout. After bathing, disconnect supply hose from tub spout or lift.**
- Please use the bathlift with cold water only.
- NB: To avoid squeezing, please do not place hands between the seat and bottom plate.
- The bathlift must only be used in conjunction with the security combination. Please ensure the distance between the water level and security combination is at least 15 cm. Please have the security combination checked annually by a specialist.
- The bathlift is connected to a suitable water supply (e.g. separate tap, tap of washing machine). Therefore, you can use your shower as usual.
- It is possible to connect the **AQUATEC** bathlifts to the shower hose supply of one handle mixers and two tap fittings. Please note that the tap or handle of your fitting, which opens the inflow of the shower must not be operated under pressure! We recommend that you have the bathlift installed by a specialist.
- Before bathing, please check that all pieces are in good working order (side flaps, remote control, switch-over, rocker switch), as well as the tightness of all hoses and valves. Please check whilst the tap is open if all supplies are fixed.
- In case of defects or signs of wear and tear, please contact your dealer, who can provide the necessary spare parts and repairs to your bathlift.
- **Disabled persons should have help while using the lift. Extremely disabled persons should bathe under supervision.**

- Always check the water temperature before bathing. If you use bathing oils or salts, pay attention to the correct dosage. The lift is harder to clean if you use higher concentrations and that can affect the sliding of the lifting scissor mechanism.

- **Read the connection instructions carefully, if necessary consult your dealer or a plumber.**

What to do in case of malfunction

- **Lift does not react to hand control or the rocker switch:** Check if water tub spout is open and the change-over switch is on liftoperation.
- **Lift moves slowly or comes to a stop:** Clean filter in quick disconnect (do not take out filter and never connect equipment without a filter).

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
(GB) (USA)

Assembly of AQUATEC bath lift

Screw the 4 suction cups in each of the 4 corners of the base plate. (When the suction cup is screwed on the release tap should point outwards.) Suction cups on the base will prevent slipping. Place the lift in the tub so that all four holes for the back support point towards the back of the tub (I).

Place back support prongs in one of the rows of holes (II).

Insert to the first bend of the prongs and tilt back to rest firmly on the seat. The back support can be removed by simply tilting it slightly towards the front (III).

AQUATEC MINOR with reclining back support: After placing the back support on the lift plug the small clear hose from the back support into the small quick disconnect on the control unit under the seat (II A).

Taking lift out of tub

If lift has to be removed:

- lower lift completely
- shut off tub spout
- disconnect supply hose from tub spout
- drain excess water
- take back support off
- take lift out of tub and store or put in carrying bag (accessory)

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GB USA

Assembly of AQUATEC MAJOR

Put the base in the tub so that the opening in the base plate lies above the drain (I).

Raise the lift up. Put the supporting plate onto the base so the black fixing balls fit in the recesses on the base (II) + (III).

Taking lift out of tub

If you would like to take the lift out of the tub:

- lower lift completely
- shut off tub spout
- disconnect supply hose from tub spout
- drain excess water
- take lift out of tub and store away or put in carrying bag (accessory).

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
(GB USA)

Connection

To connect the lift so it functions correctly, 3 elements have to be connected with each other; the tub spout, the change-over switch (optional accessory) and the lift.

Proceed as follows:
Connection to tub spout (I). Disconnect the hand shower hose from tub spout or shower pipe. The connection hose from the change-over switch will thread onto this. Connect the other end with the middle thread of the change-over switch. Screw the quick-connect with the safety combination onto the thread labeled with the lifter symbol (diag. I). Supply hose, shower and hand control should be connected as shown in diagram (I).

USA only:
Connection to shower pipe (I A)
Remove shower head and screw on the quick-disconnect and plug in the larger end of the white hose. (If it leaks, add plumbers or Teflon tape to shower pipe threads first.) Proceed as above.

Connection to washing machine faucet with Y-piece and Waterstop (optional accessory) (II)

Connection to sink faucet with reducing piece (III)

Connection to separate faucet or water line (dedicated line) (IV)

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
(GB USA)

Mounting the change-over Switch:

The change-over switch is used for the alternating use of the lift and a hand shower and is mounted on the wall within reach of the user. There are 3 mounting options:

1. Mounting with double-sided adhesive tape

Thoroughly clean the surface where you would like to install the change-over switch with the cleaning pad enclosed. Remove the protective film from one side of the tape and position the tape along the upper edge on the rear of the switch panel. Press firmly. Remove the remaining piece of film and press the change-over switch in place on the wall (diag. I).

2. Mounting with screws (accessories)

Take change-over switch in both hands and press with both thumbs onto the edges of the rocker switch so that the housing separates from the rear panel (II). Holes in back plate are used with screws and dowels for mounting onto tile joints (II A).

3. Mounting with suction cups (accessories)

Screw 2 suction cups firmly onto the threads on the backplate of the change-over switch (III). Clean tiles thoroughly and press moistened suction cups in place.

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Lift Operation

By Hand Control (diag. I)

The **AQUATEC MINOR** or **AQUATEC MAJOR** is raised by pushing the operating lever to position ①. Now the user can get on the lift. Move the lever to position ②, the lift lowers. The lift stops automatically when the lever is released. After bathing, push the lever to position ① until the lift has reached its raised position. Drain water, from the tub clean lift and wash covermats.

① = lift raises
② = lift lowers
centre = STOP

By Hand Control (diag. II)

Along with the functions described above, you can operate an optional reclining back support (accessory).

③ = back support lowers
④ = back support raises
centre = STOP

Showering

The change-over switch provides the opportunity to shower while bathing. Push the side marked with the shower symbol.

Pressure load chart

The diagram shows the relationship between body weight and the necessary water pressure for operating the lift.

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(GB) (USA)

Connection

Connections are possible at tub spout (diag. I), sink faucet (diag. II), shower head pipe (diag. III.A) or dedicated line (diag. III and IV).

(USA)

III A

To connect the lift properly, 3 elements have to be connected with each other, the faucet, the vaillant junction and the lift. Please proceed as follows:

Disconnect the shower hose or shower head from the faucet or pipe. Screw vaillant junction onto threaded end. Connect the supply hose from the lift to one side and the shower hose to the other side of the vaillant junction (diag. V).

The supply hose from the lift is attached in the rear of the lift seat. A switch located on the vaillant junction allows for alternate operation of the shower and lift (diag. VI).

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GB USA

Lift Operation

Your bath lift achieves its upward motion from the water pressure at your water supply. It is important that the water line always be open or turned on to operate the lift. The button or lever located on your shower spout is very important because the lift is attached to the shower connection. As you turn on your shower, your lift is ready for operation (diag. II). The vaillant junction at your spout allows you to alternate between lift operation and shower operation (diag. I). By pushing the rocker switch located on the seat to position ①, the lift is raised, pushing the switch to position ② slowly lowers the lift.

① = lift raises
 ② = lift lowers
 centre = STOP

After bathing, push the switch to position ① until the lift has reached the beginning position. Drain bath water, clean lift and wash covermats.

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GB USA

AQUATEC TRANS XL

Swivel seat for easy transfer onto the lift.

AQUATEC DISK XL

Rotary transfer aid for easy transfer onto the lift.

AQUATEC STEP

Step used for entering bath and bath lifts.

Special cushion

with side flaps. For use with the special back support.

Special stationary back support

with individually adjustable side flaps.

Adjustable head support

Adjustable height and side flaps. For use with the stationary or reclining back support.

Carrying bag

Waterproof for storage or easy transport when traveling.

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GB USA

Care and cleaning of AQUATEC bath lift

AQUATEC bath lifts do not need servicing. To ensure long lasting operation, we recommend the following steps:

The covermats can be laundered on normal cycle in a washing machine (60° C).

The anti-bacterial properties are not affected.

The mats and covers should be line dried. Do not put in a dryer.

The lift should be showered off frequently and cleaned with household cleaning products.

When bathing salts are used, scrub the lift with bathroom cleaning products. Grease the guiding rail with petroleum jelly.

To ensure proper operation of your **AQUATEC** bath lift, it should be inspected every 2 years by your local specialist.

All **AQUATEC** bath lifts come with a 3 year guarantee effective from date of production. Covermats are excluded.

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