Aquatec® Sansibar

| | | |
|----|--|----|
| en | Wall Mounted Shower Seat User Manual | 3 |
| de | Duschsitz für die Wandmontage Gebrauchsanweisung | 5 |
| fr | Siège de douche à montage mural Manuel d'utilisation | 8 |
| it | Seggiolino da doccia a parete Manuale d'uso | 11 |
| es | Asiento de ducha plegable de pared Manual del usuario | 14 |
| pt | Assento do duche de montagem na parede Manual de utilização | 17 |
| nl | Douchezitje met wandmontage Gebruiksaanwijzing | 20 |
| sv | Väggmonterad duschsits Bruksanvisning | 23 |
| fi | Seinään kiinnitettävä suihkuistuin Käyttöohje | 25 |
| cs | Nástěnné sedadlo do sprchy Návod k obsluze | 28 |
| hr | Zidna sjedalica za tuš Upute za uporabu | 31 |


Fig. 1


Fig. 2


Fig. 3

AQUATEC®

This manual MUST be given to the user of the product.
BEFORE using this product, this manual MUST be read and saved for future reference.

Yes, you can.®


Fig. 4


Fig. 5


Fig. 6


Fig. 7

Fig. 8


Fig. 9


Fig. 10


Fig. 11


Fig. 12


Fig. 13


Fig. 14


Fig. 15

1 General

1.1 Introduction

This user manual contains important information about the handling of the product. To ensure safety when using the product, read the user manual carefully and follow the safety instructions.

Invacare reserves the right to alter product specifications without further notice.

Before reading this document, make sure you have the latest version. You find the latest version as a PDF on the Invacare website.

If you find that the font size in the printed document is difficult to read, you can download the PDF version from the website. The PDF can then be scaled on screen to a font size that is more comfortable for you.

In case of a serious incident with the product, you should inform the manufacturer and the competent authority in your country.

1.2 Warranty Information

We provide a manufacturer's warranty for the product in accordance with our General Terms and Conditions of Business in the respective countries.

Warranty claims can only be made through the provider from whom the product was obtained.

1.3 Compliance

Quality is fundamental to the company's operation, working within the disciplines of ISO 13485.

This product features the CE mark, in compliance with the Medical Device Regulation 2017/745 Class I.

This product features the UKCA mark, in compliance with Part II UK MDR 2002 (as amended) Class I.

We are continuously working towards ensuring that the company's impact on the environment, locally and globally, is reduced to a minimum.

We only use REACH compliant materials and components.

Product-specific Standards

The product has been tested and conforms to EN 12182 (Assistive products for persons with disability) and related standards. For further information about local standards and regulations, contact your local Invacare representative. See addresses at the end of this document.

1.4 Service Life

The expected service life of this product is three years when used daily and in accordance with the safety instructions, maintenance intervals and correct use, stated in this manual. The effective service life can vary according to frequency and intensity of use.

1.5 Limitation of Liability

Invacare accepts no liability for damage arising from:

- Non-compliance with the user manual
- Incorrect use
- Natural wear and tear
- Incorrect assembly or set-up by the purchaser or a third party
- Technical modifications
- Unauthorised modifications and/or use of unsuitable spare parts

2 Safety

2.1 Safety Information

WARNING!

Risk of Serious Injury or Damage

Improper use of this product may cause injury or damage.

- If you are unable to understand the warnings, cautions or instructions, contact a health care professional or provider before attempting to use this equipment.
- Do not use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as user manual, service manual or instruction sheet supplied with this product or optional equipment.

WARNING!

Risk of Falling!

- Do not use if defective.
- Do not make any unauthorized alterations or modifications to the product.
- Do not stand on this product.
- All legs must be set at the same height.
- Always use new wall plugs and bolts when repositioning the product.

WARNING!

- Before first use, clean the product thoroughly.
- Before every use, check the product for damage and ensure that it is securely assembled.

CAUTION!

Risk of Injury

Non-original or wrong parts may affect the function and safety of the product.

- Only use original parts for the product in use.
- Due to regional differences, refer to your local Invacare catalog or website for available options or contact your Invacare distributor. See addresses at the end of this document.

2.2 Labels and Symbols on the Product (Fig. 15)

| | | | |
|---|------------------------|---|--------------------------|
| Ⓐ | Manufacturer | Ⓒ | Production date |
| Ⓑ | UK Conformity Assessed | Ⓓ | Unique Device Identifier |
| Ⓒ | European Conformity | Ⓔ | Medical Device |
| Ⓓ | Lot number | Ⓕ | Reference number |
| Ⓔ | Read user manual | Ⓖ | Product designation |
| Ⓕ | Maximum user weight | | |

The identification label is attached to the underside of the seatplate.

3 Product Overview

3.1 Intended Use

Aquatec Sansibar is a wall-mounted shower seat composed of a seat with hygiene cut-out and supportive legs. This product is intended to be used as a shower aid.

Intended Users

Users are any adults and adolescents with restricted or limited mobility or limited muscle strength. If used without an assistant, the users must be able to keep their own balance.

Indication

- Restricted or limited mobility (e.g. legs and/or hips) or
- having balance issues to stand up or walk due to transitional or permanent disabilities.

There are no contraindications known when using the product as intended. For the maximum user weight, see 8.1 *Technical Data (Fig. 14), page 5*.

3.2 Scope of Delivery (Fig. 1)

- **A** Template
- **B** Shower seat
- **C** Plug (5x)
- **D** Adapter (4x)
- **E** Adjustable foot (2x)
- **F** User manual

4 Setup (Fig. 2 – 12)

CAUTION! Risk of Injury

- Mount only to walls that have a sufficient load-bearing capacity. (Fig. 2)
- Use only suitable fastening materials.

1. Measure the individual seat height of the user. (Fig. 3)
2. Use the adapters to adjust the height of the supporting legs to the required seat height h (Fig. 4).

The total height G of the supporting legs of the folding shower seat is 291 mm (+ 23 mm per adapter) (Fig. 5). Use max. three adapters per side.

3. Measure the depth T of the shower tray. (Fig. 6)
4. Deduct the depth T of the shower tray from overall height G .
5. Measure the resulting height H on the template and cut the template horizontally. (Fig. 7)
6. Place the template on the wall so that the cut edge is flush with the upper edge of the shower tray. (Fig. 8)

NOTICE!

- Make sure that the mounting surface on the wall is even.

7. Use a scribing tool to mark the five mounting points indicated by the template. (Fig. 8)
8. Use a drill to make suitably sized holes (8 mm) at the points marked on the wall. (Fig. 9)
9. Use suitable materials (wall plug, washer and screw) to fix the folding shower seat to the wall. (Fig. 10)

Make sure that the pivot mechanism at the back sits correctly into the groove on the mounting plate. (Fig. 10)

10. Place the plugs supplied over the holes. (Fig. 11)
11. Use a spirit level to level out correctly (Fig. 11). Rotate the leg extension tip counterclockwise (for lengthening) or clockwise (for shortening) (Fig. 12).

5 Usage

WARNING! Risk of Overturning

Ensure that both supporting legs are set to the same height.

WARNING!

Risk of Trapping/Squeezing

- When folding the seat down, ensure that your feet do not get trapped under the legs.

1. Fold up the shower seat after use. (Fig. 13)

6 Maintenance

6.1 General Maintenance Information

This product is maintenance-free, provided that the cleaning and disinfection instructions are adhered to.

6.2 Cleaning and Disinfection

General Safety Information

CAUTION!

Risk of Contamination

- Take precautions for yourself and use appropriate protective equipment.

NOTICE!

Wrong fluids or methods can harm or damage the product.

- All cleaning agents and disinfectants used must be effective, compatible with one another and must protect the materials they are used to clean.
- Never use corrosive fluids (alkalines, acid etc.) or abrasive cleaning agents. We recommend an ordinary household cleaning agent such as dishwashing liquid, if not specified otherwise in the cleaning instructions.
- Never use a solvent (cellulose thinner, acetone etc.) that changes the structure of the plastic or dissolves the attached labels.
- Always make sure that the product is completely dried before taking it into use again.

For cleaning and disinfection in clinical or long-term care environments, follow your in-house procedures.

Cleaning Intervals

NOTICE!

Regular cleaning and disinfection enhances smooth operation, increases the service life and prevents contamination. Clean and disinfect the product:

- regularly while in use,
- before and after any service procedure,
- when it has been in contact with any body fluids,
- before using it for a new user.

Cleaning Instructions

NOTICE!

Maximum cleaning temperature:
– 60 °C by hand

1. Clean the product by wiping down all generally accessible surfaces.
2. Dry the parts with a cloth.

Disinfection Instructions

Method: Follow the application notes for the used disinfectant and wipe-disinfect all accessible surfaces.

Disinfectant: Ordinary household disinfectant.

Drying: Allow the product to air-dry.

7 After Use

7.1 Reconditioning

This product is suitable for reuse. To recondition the product for a new user, carry out the following actions:

- Inspection
- Cleaning and disinfection

For detailed information, see 6 Maintenance, page 4.

Make sure that the user manual is handed over with the product.

If any damage or malfunction is detected, do not reuse the product.

7.2 Disposal

Be environmentally responsible and recycle this product through your recycling facility at its end of life.

Disassemble the product and its components, so the different materials can be separated and recycled individually.

The disposal and recycling of used products and packaging must comply with the laws and regulations for waste handling in each country. Contact your local waste management company for information.

8 Technical Data

8.1 Technical Data (Fig. 14)

| Dimensions | |
|-------------------------------|---------------|
| Ⓐ | 400 mm |
| Ⓑ | 410 mm |
| Ⓒ | 460 — 530* mm |
| Ⓔ | 460 mm |
| Ⓕ | 530 — 600 mm |
| Width, hygiene recess | 165 mm |
| Depth, hygiene recess | 155 mm |
| Depth, folded | 110 mm |
| Height, folded | 540 — 610 mm |
| Graduation, height adjustment | 23 mm |

| Weights | |
|------------------|--------|
| Product weight | 2.6 kg |
| Max. user weight | 130 kg |

| Materials | | | |
|-----------|----------------------------|---|----|
| ① | PP, glass fiber reinforced | ② | PP |

* Third height adapter optional

All product components are corrosion resistant and do not contain natural rubber latex.

8.2 Environmental Parameters

| Operating Conditions | |
|---------------------------------|---------------------------|
| Ambient temperature | 10 - 40 °C |
| Relative humidity | 15% - 93%, non-condensing |
| Storage and Shipping Conditions | |

| | |
|---------------------|---------------------------|
| Ambient temperature | 0 - 40 °C |
| Relative humidity | 30% - 75%, non-condensing |

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Invacare AG
Invacare Poirier SAS
Invacare GmbH
Invacare UK Operations Limited