Yes, you can.[®]

Invacare[®] ScanBed 750TM

User Manual (GB)

Brugermanual (DK)

Bruksanvisning (SE)

Brukerveiledning (NO)

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Manuale d'uso (IT)

Table of contents

User part	6
1. General information	6
2. Operation of the SB750™	7
Technical part	10
3. Information	10
4. Receiving the SB750™	11
5. Fitting the accessories	13
6. Disassembly/assembly of the SB750™	15
7. Wiring	17
8. Order numbers for accessories	19
9. Maintenance and check-ups	21
10. Beds equipped with accumulator back up	21
11. Maintenance chart	22
12. Trouble-shooting the electrical system	23
13. Lubrication plan	23
14. Cleaning	24
15. Technical specifications	24
16. Electrical data	25
17. Weight	26
18. Disposal	26

User part

Congratulations on your choice of *Invacare*® EC-Høng's **SB750™** nursing bed.

- The bed has been developed for patients over 12 years of age requiring care in a nursing home or home care nursing.
- The bed combines excellent stability and ergonomic design with easy disassembly and operation.
- The unique backrest return of the bed provides the patient with excellent comfort and the carer with ideal conditions for an ergonomically correct working position.

In order to optimise the patient's comfort, *Invacare*® recommends using a 12-cm mattress.

I. General information

The **SB750™** bed fulfils all requirements regarding maximum distances. However, if the bed is used for the care of patients with small body dimensions, it must be especially noted that there is a risk of such a patient slipping through the openings between the side rails or through the opening between the side rail and the mattress support.

The bed must not be used by patients under 12 years of age, or by patients with body size equivalent to an average 12 year old or smaller.

The bed, in combination with side rails, must not be used for persons under 45 kg or under 150 cm's height, or for persons that are restless (spasms) or confused, unless

- a professional risk assessment has taken place and been accepted, or
- a correctly installed safety cover on the side rails are used.

When using side rails and the according safety cover, it is essential to ensure correct fitting - **otherwise there is danger of entrapment / suffocation between mattress support, side rail and bed end!**

If the bed is used by restless (spasms) or confused persons, an ACP box can be used for blocking the hand control functions, and / or a limitation cross bar installed under the leg section of the mattress support. This cross bar will prevent the leg section from being lowered below the level of the mattress support, thereby reducing the opening between leg section and bed end.

Ensure the bed is adjusted to it's lowest position before leaving the bed unattended - thereby the effect of fall-down / entrapment accidents is reduced.

Invacare® accepts no liability for any use, change or assembly of the product other than as stated in this User's Manual.


The bed must be placed so that the height adjustment (up/down) is not obstructed by, for example, lifts or furniture. Otherwise there is a risk for injury on user and/or equipment.


When entering or exiting the bed, the height adjustment should be used. Furthermore, the backrest can be used as support (yellow button on the hand control). However, the thigh- and leg section must be horizontal, otherwise, there is risk of overloading of the mattress support.

2. Operation of the SB750™

The bed can be equipped with this type of hand control (Soft Control):


Or with this type of hand control (HB 70):

Operating the leg section on beds with gas springs

Lift the handle with one hand, while adjusting the leg section to the desired position with the other hand.
Release the handle.

Operating the leg section on SB750, 120 cm width, with Rastofix

To raise the leg section: Lift the leg section to desired height.

NB: Please ensure that both Rastofix brackets are correctly locked in place.

To lower the leg section: Lift the leg section slightly and lower to horizontal.


Always leave the bed in the lowest position. Otherwise there is a risk of squeezing due to accidental lowering of the mattress support. A person can be seriously injured under the bed during height adjustment.

Operating the braking castors

Operating castors without central braking system

When the bed is positioned correctly, at least one castor at the head end and one castor at the foot end must be locked.

- 1) **Braking:** Step on the pedal.
- 2) **Releasing the brake:** Step on the release pedal.

Operating castors with central braking system

When the bed is positioned correctly, it must be locked.

- 1) **Braking:** If the brake is in neutral, step on the red pedal.
- 2) **Releasing the brake:** Step on the green pedal, until the brake is in neutral position.

Operating the steerable castor

The **SB750TM** with central braking system may be equipped with a steerable castor. The steerable castor is operated by means of the central braking pedal.

- 1) **Activating the steering:** If the brake is in neutral, step on the green pedal.
- 2) **Deactivating the steering:** Step on the red pedal, until the brake is in neutral position.



Set-off from the castors might under special conditions appear at different types of absorbing floor covering - including untreated or badly treated floors. In matters of doubt, *Invacare*[®] recommends to place a suitable kind of protection between the castors and the floor.

Operating the BRITT II/III/IV wooden side rail

The bed end may be equipped to lock the side rail at half height as well as in the top position.
Normally the side rail can only be locked in top position.

Up: Pull up the top wooden side rail bar, until the locking pin locks with an audible click.


Down: Lift the top wooden side rail bar and press the two locking pins together. Lower the side rail.

Operating the steel side rail

Up: Lift the top steel bar of the side rail towards the end with the locking mechanism.

Down: Press the release button and pull the top bar of the steel side rail away from the locking mechanism.

There is risk of entrapment of fingers on assembly and operation of the side rails.
Please ensure that the side rail locking system has engaged properly by pulling / pushing the upper bar.


Adjusting the height of the lifting pole handle

Loosen the cord as shown in fig. A. The lifting handle can now be adjusted to the desired height.

Press the cord together as shown in fig. B and check that the cord is locked in the cord lock by pulling the handle.

Position the lifting pole in such a way that the handle is above the bed. If the lifting pole is used while the handle has been turned away from the bed - the bed can tip when the handle is used.

A

B

Technical part

Invacare[®] is certified according to DS/EN ISO 9001 and ISO 13485 which ensures that our customers are always supplied with *Invacare*[®] products of uniform quality.

Throughout the entire production process our materials and products are quality controlled by the operators. Furthermore there is a final test when the product is fully assembled.

The operator carrying out the final test, which comprises of checking all moveable parts, motors/gas springs and castors, issues a QA number on the product thereby confirming its quality.

QA XXX

If the product does not correspond to the quality demands of *Invacare*[®] it will be discarded.

If, contrary to our expectations, a problem should arise in connection with the delivered product, please contact your *Invacare*[®] supplier.

3. Information

Please read the entire technical information section carefully before using or servicing the bed. All indications of right and left are based on a patient lying on his back in the bed. Please note that there may be sections in this User's Manual, which are not relevant for your bed.

- The **SB750TM** is CE-marked in accordance with directive 93/42/EEC concerning Medical Apparatus.
- The **SB750TM** is tested without remarks by HMI according to EN 1970.
- **SB750TM** is tested and approved in accordance with NF (Norme Francaise).
- **SB750TM** is EMC tested and approved in accordance with EN 60601-1-2:93 by UL Demko.
- The motors and control unit of the **SB750TM** are approved according to EN 60601:1996-03.
- **SB750TM** has undergone a risk analysis according to EN 14971.

Max. patient weight: 185 kg.

The hand control, control unit and motors are protected according to IP 66.

A lock cam must be used on the control unit - if not, *Invacare*[®] cannot guarantee the IP protection.

Electromagnetic interference between the bed and other electrical products can occur. To reduce or eliminate such electromagnetic interference, increase the distance between the bed and the products or switch them off.

This medical bed can be used together with medical electrical equipment connected to the heart (intracardially) or the blood vessels (intravascularly), provided that following points are respected:

- The medical bed should be equipped with means for potential equalization connection marked out by a symbol shown in the back of this manual.
- Medical electrical equipment should not be fixed on the bed's metallic accessories such as side rails, lifting pole, drip rod, bed ends ect.

In addition, the medical electrical equipment power supply cord should be kept clear of the accessories or any other moving part of the bed.

Invacare[®] accepts no liability for any use, change or assembly of the product other than as stated in this User's Manual.

If the functions of the bed change, please send the bed immediately for a check-up according to the maintenance chart, chapter 11.

4. Receiving the SB750™

To avoid condensation, the bed should not be used until it has reached a temperature of 10-50° C.


Check whether the bed shows any signs of damage. If the bed is damaged, see terms of delivery.

If the bed has been delivered in parts, it must be assembled.
Cf. chapter 6 “Disassembly/Assembly of the SB750™”.

Check, that the connection between mattress support and base frame is properly locked (cf. the picture).

Check that all plugs of the motors and hand control are correctly connected to the control unit according to the printed symbols.

**Locked pin/
snap lock**

**Unlocked pin/
snap lock**

The bed must have the brakes engaged during nursing of a patient in the bed.

Disconnect the plug from the mains before moving the bed. The cable must be kept clear off the floor and the castors during transport.


Connect the SB750™ to the mains.

For ergonomic working position: Run the bed to approx. ½ height level.
Cf. chapter 2 “Operation of the SB750™”.

Mount the accessories: Bed ends, side rails, lifting pole and bumper wheel.
Cf. chapter 5 “Fitting the accessories”.

Parts of a standard bed

1. Allen key
2. Backrest
3. Seat section
4. 24 V motor, leg section or gas spring, leg section
5. Thigh section
6. 24 V motor, thigh section
7. Leg section
8. Top frame, foot end
9. Brake pedal for beds provided with central brake
10. Brake for beds provided with brake castors
11. Castor
12. Base frame
13. 24 V motor, backrest
14. Shear arms
15. Top frame, head end
16. Control unit

Accessories

- Battery box
- Fittings
- Bed end
- Wooden side rail
- Steel side rail
- Finger screw
- Lifting pole
- Bumper wheel

Emergency release of the mattress support parts

An emergency release of the mattress support is necessary in the case of e.g. a power- or motor failure. An emergency release of the height adjustment is NOT possible!

- Remove the plug from the mains before emergency release of the mattress support.
- In an emergency the mattress sections are released by pulling out the cotter pin from the motor/gas spring in question.


A minimum of 2 persons are required to release a mattress section. Both persons hold the mattress section in locked position. One of them pulls out the cotter pin. Both slowly lower the mattress section until it is completely down.


5. Fitting the accessories

The tools necessary for mounting/dismounting the accessories of the **SB750™** are an Allen key and an adjustable wrench. The Allen key is kept on the inside of the cross bar at the head end of the top frame. The Allen key is used for mounting/dismounting all fittings.


Mounting the mattress retainer

Mattress retainer is squeezed down over the mattress frame.

To be placed between third and fourth lamella on backrest.


To be placed between third and fourth lamella on leg section.

Mounting the bed end fittings

(If fittings have not been mounted in the factory)

1. Mount the 2 small fittings at the head end and tighten them with the Allen key.
2. Mount the large fittings and the welded fitting at the foot end and tighten them with the Allen key.

Mounting the PIGGY and ODA bed ends

Press the bed end into the U-profiles.

Mounting the FLEX III/KARIN III bed ends

Fasten the bed end to the bed end fittings with the bolts and tighten them with the Allen key.

Mounting the SANNE III bed end

Lift the bed end and place on the fittings at both sides. Push down to click into place.

The locking pin must then be rotated and locked into position. Make sure that the bed end is locked.

Dismounting the SANNE III bed end


Pull the locking pin out, and turn it halfway.

Repeat at the other end. Pull the bed end up and out.

Mounting/dismounting the steel side rail


Mounting: Referring to labelling on the side rail Tighten the side rail with 2 finger screws.

Dismounting: Loosen the finger screws, and remove the steel side rail.




Mounting the BRITT II/III/IV side rail

1. Raise the bed about 1/3.
 2. Remove the protective tape from the locking pins at the ends of the side rail.
 3. Mount one end at a time. Lift the lower wooden bar and guide the metal bow to the bed end rail (A).
 4. Press the locking pin of the side rail with a finger (B).
 5. Guide the side rail to the bed end until the spring support snaps into place and prevents the side rail from falling out of the bed end.
- Repeat points 3 to 5 when mounting the other end of the side rail.

Dismounting the BRITT II/III/IV side rail

1. Raise the bed about 1/3. Lower the side rail.
 2. Dismount one end at a time. Lift the lower wooden bar so that the spring support is visible.
 3. Press the spring support down with a wrench or a screwdriver (C).
 4. Hold the lower wooden bar. Lower the side rail, until the metal bow is free of the bed end rail.
- Repeat points 2 to 4 when dismounting the other end of the side rail.

Fitting the lifting pole

Remove the plastic stopper in the lifting pole tube in which the lifting pole is to be mounted.
Mount the lifting pole - the finger screw must be used.


Lifting pole tube with plastic stopper

Finger screw

Fitting the bumper wheel

Mount the bushing in the bumper wheel - and tighten it with a screw under the fitting.
Bumper wheels may not be mounted on SANNE bed ends.


Mounting the limitation crossbar under the leg section


If the bed is used by restless (spasms) or confused persons a limitation cross bar installed under the leg section of the mattress support. This cross bar will prevent the leg section from being lowered below the level of the mattress support, thereby reducing the opening between leg section and bed end.

Raise the leg section of the mattress support and the lifter arm, and mount the limitation crossbar.
Tighten both screws.

The limitation crossbar is placed here

The screws are tightened

Correctly mounted crossbar


Operating the built-in mattress support extension

Warning! In order to avoid entrapment / suffocation, it is essential to avoid extending the bed ends without extending the mattress support accordingly.
See page 14: Mounting the limitation crossbar under the leg section

The mattress support extension is mounted using the Allen key.
The Allen key is kept on the inside of the cross bar at the head end.

1. Remove the wooden side rails (cf. chapter 5).
2. Remove the 2 screws for fastening the fitting below the side bar by means of the Allen key.
3. Pull out the fitting so that it extends 16 cm, insert the 2 screws again and retighten them.
4. Loosen the 2 screws below the mattress handle on the leg section, pull the handle completely up, and retighten the screws.
5. Insert the mattress support extension.

6. Disassembly/assembly of the SB750™


Disassembly

1. Make sure that the castors have the brakes on.
2. Dismount all accessories.
3. Run the bed to the marked height level so that the cams on the base frame's guide way and the shear arms are opposite each other.
4. Disconnect the power supply.
5. Remove the locking cam from the control unit.
6. Disconnect the plugs of backrest, thigh and leg motors and check that the cables are not caught in the base frame.
7. Remove the 2x2 finger screws in the side bars.
8. Pull the head top frame out of the inserts in the side bars to stop and lower the section onto the floor.
9. Release the 2 lock fittings/snap locks on the foot top frame by pulling out the locking pin and turning it to the right.
10. Lower the foot-section of the top frame and lift it off the base frame.
11. Connect the base frame to the mains again and lower it completely.


3.


7.

8.


10.

Assembly


1. Make sure that the castors have their brakes on.
2. Connect the base frame to the mains.
3. Run the bed to the marked height level so that the cams (see picture 3, page 15) on the base frame's guide way and the shear arms are opposite each other.
4. Lift the foot end of the top frame so that the two 2 lock fittings/snap locks fit onto the pins on the shearing arms. (4.1)


Make sure that the locking pin has clicked into place and is locked (4.2). Then move the foot end of the top frame into horizontal position.

5. Place the head end of the top frame in front of the base frame and lower it down so that the plastic gliders slide into the guide way.


6. Lift the head end of the top frame and push it into the inserts until it stops.



7. Screw the 2x2 finger screws into the side bars.
8. Connect the plugs of the backrest, thigh and leg motors to the control unit, see section 7 regarding wiring.
9. Mount the accessories.


7. Wiring

The wires are mounted from the bottom and through the case of the shaft end of the motor.

Guide the wires under the lamellas.


Mount the wires through the cable bow and further to the control unit.

These photos show a bed with 2 motors for the mattress support.


Control the wiring by operating all the motors of the bed to their outer positions.

When all wires have been correctly mounted, there is no squeezing risk of the wires.

It is normal that the wires after short usage loosen a bit.

ACP Box (optional)

The hand control can be locked by activating the switch on the ACP box located in the head end of the bed.
To check if the locking function has been activated, depress a button on the hand control.

Open

Locked

**Removal of
safety lock**

Replacement of hand control fitted to the ACP-Box

The white locking ring on the ACP box is depressed, after which the hand control cable can be removed.

8. Order numbers for accessories

Article	Dimensions	Order no.:
WOODEN BED END		
SANNE III (beech)	85 cm	020429.01
SANNE III (beech)	90 cm	020740.01
SANNE III (beech)	105 cm	021953.01
SANNE III (beech)	120 cm	021954.01
SANNE III (cherry / beech)	90 cm	021956.07
SANNE III, increased and lowered (beech,) *	85 cm	1451258-0101
SANNE III, increased and lowered (beech,) *	90 cm	1451257-0101
SANNE III, increased and lowered (beech,) *	105 cm	1484403-0101
SANNE III, increased and lowered (beech,) *	120 cm	1484404-0101
*) To be used with the LILLI III wooden side rail		
FLEX III (beech)	85 cm	020417.01
FLEX III (beech)	90 cm	020737.01
PIGGY (beech / U-profiles)	85 cm	020431.01
PIGGY (beech / U-profiles)	90 cm	018940.01
ODA (beech / U-profiles)	85 cm	020430.01
ODA (beech / U-profiles)	90 cm	016609.01
Bed ends are available in other versions than the above mentioned - contact <i>Invacare</i> ®		
Side rails		
Scala Basic Metal side rail (foldable)		1432784-1015
Scala Basic Plus Metal side rail (foldable, with plastic insert)		1432787-1015
Scala Medium Metal side rail (foldable)		1432781-1015
Scala Decubi Metal side rail (foldable)		1432793-1015
Metal side rail Model 5539, extended (foldable, with loose finger screws)		50.55390.D0/L+R
BRITT IV wooden side rail (high + 2.5 cm, beech)		1427875-0101
BRITT IV wooden side rail (high + 2.5 cm, cherry)		1427875-0102
BRITT IV wooden side rail (high + 2.5 cm, white pigmented)		1427875-0105
BRITT III wooden side rail (beech)		020434.01
BRITT III wooden side rail, short (beech)		021342.01
BRITT III wooden side rail (cherry)		020434.02
BRITT III wooden side rail, short (white pigmented)		021342.05
BRITT III wooden side rail (white pigmented)		020434.05
BRITT II wooden side rail, extended (beech)		019372.01
BRITT II wooden side rail, extended (cherry)		019372.02
BRITT II wooden side rail, extended (white pigmented)		019372.05
LILLI III wooden side rail (beech)		1433812-0101
LILLI III wooden side rail (cherry)		1433812-0102
LILLI III wooden side rail, extended (beech)		1434044-0101
Covering tubes for mattress support extension		017639.D0
SUPPORT HANDLES		
Support handle	25x30 cm	021964.D0
Support handle	40x50 cm	021963.D0
Support handle	25x80 cm	1417510-1015
Support handle	40x30 cm	1417511-1015
Support handle	40x95 cm	1417512-1015
LIFTING POLE		
Lifting pole		50.57600.D0
LIMITATION CROSSBAR		
Limitation crossbar (prevents the leg section from being lowered below the level of the top frame)		1510820-1015

Article	Dimensions	Order no.:
DRIP		
Drip rod		50.60030.D0
Drip rod (for mounting in lifting pole tube)		50.60910.00
BUMPER WHEELS		
Bumper wheels (one)	75 mm.....	50.59600.00
Bumper wheels (one)	100 mm.....	50.59610.00
Bumper wheels (one)	150 mm.....	50.59640.00
Bumper wheels (one, vertical)	100 mm.....	50.59630.D0
MATTRESSES		
Dacapo Standard (38 kg/m ³ cold foam, cotton cover)	85x200x10.....	1421530
Dacapo Standard (38 kg/m ³ cold foam, cotton cover)	90x200x10.....	1421531
Dacapo Standard (38 kg/m ³ cold foam, cotton cover)	105x200x10.....	1421532
Dacapo Standard (38 kg/m ³ cold foam, cotton cover)	120x200x10.....	1421533
Dacapo Basic (38 kg/m ³ cold foam, incontinence cover)	88x200x12.....	1421551
Dacapo Basic (38 kg/m ³ cold foam, incontinence cover)	83x200x12.....	1421552
Dacapo Basic (38 kg/m ³ cold foam, incontinence cover)	88x183x12.....	1421555
Dacapo Basic (38 kg/m ³ cold foam, incontinence cover)	83x183x12.....	1421556
Dacapo Basic Heavy user (3 coated 50/50/38 kg/m ³ cold foam, incontinence cover) ..	105x200x15.....	1484429
Dacapo Basic Heavy user (3 coated 50/50/38 kg/m ³ cold foam, incontinence cover) ..	120x200x15.....	1451262
Dacapo Organic (50 kg/m ³ Organic foam, egg tray shape, incontinence cover)	85x200x12.....	1421572
Dacapo Organic (50 kg/m ³ Organic foam, egg tray shape, incontinence cover)	90x200x12.....	1421573
Dacapo Combi (2 coated 50/38 kg/m ³ , Visco-elastic/cold foam, incontinence cover) ...	85x200x14.....	1421489
Dacapo Combi (2 coated 50/38 kg/m ³ , Visco-elastic/cold foam, incontinence cover) ...	90x200x14.....	1421490
Dacapo Square (2 coated 50/38 kg/m ³ , Organic foam/cold foam, incontinence cover) ..	83x200x14.....	1421515
Dacapo Square (2 coated 50/38 kg/m ³ , Organic foam/cold foam, incontinence cover) ..	88x200x14.....	1421516
Dacapo Top (50 kg/m ³ , Visco-elastic foam, incontinence cover)	85x200x6.....	1421587
Dacapo Top (50 kg/m ³ , Visco-elastic foam, incontinence cover)	90x200x6.....	1421588
It is possible to buy custom-fit insertion mattresses for extended beds.		
LOW - SLUNG SET		
Low-slung set with Ø75 mm braking castor		021343.D0
SIDE RAIL COVERS		
Scala Cover for Britt III, IV and Lilli III wooden side rail		1449412
Scala Cover for Britt III wooden side rail, short		1449411
Scala Cover for Britt II and Lilli III wooden side rail, extended		1493565
Scala Cover for Britt III, IV and Lilli III wooden side rail, net		1449413
Scala Cover for Britt II and Lilli III wooden side rail, net, extended		1449415
Scala Cover for metal side rail, Basic		1449423
Scala Cover for metal side rail, increased		1449419
Scala Cover for metal side rail, extended		1449421
Safety Cover for Britt III, IV and Lilli III wooden side rail		1511890
Safety Cover for Britt III wooden side rail, short		1513592
Safety Cover for Britt II and Lilli III wooden side rail, extended		1511892

9. Maintenance and check-ups

Service and maintenance of the **SB750™** must only be performed by personnel who have received the necessary instruction or training.

With normal daily use, service must be carried out according to the service schedule after 2 years use and thereafter every second year.

When moving the bed, service must be carried out according to the service schedule.

When using Flex or Karin bed ends, which are made of a flexible type of wood, it will be necessary to tighten the hexagon head screws at intervals other than described in the service schedule.

Note! The top frame must be supported during service inspections to prevent accidental lowering.

We recommend a safety test comprising the motors' performance and mechanical state once a year.

Motors, hand control and control unit

These parts are serviced by exchanging the faulty part.

10. Beds equipped with accumulator back up

A bed with accumulator back up is equipped with an additional box next to the control unit.


Preventive maintenance

The accumulators have to be exchanged after 4 years. Depending on how the back up is used the exchange may have to take place earlier.

Frequent, sudden discharges reduce the back-up's life.

We recommend the accumulator be tested at least once a year.

Accumulators are not damaged by continuous connection to the mains.

Exchanging the accumulators

Accumulators must be exchanged as sets and replaced with accumulators of the same type or mechanically and electrically compatible accumulators:

The accumulators are delivered as a pair under order no.: 818323 (12V - 1,2 Ah).

Accumulators must be new and charged at least every 6 months for maintenance purposes.

Accumulators in a set must have identical production codes.

Prior to insertion make sure that the accumulator set is correctly connected cf. the drawing (in the accumulator box), and that none of the connections are loose.

Warning

Old or faulty accumulators may generate an explosive gas mixture during charging.

The accumulator box is provided with vents to ensure adequate ventilation of the box.

The vents must not be blocked or covered, as this may result in pressure built-up and the risk of explosion.

Waste disposal

Old accumulators may be returned to *Invacare®*, alternatively they are disposed of like car batteries.

II. Maintenance chart

Service and maintenance of the SB750™ must only be performed by personnel who have received the necessary instruction or training

S/N (located on mattress support): _____

Date:	Initials:					
Circlips, cotter pins and plastic fixing ring checked.						
Screws tightened.						
Weldings checked.						
Side rail locking and moving system tested.						
Castor fittings tightened.						
Castor brakes checked.						
Height adjustment motor checked.						
Back rest motor checked.						
Thigh section motor checked.						
Leg section gas spring/motor checked.						
Lead-in connections correctly suspended and undamaged. Plugs undamaged.						
Damaged coating repaired.						
Lubrication according to the lubrication plan. NB! The wooden side rails' gliding system must not be lubricated with oil - otherwise the wooden bars will move sluggishly.						
Accessories checked.						


A service contract can be made in the countries, where *Invacare*® has its own sales company.

In certain countries *Invacare*® offers courses in service and maintenance of the **SB750™**.


Spare parts lists and additional user manuals are available from *Invacare*®.

12. Trouble-shooting the electrical system

Note! Make sure that the mains are connected.

13. Lubrication plan

We recommend lubricating the **SB750™** according to the following instructions:

1. Lifting arm bearings - lubricate with oil
2. Motor/gas spring bearings - lubricate with oil
3. Gliders and guide ways - lubricate with grease
4. Shear arms and rollers shaft - lubricate with oil/grease

Lubricate with medically clean oil, e.g. KEW-WO 50, order no.: 813239 and grease, order no.: 1497607.


14. Cleaning

The **SB750™** cannot withstand cleaning in washing machines nor flushing with a water jet. Use ordinary disinfecting detergents. Dry the bed after cleaning. Never use acids, alkalines or solvents. Please make sure the power plug is removed from the socket outlet.

The back rest and the leg/thigh sections can be lowered onto the mattress support, which facilitates cleaning and access to the control unit.

Important details for cleaning hand control, control unit and motors IP 66: Always clean with water and brush. Water may be pressurised, but never use high-pressure cleaning or steam cleaning directly on the electrical parts.

15. Technical specifications

WxL	A	B	C	D	E
85x200	200	213	85	95	60
90x200	200	213	90	101	60
85x183	183	196	85	95	43
90x183	183	196	90	101	43
105x200	200	213	105	116	60
120x200	200	213	120	131	60


All measurements are stated in cm. All angles are stated in degrees. All measurements and angles are stated without tolerances. *Invacare®* reserves the right to change the stated measurements and angles without warning.

Transportation and storage conditions:

Temperature between -10° and 50° C

Humidity between 10 and 80 % RH

Pressure between 700 and 1.100 hPa

16. Electrical data


Voltage supply:	230 V AC/50 Hz.
Max. current input:	1,5 A.
Intermittent (periodic use of motors):	10% max. 6 min/h.
Max. accumulator capacity:	1,2 Ah.
Protection class:	IP 66.
Output:	70 VA.
Sound level:	48 dB (A).

Double insulated, class II, type B.

Alternating current: 

Direct current: 

Max. load (SWL)
(Patient + mattress + siderail + lifting pole + other equipment).

The patient is not separated from the ground and the chassis.

Potential equalization.

The control unit has no contact breaker, disconnecting the plug is the only means of interrupting the power supply.

17. Weight

Measurement / Weight	85x200 90x200	105x200	120x200
Top frame, head, with mattress support	25 kg	34 kg	36 kg
Top frame, foot, with mattress support	25 kg	36 kg	39 kg
Base frame	40 kg	40 kg	42 kg
Complete bed, excluding accessories	90 kg	110 kg	117 kg
Bed end (SANNE III/FLEX III/KARIN III)	6 kg/pcs.	7 kg/pcs.	8 kg/pcs.
Wooden side rail (BRITT II/BRITT III)	7 kg/pcs.	7 kg/pcs.	7 kg/pcs.
Wooden side rail (BRITT IV)	8 kg/pcs.	8 kg/pcs.	8 kg/pcs.
Bed end (ODA/PIGGY)	6 kg/pcs.	-	-
Steel side rail	7 kg/pcs.	7 kg/pcs.	7 kg/pcs.
Lifting pole	6,5 kg	6,5 kg	6,5 kg

18. Disposal

This product has been supplied from an environmentally aware manufacturer that complies with the Waste Electrical and Electronic Equipment (WEEE) Directive 2002/96/CE.

This product may contain substances that could be harmful to the environment if disposed of in places (land fills) that are not appropriate according to legislation.

The »crossed out wheelie bin« symbol is placed on this product to encourage you to recycle wherever possible.

Please be environmentally responsible and recycle this product through your recycling facility at its end of life.

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Invacare AG, Benkenstrasse 260, CH-4108
Invacare Limited, Pencoed Technology Park,
INVACARE Rea AB