Yes, you can.®

Invacare® **Hand control**

Bed Accessories

User Manual

EN

Brugsanvisning

DA

Bruksanvisning

SV

Brukerveiledning

NO

Manuale d'uso

IT

This manual must be given to the user of the the product.

Before using this product, read this manual and safe for future reference

Hand control

type RA no. 1426936

Yes, you can.®

Description:

The hand control is intended for persons with limited use of their hand or fingers. It can be used alone or together with the hand control there was delivered together with the bed. In that case, use the Y-cable (order no. 819465)

The hand control can be used on SB400, SB750 and SB755 beds.

The hand control may be used only as directed by a doctor or therapist, and after it has assessed that the user is able to use the hand control safely.

There is an increased risk of unintended moves of the bed, because of the large buttons.

The hand control should not lie in the bed, but after use, placed on the side rail or on the bed frame.

The handcontrol is not water resistant and must be removed while the bed cleaned.

Approvals:

The hand control is not CE-marked and may only used after the prescription form is filled out and signed.

If the hand control is transferred to another user/patients a new form should be filled out. New form are available from *Invacare*®.

NB! The TÜV approval of the bed is invalid, when the hand control is connected.

Prescription form for use of hand control RA

This form must be filled in, signed and archived with the users/patient medical record before the hand control is used.

User/patient name: _____ User ID: _____

I hereby declare that :

- User/patient assess the ability to use the hand control type RA, especially with respect that the user/patient does not activate the hand control accidental.
- User/patient and attendants have been instructed in safe use of the hand control according to the user manual.

If the condition of the user/patient is changes, a new assessment repeated.
The hand control must be removed from the bed and not transferable to another user, before a new assessment is done and documented.

If it is necessary the hand control can be retracted.

Date Signature - Doctor / Therapist

Recall

The bed of the above mentioned user /patient is no longer equipped with hand control RA.

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INVACARE AUSTRIA GmbH
INVACARE EC-Høng A/S