# Invacare® Groin Band for Standing Transfer Vest

EN	<b>Slings</b>	User Manual . . . . .	1
DA	<b>Sejl</b>	Brugsanvisning . . . . .	3
DE	<b>Patientengurte</b>	Gebrauchsanweisung . . . . .	4
ES	<b>Eslinga</b>	Manual del usuario . . . . .	5
FR	<b>Sangles</b>	Manuel d'utilisation . . . . .	6
IT	<b>Imbracatura</b>	Manuale d'uso . . . . .	7
NL	<b>Liftmat</b>	Gebruiksaanwijzing . . . . .	9
NO	<b>Seil</b>	Bruksanvisning . . . . .	10
PT	<b>Cesta de transferência</b>	Manual de utilização . . . . .	11
SV	<b>Lyftselar</b>	Bruksanvisning . . . . .	12

**EN**

## I General

### I.1 Figures


Figur

Abbildung

Afbeelding

Figures

Figura

This manual **MUST** be given to the user of the product.  
BEFORE using this product, read this manual and save for future reference.

Yes, you can.®

## 1.2 Symbols

	Signal words are used in this manual and apply to hazards or unsafe practices which could result in personal injury or property damage. See the information below for definitions of the signal words. <b>WARNING!</b> Warning indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury. <b>CAUTION!</b> Caution indicates a potentially hazardous situation which, if not avoided, may result in property damage or minor injury or both.
	Gives useful tips, recommendations and information for efficient, trouble-free use.
	<b>IMPORTANT</b> indicates a hazardous situation that could result in damage to property if it is not avoided.
	This product complies with Directive 93/42/EEC concerning medical devices. The launch date of this product is stated in the CE declaration of conformity.
	Safe Working Load
	Consult the user manual.
	DO NOT exceed 95°C washing temperature.
	DO NOT bleach.
	Tumble Dry, Normal, Low Heat
	DO NOT iron.

## 1.3 Intended Use

The Invacare® Standing Transfer Vest sling is available with or without the groin band, the groin band can be attached to the transfer standing vest very easily.

The Standing Transfer Vest with groin band is suitable for clients with reduced standing/walking abilities and/or reduced balance. The sling provides maximum personal access for dressing and hygiene.

-  The Standing Transfer Vest lifts is applied around the upper part of the body, located in the under arm area. It is essential that the user has the ability to stand up, and has some degree of weight bearing ability to ensure that he feels contact with the floor.

## 1.4 Service Life

All Invacare slings are guaranteed for one year. Expected service life of the product is 1-5 years. Service life time varies depending on the fabric, frequency of use, laundering practices, and weight carried.

## 2 Safety

### 2.1 General Guidelines

#### WARNING!

- DO NOT use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as user manuals, service manuals or instruction sheets supplied with this product or optional equipment. If you are unable to understand the warnings, cautions or instructions, contact a healthcare professional, dealer or technical personnel before attempting to use this equipment - otherwise, injury or damage may occur.
- Please check your sling on delivery. If the sling is not in perfect condition, please contact your supplier.

### 2.2 Safe Working Load

		
<b>S</b>	1492662	159 kg
<b>L</b>	1492663	159 kg

## 3 Usage

### 3.1 Sling application

#### WARNING!

- Slings should be checked regularly for wear and damage (in accordance with each country's regulations). NEVER use a damaged sling.
- If you are in any doubt, please do not attempt to use the sling and seek the advice of a professional or your supplier.
- Invacare recommends that a risk assessment is always carried out by a professional prior to issuing lifting equipment. Always choose the sling design and size according to the client's weight, size and physical ability.
- If you are in doubt about the use of Invacare's slings with other manufacturers' equipment, contact Invacare or your supplier for advice.

### Sling application — Sitting in chair or wheelchair

1. Lean the client forward ensuring that he is well supported to reduce risk of falling and to promote feelings of security.
2. Slide the vest down behind the back of the client in a central position aligned with the waist.
3. Sit the client back against the sling. **Figure 2.**
4. Take the arms through the straps and fasten the vest around the client's stomach by using the adjustable buckle. **Figure 3.**
5. Place the groin band round the thigh and close the buckle (an audible click should be heard). The groin band helps to keep the Standing Transfer Vest in place.

### 3.2 How to fix the groin band to the Standing Transfer Vest

1. Attach the groin band to the loops placed at the back of the Standing Transfer Vest.
2. Attach the buckles to the loops placed at the front of the Standing Transfer Vest. **Figure 1.**

### 3.3 Attach loops to the spreader bar or lifting strap

#### WARNING!

- Use only CE marked spreader bars and hoists, preferably with a rated Safe Working Load equal to or higher than that of the sling. If the SWL is lower than that of the sling, the SWL of the spreader bar or hoist is applicable.
- Always check that all the straps are properly attached to the spreader bar or lifting strap before attempting to lift the client.

1. Attach the loops to the spreader bar or lifting straps of the lifting equipment.

The sling can be used with spreader bars with 2 attachment points.

For security and comfort reasons, the size (width) of the spreader bar must correspond with the size of the sling.

2. Before lifting the client, always check that both feet are placed on the floor for weight bearing.
3. Help the client to lean forward (in his natural way of movement) before beginning the lift. **Figure 4.**

Ensure that the client is lifted only high enough so that both feet do not lose contact with the floor.

4. When the client is in a standing position, the posture of the client can be adjusted by altering the shoulder straps. **Figure 5.**

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© 2013 Invacare® Corporation. All rights reserved. Reproduction, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.
#### Invacare Corporation