# Invacare® Sonata

en	Bed	
	User Manual .....	3
es	Cama	
	Manual del usuario .....	27
fr	Lit	
	Manuel d'utilisation .....	53
it	Letto ortopedico	
	Manuale d'uso .....	79
pt	Cama	
	Manual de utilização .....	105

This manual MUST be given to the user of the product.  
BEFORE using this product, read this manual and save for future  
reference.

**Yes, you can.**

©2017 Invacare Corporation

All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

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# Contents

---

This manual MUST be given to the user of the product.  
BEFORE using this product, read this manual and save for future reference.

<b>1 General</b>	<b>4</b>
1.1 Introduction	4
1.1.1 Symbols in this manual	4
1.2 Intended Use	4
1.3 Service life	4
1.4 Compliance	4
1.5 Warranty	4
1.6 Limitation of liability	4
<b>2 Safety</b>	<b>6</b>
2.1 General safety information	6
2.2 Mattresses	6
2.3 Labels and symbols on the product	7
2.3.1 Product label	7
2.3.2 Other labels and symbols	7
<b>3 Setup</b>	<b>8</b>
3.1 General Safety Information	8
3.2 Main parts of bed	8
3.3 Assembly of mattress support	8
3.4 Mounting cranks/motors	9
3.4.1 Manually operated mattress support	9
3.4.2 Electrically operated mattress support	9
3.5 Mounting of bed ends and legs	9
3.6 Installing side rails	11
3.6.1 Verso II side rail	11
3.6.2 Diana side rail	11
3.6.3 Flora side rail	12
3.7 Wiring	13
<b>4 Operating the bed</b>	<b>14</b>
4.1 General Safety Information	14
4.2 Operating the hand control unit	14
4.3 Operating the bed functions manually	14
4.4 Operating side rails	15
4.4.1 Verso II side rails	15
4.4.2 Diana side rail	15
4.4.3 Flora side rail	15
4.5 Castors and brakes	16
4.5.1 Castor brake	16
4.6 Emergency release of a mattress support section	16
4.7 Mattress retainers	16
4.8 Lifting pole	16
4.8.1 Place the lifting pole	16
4.8.2 Adjusting the handle height	17
<b>5 Accessories</b>	<b>18</b>
5.1 List of available accessories	18
<b>6 Maintenance</b>	<b>19</b>
6.1 Maintenance	19
6.1.1 Checklist maintenance	19
6.2 Cleaning and disinfection	19
6.2.1 Cleaning methods	19
6.3 Lubrication	19
<b>7 After Use</b>	<b>21</b>
7.1 Waste disposal	21
7.2 Reconditioning	21
<b>8 Troubleshooting</b>	<b>22</b>
8.1 Troubleshooting electrical system	22
<b>9 Technical data</b>	<b>23</b>
9.1 Dimensions	23
9.2 Mattress dimensions	23
9.3 Weights	23
9.4 Environmental conditions	23
9.5 Electrical system	23
9.6 Electromagnetic compliance (EMC)	24

# 1 General

## 1.1 Introduction

This user manual contains important information about the handling of the product. In order to ensure safety when using the product, read the user manual carefully and follow the safety instructions.

To ensure correct use, the bed must be tested and adjusted by qualified personnel.

All references to left and right are based on a person lying on his back in the bed, with his head in the head end.

If a problem should arise in connection with the delivered product, please contact your *Invacare*® dealer. An address list is shown on the back side at this manual.

### 1.1.1 Symbols in this manual

Symbols and signal words are used in this manual and apply to hazards or unsafe practices which could result in personal injury or property damage. See the information below for definitions of the signal words.

#### WARNING

Indicates a hazardous situation that could result in serious injury or death if it is not avoided.

#### CAUTION

Indicates a hazardous situation that could result in minor or slight injury if it is not avoided.

#### IMPORTANT

Indicates a hazardous situation that could result in damage to property if it is not avoided.

#### Tips and Recommendations

Gives useful tips, recommendations and information for efficient, trouble-free use.

This product complies with Directive 93/42/EEC concerning medical devices. The launch date of this product is stated in the CE declaration of conformity.

Manufacturer of the product.

## 1.2 Intended Use

The bed has been developed for domestic care and long term care:

- *Application environment 3*; Long-term care in a medical area where medical supervision is required and monitoring is provided if necessary and medical electrical equipment used in medical procedures may be provided to help maintain or improve the condition of the patient.
- *Application environment 4*; Care provided in a domestic area where medical electrical equipment is used to alleviate or compensate for an injury, disability or disease.
- The bed is intended for indoor use only.
- The bed is intended for adult users, having a physical size equal to or more than 146 cm, a weight equal to or more than 40 kg and a body mass index (BMI) equal to or more than 17.

- The bed is not intended for transportation of users. It is mobile within a room with a user in it. Castors are lockable.
- The bed is not intended for psychiatric patients.
- Maximum user weight and safe working load are listed on the product label and in the technical data section in this manual.

#### WARNING!

Any other or incorrect use could lead to hazardous situations.

*Invacare* accepts no liability for any use, change or assembly of the product, other than stated in this user manual.

## 1.3 Service life

The expected service life of this product is seven years when used daily and in accordance with the safety instructions, maintenance intervals and correct use, stated in this manual. The effective service life can vary according to frequency and intensity of use.

## 1.4 Compliance

Quality is fundamental to the company's operation, working within the disciplines of ISO 13485.

We are continuously working towards ensuring that the company's impact on the environment, locally and globally, is reduced to a minimum. We use only REACH and RoHS compliant materials and components.

The product is compliant with the European Directive 93/42/EEC concerning Class 1 medical devices.

The product has been tested and conforms to IEC 60601-2-52 – Medical Beds and all related standards. This includes tests regarding flammability and biocompatibility.

## 1.5 Warranty

The warranty covers all material and production defects for two years from the date of delivery, provided it can be demonstrated that such defects were present before delivery. All manufacturing faults or defects must be promptly reported.

*Invacare*® may repair the fault or replace the component. The warranty provided by *Invacare*® does not cover additional costs (transport, packaging, labour, sundry expenses, etc. are for the customer's account).

The warranty does not cover:

- Damage caused during transport that is not directly reported to the forwarder at the moment of delivery.
- Repairs performed by unauthorized centers and personnel.
- Parts subject to normal wear.
- Malicious damages or damaged caused by improper use of the bed.

## 1.6 Limitation of liability

*Invacare* accepts no liability for damage arising from:

- Non-compliance with the user manual
- Incorrect use
- Natural wear and tear

- Incorrect assembly or set-up by the purchaser or a third party
- Technical modifications
- Unauthorized modifications and/or use of unsuitable spare parts

## 2 Safety

### 2.1 General safety information

#### **WARNING!**

##### **Risk of entrapment / suffocation**

There's a risk of entrapment / suffocation between mattress support, side rail and bed end or between moving parts and objects placed nearby the bed.

- The bed must not be used by persons under 12 years of age, or by persons with a body size equivalent to an average 12 years old or smaller.
- The bed, in combination with side rails must not be used by persons having a physical size less than 146 cm, a weight less than 40 kg or a body mass index (BMI) less than 17.
- Due to mattress compression, an increased risk may occur over time. Periodically monitor gaps between the bed, mattress and/or side rail. Replace mattress if the gaps may lead to entrapment.

#### **WARNING!**

##### **Risk of slipping through the openings**

The bed fulfils all requirements regarding maximum distances. However, it is possible that persons with small body dimensions slip through the openings between the side rails or through the opening between the side rail and the mattress support.

- Pay special attention, if the bed is used for the care of persons with small body dimensions.

#### **WARNING!**

##### **Risk due to electromagnetic interference**

Electromagnetic interference between the bed and other electrical products can occur.

- To reduce or eliminate such electromagnetic interference, increase the distance between the bed and the products or switch them off.

This medical bed can be used together with medical electrical equipment connected to the heart (intracardially) or blood vessels (intravascularly) provided that following points are respected:

- Medical electrical equipment should not be fixed on the bed's metallic accessories such as side rails, lifting pole, drip rod, bed ends, etc.
- The medical electrical equipment power supply cord should be kept clear of the accessories or any other moving parts of the bed.

#### **WARNING!**

##### **Risk of injury or damage to property**

- Do not roll the castors over the main power cord.
- Do not bring main power cord into moving parts.
- Disconnect the plug from the mains before moving the bed.
- Make sure that no wires (mains or from other equipment) are jammed or damaged, when the bed is used.
- Keep bed components and accessories at least 30 cm away from a heated surface and not in direct sunlight.

#### **CAUTION!**

There is a risk of entrapment of fingers in the bed moving parts.

- Pay attention to your fingers.

#### **CAUTION!**

- For a user entering or exiting the bed, always lower the bed to an appropriate height. The backrest can be used as a support. Make sure, the thigh- and leg section is horizontal to avoid overloading the mattress support.
- Always lower the bed to the lowest position before leaving the patient in the bed unattended.
- Make sure that there is nothing under, over, or near the bed that can obstruct the height adjustment, like for example furniture, lifts or window frames.

### 2.2 Mattresses

#### **WARNING!**

##### **Safety aspects regarding combination of side rails and mattresses:**

To get the highest possible safety level, when using side rails on the bed, the minimum and maximum measures for mattresses, must be respected.


- For correct mattress measures see mattress table in chapter 9 *Technical data*, page 23.

#### **WARNING!**

##### **Risk of entrapment and/or suffocation**

- The user could get trapped and/or suffocate, if the horizontal space, between the mattress side and the inside of the side rail, is too big. Follow the minimum width (and length) of mattresses in combination with a side rail, as stated in the mattress table in chapter 9 *Technical data*, page 23 .
- Be aware that using very thick or soft mattresses (low density), or a combination of these, increases the risk.

**WARNING!**  
**Risk of falling**

The user can fall over the edge and get seriously injured, if the vertical distance  $\text{\textcircled{A}}$  between the top of the mattress and the edge of the side rail/bed end, is too short. See image above.  
 – Always keep a minimum distance  $\text{\textcircled{A}}$  of 22 cm.  
 – Follow the maximum mattress height in combination with the side rail as stated in the mattress table in chapter 9 *Technical data, page 23.*

	Max. Safe Working Load
	CLASS II equipment
	Type B Applied Part
	WEEE conform
	This product complies with Directive 93/42/EEC concerning medical devices.


Abbreviations for technical data:

- I<sub>in</sub> = Incoming Current
- U<sub>in</sub> = Incoming Voltage
- Int. = Intermittence
- AC = Alternating Current
- Max = maximum
- min = minutes

For more information about technical data, refer to 9 Technical data, page 23.

**2.3 Labels and symbols on the product**

**2.3.1 Product label**

The product label is placed on the frame of the bed and contains the main product information, including technical data.

**Symbols on the product label**

	Serial Number
	Reference Number
	Manufacturer Address
	Manufacturing Date
	Max. User Weight

**2.3.2 Other labels and symbols**

	Refer to User Manual
--	----------------------

**Label – user and mattress sizes**

	Definition of min. weight, min. height and min. body mass index of an adult user  Refer to user documentation for the correct mattress measures.
--	--------------------------------------------------------------------------------------------------------------------------------------------------------

**Label on side rails**

	Indicates the maximum space between side rail and the head end of the bed. (removable side rails only).
--	---------------------------------------------------------------------------------------------------------

## 3 Setup

### 3.1 General Safety Information

When you receive the bed, check the packaging. If the packaging shows any signs of damage upon delivery, contact the shipping company.

**WARNING!**

**Risk of injury or damage to property**

The assembly of the bed and installation of accessories must be done by authorized or trained personnel.

- Follow instructions carefully. If you have any question regarding the assembly, contact your local dealer or Invacare representative.
- If the bed shows any sign of damages, do not use the bed. Contact your local dealer.
- The electrical equipment of the bed must not be dismantled or combined with other electrical equipment.
- After each assembly, check that all fittings are properly tightened and that all parts have the correct function.
- 


**IMPORTANT!**

- To prevent condensation, the bed should not be used until it has reached a temperature of 10-40° C.

### 3.2 Main parts of bed

The bed can be configured either for manual or electrical operation.

Scope of delivery:

- Ⓐ Mattress support, head section
- Ⓑ Mattress support, foot section
- Ⓓ Pair of side rails
- Ⓔ Motor / control box or crank for head section
- Ⓕ Motor or crank for thigh section
- Ⓖ Pair of bed ends

- ⓓ Lifting pole (optional)
- Hand control
- 2 pcs. of mattress support retainers (2-section beds)
- 4 pcs. of mattress support retainers (4-section beds)

### 3.3 Assembly of mattress support

**CAUTION!**

**Risk of injuries or finger pinch**

- Keep your fingers away from movable parts during assembling, disassembling or adjusting.

**IMPORTANT!**

Before any use and before operating the bed:

- Check that the bed is correctly mounted according to instructions.
- Check that all screws are properly tightened.
- Ensure that the pipe pins are fully inserted through both holes of the fixation and that the spring clip is completely locked.

**After assembly of the mattress support; check:**


- that all plugs of the motors and hand control are correctly connected to the control unit according to the printed symbols.
- that top frame, head- and leg end are completely pushed together.
- that the lock pin is engaged.
- that bed ends are correctly attached.

Sonata can have either a complete or a separated mattress support.

If your bed is equipped with a complete mattress support, proceed to 3.5 Mounting of bed ends and legs, page 9 .

#### Assemble the separated mattress support

1.


Turn the mattress support up-side-down on a flat, clean surface.

2.

Place the two inserts in the head end of the mattress support and push the ends together.

3.

Tighten with the two thumb screws.

### 3.4 Mounting cranks/motors

- ! IMPORTANT!**  
If the pipe pin is not properly mounted the bed will not function as intended.
- Make sure the pipe pin is correctly mounted through all holes in each joint and properly locked with the brace.
  - Make sure that the pipe pin is mounted in the correct hole in the bed bracket.

#### 3.4.1 Manually operated mattress support

##### Mattress support; non-detachable (two or four sections)

No motors used. Crank/tubes already mounted.

##### Mattress support; detachable (four sections)

No motors used. Crank/tubes for leg section already mounted. Backrest crank/tube needs mounting.

##### Mounting of backrest crank/tube

1.

Join the backrest crank/tube with the backrest bracket, using a pipe pin in the hole of the bracket.

#### 3.4.2 Electrically operated mattress support

The longer motor (assembled on the control box) operates the backrest.

The shorter motor operates the leg section.

1.

Join the leg section tube with the leg section bracket, using a pipe pin in the upper hole of the bracket.

2.

Join the backrest tube with the backrest bracket, using a pipe pin in the hole of the bracket.


3.

Connect the section motor wires to the control box.

- Backrest motor to socket 1.
- Leg section motor to socket 2.

4.

Attach the locking cam over all motor plugs.

### 3.5 Mounting of bed ends and legs

- i** Depending on the configuration of the bed, the mattress support is either equipped with or without a bayonet system to mount the bed ends, which results in different options for bed ends and legs.

- Bed ends with integrated legs with bayonet system.
- Bed ends with integrated legs and brackets for mounting on the side tubes.
- Bed ends with brackets for mounting on the side tubes in combination with separate legs mounted on the mattress support platform.

#### CAUTION!

**Different height fixation of legs can cause instability in the bed**

On the bayonet system the mattress support platform can be fixed to the bed ends in two different heights.

Legs with adjustable height can be fixed in three different heights.

- Before using the bed, make sure that all fixations are set to the same height on all four legs.

#### Bed end with integrated legs and bayonet system

Latch the mattress support to the bed ends and press firmly into position. Do this one end at a time.

(1) High position

(2) Low position

Make sure that the bed ends are on the same level. It is possible to level the bed end in the high position (1) or the low position (2).

**Bed end with integrated legs without bayonet system**

Insert the integrated brackets **A** of the bed end into the side tube of the mattress platform and tighten the screws firmly with a 5 mm Allen Key. Do this one end at a time.

**Bed ends in combination with separate legs**

1.


For Sabrina, Camila bed ends, insert the brackets for the bed ends into the side tube of the mattress platform. Tighten the screws firmly with a 5 mm Allen Key.

2.

For the Vibeke bed end, insert the separate brackets for the bed ends into the side tube of the mattress platform and tighten the screws firmly with a 5 mm Allen Key. Lower the bed end board in the U-bracket fitting and press down firmly to ensure that bed end is fully seated at bottom of U-bracket fitting.

3.


For Nadine bed ends, insert the integrated brackets **A** of the bed end into the side tube of the mattress platform and tighten the screws firmly with a 5 mm Allen Key. Lower the bed end board in the U-bracket fitting **B** and press down firmly to ensure that bed end is fully seated at bottom of U-bracket fitting.

**Separate legs mounted on mattress support**

For bed ends without integrated legs, separate legs need to be mounted to the cross tubes of the mattress support.

Separate legs are available as fixed and adjustable height versions, with and without castors.

1.

Place the bracket on the cross tube.

2.

Insert the bolts and tighten the nuts.

### Height adjustable legs

The legs can be adjusted in 3 different heights with a pipe pin. Ensure the pipe pin is fully inserted through both holes and that the pipe pin is completely locked. Make sure that all four legs are adjusted to the same height.

## 3.6 Installing side rails

### WARNING!

#### Risk of personal damage

There is a risk of entrapment or squeezing, while assembling or disassembling the side rail.

- Pay attention to your fingers.
- Follow instructions carefully.
- After each assembly, check that all fittings are properly tightened and that the side rail has the correct function.

### Removable 3/4 length side rails

### WARNING!


#### Risk of entrapment

There's a risk of entrapment or suffocation between mattress support, side rail and bed end.

- When using removable side rails, always make sure that the distance between the bed end and the side rail's handle upper edge is less than 6 cm in the head end and more than 32 cm in the foot end of the bed.

### 3.6.1 Verso II side rail

#### Mounting

The mounting of the side rail must be done with the side rail in locked position.

1. Place both brackets of the side rail on the frame of the bed, with the locking button **A** facing towards the foot end of the bed.
2. Ensure that the bracket near the head end of the bed is placed encompassing the lifting pole support as shown in detail **B**.
3. For beds equipped with bed ends with Sanne style fitting, it is beforehand necessary to temporarily disengage the locking pin **C** from the head bed end to place the bracket.
  - a. Pull the locking pin **C** and rotate 1/4 turn to lock it in the disengaged position.
  - b. Position the bracket under the locking pin as described in step 2.
  - c. Rotate the locking pin and allow it to engage with the bed end.
4. After the brackets are positioned correctly on the frame of the bed, fix the side rail by firmly tightening the finger screws **D**.

#### Dismounting

1. Untighten the finger screws and lift the side rail of the bed frame.

### 3.6.2 Diana side rail

#### Compatibility


### WARNING!

- The Diana side rail cannot be used on Sonata beds, when the bed is configured with adjustable legs without castors and when the legs are mounted in the lowest height.

### Installing the side rail


1.

Pull out the pin and mount the bracket on the side rail. Make sure the pin locks the bracket and that the finger screw is on the outside of the side rail.


2. Do the same thing on the other side.
3. The installation on the bed must be done with the side rail in locked position.

4.

Loosen the finger screws Ⓐ on both brackets of the side rail.

5.

Place the side rail on the frame as close as possible to the head end of the bed and tighten the finger screws on both brackets.

#### ! IMPORTANT!

- Risk of squeezing or entrapment
  - The distance between the side rail and the head bed end must be smaller than 6 cm.

#### ! IMPORTANT!

- Make sure the side rail is properly locked.
  - Verify that the locking pins are in place and that the finger screws are properly tightened.

### 3.6.3 Flora side rail

The Flora side rail is only compatible with Camila bed ends.

#### Installing the side rail


#### WARNING! Risk of entrapment / suffocation

The user can get trapped or fall out of the bed if the side rail is not correctly mounted or damaged.

- Make sure that all gliding shoes are correctly guided into the bed end guideways.
- Make sure that the locking screw and tab are properly installed at the bottom of the guideways.

#### Assembling the gliding system

1.


Push the pawl Ⓐ from the backside, through the glider Ⓑ and make sure it locks with an audible click.

2. Do the same for all four pawls in each glider.


#### Attaching the gliding system and installing the side rail

1. Raise the bed to 1/3 of full height.
- 2.

Loosen the screw and tab at the bottom of the guideways on one of the bed ends and push the glider halfway up into the guideway until it attaches with an audible click.


3.

Push the upper side rail onto the two upper pawls.


4.

In the other end of the side rail bar; push the two upper pawls into the end of the bar.

5.


Push the attachment into the guideway until it attaches with an audible click (=properly engaged in the lower position).

6. Mount the lower side rail bar in both ends as in step 4.

7. Push the side rail in both ends, all the way up until the upper bar locks in the upper position.

8.

Retighten the screw and tab under the lower side rail bar, in each corner of the bed ends.

### 3.7 Wiring

#### **IMPORTANT!**

Not correctly drawn wires can be damaged when operating the sections with the motors.

- All wires must be kept clear of the floor and must not block the castors.
- All wires (motors, hand control and mains) must be connected to the control box
- The locking cam must be attached over the plugs in the control box.

Connect the mains wire plug to the 230 V socket.

## 4 Operating the bed

### 4.1 General Safety Information

**WARNING!**

**Risk of personal injury and damage to property.**

- The bed must be placed so that the height adjustment is not obstructed by, for example, lifts or furniture.
- Take care that no body parts are being squeezed between fixed parts (such as side rails, bed ends etc) and moving parts.
- The hand control must not be used by children.


**WARNING!**

**Risk of injury**

If the bed is used by confused or restless users or users with spasms:

- either, when the bed is equipped with a lockable hand control, lock the hand control functions,
- or make sure that the hand control is out of reach for the user.

### 4.2 Operating the hand control unit


- ▲ Move up (press left side of corresponding button)
- ▼ Move down (press right side of corresponding button)
- Ⓐ Backrest adjustment
- Ⓑ Thigh section adjustment
- Ⓒ Sitting position adjustment

1. Use the buttons on the hand control to operate the electrical functions of the bed.


### 4.3 Operating the bed functions manually

If the bed is configured for manual operation, it's not equipped with a hand control unit, but with a crank, positioned at the foot end of the bed.


The two-section bed has one crank to operate the head section. The four-section bed has one crank to operate the head section and one crank to operate the thigh section.

- Ⓐ for lower half of the mattress support.
- Ⓑ for upper half of the mattress support

**UP** Turn the crank handle clockwise

**DOWN** Turn the crank handle counterclockwise

**WARNING!**

- Do not use the crank if an extra weight, apart from the user, is placed on the bed

**IMPORTANT**

- The handle has to be locked/secured, when the crank is not used.


#### Lock/secure the crank

**IMPORTANT!**

- Twisting or rotating the white sleeve will damage its internal components and the locking function will cease to work.
- Only push the white sleeve inwards without twisting.

- 1.

2. Push the white sleeve Ⓒ inwards
2. Bend the crank by the knuckle Ⓓ and move it to the right.

## 4.4 Operating side rails

### WARNING!

#### Risk of entrapment or suffocation

There's a risk of entrapment or suffocation between mattress support, side rail and bed end.  
– Always ensure correct fitting and tightening of the side rails.

### WARNING!

#### Risk of falling

– Never leave the user unattended in the bed with the side rail down.  
– Make sure the side rail is in its highest position and properly locked while leaving the user unattended.

### CAUTION!

If the side rail is not locked properly it can fall down.

– Pull/push the top bar of the side rail to ensure that the locking system is properly engaged.

### CAUTION!

#### Risk of Injury

There is a risk of entrapment or squeezing while operating the side rail.

– Pay attention to your fingers and body parts of the patient.  
– Never force or drop the side rail while handling it.

### Removable 3/4 length side rails


### WARNING!

#### Risk of entrapment

There's a risk of entrapment or suffocation between mattress support, side rail and bed end.  
– When using removable side rails, always make sure that the distance between the bed end and the side rail's handle upper edge is less than 6 cm in the head end and more than 32 cm in the foot end of the bed.

#### 4.4.1 Verso II side rails

#### Folding down to open the side rail

1. Hold the top bar **B** with one hand and pull the locking button **A** with the other hand.
2. Pull the top bar sideways and release the locking button.
3. Fold the side rail down to its lowest position.

#### Folding up to close the side rail

1. Pull the top bar **B** and raise the side rail until it is properly engaged in the locking system.

### CAUTION!

If the side rail is not properly locked it can fall down.

– Pull/push the top bar of the side rail to ensure that the locking system is properly engaged.

#### 4.4.2 Diana side rail

#### Operating the side rail


Always stand in the middle of the side rail, when moving it up and down, to ensure that both sides move evenly. If not, the side rail can get stuck. To unlock it, pull/push it slightly up-/downwards to level the side rail again.

#### Move the side rail up


1. Lift the side rail up until the locking pin is securely pushed into a housing with an audible click.
2. Make sure the locking pin is securely pushed into a housing.
3. Try moving the side rail up-/downwards to make sure it is properly locked.

#### Move the side rail down

- 1.

Pull out the locking pins **A** on both sides and push the side rail down, slowly.

### WARNING!

#### Risk of squeezing

There is a risk of getting trapped between the extended tubes on the side rail **B** and the floor or getting squeezed while moving the side rail up and/or down.

– Be careful while lowering the side rail.

#### 4.4.3 Flora side rail

The Flora side rail is only compatible with Camila bed ends.


## Operating the side rail

### Raise the side rail

1. Use the two handles on the upper bar to raise the side rail until it locks in the upper position with an audible click.

### Lower the side rail

- 1.

Press the gliding system latch **A** on one end of the bed and lower the side rail to the lowest position.

2. Repeat the procedure at the other end of the bed.

## 4.5 Castors and brakes

### CAUTION!

#### Risk of trapping/squeezing

All brakes are foot-operated.  
– Do not release the brake with the fingers.

### CAUTION!

#### Risk of injury


User can fall, while entering or exiting the bed, if brakes are not locked.

- Always lock the brakes before the user is moving in or out of bed or when attending to the user.
- At least one castor at the head end and one castor in the foot end must be locked.

Castors might make marks on different types of absorbing floor surfaces - such as untreated or badly treated floors. To prevent marks, *Invacare®* recommends to place a suitable kind of protection between the castors and the floor.

### 4.5.1 Castor brake

**Locking the brake** - Step on the outer part **B** of the brake pedal

**Unlocking the brake** - Step on the inner part **A** of the brake pedal

## 4.6 Emergency release of a mattress support section

In case of a power- or motor failure, an emergency release of the back, thigh or leg section could be necessary. An emergency release of the height adjustment is NOT possible.

### CAUTION!

#### Risk of injury


- A minimum of two persons is required for an emergency release of a mattress support section.
- When releasing a mattress support section, it might lower fast. Do not reach under the mattress support while lowering it.

### IMPORTANT!

- Before an emergency release of the mattress support, remove the plug from the mains socket.

1. Both persons hold the mattress section.
2. One of them locates the motor in question and pulls out the safety pin.

3. Both persons slowly lower the mattress section until it is completely down.

## 4.7 Mattress retainers

### WARNING!

#### Risk of injury

If the mattress retainers are removed the mattress will not follow the movement, when adjusting the head-/legsections. The mattress can slip sideways and cause the user to fall out of the bed or get trapped in the bed.

- Always use mattress retainers and be careful to put them back after any adjustments.

## 4.8 Lifting pole

### WARNING!

#### Risk of injury

The bed can tip if the handle is used, while the lifting pole is turned away from the bed.

- The lifting pole always has to be positioned with the handle hanging over the bed area.
- Do not exceed the maximum load of the lifting pole; 80 kg.

### 4.8.1 Place the lifting pole

The lifting pole can be placed either on the left or the right side of the head end of the bed.

1. Insert the lifting pole into the lifting pole tube and fix it with the finger screw.


It's not required to tighten the finger screw, in case you want the lifting pole to swing away to the side of the bed.


### 4.8.2 Adjusting the handle height

The handle height should always be adjusted to the user's need.

1.

2.

1. Hold the cord out while sliding the plastic cord lock (A) upwards or downwards until the handle reaches the preferred height.
2. Lock the handle by pressing back the cord in the lock (A) and pull the handle downwards.

#### IMPORTANT!

After adjusting the handle height:

- Check that the two cords above the cord lock are parallel and inside the cord lock.
- Make sure the cord is properly locked by pulling the handle hard.

---

## 5 Accessories

---

### 5.1 List of available accessories

- Side rails
  - Side rail Verso II  
3/4 length collapsible steel side rail (156 cm x 40 cm), mounted on the side tube. Release button in the foot section
  - Side rail Diana  
3/4 length lowering side rail (153cm x 37 cm). Mounted on the side tube.
  - Side rail Flora  
Full length wooden side rail, in cherry.

- Side rail cover for Verso II
- Lifting pole
- Holder for hand control
- Drip rod bar
- Electrical up-grade kit for the mattress support

**Important**

- – Use only original accessories and spare parts. Spare parts lists and extra user manuals can be ordered from *Invacare* or via homepage.

## 6 Maintenance

### 6.1 Maintenance

A service contract can be made in the countries, where *Invacare*® has its own sales company. In certain countries *Invacare*® offers courses in service and maintenance of the bed. Spare parts lists and additional user manuals are available from *Invacare*®.

With normal daily use, service must be carried out according to the check list after two years use and thereafter every second year.

- !** **IMPORTANT!**
- The mattress support must be supported during service inspections, to prevent accidental lowering.
  - Only personnel who have received the necessary instructions or training may perform service and maintenance of the bed.
  - After reconditioning the bed, or if bed functions change, service must be carried out according to the check list.

#### Before Use

- Ensure that all manual and electrical parts functions correctly and are in a secure state.

#### After three months

- Ensure that all manual and electrical parts are functioning, and tighten bolts, screws, nuts, etc.

#### Every year

- We recommend a safety test comprising the motors' performance and mechanical state.

#### Every second year

- Service should be performed according to the check list.

 **Motors, hand control and control units** are serviced by exchanging the faulty part.

#### 6.1.1 Checklist maintenance

##### Checkpoints

- Visual inspection of all parts of the bed (no deformations)
- Circlips, cotter pins and plastic fixing ring - properly locked and intact.
- Screws - tightened.
- Weldings - intact (no wear and tear)
- Castors (rolling smoothly)
- Castor fittings - tightened.
- Castor brakes - locking properly.
- All motors - running properly (with regular speed and at low noise).
- Wires - correctly wired and undamaged.
- Electric plugs - undamaged.
- Side rail locking and moving system - properly locking and running smoothly.
- Accessories - correctly assembled and correct function.
- Damaged coating - repaired.
- Leakage of grease - check for grease.

#### CAUTION!

##### Risk of injury or damage to property

Grease leaks can lead to accidental falls or fire.

- Contact service provider if leaks of grease or other residue is detected.

### 6.2 Cleaning and disinfection

- !** **IMPORTANT!**
- Wrong fluids or methods can harm or damage your product.
- Follow instructions carefully for either Non-washable or Washable components.
  - Never use corrosive fluids (alkalines, acid, cellulose thinner, acetone etc)
  - Never use a solvent that changes the structure of the plastic or dissolves the attached labels.
  - Always make sure that the bed is carefully dried before taking it into use again.

#### 6.2.1 Cleaning methods

##### Electrical components

- !** **IMPORTANT!**
- Non-washable electronics can not withstand high temperatures.
- Do not wash or dry in higher temperatures than 40 °C.

Method: Wipe off with a wet cloth or soft brush.

Max. temperature: 40 °C

Solvent/chemicals: Water

##### Metal components

Method: Wipe off with a wet cloth or soft brush. Water may be pressurized, but not high pressure or steam.

Max. temperature: 40 °C

Solvent/chemicals: Household detergent or soap and water, 6-8 pH

##### Wood (including textile straps on side rails, if existing)

Method: Wipe off with a wet cloth or soft brush.

Max. temperature: 40 °C

Solvent/chemicals: Household detergent or soap and water, 6-8 pH

##### Textiles (including upholstery and mattresses)

See attached label on each product.

### 6.3 Lubrication

#### Lubrication plan

We recommend lubricating the bed according to the following instructions:

- Points of rotation in mattress support and base frame - lubricate with medically clean oil.
- Motor attachment points to mattress platform - lubricate with medically clean oil.

- All of the motors' tension rod - lubricate with medically clean oil.
- All motor bearings - lubricate with medically clean oil.

The wooden side rails gliding system must not be lubricated with oil - otherwise the wooden bars will move sluggishly.

## 7 After Use

### 7.1 Waste disposal

 Waste disposal/recycle must comply with the laws and regulations for waste handling in each country.

Invacare® is continuously working towards ensuring that the company's impact on environment, locally and globally, is reduced to a minimum.

We comply with the current environment legislation (e.g. WEEE and RoHS directive).

We only use REACH compliant materials and components.

- All wooden parts must be dismantled and sent for incineration.

- All electric parts must be dismantled and be disposed of as electric components.
- Plastic parts must be sent for incineration or recycling.
- Steel parts and castors must be disposed of as waste metals.

#### **IMPORTANT!**

- Accumulator back-up
  - Old accumulators must be returned to *Invacare®* or recycled as car batteries

### 7.2 Reconditioning

This product is suitable for reuse. To recondition the product for a new user, carry out the following actions:

- Inspection according to service plan
- Cleaning and disinfection

For detailed information about inspection, cleaning and disinfection, refer to 6 Maintenance, page 19.

## 8 Troubleshooting

### 8.1 Troubleshooting electrical system

Symptom	Possible cause	Remedy
Mains indicator does not light up	Mains are not connected	Connect mains
	Fuse in the control unit is blown	* Replace the control unit
	Control unit is defective	* Replace the control unit
Mains indicator lights up, but the motor is not running. The relay in the control unit makes a clicking noise.	Motor plug is not fully inserted into the control unit.	Insert the motor plug properly into the control unit
	The motor is defective.	* Replace the motor
	Motor cable is damaged.	* Replace the cable
	Control unit is defective	* Replace the control unit
Mains indicator lights up, but the motor is not running. No relay sound is heard from the control unit.	Control unit is defective	* Replace the control unit
	Hand control is defective	* Replace the hand control
Control unit is in order except for one direction on one channel.	Control unit is defective	* Replace the control unit
	Hand control is defective	* Replace the hand control
Motor is running, but the piston rod does not move.	Motor is damaged	* Replace the motor
The motor cannot lift full load.		
Motor noise, but no movement of piston rod.		
Piston rod operates inwards and not outwards.		

\* Service and maintenance of the bed must only be performed by personnel who have received the necessary instruction or training.

**Risk of personal injury and damage to the product.**

– The bed must be unplugged from the main power source before opening or repairing electrical parts.

**WARNING!**

**Risk of personal injury and damage to the product.**


– The bed must be unplugged from the main power source before opening or repairing electrical parts.

## 9 Technical data

### 9.1 Dimensions

All measurements are stated in mm. All angles are stated in degrees. All measurements and angles are stated without tolerances.

Invacare® reserves the right to change the stated measurements and angles.

A (outer)	890 – 1025 mm
B (inner)	2000 mm
C (outer)	2080 – 2130 mm
D	800 mm
E	230 mm
F	300 mm
G	610 mm
H	1230 mm
I	365 – 635 mm
J	0 – 72°
K	0 – 32°

### 9.2 Mattress dimensions

Allowed mattress sizes, depending on the side rail  
(Two-sectioned and four-sectioned bed)

Side rail:	Mattress measures (in cm):			
	Min. height	Max. height	Min. width	Min. length
Diana	13	15	85	190
Flora	12	17	85	190
Verso II	8	18	85	190

Min. density: 38kg/m<sup>3</sup>

### 9.3 Weights

	Max. patient weight (provided that the weight of the mattress and the accessories do not exceed 35 kg)	135 kg
	Max. safe working load (patient + accessories)	170 kg

Sonata standard bed, complete, excluding accessories	
Mattress support, head section four sections (2 pcs.)	15,5 kg
Mattress support, leg section four sections (2 pcs.)	14.0 kg
Mattress support, four sections (1 pc)	28,5 kg
Mattress support, four sections (1 pc)	22,5 kg
Side rail Verso II (1 pc.)	7.0 kg
Side rail Diana (1 pc.)	6.0 kg
Side rail Flora (1 pc.)	8,5 kg
Bed end Valeria (1 pc.)	4.9 kg
Bed end Camila (1 pc.)	11,0 kg
Bed end Sabrina (1 pc.)	6,8 kg
Bed end Nadine (1 pc.)	5,0 kg
Bed end Vibeke (1 pc.)	6,5 kg
Lifting pole	5,0 kg

### 9.4 Environmental conditions

	Storage and transportation	Operation
Temperature	-10°C to +50°C	+5°C to +40°C
Relative humidity	20% to 75%	
Atmospheric pressure	800 hPa to 1060 hPa	

Be aware that when a bed has been stored under low temperatures, it must be adjusted to operating conditions before use.

### 9.5 Electrical system

Voltage supply: U <sub>in</sub> = 230 Volt, AC, 50/60 Hz (AC = Alternating current)
Maximum current input: I <sub>in</sub> max. = 1 A
Intermittent (periodic motor operation): I <sub>nt</sub> = max. 10%, 2 min / 18 min
Degree of protection: IPX6*
Insulation class: Class II

## Type B Applied Part

Applied Part complying with the specified requirements for protection against electrical shock according to IEC60601-1.

(An applied parts is a part of the medical equipment which is designed to come into physical contact with the patient or parts that are likely to be brought into contact with the patient.)

Sound level: 45 to 50 dB (A)

\* IPX6 classification means that the electrical system is protected against water projected in powerful jets from any direction.

The bed is not provided with a mains switch, so the mains plug is the only separation from the mains.

## 9.6 Electromagnetic compliance (EMC)

### Guidance and manufacturer's declaration – electromagnetic emission

The medical bed is intended for use in the electromagnetic environment specified below. The customer or the user of the bed should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11(partly)	Group I	The medical bed uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11(partly)	Class B	The medical bed is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

### Guidance and manufacturer's declaration – electromagnetic immunity

The medical bed is intended for use in the electromagnetic environment specified below. The customer or the user of the bed should assure that it is used in such an environment.

Immunity test	IEC 60601–1–2 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% $U_T$ (>95% dip in $U_T$ ) for 0,5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles < 5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	< 5% $U_T$ (>95% dip in $U_T$ ) for 0,5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles < 5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the medical bed requires continued operation during power mains interruptions, it is recommended that the medical bed be powered from an un-interruptible power supply or a battery.  $U_T$ is the a. c. mains voltage prior to application of the test level.

Immunity test	IEC 60601–1–2 test level	Compliance level	Electromagnetic environment – guidance
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 V  3 V/m	3 V  3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the medical bed including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> $d = \left[ \frac{3,5}{V_1} \right] \sqrt{P}$ $d = \left[ \frac{3,5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[ \frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in metres (m).<sup>b</sup></p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the medical bed is used exceeds the applicable RF compliance level above, the medical bed should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the medical bed.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

At 80 MHz and 800 MHz, the higher frequency range applies.

**Recommended separation distances between portable and mobile RF communications equipment and the medical bed**

The medical bed is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the medical bed can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the medical bed as recommended below, according to the maximum output power of the communications equipment

Rated maximum output of transmitter [W]	Separation distance according to frequency of transmitter [m]		
	150 kHz to 80 MHz $d = \left[ \frac{3,5}{V_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[ \frac{3,5}{E_1} \right] \sqrt{P}$	800 MHz to 2,5 GHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33

10	3.69	3.69	7.38
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above the recommended separation, distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Invacare Limited