Invacare® LiNX

DLX-REM400

en	Joystick Gaiter	
	Service Manual	1
de	Joystick-Manschette	
	Servicehandbuch	5
es	Pieza de joystick	
	Manual de servicio	10
fr	Soufflet d'étanchéité du joystick	
	Manuel de maintenance	15
it	Cuffia per joystick	
	Manuale per la manutenzione	19
pt	Fole do joystick	
	Manual de assistência	24
sv	Joystickdamask	
	Service manual	28

1 Before You Start en

Distribution

These instructions are intended for use by experienced service technicians only. Gaiters and associated parts are not supplied for retail sale or installation by end users.

Warranty

This procedure is intended for the replacement of gaiters on units not covered by warranty. Undertaking this procedure on a product within warranty will void the warranty. All units under warranty should be returned to an authorised Dynamic Controls service centre for repair, refund or replacement under warranty.

Liability

Dynamic Controls accepts no responsibility for any liability, whether expressed or implied, resulting from use of unauthorized parts or servicing that deviates from this procedure.

Limitations

This procedure covers the replacement of gaiters only. If any secondary fault, including, but not limited to, the presence of internal contamination, is discovered during the process of the gaiter replacement, the unit must be withdrawn from service and either exchanged or returned to a Dynamic Controls service centre for further analysis and where appropriate repair.

Disclaimer

The remote modules covered under this procedure are supplied with tamper evident warranty void seals fitted. If a remote module is received for gaiter replacement without the seal in an intact condition, Dynamic Controls recommends that the unit be returned 'as is' without

undertaking any service or repair. This recommendation is intended to safeguard the healthcare professional or service technician from the consequences of damage or defects resulting from user tampering, whether intentional or inadvertent. Following repairs, the technician undertaking gaiter replacement under this procedure may choose to apply their own tamper evident labels, to protect against subsequent opening or tampering of the unit post-service.

ESD precautions

Precautions must be taken to protect the remote from damage by ESD, once the case is opened. The following procedure should be performed in an ESD-protected area using appropriate protection devices, such as ESD-safe matting, grounding cords and wrist straps.

2 Equipment

Parts List

 Different variants of the gaiter and gaiter retainer are available on the market, which differ from each other. Ensure that you use the enclosed part for repair only!

Quantity	Description
1	Gaiter (required)
1	DOWSIL (formerly Dow Corning) DC4 Electrical Insulating Compound (required)

 DC4 electrical insulating compound can be sourced from: Element14 / Farnell — Order Code 537019

DEALER: Keep this manual.
The procedures in this manual MUST be performed by a qualified technician.

Yes, you can.®

Tools required

- 1 x TORX screwdriver bit, size T8
- 1 x TORX screwdriver bit, size T10
- 1 x torque driver capable of taking the T8 and T10 TORX bits and reaching at least 120 cNm torque
- Cotton buds / swabs

3 Initial Checks

Before removing, the gaiter should be checked carefully to confirm the extent of wear or damage. If the gaiter is damaged to an extent that may have allowed foreign objects or significant ingress to pass into the joystick housing, the remote module should be returned for thorough internal assessment by a Dynamic Controls service centre.

If the remote module shows any signs of significant mechanical damage or impact, the unit should be returned to a Dynamic Controls service centre for assessment. Significant damage would typically include loss of function or openings in the unit chassis.

4 Disassembly

WARNING!

– Do not disassemble beyond the point described in this procedure. Further disassembly may result in product calibration moving out of specification, compromising the product’s safe operation.

1.


Fig. 4-1 Remove joystick knob

Remove the joystick knob: pull it directly away from the remote module's body but **DO NOT** twist it. Some force may be required if the knob is tight.

WARNING!

– Twisting the knob when it is still attached to the shaft may result in damage to the joystick mechanism.

2.


Fig. 4-2 Remove the decorative cover

Carefully remove the square, decorative cover **A** by inserting a thin-bladed tool (e.g. a screwdriver) into the slot **B** and levering the cover up and over the joystick shaft.

3.


Fig. 4-3 Remove four screws

Using the T10 driver, remove the four screws **A** that were exposed after removing the decorative cover in step 2. The screws can be removed in any order. One screw is located underneath the warranty label **B**.

4.


Fig. 4-4 Remove rear cover

Remove the rear cover (A) by inserting and sliding a thin object, or fingernail, along one of the cover's edges, whilst gently prising it away from the remote module's body. The areas marked (B) are a good place to start.

5.


Fig. 4-5 Remove screws (A) and (B)

Using the T10 driver, remove the four screws (A). These screws can be removed in any order. Using the T8 driver, remove the four screws (B). These screws can be removed in any order.

6.


Fig. 4-6 Separate upper and lower case and disconnect connector

Separate the upper and lower cases (A), (B). Disconnect the interconnecting internal looms from connectors (C) and (D).

7.


Fig. 4-7 Remove the joystick sub-assembly and gaiter from the chassis

Carefully remove the joystick sub-assembly and gaiter (A) from the chassis (B).

8. Remove and discard the gaiter. Note that the gaiter may be attached to the joystick sub-assembly or the chassis.

5 Component Inspection

Contamination

With the gaiter and retainer removed, inspect the visible joystick mechanism for signs of ingress, including signs of liquid residue and solid debris.

Deflect the joystick shaft in all directions so that the inside of the joystick mechanism can also be checked for signs of debris or ingress. If there is any sign of contamination, the complete remote module should be replaced or, where appropriate, returned for service by a Dynamic Controls service centre.

WARNING!

– Unauthorised cleaning of the joystick mechanism or application of unsuitable lubricants may adversely affect the performance of the joystick mechanism, creating a safety risk for the user.

Spring Retainer And Shaft


Fig. 5-1 Spring retainer and shaft

While the gaiter is removed, check the condition of the spring retainer (A) and shaft. Check for signs of physical damage, including deformed shafts and cracks in the spring retainer. Ensure the spring retainer's locating pin (B) is secure and not protruding from the retainer. Note that the mould-parting line on the spring retainer is normal and not a sign of damage.

If any part is damaged, loose or missing, the complete remote module should be replaced or, where appropriate, returned for service by a Dynamic Controls service centre.

6 Reassembly

1.


Fig. 6-1 Apply electrical insulating compound

Using a cotton bud, apply a thin coating¹ of DC4 Electrical Insulating Compound to the middle and lower layers of the gaiter, as highlighted with the grid in figure above. Ensure that the grease is smeared evenly all the way around the highlighted areas and not placed in just one spot.

2.


Fig. 6-2 Push the gaiter through the chassis

From the outer side of the upper case, push the gaiter (A) through the chassis (B), whilst rotating the gaiter within the chassis to ensure that it can move freely.

3.


Fig. 6-3 Push the joystick assembly through the chassis and gaiter

From the inner side of the upper case, push the joystick sub-assembly (A) through the chassis (B) and gaiter (C), ensuring that the gaiter sits down fully and does not have any folds or wrinkles. Check that it sits correctly on the spring retainer as shown (D).


Fig. 6-4 Top of gaiter sits below the top shoulder of the spring retainer

Ensure that the top of the gaiter sits below the top shoulder of the spring retainer (6-4 Top of gaiter sits below the top shoulder of the spring retainer, page 4, right), not on top of it (6-4 Top of gaiter sits below the top shoulder of the spring retainer, page 4, left).

4.


Fig. 6-5 Reconnect connector and close upper and lower case

Reconnect the interconnecting internal looms to connectors (C) and (D). Check that the speaker (E) is still fitted and correctly positioned. If it has been dislodged, refit it with its connectors (F) facing outwards and towards the left-hand side of the lower case when viewing from the inside. Close the upper and lower cases (A), (B).

1. The amount of insulating compound that you use is not critical. Use just enough so that the surfaces, as indicated in the images, are lightly covered.

5.


Fig. 6-6 Replace and tighten screws **A**

Using the T10 driver, replace and tighten the four screws **A** to 120 cNm, the screws can be replaced in any order. For security, dealers may wish to attach their own anti-tamper labels so that they can determine if the unit is subsequently opened.

6.


Fig. 6-7 Replace and tighten screws **A** and **B**

Using the T10 driver, replace and tighten the four screws **A** to 50 cNm, the screws can be tightened in any order. Using the T8 driver, replace and tighten the four screws **B** to 40 cNm, the screws can be tightened in any order.

7.


Fig. 6-8 Re-attach the joystick decorative cover and back cover

Re-attach the joystick decorative cover **A** and back cover **B**, ensuring all clips are fully engaged.

8.


Fig. 6-9 Replace the joystick

Replace the joystick knob **A**. If necessary, and before placing it on the shaft, rotate the knob to align the ribs **B** with the notches **C**. Push the knob firmly towards the body of the remote module – when fitted correctly, the knob should not be able to be rotated about the shaft.

7 Testing

Visually inspect the gaiter's position, ensuring that it sits down evenly, with no folds or wrinkles in it.

Test the joystick movement by moving the joystick in a circular motion about its full extent. It must move freely in all directions and consistently return to the centre, neutral position when released from any position.

Connect to a system. Check for any out-of-neutral at power up or other faults.

Use the remote module to drive a wheelchair, ensuring that the wheelchair drives and stops correctly, the joystick returns to neutral and that, after a short while, the park brakes engage (usually indicated by an audible click).

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Invacare Ireland Ltd,
Invacare AB
Invacare Deutschland GmbH