Aquatec® KOGIA

en	Bath lift Service Manual	3
de	Badewannenlifter Servicehandbuch	9

PROVIDER: Keep this manual.
The procedures in this manual MUST be performed by a qualified technician.

AQUATEC®

Yes, you can.®

©2022 Invacare® Corporation

All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

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Contents

1 General Information	4
1.1 Safety Advice	4
1.2 Safety and Fixation Instructions	4
2 Product Safety	5
2.1 Periodic Checks	5
2.1.1 Visual Check	5
2.1.2 Function Check	5
2.1.3 Checking Electrical Safety	5
2.1.4 Check List	5
2.2 Servicing	5
3 Reconditioning	6
3.1 Cleaning and Disinfection	6
4 Repair	7
4.1 Spare Parts	7
4.2 Replacing Side Flaps	7
4.3 Replacing Hand Control	7
4.4 Replacing Suction Cups	8

1 General Information

This service manual contains the technical information for maintaining product safety, reconditioning and repairing Aquatec® bath lifts.

1.1 Safety Advice

WARNING!

Danger of accident and severe injuries.

If a bath lift is faulty it can cause accidents and severe injuries. The product requires regular checks and servicing is required for this reason.

IMPORTANT!

- The product must be cleaned, disinfected, the suction cups must be replaced and a check (see *2.1 Periodic Checks, page 5*) is required before re-use.
- Checks and servicing, reconditioning and repairs for bath lifts may only be carried out by authorized specialist dealers.

1.2 Safety and Fixation Instructions

These safety instructions are intended to help avoid accidents during work and must be observed under all circumstances.

All employees coming into contact with contaminated products must regularly consult a company doctor. Work clothing and personal protective equipment must be available in necessary quantities and be in proper condition. Reliable hand and surface disinfection must be ensured.

Prior to all Inspections or Repairs

- Read this repair manual and the associated user manual carefully and follow the instructions.

Personal Safety Equipment

Safety Shoes

The bath lift and some of its components are very heavy. If these parts are dropped, they can cause injury to the feet.

- Wear safety shoes at all times when working.

General Safety Information and Installation and Removal Information

CAUTION!

Risk of crushing

Different components, such as the drive unit, seat, etc. are very heavy. They can therefore pose a risk of injury for the hands.

- Please note the heavy weight of certain components. This particularly applies to the removal of the drive unit, seat or backrest.

CAUTION!

There is a risk of injury and risk of damage to the bath lift if service work is carried out incorrectly or incompletely

- Only use undamaged tools in good condition.
- After the conclusion of the work and before restarting the bath lift, check that all connections are secure.
- After the conclusion of the work and before restarting the bath lift, check that all parts are correctly engaged.
- When the work is complete, always carry out a trial run.

2 Product Safety

2.1 Periodic Checks

IMPORTANT!

The product requires regular checks to ensure product safety. The required intervals are to be observed.

WARNING!

- It is recommended to check the bath lift after 24 months at the latest (in accordance with EN ISO 10535) and the suction cups must be replaced.
- Check the electrical safety of the product in line with national requirements.

2.1.1 Visual Check

1. Check the seat plate, backrest plate and base plate for damage and ensure they are complete.
2. Check the shear arms and glide surfaces by lifting them and pushing them together again.
3. Check the bearings of the seat plate / shear arms and base plate / shear arms.
4. Examine the hand control and battery box for any damage. Replace the hand control and battery box if the housing is damaged or cracked.
5. Examine the plug-in charger for any damage.

2.1.2 Function Check

1. Use a voltmeter to check the charger for open-circuit voltage of 14 ± 0.2 V at the charging poles of the jack plug.
2. Raise and lower the bath lift under maximum load to check its functions and listen for any unusual noises. If the battery is fully charged, the product can be lifted up to 22 times under nominal load. The age of the battery can decrease the number of lifts.

2.1.3 Checking Electrical Safety

IMPORTANT!

Regular electrical checks according to EN 62353 (VDE0751) must be carried out in line with national requirements.

1. Inspect the charger together with the battery for damage to the casing or cable, test functions and check for stray current.
2. Check the drive for damage to the casing or cable.

2.1.4 Check List

General Inspection	☺	☹
Is the device in good condition and is it complete (device and accessories)?	<input type="checkbox"/>	<input type="checkbox"/>
Is the identification label easily readable and is it securely mounted on the product?	<input type="checkbox"/>	<input type="checkbox"/>
Have subsequent changes been made to the product?	<input type="checkbox"/>	<input type="checkbox"/>
Is the device free from damage or weaknesses of any kind? Does the device operate correctly under nominal load?	<input type="checkbox"/>	<input type="checkbox"/>
Is the device fully functional in accordance with the user manual? Is the product accompanied by a user manual?	<input type="checkbox"/>	<input type="checkbox"/>
Is all safety equipment working as intended?	<input type="checkbox"/>	<input type="checkbox"/>
Is the device electrically safe, where relevant?	<input type="checkbox"/>	<input type="checkbox"/>
Are all safety-related notices sufficiently legible? Is there an identification label and is this securely attached to the supporting surface?	<input type="checkbox"/>	<input type="checkbox"/>
Eradication of Faults	☺	☹
Have all the faults found been eliminated and have faulty components been replaced?	<input type="checkbox"/>	<input type="checkbox"/>
Does the device now operate correctly?	<input type="checkbox"/>	<input type="checkbox"/>
Is all safety equipment now working as intended?	<input type="checkbox"/>	<input type="checkbox"/>
Are all the necessary safety-related notices now in place?	<input type="checkbox"/>	<input type="checkbox"/>
Completion of Checks	☺	☹
Has the device been cleaned and disinfected?	<input type="checkbox"/>	<input type="checkbox"/>
Is the device technically and functionally safe?	<input type="checkbox"/>	<input type="checkbox"/>
Is the device electrically safe, where relevant?	<input type="checkbox"/>	<input type="checkbox"/>
Is the identification label easily readable and is it securely mounted on the product?	<input type="checkbox"/>	<input type="checkbox"/>
Has the device been packed correctly?	<input type="checkbox"/>	<input type="checkbox"/>

2.2 Servicing

If the cleaning and safety notices are observed, the product will be maintenance free.

3 Reconditioning

Before the conditioning decision is made, it should also be verified beyond a critical evaluation of objective feasibility that this is sensible from an economic and ecological perspective, taking into consideration the risk associated with the conditioning and use of the medical device and the expense and quality assurance, or whether it would be better to dispose of the medical device.

The reconditioning of the product encompasses the following points:

- 3.1 *Cleaning and Disinfection, page 6*
- 2.1 *Periodic Checks, page 5*
- 2.2 *Servicing, page 5*

IMPORTANT!

- If possible, continue to use the old identification label; if this is not possible, the new identification label must contain the same information and the old serial number. (Replacement of spare parts with serial numbers).
- If screws with thread locker are loosened, these must be replaced with new ones.

3.1 Cleaning and Disinfection

Manual cleaning and disinfection procedures must always be conducted in accordance with documented work instructions and with means and processes checked for effectiveness and tailored to the medical device. For mechanical cleaning and disinfection procedures it can be procedurally ensured that the parameters necessary to attain quantifiable cleaning and disinfection performance are complied with. In order to guarantee the continuous quality of manual and/or mechanical conditioning procedures it must be proven by at least periodical checks that the procedure is still effective and no unintentional changes have arisen.

Conditioning ends with the documented approval for use of the medical device. The chain of actions must be optimized as weaknesses in one of the necessary individual steps threaten the overall success. Transport and storage must not have a negative impact on the characteristics of the conditioned medical device. During transportation, it must be ensured that contamination of a conditioned or new medical device is precluded until its use by means of strict separation of clean/unclean devices.

All detergents and disinfectants used must be effective, compatible with one another and be gentle on the

materials used. Relevant cleaning and disinfection baths must be regularly replaced at defined intervals in order to guarantee their effectiveness.

With respect to contact time and concentration, the disinfectants used must comply with the list of disinfectants issued by 'The National Institute for Clinical Excellence' guidelines on Infection Control www.nice.org.uk/CG139 and your local infection control policy.

Detergents and disinfectants must be selected in accordance with the equipment and facilities available in each instance in order to guarantee a reliable cleaning result.

Exchanged products and product components must be disposed of in accordance with their respective material composition.

Preparation for Cleaning

1. Remove battery box.
2. Detach the backrest from the seat plate.
3. Cover electronic connections with protective caps to ensure that water cannot get in.

Cleaning

1. Remove any superficial contamination.
2. Treat all components of the bath lifter except battery box and charger with detergent and rinse down with clean water.
3. Replace the covers or wash them in the washing machine at a maximum temperature of 60 °C.

IMPORTANT!

Risk of damage

- Washing at high temperatures can cause shrinkage.
- Do not dry covers in the tumble dryer.

4. Use a universal detergent to clean the battery box, covers and charger by hand.
5. Allow all components of the bath lifter to dry in the ambient air or dry off with hand towels.
6. Check the suction cups and replace if necessary.

Disinfection

1. Disinfect the seat surface, the front of the backrest and the hand control by wiping them down with disinfectant.
2. Dry all components of the bath lifter in the ambient air or allow them to dry in the air until there is no more superficial moisture on the product.

4 Repair

Repairs require a high level of expertise. These assembly instructions therefore break down the various tasks into 3 categories:

Requirement	Symbol
Easy – technical understanding required	●○○
Medium – technical knowledge required	●●○
Difficult – technical knowledge and expertise in assembling bath lifts required	●●●

The required tools and their sizes are listed before the instructions.

Tool	Symbol
Phillips head screwdriver	⊕
Slot head screwdriver	●
Torx screwdriver	⊛
Snap ring pliers	○

WARNING!

– To avoid the risk of contamination, clean and disinfect the bath lift before carrying out repairs.

IMPORTANT!

– Repairs may only be carried out by an authorized dealer.
 – If possible, continue to use the old identification label; if this is not possible, the new identification label must contain the same information and the old serial number. (Replacement of spare parts with serial numbers).
 – When components are replaced it is necessary to ensure the traceability of the components replaced.
 – If screws with thread locker are loosened, these must be replaced with new ones.
 – If screws with circlip rings are loosened, these must be replaced with new ones.
 – Parts that become damaged during removal must be replaced with new ones.

4.1 Spare Parts

WARNING!

Original spare parts must be used in all repairs. Otherwise the warranty and product declaration of conformity shall be rendered invalid.

All spare parts must be obtained from the Invacare customer service department. An electronic spare parts catalogue can be found at www.invacare.eu.com.

4.2 Replacing Side Flaps

Difficulty: ●○○

Tool: —

Pos.	Designation	Article number	Quantity
Ⓐ	Side flap	1563646	2


Fig. 4-1

Removal

1. Press in the clips of side flap Ⓐ and slowly draw out the side flap.

Fitting

1. Slide the clips of side flap Ⓐ evenly into the slots of the seat plate Ⓑ until they engage.

4.3 Replacing Hand Control

Difficulty: ●●●

Tool: ⊕

Pos.	Designation	Article number	Quantity
Ⓐ	Hand control incl. sealing	1590349	1

Removal

1. Remove battery.
2. Remove backrest from seat and put backrest on worktop.
3. Remove screws of electronic cover.
4. Carefully lift up electronic cover from electronic box.
5. Remove sealing.
6. Disconnect and remove old hand control from circuit board.

Fitting


Fig. 4-2

1. Put black grommet Ⓐ of hand control into intended gap Ⓑ of electronic box.
2. Insert plug Ⓒ of new hand control into connector Ⓓ on circuit board.
3. Insert new sealing.

4. Replace electronic cover.

 Make sure that no cables get pinched.

5. Tighten screws of electronic cover with 1.2 Nm (\pm 0.5 Nm).

6. Set up backrest to seat.

7. Check all the functions of the bath lift.

WARNING!

- The hand control should only be replaced as a complete component.
- Do not open or attempt to repair the hand control.

Please contact Invacare customer service by phone or order a new hand control using the electronic spare parts catalog at www.invacare.eu.com.

4.4 Replacing Suction Cups

Difficulty: ●○○

Tool: —

Pos.	Designation	Article number	Quantity
Ⓐ	Pair of suction cups	15071	2

Fitting

1.


Fig. 4-3

Put rear suction cup with threaded bolt Ⓐ on sleeve nut Ⓑ and screw it in. Screw as far as the suction cup fits the orientation on the graphic.

2.


Fig. 4-4

Put front suction cup with threaded bolt Ⓐ on sleeve nut Ⓑ and screw it in. Screw as far as the suction cup fits the orientation on the graphic.

3. Insert toggle Ⓒ in corresponding holding fixture Ⓓ.

Removal

1. Removal takes place in reverse order.

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Invacare Poirier SAS
Invacare Austria GmbH
Invacare AG
Invacare Limited