Invacare® AVIVA™ RX10

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# Contents

This manual MUST be given to the user of the product.  
BEFORE using this product, this manual MUST be read and saved for future reference.

<b>1 General</b>	<b>73</b>
1.1 Introduction	73
1.2 Symbols in This Manual	73
1.3 Compliance	73
1.3.1 Product-Specific Standards	73
1.4 Usability	73
1.5 Warranty Information	74
1.6 Service Life	74
1.7 Limitation of Liability	74
<b>2 Safety</b>	<b>74</b>
2.1 General Safety Notes	74
2.2 Safety Information for the Electrical System	75
2.3 Safety Information for Electromagnetic Interference	76
2.4 Safety Information for Driving and Freewheel Mode	76
2.5 Safety Information for Care and Maintenance	77
2.6 Safety Information for Changes and Modifications to the Power Wheelchair	77
2.7 General Safety Notes for the Remote	78
<b>3 Product Overview</b>	<b>79</b>
3.1 Intended Use	79
3.1.1 Product Description	79
3.1.2 Intended User	79
3.1.3 Indications	79
3.2 Type Classification	79
3.3 Main Parts of Power Wheelchair	79
3.4 Main Parts of Remotes	79
3.4.1 The Status Indicator	79
3.4.2 Battery Gauge	79
3.4.3 User Interface DLX-REM060	80
3.4.4 User Interface DLX-REM150	80
3.4.5 User Interface DLX-REM210	80
3.4.6 User Interface DLX-REM215	80
3.5 Labels on the Power Wheelchair	81
3.6 Labels on the Remote	82
<b>4 Setup</b>	<b>83</b>
4.1 General Setup Information	83
4.2 Setting Up the Remote	83
4.2.1 Wiring	83
4.2.2 Connecting	84
4.3 Adjusting the Seat	84
4.3.1 Installing the Backrest	84
4.3.2 Adjusting Manual Backrest	84
4.3.3 Adjusting Seat Height and Manually Seat Tilt	84
4.3.4 Adjusting Seat Depth	85
4.4 Adjusting Armrest	85
4.4.1 Adjusting Armrest Height	85
4.4.2 Adjusting Armrest Depth	85
4.4.3 Adjusting Armrest Width	85
4.5 Adjusting Remote for Length of User's Arm	85
4.6 Adjusting the Headrest	85
4.6.1 Adjusting Position of Headrest	86
4.6.2 Adjusting Height of Headrest	86
4.7 Adjusting Standard 80° Legrest	86
4.7.1 Adjusting Width of Legrest	86
4.7.2 Adjusting Length of Legrest	86
4.7.3 Swivelling and/or Removing Legrest	86
<b>5 Using Power Wheelchair</b>	<b>86</b>
5.1 Getting in and out of Power Wheelchair	86
5.1.1 Removing Standard Armrest for Side Transfer	86
5.1.2 Swivelling Remote to Side	86
5.1.3 Getting in and out of Power Wheelchair	86
5.2 Before Driving	87
5.2.1 Operating the Remote	87
5.2.2 Activating the Drive Function	88
5.2.3 Operating Powered Seating Functions	88
5.2.4 Operating the Light System	88
5.2.5 Operating the Horn	88
5.2.6 Operating an Emergency Stop	89
5.2.7 The Sleep Mode	89
5.2.8 Locking / Unlocking the Remote	89
5.2.9 Audible Cues	89
5.3 Taking Obstacles	90
5.3.1 Maximum Obstacle Height	90
5.3.2 Correct Way to Take Obstacles	90
5.4 Driving up and down Gradients	90
5.5 Use on Public Roads	91
5.6 Parking and Stationary	91
5.7 Pushing the Power Wheelchair in Freewheel Mode	91
5.7.1 Disengaging / Engaging Motors	91
<b>6 Control System</b>	<b>91</b>
6.1 Control Protection System	91
6.1.1 Using Main Fuse	91
6.2 Batteries	91
6.2.1 Connecting / Disconnecting Batteries	92
6.2.2 General Instructions on Charging	92
6.2.3 Charging Batteries	92
6.2.4 Battery Alarms	93
6.2.5 Disconnecting Power Wheelchair After Charging	93
6.2.6 Storage and Maintenance	93
6.2.7 Instructions on Using Batteries	93
6.2.8 Cleaning Battery Terminals	94
6.2.9 Transporting Batteries	94
6.2.10 General Instructions on Handling Batteries	94
6.2.11 Handling Damaged Batteries Correctly	94
6.3 Battery Charger	94
6.3.1 Symbols on the Charger	94
6.3.2 Product Overview	95
6.3.3 Caution	95
6.3.4 Operational Instructions	95
6.3.5 Troubleshooting	95
6.3.6 Technical Specifications	95
<b>7 Transport</b>	<b>96</b>
7.1 Transport — General Information	96
7.2 Transferring Power Wheelchair to Vehicle	96
7.3 Use Power Wheelchair as Vehicle Seat	96
7.3.1 Anchoring Power Wheelchair Used as Vehicle Seat	97
7.3.2 Securing User in Power Wheelchair	97
7.4 Transporting Power Wheelchair Without Occupant	98
<b>8 Maintenance</b>	<b>98</b>
8.1 Maintenance Introduction	98
8.2 Inspection Checks	98
8.2.1 Before Each Use of Power Wheelchair	99
8.2.2 Weekly	99
8.2.3 Monthly	99
8.3 Troubleshooting Remote	100
8.3.1 Fault Diagnosis	100

8.3.2 Fault Codes and Diagnosis Codes .....	100
8.3.3 OON ("Out Of Neutral") .....	101
8.3.4 Drive Inhibit Indication .....	101
8.3.5 Cut-off Voltage .....	101
8.4 Wheels and Tyres .....	101
8.5 Short-term Storage .....	101
8.6 Long-term Storage .....	101
8.7 Cleaning and Disinfection .....	102
8.7.1 General Safety Information .....	102
8.7.2 Cleaning Intervals .....	102
8.7.3 Cleaning .....	102
8.7.4 Disinfection Instructions .....	102
<b>9 After Use .....</b>	<b>102</b>
9.1 Reconditioning .....	102
9.2 Disposal .....	102
<b>10 Technical Data .....</b>	<b>103</b>
10.1 Technical Specifications .....	103
<b>11 Service .....</b>	<b>106</b>
11.1 Inspections Performed .....	106

# 1 General

## 1.1 Introduction

This user manual contains important information about the handling of the product. To ensure safety when using the product, read the user manual carefully and follow the safety instructions.

Only use this product if you have read and understood this manual. Seek additional advice from a healthcare professional who is familiar with your medical condition and clarify any questions regarding the correct use and necessary adjustment with the healthcare professional.

Note that there may be sections in this document, which are not relevant to your product, since this document applies to all available models (on the date of printing). If not otherwise stated, each section in this document refers to all models of the product.

The models and configurations available in your country can be found in the country-specific sales documents.

Invacare reserves the right to alter product specifications without further notice.

Before reading this document, make sure you have the latest version. You find the latest version as a PDF on the Invacare website. Previous product versions may not be described in this Manual's current revision. If you require assistance, please contact Invacare.

If you find that the font size in the printed document is difficult to read, you can download the PDF version from the website. The PDF can then be scaled on screen to a font size that is more comfortable for you.

For more information about the product, for example product safety notices and product recalls, contact your Invacare distributor. See addresses at the end of this document.

In case of a serious incident with the product, you should inform the manufacturer and the competent authority in your country.

## 1.2 Symbols in This Manual

Symbols and signal words are used in this manual and apply to hazards or unsafe practices which could result in personal injury or property damage. This document is printed in greyscale. For your information, the safety messages have the following colour coding according to ANSI Z535.6: Danger (Red), Warning (Orange), Caution (Yellow) and Notice (Blue). See the information below for definitions of the signal words.

### DANGER!

Indicates a hazardous situation that will result in serious injury or death if it is not avoided.

### WARNING!

Indicates a hazardous situation that could result in serious injury or death if it is not avoided.

### CAUTION!

Indicates a hazardous situation that could result in minor or slight injury if it is not avoided.

### NOTICE!

Indicates a hazardous situation that could result in damage to property if it is not avoided.

### Tips and Recommendations

Gives useful tips, recommendations, and information for efficient, trouble-free use.

### Tools

Identifies required tools, components and items which are needed to carry out certain work.

### Other Symbols

### UK Responsible Person

Indicates if a product is not manufactured in the UK.

## 1.3 Compliance

Quality is fundamental to the company's operation, working within the disciplines of ISO 13485.

This product features the CE mark, in compliance with the Medical Device Regulation 2017/745 Class I.

This product features the UKCA mark, in compliance with Part II UK MDR 2002 (as amended) Class I.

We are continuously working towards ensuring that the company's impact on the environment, locally and globally, is reduced to a minimum.

We only use REACH compliant materials and components.

We comply with the current environmental legislations WEEE and RoHS.

### 1.3.1 Product-Specific Standards

The product has been tested and conforms to EN 12184 (Electrically powered wheelchairs, scooters and their chargers) and all related standards (e.g. ISO 7176).

Hereby, Invacare declares that the radio equipment type DLX-REM150 is in compliance with Directive 2014/53/EU. The full text of the EU declaration of conformity is available at the following internet address: [www.invacare.eu.com](http://www.invacare.eu.com).

When equipped with an appropriate lighting system, the product is suitable for use on public roads.

For further information about local standards and regulations, contact your local Invacare distributor. See addresses at the end of this document.

## 1.4 Usability

Only use a power wheelchair when it is in perfect working order. Otherwise, you might put yourself and others at risk.

The following list does not claim to be exhaustive. It is only intended to show some of the situations that could affect the usability of your power wheelchair.

In certain situations, you should immediately stop using your power wheelchair. Other situations allow you to use the power wheelchair to get to your provider.

**You should immediately stop using your power wheelchair if its usability is restricted due to:**

- Unexpected driving behaviour
- brake failure

**You should immediately contact an authorised Invacare provider if the usability of your power wheelchair is restricted due to:**

- the lighting system (if fitted) failing or being defective
- reflectors falling off
- worn thread or insufficient tire pressure
- damage to the armrests (e.g. torn armrest padding)
- damage to the legrest hangers (e.g. missing or torn heel straps)
- damage to the postural belt
- damage to the joystick (joystick cannot be moved into the neutral position)
- cables that are damaged, kinked, pinched or have come loose from the fixation

- the power wheelchair drifting when braking
- the power wheelchair pulling to one side when moving
- unusual sounds developing or occurring

Or if you have the feeling that something is wrong with your power wheelchair.

## 1.5 Warranty Information

We provide a manufacturer's warranty for the product in accordance with our General Terms and Conditions of Business in the respective countries.

Warranty claims can only be made through the provider from whom the product was obtained.

## 1.6 Service Life

We estimate a service life of five years for this product, provided it is used in strict accordance with the intended use as set out in this document and all maintenance and service requirements are met. The estimated service life can be exceeded if the product is carefully used and properly maintained, and provided technical and scientific advances do not result in technical limitations. The service life can also be considerably reduced by extreme or incorrect usage. The fact that we estimate a service life for this product does not constitute an additional warranty.

## 1.7 Limitation of Liability

Invacare accepts no liability for damage arising from:

- Non-compliance with the user manual
- Incorrect use
- Natural wear and tear
- Incorrect assembly or set-up by the purchaser or a third party
- Technical modifications
- Unauthorised modifications and/or use of unsuitable spare parts

# 2 Safety

## 2.1 General Safety Notes

### WARNING!

#### Risk of serious injury or damage

- Improper use of this product may cause injury or damage.
- If you are unable to understand the warnings, cautions or instructions, contact a health care professional or provider before attempting to use this equipment.
  - Do not use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as user manual, service manual or instruction sheet supplied with this product or optional equipment.

### DANGER!

#### Risk of death, serious injury or damage

- Lighted cigarettes dropped onto an upholstered seating system can cause a fire resulting in death, serious injury, or damage. Power wheelchair occupants are at particular risk of death or serious injury from these fires and resulting fumes because they may not have the ability to move away from the power wheelchair.
- DO NOT smoke while using this power wheelchair.

### WARNING!

#### Risk of serious injury or damage

- Storing or using the power wheelchair near open flame or combustible products can result in serious injury or damage.
- Avoid storing or using the power wheelchair near open flame or combustible products.

### WARNING!

#### Risk of death, serious injury or damage

- The power wheelchair may tip over or collide with surroundings when you change its stability characteristics by changing your seating position.
- It is highly recommended not to drive on a slope with a tilted or reclined seat. If it is not possible to avoid such situations please consult your healthcare professional to determine safe conditions. Body weight and height as well as slope surface conditions and the chair setup have an impact on real world stability of the chair.

### WARNING!

#### Risk of damage or injury if power wheelchair is accidentally set into motion

- Switch the power wheelchair off before you get in, get out or handle unwieldy objects.
- When the drive is disengaged, the brake inside the drive is deactivated. For this reason, pushing the power wheelchair by an attendant is only recommended on flat surfaces, never on gradients. Never leave your power wheelchair on a gradient with its motors disengaged. Always re-engage the motors immediately after pushing the power wheelchair (refer to 5.7 *Pushing the Power Wheelchair in Freewheel Mode*, page 91).

### WARNING!

#### Risk of death, serious injury or damage

- Improper monitoring or maintenance may cause injury, damage or death due to ingestion or choking on parts or materials.
- Closely supervise children, pets, or people with physical or mental disabilities.

### WARNING!

#### Risk of death, serious injury or damage

- Risk of entrapment and strangulation when loose personal belongings (e.g. jewellery, scarfs) get caught by moving or protruding parts.
- Make sure that any loose items are clear of moving parts of the power wheelchair, e.g. wheels or powered seating components.
  - Keep your hands, clothing and all other objects away from wheels or powered seating components when they are in operation.
  - Power off power wheelchair immediately to stop any movement.

### WARNING!

#### Risk of death, serious injury or damage

- Improper routing of cables may cause a tripping entanglement or strangulation hazard that may result in death, serious injury or damage.
- Ensure all cables are routed and secured properly.
  - Ensure there are no loops of excess cable extending away from the wheelchair.

**WARNING!**

**Risk of injury if the power wheelchair is driven when ability to operate a vehicle is impaired by medication or alcohol**

- Never drive the power wheelchair under the influence of medication or alcohol. If necessary, the power wheelchair must be operated by an attendant who is physically and mentally able.

**WARNING!**

**Risk of injury if the power wheelchair is switched off while driving, for example by pressing the power button or disconnecting a cable, due to it coming to an abrupt, sharp stop**

- If you have to brake in an emergency, simply release the joystick which will bring you to a halt (refer to the remote user manual for more information).

**WARNING!**

**Risk of injury when transferring power wheelchair to a vehicle for transport with the occupant seated in it**

- It is always better to transfer the power wheelchair to a vehicle without the occupant seated in it.
- If the power wheelchair needs to be loaded up a ramp together with its driver, ensure that the ramp does not exceed the rated slope (refer to *10.1 Technical Specifications, page 103*).
- If the power wheelchair does need to be loaded using a ramp which exceeds the rated slope (refer to *10.1 Technical Specifications, page 103*), then you must use a winch. An attendant can safely monitor and assist the loading process.
- As an alternative you can use a platform lift. Ensure that the total weight of the power wheelchair including the user does not exceed the maximum permissible weight for the platform lift or winch if you are using.

**WARNING!**

**Risk of falling out of the power wheelchair**

- Do not slide forward on the seat, do not lean forward between your knees, do not lean backwards out over the top of the backrest, for example to reach an object.
- If a posture belt is installed, it should be correctly adjusted and used each time you use the power wheelchair.
- When transferring to a different seat, position the power wheelchair as close as possible to the new seat.

**CAUTION!**

**Risk of injury if maximum permissible load is exceeded**

- Do not exceed the maximum permissible load (refer to *10.1 Technical Specifications, page 103*).
- The power wheelchair is only designed for use by a single occupant whose maximum weight does not exceed the maximum permissible load of the device. Never use the power wheelchair to transport more than one person.

**CAUTION!**

**Risk of injury by moving parts**

- Make sure that no injury is incurred by moving parts of the power wheelchair, like wheels or one of the actuator modules (if fitted), especially when children are around.

**CAUTION!**

**Risk of injury due to wrong lifting or dropping of heavy components**

- When maintaining, servicing or lifting any part of your power wheelchair, take into account the weight of the individual components especially the batteries. Be sure at all times to adopt the correct lifting posture and ask for assistance if necessary.

**CAUTION!**

**Risk of injury from hot surfaces**

- Do not leave the power wheelchair in direct sunlight for prolonged periods. Metal parts and surfaces such as the seat and armrests can become very hot.

**CAUTION!**

**Risk of fire or breaking down due to electric devices being connected**

- Do not connect any electric devices to your power wheelchair that are not expressly certified by Invacare for this purpose. Have all electrical installations done by your authorised Invacare provider.

## 2.2 Safety Information for the Electrical System

**WARNING!**

**Risk of death, serious injury or damage**

- Misuse of the power wheelchair may cause the power wheelchair to start smoking, sparking, or burning. Death, serious injury, or damage may occur due to fire.
- DO NOT use the power wheelchair other than its intended purpose.
  - If the power wheelchair starts smoking, sparking, or burning, discontinue using the power wheelchair and seek service IMMEDIATELY.

**WARNING!**

**Risk of death or serious injury**

- Electric shock can cause death or serious injury
- To avoid electric shock, inspect plug and cord for cuts and / or frayed wires. Replace cut cords or frayed wires immediately.

**WARNING!**

**Risk of death or serious injury**

- Failure to observe these warnings can cause an electrical short resulting in death, serious injury, or damage to the electrical system.
- The POSITIVE (+) RED battery cable MUST connect to the POSITIVE (+) battery terminal(s) / post(s).
  - The NEGATIVE (-) BLACK battery cable MUST connect to the NEGATIVE (-) battery terminal(s) / post(s).
  - NEVER allow any of your tools and / or battery cable(s) to contact BOTH battery post(s) at the same time. An electrical short may occur and serious injury or damage may occur.
  - Install protective caps on positive and negative battery terminals.
  - Replace cable(s) immediately if cable(s) insulation becomes damaged.
  - DO NOT remove fuse or mounting hardware from POSITIVE (+) red battery cable mounting screw.

**WARNING!****Risk of death, serious injury, or damage**

Corroded electrical components due to water or liquid exposure can result in death, serious injury, or damage.

- Minimize exposure of electrical components to water and / or liquids.
- Electrical components damaged by corrosion **MUST** be replaced immediately.
- Power wheelchairs that are frequently exposed to water / liquids may require replacement of electrical components more frequently.

**WARNING!****Risk of fire**

Switched on lamps produce heat. If you cover the lamps with fabrics such as clothes, there is a risk that the fabric may catch fire.

- NEVER cover the light system with fabric.

**WARNING!****Risk of death, serious injury or damage when carrying along oxygen systems**

Textiles and other materials that normally would not burn are easily ignited and burn with great intensity in oxygen enriched air.

- Check the oxygen tubing daily, from the cylinder to the delivery site, for leaks and hold away from electrical sparks and any source of ignition.

**WARNING!****Risk of injury or damage due to electrical shorts**

Connector pins on cables connected to the power module can still be live even when the system is off.

- Cables with live pins should be connected, restrained or covered (with non-conductive materials) so that they are not exposed to human contact or materials that could cause electrical shorts.
- When cables with live pins have to be disconnected, for example, when removing the bus cable from the remote for safety reasons, make sure to restrain or cover the pins (with non-conductive materials).

**NOTICE!**

A failure in the electric system can lead to unusual behaviour such as continuous light, no light, or noises from the magnetic brakes.

- If a failure exists, switch off the remote and switch it on again.
- If a failure still exists, then disconnect or remove the power source. Depending on the power wheelchair model, you can either remove the battery packs or disconnect the batteries from the power module. If in doubt which cable to disconnect, contact your provider.
- In any case, contact your provider.

## 2.3 Safety Information for Electromagnetic Interference

This powered power wheelchair is meant to operate without introducing significant electromagnetic disturbances into the environment and without significant degradation of operational performance in the presence of electromagnetic disturbances expected in normal use. Therefore, the wheelchair was successfully tested in accordance with International standards as to its compliance with Elettromagnetic Interference (EMI) regulations.

It is suitable for use in locations in residential environments and in establishments directly connected to a low voltage power supply network which supplies buildings used for domestic purposes. However, electromagnetic fields, such as those generated by radio and television transmitters, and cellular phones can influence the functions of powered power wheelchairs.

Also, the electronics used in our power wheelchairs can generate a low level of electromagnetic interference, which however will remain within the tolerance permitted by law. For these reasons we ask you to please observe the following precautions:

**WARNING!****Risk of malfunction due to electromagnetic interference**

- Do not switch on or operate portable transceivers or communication devices (such as radio transceivers or cellular phones) when the power wheelchair is switched on.
- Avoid getting near strong radio and television transmitters.
- In case the power wheelchair should be set in motion unintentionally or the brakes are released, switch it off immediately.
- Adding electrical accessories / options and other components or modifying the power wheelchair in any way can make it susceptible to electromagnetic interference. Keep in mind that there is no sure way to determine the effect such modifications will have on the overall immunity of the electronic system.
- Report all occurrences of unintentional movement of the power wheelchair, or release of the electric brakes to the manufacturer.

**NOTICE!**

Electromagnetic Interference can result in unintentional movement of the power wheelchair.

- Switch off the remote and, if possible, nearby electronic equipment before switching the wheelchair on again.
- Move away from the source of Electromagnetic interference.

## 2.4 Safety Information for Driving and Freewheel Mode

**DANGER!****Risk of death, serious injury, or damage**

Malfunctioning joystick could cause unintended / erratic movement resulting in death, serious injury, or damage

- If unintended / erratic movement occurs, stop using the wheelchair immediately and contact a qualified technician.

**WARNING!****Risk of serious injury or damage**

Improper positioning while leaning or bending could cause the wheelchair to tip forward resulting in serious injury or damage

- To assure stability and proper operation of your power wheelchair, you must at all times maintain proper balance. Your power wheelchair has been designed to remain upright and stable during normal daily activities as long as you DO NOT move beyond the centre of gravity.
- DO NOT lean forward out of the power wheelchair any further than the length of the armrests.

- DO NOT attempt to reach objects if you have to move forward in the seat or pick them up from the floor by reaching down between your knees.

**WARNING!****Risk of breaking down in adverse weather conditions, i.e. extreme cold, in an isolated area**

- If you are a user with severely limited mobility, we advise that in the case of adverse weather conditions DO NOT attempt a journey without an accompanying attendant.

**WARNING!****Risk of injury if the power wheelchair tips over**

- Inclines and declines can only be travelled up to the maximum safe slope (refer to *10.1 Technical Specifications, page 103*).
- Always return the backrest of your seat or the seat tilt to an upright position before ascending slopes. We recommend that you position the seat backrest and the seat tilt (if fitted) slightly to the rear before descending slopes.
- Only ever drive downhill at a maximum of 50% of the top speed.
- Avoid abrupt braking or accelerating on gradients.
- If at all possible, avoid driving on wet, slippery, icy, or oily surfaces (such as snow, gravel, ice etc.) where there is a risk of you losing control over the vehicle, especially on a gradient. This may include certain painted or otherwise treated wood surfaces. If driving on such a surface is inevitable, then always drive slowly and with the utmost caution.
- Never attempt to overcome an obstacle when on an uphill or downhill gradient.
- Never attempt to drive up or down a flight of steps with your power wheelchair.
- When overcoming obstacles, always observe the maximum obstacle height and information about overcoming obstacles (refer to *5.3.2 Correct Way to Take Obstacles, page 90*).
- Avoid shifting your centre of gravity as well as abrupt joystick movements and changes of direction when the power wheelchair is in motion.
- Never use the power wheelchair to transport more than one person.
- Do not exceed the overall maximum permissible load (refer to *10.1 Technical Specifications, page 103*).
- Note that the power wheelchair will brake or accelerate if you change the driving mode whilst the power wheelchair is in motion.

**WARNING!****Risk of injury if your foot slides off the footrest and gets caught underneath the power wheelchair when it is in motion**

- Make sure each time before you drive the power wheelchair that your feet are squarely and securely in place on the footplates, and that both legrests are properly locked into place.

**WARNING!****Risk of injury if you collide with an obstacle when driving through narrow passages such as doorways and entrances**

- Drive through narrow passages in the lowest driving mode and with due caution.

**WARNING!****Tippling hazard if antitippers are removed, damaged or changed to a position different to the factory settings**

- Antitippers should only ever be removed for dismantling the power wheelchair for transport in a vehicle or for storage.
- The antitippers must always be fitted if the power wheelchair is being used.

**WARNING!****Risk of tipping**

Antitippers (stabilizers) are only effective on firm ground. They sink in on soft ground such as grass, snow or mud if the power wheelchair rests itself on them. They lose their effect and the power wheelchair can tip over.

- Only drive with extreme care on soft ground, especially during uphill and downhill journeys. In the process pay increased attention to the tip stability of the power wheelchair.

## 2.5 Safety Information for Care and Maintenance

**WARNING!****Risk of death, serious injury, or damage**

Incorrect repair and/or servicing of this power wheelchair performed by users/caregivers or unqualified technicians can result in death, serious injury, or damage.

- DO NOT attempt to carry out maintenance work that is not described in this user manual. Such repair and/or service MUST be performed by a qualified technician. Contact a provider or Invacare technician.

**CAUTION!****Risk of accident and loss of warranty if maintenance is insufficient**

- For reasons of safety and in order to avoid accidents which result from unnoticed wear, it is important that this power wheelchair undergoes an inspection once every year under normal operating conditions (see inspection plan contained in service instructions).
- Under difficult operating conditions such as daily travel on steep slopes, or in the case of use in medical care cases with frequently changing power wheelchair users, it would be expedient to carry out intermediate checks on the brakes, accessories / options and running gear.
- If the power wheelchair is to be operated on public roads, the vehicle driver is responsible for ensuring that it is in an operationally reliable condition. Inadequate or neglected care and maintenance of the power wheelchair will result in a limitation of the manufacturer's liability.

## 2.6 Safety Information for Changes and Modifications to the Power Wheelchair

**CE marking of the power wheelchair:**

- The conformity assessment / CE marking was carried out according to the respective valid regulations and only applies to the complete product.
- The CE marking is invalidated if components or accessories / options are replaced or added that have not been approved for this product by Invacare.

- In this case, the company that adds or replaces the components or accessories / options is responsible for the conformity assessment / CE marking or for registering the power wheelchair as a special design and for the relevant documentation.

**WARNING!****Risk of serious injury or damage**

Use of incorrect or improper replacement (service) parts may cause injury or damage

- Replacement parts **MUST** match original Invacare parts.
- Always provide the wheelchair serial number to assist in ordering the correct replacement parts.

**CAUTION!****Risk of injuries and damage to power wheelchair due to unapproved components and accessories / options**

Seating systems, additions and accessories / options which have not been approved by Invacare for use with this power wheelchair can affect the tipping stability and increase tipping hazards.

- Only ever use seating systems, additions and accessories / options which have been approved by Invacare for this power wheelchair.

Seating systems which are not approved by Invacare for use with this power wheelchair do not, under certain circumstances, comply with the valid standards and could increase the flammability and the risk of skin irritation.

- Only use seating systems that have been approved by Invacare for this power wheelchair.

Electrical and electronic components which have not been approved by Invacare for use with this power wheelchair can cause fire hazards and lead to electromagnetic damage.

- Only ever use electrical and electronic components which have been approved by Invacare for this power wheelchair.

Batteries which have not been approved by Invacare for use with this power wheelchair can cause chemical burns.

- Only ever use batteries which have been approved by Invacare for this power wheelchair.

**CAUTION!****Risk of injuries, and damage to the power wheelchair, if unapproved backrests are used**

A retrofitted backrest which is not approved by Invacare for use with this power wheelchair may overload the backrest tube and thus increase the risk of injuries and of damage to the power wheelchair.

- Please contact your Invacare specialist provider who will perform risk analyses, calculations, stability checks etc. to ensure that the backrest can be used safely.

**Important information about maintenance work tools**

Some maintenance work which is described in this manual and can be carried out by the user without problems require the correct tools for proper work. If you do not have the correct tool available we do not recommend that you try to carry out the relevant work. In this case, we urgently recommend that you contact an authorised specialist workshop.

**2.7 General Safety Notes for the Remote****WARNING!****Risk of injury or damage to the power wheelchair**

Do not install, maintain or operate this equipment before you have read and understood all the instructions and all the manuals for this product and all other products that you use or install together with this product.

- Follow the instructions in the user manuals.

**WARNING!****Risk of serious injury or damage to the power wheelchair or surrounding property**

Wrong settings can make the Power wheelchair uncontrollable or unstable. An uncontrolled or unstable wheelchair can cause an unsafe situation such as a crash.

- Performance adjustments must only be made by qualified technicians or by persons who completely understand the programming parameters, the adjustment process, the configuration of the Power wheelchair and the capabilities of the driver.
- Performance adjustments must only be made in dry conditions.

**WARNING!****Risk of injury or damage due to electrical shorts**

Connector pins on cables connected to the power module can still be live even when the system is off.

- Cables with live pins should be connected, restrained or covered (with non-conductive materials) so that they are not exposed to human contact or materials that could cause electrical shorts.
- When cables with live pins have to be disconnected, for example, when removing the bus cable from the remote for safety reasons, make sure to restrain or cover the pins (with non-conductive materials).

**WARNING!****Risk of injury or damage to the power wheelchair**

Risk of unintended movement of the Power wheelchair or seating system when loose personal belongings (e.g. jewellery, scarfs) become entangled around the joystick.

- Make sure that any loose items are clear of the joystick when your power wheelchair is powered up.
- Power off your power wheelchair immediately to stop any movement.

**CAUTION!****Risk of injury from hot surfaces**

Remote module can get hot when exposed to strong sunlight for long periods.

- Do not leave power wheelchair in direct sunlight for long periods.

**NOTICE!**

If you touch the connector pins, they can become dirty or they can be damaged by electrostatic discharge.

- Do not touch the connector pins.

**NOTICE!**

There are no user-serviceable parts inside any case.

- Do not open or disassemble any case.

The Aviva RX10 contains a Bluetooth interface depending on the configuration. Smartphones can connect over the MyLiNX App. This interface can only be used to transfer statistical data, e.g. the battery status. The power wheelchair cannot receive data over this interface. The power wheelchair operates in every environment without any limitations if Bluetooth is enabled. No specific safety precautions are required.

## 3 Product Overview

### 3.1 Intended Use

#### 3.1.1 Product Description

The AVIVA RX10 is a rear wheel drive power wheelchair. It offers a compact design and is easy to control and manoeuvre.

#### 3.1.2 Intended User

This power wheelchair was designed for adults and adolescents whose ability to walk is impaired, but who are still in terms of their eyesight and physically and mentally able to operate an electric power wheelchair.

The maximum user weight for the AVIVA RX10 is 136 kg.

#### 3.1.3 Indications

The use of this power wheelchair is recommended for the following indications:

- The inability or a greatly restricted ability to walk within the scope of the basic requirement to be able to move within one's own four walls.
- The need to leave the dwelling place in order to get some fresh air during a short walk or to reach those places generally to be found at close distance to the dwelling and where everyday business is carried out.

Provision of power wheelchairs for interior and exterior areas is advisable if the use of hand-operated wheelchairs is no longer possible on account of the disability, yet proper operation of an electromotive drive unit is still practicable.


#### Contraindications

There are no contraindications known.

### 3.2 Type Classification

This vehicle has been classified according to EN 12184 as a **class B mobility product** (for indoor and outdoor areas). It is therefore compact and agile enough for indoor areas, but also able to overcome many obstacles in outdoor areas.

### 3.3 Main Parts of Power Wheelchair

### 3.4 Main Parts of Remotes

#### 3.4.1 The Status Indicator

The status indicator is located inside the power button. When the LiNX remote is not powered up, the status indicator is not lit.

When the LiNX remote is powered up and there are no faults with the system, the status indicator lights green.

If there is a fault with the system when powered up, the status indicator flashes red. The number of flashes indicates the type of fault, refer to 8.3.2 *Fault Codes and Diagnosis Codes*, page 100.

#### 3.4.2 Battery Gauge

The driving range status is shown in the battery gauge. The driving range status can vary depending on the driving behaviour, for example when driving up or down hill.


#### CAUTION!

##### Risk of injury or damage due to depleted batteries


If the power wheelchair is used with depleted batteries, you may get stranded in a hazardous situation that can lead to injury or damage.

- Make sure the driving range status is sufficient for the distance you intend to travel.
- If the driving range status is low or very low, it is recommended to charge the batteries before driving.


### 3.4.3 User Interface DLX-REM060

- Drive function


### 3.4.4 User Interface DLX-REM150

- Drive function
- Light system


### 3.4.5 User Interface DLX-REM210

- Drive function
- Seating function



### 3.4.6 User Interface DLX-REM215






- Drive function
- Light system
- Seating function

## 3.5 Labels on the Power Wheelchair

	(A)		<p>Identification of the position of the coupling lever for driving and push operation (only the right side visible in the picture). For details see below.</p>
	(B)	 <p>The colour of the symbol background is yellow on product labels.</p>	<p>Indication of pinch points that could occur on the power wheelchair.</p>
	(C)	 <p>ISO 7176-19</p>	<p>Identification of the tie-down points at the front and back: If the symbol appears on a bright yellow sticker, the anchoring point is suitable for fixation of the power wheelchair in a vehicle for use as a vehicle seat.</p>
	(D)		<p>This symbol indicates the circuit breaker.</p>
	(E)	 <p>ISO 7176-19</p>	<p>The label indicates, that the power wheelchair must be facing forward when used as vehicle seat regarding the requirements of ISO 7176-19.</p>
	(F)		<p>Identification label sticker on the chassis at the right. For details see below.</p>

## Explanation of Symbols on Labels

	Manufacturer		Unique Device Identification
	Date of manufacture		Battery type
	Medical device		Factory setting
	European Conformity		Serial number
	UK Conformity Assessed		Maximum speed
	QR code contains link to user manual		Rated slope
	Read the user manual		Unladen weight
	WEEE Conformity		Maximum user weight

	<p>This symbol indicates the “Drive” position of the coupling lever. In this position the motor is engaged and the motor brakes are operational. You can drive the power wheelchair.</p> <ul style="list-style-type: none"> <li>Note that for driving purposes both motors must always be engaged.</li> </ul>
	<p>This symbol indicates the “Push” position of the coupling lever. In this position the motor is disengaged and the motor brakes are not operational. The power wheelchair can be pushed by an attendant and the wheels turn freely.</p> <ul style="list-style-type: none"> <li>Note that the remote must be switched off.</li> <li>Also refer to 5.7 Pushing the Power Wheelchair in Freewheel Mode , page 91.</li> </ul>
	<p>Read the user manual. This symbol appears on different labels and positions.</p>

### 3.6 Labels on the Remote

	A		Recommendation to read the instruction manual before using the module.					
	B	IPx4	This is the enclosure's ingress protection rating.					
	C		This is the WEEE symbol (Waste Electrical and Electronic Equipment Directive) <sup>1</sup> .					
	D		Tamper evident seal.					
	E		Product label containing: <ul style="list-style-type: none"> <li>Dynamic Controls' website</li> <li>Dynamic Controls' part description</li> </ul>					
	F		Product label containing: <ul style="list-style-type: none"> <li>The serial number</li> <li>The part number</li> </ul>					
	G		The petrol pump indicates the battery charger input.					
	H		<table border="0"> <tr> <td>1. Hardware version</td> <td>4. Application version</td> </tr> <tr> <td>2. Hardware major version</td> <td>5. Application major version</td> </tr> <tr> <td>3. Hardware minor version</td> <td>6. Application minor version</td> </tr> </table>	1. Hardware version	4. Application version	2. Hardware major version	5. Application major version	3. Hardware minor version
1. Hardware version	4. Application version							
2. Hardware major version	5. Application major version							
3. Hardware minor version	6. Application minor version							

<sup>1</sup> This product has been supplied from an environmentally aware manufacturer. This product may contain substances that could be harmful to the environment if disposed of in places (landfills) that are not appropriate according to legislation.

- The 'crossed out wheelie bin' symbol is placed on this product to encourage you to recycle wherever possible.
- Please be environmentally responsible and recycle this product through your recycling facility at its end of life.

## 4 Setup

### 4.1 General Setup Information

 Because of safety reasons the battery is disconnected from the power module ex works. To connect the battery with the power module, refer to 6.2.1 *Connecting / Disconnecting Batteries*, page 92.

#### WARNING!

##### Risk of death, serious injury or damage

Continued use of the power wheelchair that is not set to the correct specifications may cause erratic behaviour of the power wheelchair resulting in death, serious injury, or damage.

- Performance adjustments should only be made by healthcare professionals or persons fully conversant with this process and the driver's capabilities.
- After the power wheelchair has been set-up / adjusted, check to make sure that the power wheelchair performs to the specifications entered during the set-up procedure. If the power wheelchair does not perform to specifications, IMMEDIATELY turn the power wheelchair Off and re-enter set-up specifications. Contact Invacare, if power wheelchair still does not perform to correct specifications.

#### WARNING!

##### Risk of serious injury or damage

Incorrect adjustment to the center of gravity may lead to serious lack of product stability and a risk of tipping or sliding. This can result in serious injury.

- Adjustments to the centre of gravity must ONLY be carried out by qualified technicians. Instructions of how to make such adjustments are available to them and not listed here.

#### WARNING!

##### Risk of death, serious injury or damage

Attaching hardware that is loosely secured or missing could cause instability resulting in death, serious personal injury, or property damage.

- After ANY adjustments, repair or service and before use, make sure that all attaching hardware is present and tightened securely.

#### WARNING!

##### Risk of injury or damage

Incorrect set up of this power wheelchair performed by users / caregivers or unqualified technicians can result in injury or damage.

- DO NOT attempt to set up this power wheelchair. Initial set up of this power wheelchair MUST be performed by a qualified technician.
- Adjustment by the user is only recommended after they have been given appropriate guidance by the healthcare professional.
- DO NOT attempt to carry out the work if you do not have the listed tools available.

#### CAUTION!

##### Risk of injury or damage

The power wheelchair is fitted with an individual, multiply adjustable seating system including adjustable legrests, armrests, a headrest or other options which are used to adapt the seat to the physical requirements and the condition of the user. It is possible that collisions or pinch points can occur between power wheelchair components due to various combinations of adjustment options and their individual settings. When adapting the seating system and the seat functions to the user:

- Beware of pinch points when adjusting the power wheelchair components and
- ensure that no power wheelchair components collide.

#### NOTICE!

The power wheelchair is manufactured and configured individually to the specifications of the order. The assessment must be performed by a healthcare professional according to the user's requirements and health conditions.

- Consult a healthcare professional if you intend to adapt the power wheelchair configuration.
- Any adaptation should be performed by a qualified technician.

Initial setup should always be done by a healthcare professional. Adjustment by the user is only recommended after they have been given appropriate guidance by the healthcare professional.

#### Electrical Adjustment Options

Refer to 5.2.3 *Operating Powered Seating Functions*, page 88 for more information on operating electrical adjustment options.

#### Footplates

All footplates offered by Invacare can be folded upwards.

### 4.2 Setting Up the Remote

The tasks described in this chapter are intended to be performed by trained and authorised service technicians for initial setup. They are not intended to be performed by the user.

#### 4.2.1 Wiring

For safe and reliable operation, the installation of looms and cables must follow the basic principles of power wiring.

Cables must be secured between their connectors and any point of flexing so that flexing forces are not transferred to the connectors.

#### CAUTION!

##### Risk of injury and damage to the remote

Damage to cables increases wiring impedance. A damaged cable can potentially produce localised heat, sparks or arcing and become a source of ignition to surrounding flammable material.

- The installation must ensure that all power cables, including the bus cable, are protected against damage and potential contact with flammable materials.

#### NOTICE!

Cables and remote modules can get damaged if not positioned properly.

- Route and position cables and remote modules so that they are free from physical strain, abuse or damage, such as snagging, crushing, impacts from external objects, pinching or abrasion.

Adequate strain relief must be provided for all cables, and the mechanical limits of the cables and looms must not be exceeded.

Ensure that connectors and connector sockets are shielded from water splashes and water ingress. Cables with female connectors should face horizontally or downwards. Ensure all connectors are fully mated.

**CAUTION!**

**Risk of injury and damage to the remote**

Connector pins on cables connected to the power module can still be live even when the system is off.

- Cables with live pins should be connected, restrained or covered so that they are not exposed to human contact or materials that could cause electrical shorts.

Make sure that the cables do not extend beyond the wheelchair to prevent them from being caught or damaged by external objects. Take particular care on wheelchairs with movable structures such as a powered tilt.

**WARNING!**

**Risk of injury or damage due to electrical shorts**

Continuous contact between user and cable can result in frayed cable jacket. This increases risk of electric shorts.

- Avoid routing the cable where it will come into continuous contact with the end user.

When installing the bus cable, avoid undue straining of the cable and connection points. Flexing of the cable should be minimised wherever possible, to extend service life and minimize the risk of accidental damage.

**NOTICE!**

Regular bending can damage bus cable

- It is recommended to use a cable chain to support the bus cable, where the cable is subject to regular cyclic bending. The maximum stretch of the chain should be less than the length of the bus cable. The force applied to flex the cable should never exceed 10 N.

Appropriate life testing should be carried out to determine / confirm the expected service life and inspection and maintenance schedule.

**4.2.2 Connecting**

**CAUTION!**

**Risk of unintended stops**

If the plug of the remote cable is damaged, the remote cable may come loose while driving. The remote may lose power, could suddenly power down and force an unintended stop.

- Always check the plug of the remote for damage. Contact your provider immediately in case of a damage.

**NOTICE!**

The remote plug and connector socket fit together in one way only.


- Do not force them together.

1. Lightly push to connect the plug of the remote cable and the connector socket with an audible click.

**4.3 Adjusting the Seat**

**4.3.1 Installing the Backrest**

**Installing Manual Backrest**

1. Loosen clamping lever **A** on both sides of the seat.
2. Insert backrest into backrest holders **B**.
3. Tighten clamping lever on both sides of the seat.


**Installing Powered Backrest**

1. Fold up backrest actuator and secure backrest actuator with safety pin to backrest

**4.3.2 Adjusting Manual Backrest**

- 5 mm Allen key


1. Remove screw **A** on both sides of the backrest support **B**.
2. Adjust backrest **C** to desired angle for the user.
3. Tighten screw on both sides of the backrest support.

**4.3.3 Adjusting Seat Height and Manually Seat Tilt**


- 5 mm Allen key
- 13 mm wrench

**Adjusting Seat Height**

1. Remove all seat plates.
2. Remove screws **A**, nuts and washers **B** and washers **C**.
3. Adjust seat height by installing the seat frame with one of the two available drilling holes **D** to the base.
4. Tighten screws.

## Adjusting Manually Seat Tilt

1. Remove front cover.
2. Remove screw, nut and washer (E).
3. Adjust seat tilt by moving the seat frame up or down until the required angle is reached (see following table for depending seat tilt angle and drilling hole position).
4. Install the seat tilt tube (F) to the depending drilling hole and tighten the screw.

Depending on the seat height some combinations of manual adjustable seat tilt angles are not valid, refer to the following table:


Pos. 1	Pos. 2	Pos. 3	Pos. 4	Pos. 5	Pos. 6
Seat-to-floor Height 450 mm					
0°	2.9°	6.9°	10.9°	15.1°	23° <sup>1</sup>
Seat-to-floor Height 470 mm					
-3.2° <sup>1</sup>	0.3°	3.7°	7.8°	11.9°	19.9°

<sup>1</sup> Invalid setting

### 4.3.4 Adjusting Seat Depth

- 5 mm Allen key


1. Remove screws (A) on both sides of the seat.
2. Adjust backrest to desired seat depth for the user.
  - Make sure the backrest support (B) is adjusted to the same drilling holes (C) on both sides equally.
3. Tighten screws on both sides of the seat.

## 4.4 Adjusting Armrest

### 4.4.1 Adjusting Armrest Height

#### Standard Armrest


1. Loosen hand screw (A).
2. Adjust component to desired position.
3. Tighten wing screw.

### 4.4.2 Adjusting Armrest Depth

- 5 mm Allen key


1. Loosen screws (A).
2. Adjust component to desired position.
3. Tighten screws.

### 4.4.3 Adjusting Armrest Width

- 5 mm Allen key


1. Loosen screws (A).
2. Adjust component to desired position.
3. Tighten screw.

## 4.5 Adjusting Remote for Length of User's Arm

Before adjusting the remote for the length of the user's arm, adjust the armrest depth to the desired position, refer to 4.4.2 *Adjusting Armrest Depth*, page 85.

1. Loosen hand screw (A).
2. Adjust component to desired position.
3. Tighten set screw.

## 4.6 Adjusting the Headrest

### CAUTION!

**Risk of injury during use of the power wheelchair as a vehicle seat if a headrest is wrongly adjusted or not installed**

This can cause the neck to be hyperextended during collisions.

- A headrest must be installed. The headrest optionally supplied for this power wheelchair by Invacare is perfectly suitable for use during transport.
- The headrest must be adjusted to the user's ear height.

- It may be necessary to remove and modify the back cushion's cover in order to access the headrest mounting holes on the back pan.
- An optional shim plate is available. It may be installed between the clamp assembly and the back pan to provide additional spacing / clearance on Posture Back and Deep Back.

The headrest clamp hardware is designed to install into existing mounting holes in the backrest pan.

#### 4.6.1 Adjusting Position of Headrest

- 5 mm Allen key

1. Loosen the screws **A**, **B** or the clamping lever **C**.
2. Adjust the headrest to the required position.
3. Re-tighten screws and clamping lever.


#### 4.6.2 Adjusting Height of Headrest

1. Loosen hand screw **A**.
2. Adjust component to desired position.
3. Tighten hand screw.


### 4.7 Adjusting Standard 80° Legrest

#### 4.7.1 Adjusting Width of Legrest


1. Remove screws **A**.
2. Adjust component to desired position.
3. Tighten screw.

#### 4.7.2 Adjusting Length of Legrest

1. Loosen hand lever **A**.
2. Adjust component to desired position.
3. Tighten hand lever.

#### 4.7.3 Swivelling and/or Removing Legrest

1. Push release lever **A** inward or outward.
2. Adjust component to desired position.
3. Pull upward to remove component.

## 5 Using Power Wheelchair

### 5.1 Getting in and out of Power Wheelchair


#### NOTICE!

- The armrest must be removed or swiveled up in order to get into or out of the power wheelchair from the side.

#### 5.1.1 Removing Standard Armrest for Side Transfer

Depending on which side the remote is installed on, you need to disconnect the remote cable before removing the armrest.


1. Pull plug **A** of remote cable to disconnect remote.
2. Loosen wing screw **B**.
3. Remove armrest from holder.

*The graphic serves as example.*

#### 5.1.2 Swivelling Remote to Side

If your power wheelchair is fitted with a swivelling remote holder, then the remote can be moved away to the side, for example, to drive up close to a table.

#### Swing-Away Remote Holder

1. Push remote to swivel remote holder to the side.

#### 5.1.3 Getting in and out of Power Wheelchair

#### WARNING!


#### Risk of serious injury or damage

Improper transfer techniques may cause serious injury or damage

- Before attempting transfers, consult a healthcare professional to determine proper transfer techniques for the user and type of wheelchair.
- Follow the instructions below.

If you do not have sufficient muscle strength, you should ask other persons for help. Use a sliding board, if possible.

1. Reduce the gap between transfer surface (B) and power wheelchair seat (C) to the minimum distance (A) necessary to perform transfer. This might have to be done by an attendant.
2. Align castors parallel to drive wheels to improve stability during transfer.
3. Always switch your power wheelchair off.
4. Always engage both motor locks / clutches and free wheel hubs (if fitted) to prevent wheels from moving.
5. Depending on armrest type of your power wheelchair, detach armrest or swivel it up.
6. Now slide in or out of your power wheelchair.

## 5.2 Before Driving

Before you take your **first trip**, you should familiarise yourself well with the operation of the power wheelchair and with all operating elements. Take your time to test all functions and driving modes.

 If installed, make sure to properly adjust and use the posture belt each time you use the power wheelchair.

**Sitting comfortably = Driving safely**

### WARNING!

#### Risk of death, serious injury or damage

The power wheelchair may tip over or collide with surroundings when you change its stability characteristics by changing your seating position.

- It is highly recommended not to drive on a slope with a tilted or reclined seat. If it is not possible to avoid such situations please consult your healthcare professional to determine safe conditions. Body weight and height as well as slope surface conditions and the chair setup have an impact on real world stability of the chair.


Before **each trip**, make sure that:

- You are within easy reach of all operating controls.
- The battery charge is sufficient for the distance you intend to travel.
- The posture belt is in perfect order.

 The maximum load capacity that is stated in the technical data only states that the system is designed for this mass in total. However, this does not mean that one can sit a person with this body weight in the power wheelchair without restrictions. Attention must be paid to the body proportions, such as height, weight distribution, abdominal belt, leg and calf strap and seat depth. These factors have a strong influence on driving features such as tilt stability and traction. It may possibly be necessary to carry out adaptations to the seat system.

### 5.2.1 Operating the Remote

Your wheelchair always powers up in drive function 1 and is ready to drive. For remotes that provide multiple drive functions (DLX-REM210 or DLX-REM215) the drive function can be changed. For details about changing the drive function, refer to 5.2.2 *Activating the Drive Function*, page 88.

### Powering up the Remote

1. Press power button (A).

If there is no fault with the system, the status indicator lights up green and the battery gauge displays the current battery status, refer to 3.4.2 *Battery Gauge*, page 79.

If there is a fault with the system when powering up, the status indicator indicates the fault with a series of red flashes, refer to 8.3.2 *Fault Codes and Diagnosis Codes*, page 100. If the fault is one that prevents the system from driving, the battery gauge flashes continuously.

If you need to disable the Bluetooth interface this can be done by pressing the power button for more than 3 seconds when powering up the system. It is indicated by the status indicator pulsing for the duration of 6 seconds.

### Powering down the Remote


1. Press power button (A).
2. System powers down and status indicator switches off.

The Power button can also be used to perform an emergency stop, refer to 5.2.6 *Operating an Emergency Stop*, page 89 .


The Power button is also used to lock the system, refer to 5.2.8 *Locking / Unlocking the Remote*, page 89.

### Using the Joystick

The joystick controls the direction and speed of the wheelchair.

When the joystick is deflected from the centre (neutral) position, the wheelchair moves in the direction of the joystick movement.

The speed of the wheelchair is proportional to the joystick deflections, so that the further the joystick is moved from the neutral position, the faster the wheelchair travels.


If the user moves the joystick back to the neutral position, the wheelchair slows down and stops.

If the user releases the joystick from any position other than the neutral position, the joystick returns to the neutral position and the wheelchair slows down and stops.

The joystick can also be used to wake up the system when in sleep mode, if this parameter has been enabled by the provider, refer to 5.2.7 *The Sleep Mode*, page 89.

### Controlling the Maximum Speed


The speed dial allows you to limit the maximum speed of the power wheelchair (that is the speed when the joystick is fully deflected) to suit your preferences and environment.

The speed dial (A) offers ten discrete steps between the lowest speed (B) and the highest speed (C).

1. Turn the speed dial to adjust the maximum speed.

### 5.2.2 Activating the Drive Function

1. Press Drive function key (A). The remote switches to drive function, the Drive function indicator (C) shows the pre-selected drive function (1,2 or 3) and the wheel in the drive status display lights up green.

2. Press Drive function selector keys (A) or (B) until desired drive function lights up.
  - The Drive function indicator (C) shows the drive function.

Drive function 1

Drive function 2


Drive function 3

**i** With the Drive function selector key you can choose between three different drive functions, that are configured by Invacare and can be customised to your needs and requests by the provider.

### 5.2.3 Operating Powered Seating Functions

#### Activate Seating Function

1. Press Seating function key (A).
  - The wheelchair changes to seating function and the Drive / actuator status display (C) lights up amber.
2. Press Seating function selector keys (A) and (B) or move joystick left or right several times until desired seating function lights up, refer to table *Seating Functions*.
3. Deflect joystick to front or rear to activate actuator.

**i** The distance you deflect the joystick determines the dynamics of the movement. If you only deflect the joystick a little, the actuator only moves slowly. If you deflect the joystick as far as you can, the actuator moves faster.

#### Seating Functions

Not every wheelchair has all options.

Powered seat tilt

None

Powered recline


Unspecified

### 5.2.4 Operating the Light System


#### Operating the Lights

**i** If you drive outside, turn on the lights under bad visibility conditions or darkness.

1. Short press Light button (A). The lights are turned on or off.

#### Operating the Direction Indicators


##### Direction Indicator Right

1. Short press right direction indicator button (A). Right direction indicator is turned on or off.

##### Direction Indicator Left


1. Short press left direction indicator button (B). Left direction indicator is turned on or off.

#### Operating the Hazard Lights

1. Short press Hazard lights button (A). The hazard lights are turned on or off.


### 5.2.5 Operating the Horn

1. Press the horn button (A) to sound the horn. The horn sounds for as long as the horn button is pressed.

The horn button is also used for unlocking a locked system, refer to 5.2.8 Locking / Unlocking the Remote, page 89.

### 5.2.6 Operating an Emergency Stop

If you press the power button (A) while driving, an emergency stop is carried out. The remote powers down after this.

### 5.2.7 The Sleep Mode

The sleep mode is no factory setting, but can be enabled by your provider. If this parameter is set ON, the system goes into sleep mode after a period of time without user activity. This period of time can be set by the provider.

The transition to sleep mode is indicated by the remote module's LEDs dimming gradually. During the transition, the joystick, horn, speed dial and power button will continue to operate.


To wake the system from sleep, either press the power button or move the joystick, if this parameter has been enabled by the provider.

### 5.2.8 Locking / Unlocking the Remote


By default, lock function is disabled. Contact your provider to change the configuration.

If function is enabled, the system can be locked / unlocked using below described sequence.

#### Locking the Remote


1. Press power button (A) for more than four seconds when remote is powered up.


When entering the locked state, the battery gauge indicates the transition by flashing LEDs red, amber and green (far left, middle and far right) three times.

#### Unlocking the Remote

1. Press power button (A).
2. Press horn (B) twice within ten seconds.

If you implement the unlock sequence incorrectly or press the power button again before the unlock sequence is complete, the system returns to the locked state.

During an unlock attempt, the battery gauge indicates the system is in a locked state by flashing LEDs red, amber and green (far left, middle and far right) until either the system is powered off, unlocked or the Sequence Timeout is reached.

### 5.2.9 Audible Cues

Next to visual cues of seating and drive functions, audible cues can be set up for REM2XX remotes. Audible Cues are sounds played through the remote's speaker in response to certain system events or navigation actions. The configuration of audible cues can only be made via LiNX Access Tool (iOS or PC) and needs to be done by a qualified technician.

Audible cues are designed to help understand where you are in the LiNX system and are especially beneficial for:

- users with impaired vision or
- users who cannot see the remote or
- users who wish to have additional feedback from their actions, so there is no need to constantly monitor the remote.

There are two types of audible cues.

- Event cues: These are cues played in response to system events.
- Navigation cues: These are cues played in response to menu navigation actions.

#### Event Cues

Not all system events have an audible cue. For example, no audible cue is played when the system moves into sleep mode.

Event cues comprise two or three notes and are played on entering a specific state.

Sound	Event Cue Condition
	Played before powering off

#### Navigation Cues

Navigation cues are played when activating a function.


Sound	Navigation Cue Condition
	Played when activating a drive function.
	Played when activating a seating function.

#### Function Identifier

A function identifier is an optional audible cue that is played directly after a navigation cue. It provides a count by repeating the same note and it is useful, for example, to identify functions of the same type within the same profile.


The function identifier can be set your provider. The number of times this note repeats can be 1 to 6. This parameter can also be set to **None** or **Reverse**. If set to **None**, no function identifier cue is played after a navigation cue. If set to **Reverse**, a single note is played with a longer duration and higher frequency than the note used for the repeating function identifier.

Function = Drive 1                      Identifier = None


Function = Drive 2

Identifier = 1


Function = Drive 3

Identifier = 2

Function = Drive 4

Identifier = Reverse


In this example four drive functions of the same profile are shown. The function identifier has been set for each drive function with the following values: **None, 1, 2 and Reverse.**

**Profile Index**

**Profile Index**

Profile 1 A profile index is played when navigating between profiles, playing one note for the first profile and two notes for the second profile.

Profile 2

In this example in the first profile is entered a drive function.

**Example**

Profile Index	Function	Identifier

If function identifiers are used, then three sound elements are played:

1. profile index (e.g. one note to indicate first profile)
2. navigation cue (e.g. drive function)
3. function identifier (e.g. function identifier is set to 2)

**5.3 Taking Obstacles**

**5.3.1 Maximum Obstacle Height**

The maximum obstacle height is:

- Forward : 50 mm
- Reverse: 50 mm

For more information, refer to *10.1 Technical Specifications, page 103.*

**CAUTION! Risk of tipping over**

- Never approach obstacles at an angle but at 90 degrees as shown below.
- Approach obstacles followed by a gradient with caution. If unsure whether the gradient is too steep or not, move away from the obstacle and if possible try to find another location.
- Never approach obstacles on an uneven and / or loose ground.
- Never drive with too low tire pressure.
- Put your backrest into an upright position before ascending an obstacle.


**CAUTION! Risk of falling out of the power wheelchair and damage to the power wheelchair such as broken castors**

- Never approach obstacles that are higher than the maximum climbable obstacle height.
- Never let the footrest / legrest touch the ground when descending an obstacle.
- If unsure whether taking an obstacle is possible or not, move away from the obstacle and if possible find another location.

**5.3.2 Correct Way to Take Obstacles**

The following instructions how to take obstacles also apply for attendants if the power wheelchair is fitted with an attendant control.

Correct

Incorrect

**Ascending**

1. Approach obstacle or kerb slowly, head-on and at a right angle.
2. Stop in the following position: approx. 10 – 15 cm in front of obstacle.
3. Check position of front wheels. They must be in driving direction and at right angles to obstacle.
4. Approach slowly and keep at consistent speed until rear wheels have also passed over obstacle.

**Descending**

The approach to descend an obstacle is the same as to ascend it with the difference that you need not to stop before descending.

1. Descend obstacle with medium speed.
  - When descending an obstacle too slowly it could happen that the antitippers get stuck and lift the drive wheels off the ground. Driving the power wheelchair is then no longer possible.

**5.4 Driving up and down Gradients**

The power wheelchair was tested in the least stable configuration. For information concerning the least stable configuration and the maximum safe slope, refer to *10.1 Technical Specifications, page 103*

**WARNING! Risk of death, serious injury or damage**

The power wheelchair may tip over or collide with surroundings when you change its stability characteristics by changing your seating position.

- It is highly recommended not to drive on a slope with a tilted or reclined seat. If it is not possible to avoid such situations please consult your healthcare professional to determine safe conditions. Body weight and height as well as slope surface conditions and the chair setup have an impact on real world stability of the chair.

**CAUTION! Risk of tipping over**

- Only ever drive downhill at a maximum of 50% of the top speed. Avoid sudden changes of direction or abrupt braking when driving on slopes.
- Always return the backrest of your seat or the seat tilt (if adjustable seat tilt is available) to an upright position before ascending slopes. We recommend that you position the seat backrest or the seat tilt slightly to the rear before descending slopes.

- Start and stop driving on a slope gently.
- Never attempt to ascend or descend a slope on slippery surfaces or where there is a risk of skidding (such as wet pavement, ice etc).
- Avoid trying to get out of the power wheelchair on an incline or a gradient.
- Always drive straight in the direction the road or path you are on goes, rather than attempting to zigzag.
- Never attempt to turn around on an incline or a slope.

### CAUTION!

#### Braking distance is much longer on a downhill slope than on even terrain

- Never drive down a slope that exceeds the rated slope, refer to *10.1 Technical Specifications, page 103*.

## 5.5 Use on Public Roads

If you wish to use your power wheelchair on public roads and lighting is required by national legislation, then your power wheelchair needs to be fitted with an appropriate lighting system. Additional modifications may be required depending on the country. Contact your Invacare provider if you have any questions.

## 5.6 Parking and Stationary

When parking your power wheelchair or if your power wheelchair is stationary for a prolonged period:

1. Switch the power wheelchair's power system off (ON-/OFF key).

## 5.7 Pushing the Power Wheelchair in Freewheel Mode

The motors of the power wheelchair are equipped with automatic brakes, preventing that the power wheelchair starts rolling out of control when the remote is switched off. When pushing the power wheelchair manually whilst freewheeling, the magnetic brakes must be disengaged.

 Pushing the power wheelchair by hand may require more physical force than expected (more than 100 N). The necessary force nevertheless complies with the requirements of ISO 7176-14.

 The intended use of the freewheel mode is to maneuver the power wheelchair over short distances. The push handles or push bars support this function, but be aware that there might be some impairment between the feet of the assistant and the rear part of the power wheelchair.

### 5.7.1 Disengaging / Engaging Motors


### CAUTION!

#### Risk of the power wheelchair running away

- When the motors are disengaged (for push operation whilst freewheeling), the electromagnetic motor brakes are deactivated. When the power wheelchair is parked, the turn knobs for engaging and disengaging the motors must without fail be locked firmly into the "DRIVE" position (electromagnetic motor brakes activated).

 The motors may only be disengaged by an attendant, not by the user. This ensures that the motors are only disengaged if an attendant is available to secure the power wheelchair and prevent unintended rolling.

The engaging handles for disengaging the motors are located on each motor. For specification of symbols, refer to *3.5 Labels on the Power Wheelchair, page 81*.

1. Switch off remote.
2. Turn engaging handle of motor (A) towards the ground. Motor is disengaged, indicated by symbol (B).
3. Turn engaging handle of motor (A) towards the power wheelchair. Motor is engaged, indicated by symbol (C).

## 6 Control System

### 6.1 Control Protection System

The power wheelchair control system is fitted with an overload protection.

If the drive is severely overloaded over a long period of time (for example, when driving up a steep hill) and especially when the ambient temperature is high, the control system could overheat. In this case, the wheelchair performance is gradually reduced until it comes to a halt. The status display shows a corresponding error code (refer to the user manual of your remote). By switching the remote off and back on again, the error code is cleared and the control system is switched back on. It can however take up to five minutes until the control system has cooled down enough for the drive to restore full performance again.

If the drive is stalled by an insurmountable obstacle, for example, a kerb or similar which is too high, and the driver attempts driving for more than 20 seconds against this obstacle, the control system automatically switches off to prevent the motors from being damaged. The status display shows a corresponding error code (refer to the user manual of your remote). By switching the remote off and back on again, the error code is cleared and the control system is switched back on.

#### 6.1.1 Using Main Fuse

 The power wheelchair cannot be charged when the system is turned off by the main fuse.

 A defective main fuse may be replaced only after checking the entire controls system. A specialised Invacare provider must perform the replacement. For more information about the main fuse type, refer to *10.1 Technical Specifications, page 103*.


The main fuse is an additional safety feature of the controls protection system. When the system is overloaded, the main fuse turns the system off automatically. By pressing the main fuse button, located at the battery box (A) of your power wheelchair, the system can be turned on again.

### 6.2 Batteries


Power is supplied by two 12 V batteries. The batteries are maintenance-free and only need regular charging.

In the following, you find information on how to charge, handle, transport, store, maintain, and use batteries.

### 6.2.1 Connecting / Disconnecting Batteries

1. Loosen handscrews **A** to remove the battery cover on the rear of the power wheelchair.

2. Connect the battery to the power module:
  - a. Connect plug **B** to the B+ socket **C** of the power module.
3. Disconnect the battery from the power module:
  - a. Disconnect plug **B** from the B+ socket **C** of the power module.

### 6.2.2 General Instructions on Charging

New batteries should always be fully charged once before their first use. New batteries will be at their full capacity after having run through approx. 10 - 20 charging cycles (break-in period). This break-in period is necessary to fully activate the battery for maximum performance and longevity. Thus, range and running time of your power wheelchair could initially increase with use.

Gel/AGM lead acid batteries do not have a memory effect as NiCd batteries.

Follow the instructions listed below to ensure safe use and longevity of the batteries:

- Charge 18 hours prior to initial usage.
- We recommend charging the batteries daily after every discharge even after partly discharge, as well as each night over night. Depending on the level of discharge, it can take up to 12 hours until the batteries are fully charged again.
- When the battery indicator reached the red LED range, charge the batteries for 16 hours minimum, neglecting the charge complete display!
- Try to provide a 24 hour charge once a week to make sure that both batteries are fully charged.
- Do not cycle your batteries at a low state of charge without regularly recharging them fully.
- Do not charge your batteries under extreme temperatures. High temperatures above 30 °C are not recommended for charging as well as low temperatures below 10 °C.
- Use only charging devices in Class 2. This class of chargers may be left unattended during charging. All charging devices which are supplied by Invacare comply with these requirements.
- You cannot overcharge the batteries when using the charger supplied with your power wheelchair, or a charger that has been approved by Invacare.
- Protect your charger from sources of heat such as heaters and direct sunlight. If the battery charger overheats, charging current will be reduced and the charging process delayed.

### 6.2.3 Charging Batteries

For further information about charging the batteries, refer to 6.3 *Battery Charger*, page 94.

**WARNING!**

**Risk of injury if using the power wheelchair during charging**

- DO NOT attempt to recharge the batteries and operate the power wheelchair at the same time.
- DO NOT sit in the power wheelchair while charging the batteries.

**WARNING!**

**Risk of fire**

- Only charge the power wheelchair in a well-ventilated environment to prevent the accumulation of flammable gas.
- During the charging process explosive gases occur. Keep the power wheelchair and the charger away from sources of ignition such as flames and sparks.

**WARNING!**

**Risk of explosion and destruction of batteries if the wrong battery charger is used**

- Only ever use the battery charger supplied with your power wheelchair.

**WARNING!**

**Risk of electric shock and damage to the battery charger if it gets wet**

- Protect the battery charger from water.
- Always charge in a dry environment.

**WARNING!**

**Risk of short circuit and electric shock if the battery charger has been damaged**

- Do not use the battery charger if it has been dropped or damaged.

**WARNING!**

**Risk of electric shock and damage to the batteries**

- NEVER attempt to recharge the batteries by attaching cables directly to the battery terminals.

**WARNING!**

**Risk of fire and electric shock if a damaged extension cable is used**

- Only ever use an extension cable if it is absolutely necessary. In case you must use one, make sure it is in good condition.

**WARNING!**


**Risk of death, serious injury or damage**


Improper routing of charger cord(s) may cause tripping, entanglement or strangulation hazard that may result in injury, damage or death.

- Ensure all charger cord(s) are routed and secured properly.
- Close supervision and attention is needed when charging the wheelchair near children, pets or people with physical / mental disabilities.

 The batteries are equipped with safety vents that allow for the evaporation of gas which is generated during the charging process. If the safety vents cannot release the gas properly, the batteries may overheat and permanently deform. An unpleasant smell and reduced function of the batteries may be noticed. However, the batteries remain safe. Stop charging immediately and let the power wheelchair cool down. Please contact your provider to exchange the batteries.

 Please cycle the power prior to charging if wheelchair has not been used within 24 hours. This will ensure the enhanced battery gauge registers the charge to give an accurate reading during use of wheelchair. If it is not powered up, the battery gauge does not indicate the charging state. For more information about the charging state, see user manual of your charger.

1. Switch off power wheelchair.
2. Connect battery charger to charger socket .
3. Connect battery charger to power supply.

If the remote is powered up, the battery gauge indicates the system is connected to the charger by cycling between a left-to-right chase sequence, and then displaying the approximate battery charge state at the end of the chase sequence.

Battery charge state 1  
Red LED on.

Battery charge state 2  
Red and one amber LED on.

Battery charge state 3  
Red and two amber LEDs on.

Battery charge state 4  
Red, amber and one green LED on.

Fully charged  
Green, green, amber, amber and red LEDs on.

For more information about the drive inhibit mode, refer to [8.3.4 Drive Inhibit Indication](#), page 101.

### Battery Synchronisation

 NEW Batteries Only

The wheelchair must be powered up during charging to ensure that accurate battery charge levels display on the remote.

New batteries must be charged fully. The Battery synchronisation procedure MUST be performed within 24 hours of powering on the wheelchair.

The Battery synchronisation procedure can be found in the LiNX service manual and must be performed by a provider or qualified technician.

### 6.2.4 Battery Alarms

#### High Voltage Warning

The batteries are overcharged.  
All LEDs on and the green LEDs flashing.

1. Disconnect battery charger.

#### Low Voltage Warning

The batteries are empty.  
Only one red LED on and flashing.

1. Power down wheelchair.
2. Charge batteries immediately.

### 6.2.5 Disconnecting Power Wheelchair After Charging

1. Once charging is complete, first disconnect battery charger from power supply, then disconnect plug from remote.

### 6.2.6 Storage and Maintenance

Follow the instructions listed below to ensure safe use and longevity of the batteries:

- Always store the power wheelchair fully charged.
- Do not leave the batteries in a low state of charge for an extended length of time. Charge a discharged battery as soon as possible.
- In case your power wheelchair is not used for a longer period of time (that is more than two weeks), the batteries must be charged at least once a month to maintain a full charge and always be charged before use.
- Avoid hot and cold extremes when storing. We recommend to store power wheelchair at a temperature of 15 °C.
- Gel and AGM batteries are maintenance-free. Any performance issues should be handled by a properly trained power wheelchair technician.

### 6.2.7 Instructions on Using Batteries

**CAUTION!**  
**Risk of damaging the batteries.**

- Avoid ultra-deep discharges and never drain your batteries completely.

- Pay attention to the Battery Charge Indicator! Charge the batteries when the Battery Charge Indicator shows that battery charge is low. How fast the batteries discharge depends on many circumstances, such as ambient temperature, condition of the surface of the road, tire pressure, weight of the driver, way of driving and utilisation of lighting, if fitted.
- Always try to charge the batteries, before the indicator on the remote is displaying the state of charge with red colour. The red colour means a remaining capacity about 20 %.
- When the red LED is flashing, the Battery Safe feature is enabled. From this time, speed and acceleration is reduced drastically. It will allow you to move the power wheelchair slowly out of a dangerous situation before the electronic finally cuts off. This is deep discharging and should be avoided.
- Driving with flashing red LED means an extreme stress for the battery and should be avoided under normal circumstances.
- Be aware that for temperatures below 20 °C, the nominal battery capacity starts to decline. For example, at -10 °C the capacity is reduced to about 50 % of the nominal battery capacity.
- To avoid damaging the batteries, never allow them to be fully discharged. Do not drive on heavily discharged batteries if it is not absolutely necessary, as this will strain the batteries unduly and shorten their life expectancy.
- The earlier you recharge the batteries, the longer they live.
- The depth of discharge affects the cycle life. The harder a battery has to work, the shorter is its life expectancy. Examples:
  - One deep discharge stresses the same as 6 normal cycles (green /orange display off).
  - The battery life is about 500 cycles at 80 % discharge (first 4 LED off), or about 5000 cycles at 10 % discharge (one LED off).

- Under normal operation, once a month the battery should be discharged until all green and orange LED are off. This should be done within one day. A 16 hour charge afterwards is necessary as reconditioning.

### 6.2.8 Cleaning Battery Terminals

**WARNING!**

- Most batteries are not sold with instructions. However, warnings are frequently noted on the cell caps. Read them carefully.
- DO NOT allow the liquid in the battery to come in contact with skin, clothes or other possessions. It is a form of acid and harmful or damaging burns may result. Should the liquid touch your skin, wash the area IMMEDIATELY and thoroughly with cool water. In serious cases or if eye contact is made, seek medical attention IMMEDIATELY.

1. Examine battery terminals for corrosion.
2. Verify the plastic caps are in place over battery cell holes.
3. Clean terminals by using a battery cleaning tool, wire brush, or medium grade sand paper.

Upon completion, areas should be shiny, not dull.

4. Carefully dust off all metal particles.

### 6.2.9 Transporting Batteries

The batteries supplied with your power wheelchair are not hazardous goods. This classification is based on the German GGVS Hazardous Goods Road Transport Ordinances, and the IATA/DGR Hazardous Goods Rail Transport / Air Transport Ordinances. Batteries may be transported without restrictions, whether by road, rail or by air. Individual transport companies have, however, guidelines which can possibly restrict or forbid certain transport procedures. Please ask the transport company regarding each individual case.

### 6.2.10 General Instructions on Handling Batteries

- The batteries reach their end of life when the drive range is significantly smaller than usual. Contact your provider or service technician for details.
- Always have your batteries installed by a properly trained power wheelchair technician or a person with adequate knowledge. They have the necessary training and tools to do the job safely and correctly.

### 6.2.11 Handling Damaged Batteries Correctly

If the batteries are defective or damaged, the power wheelchair must not be used under any circumstances. Contact your provider regarding a repair or exchanging the batteries.

Damaged batteries shall only be handles by a properly trained power wheelchair technician.

**WARNING!**

**Risk of burns**

- Never touch or remove overheating batteries. Only unplug the charger.
- Never touch leaking batteries.

**CAUTION!**

**Corrosion and burns from acid leakage if batteries are damaged**

- Remove clothes that have been soiled by acid immediately.

**After contact with skin:**

- Immediately wash affected area with lots of water.

**After contact with eyes:**

- Immediately rinse eyes under running water for several minutes; consult a physician.

### Disposing of Dead or Damaged Batteries Correctly

Batteries are following special disposal rules. Your provider has all information available to safely exchange and dispose the defect batteries.


## 6.3 Battery Charger

The battery charger is intended to be used off-board only in conjunction with the Aviva RX10. Only lead-acid gel type batteries supplied with the Aviva RX10 can be charged.

Charging is done in four stages:


1. In an initialisation only small current is supplied to the battery until the voltage reaches nominal values.
2. During the fast charge stage constant current will be supplied at a maximum.
3. In the equalisation stage a constant voltage is applied to balance the batterie cells.
4. During the floating stage the current will be reduced to a set value, the charger turns off and the battery charging is complete.

### 6.3.1 Symbols on the Charger

- (A) Read the user manual
- (B) UK Conformity Assessed
- (C) For indoor use only
- (D) WEEE conformity
- (E) UL certification
- (F) BFP symbol
- (G) Class II protection class
- (H) European Conformity
- (I) Regulatory compliance mark for Australia / New Zealand
- (J) Korea certification mark
- (K) Product Safety Electrical Appliance & Material Certification for Japan
- (L) TÜV Rheinland / GS symbol
- (M) Fuse type
- (N) Connection diagram for charger output

### 6.3.2 Product Overview

Ⓐ POWER LED indication

Ⓑ CHARGING LED indication

POWER LED	CHARGING LED	Status
Green	Green	Fully charged
OFF	Green	Connected to mains, but not to power wheelchair
OFF	OFF	Not connected to mains and power wheelchair
Red	Green	Batteries are being charged
Flashing red	Green	Output protection (e.g. short-circuit, reverse polarity, under-voltage)
Flashing red	Flashing red	Output over-voltage protection

### 6.3.3 Caution

Before using the battery charger, read all instructions and cautionary markings.

This charger is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the charger by a person responsible for their safety.

#### WARNING!

**Risk of explosion and destruction of batteries if the wrong battery charger is used**

- Only ever use the battery charger supplied with your power wheelchair.

#### WARNING!

**Risk of electric shock!**

- Do not attempt to open or modify this charger in any way. No user serviceable parts inside.

#### WARNING!

**Explosive gases!**

- Prevent flames and sparks.
- Provide adequate ventilation during charging.

- Failure to follow the warnings and instructions may result in electric shock, fire and/or serious injury. Save all safety warnings and instructions for future reference.
- Only use the certified power cord, supplied with the battery charger.
- For indoor use only. Do not expose to rain, snow or damp conditions.
- The charger must be placed in a well ventilated area, protected from direct sunlight and inflammable surfaces, paper textiles etc.
- To reduce the risk of fire hazard, do not cover or obstruct the ventilation holes on the charger.

- Before each use of the charger, check that there is no sign of damage to the case, cables and socket. Contact your provider in the event of any damage. Do not open the charger yourself. Have repairs performed only by a qualified technician.
- Contact your provider if you are unsure of how to operate the charger.
- Keep the charger out of children's reach. Children should be supervised to ensure that they do not play with the charger.
- Only use the charger for the specified type and capacity.
- Disconnect the AC Input supply before plugging or unplugging the connections to the battery.
- Unplug both of the DC output and the AC Input power cord of the charger when the batteries are fully charged.
- The temperature of the case may exceed 41° C in certain circumstances when charging. Avoid to touch the case directly.
- For normal use, the charger requires no maintenance over and above general cleaning. Use a soft cloth, which may be damp if necessary, to clean the charger.
- The charger is supplied with a polarized plug. The charging plug must not be dismantled or modified.
- The off-board battery charger is not to be carried on the power wheelchair.
- Occasional use of the power wheelchair prior to charging complete indication is acceptable if the need of use is urgent.

### 6.3.4 Operational Instructions

Make sure the battery charger output voltage is the same as the output voltage of the connecting battery.

1. Connect the battery charger to the charging socket of the power wheelchair.
2. Connect the power cord. The red POWER LED and the green CHARGING LED indicate that the batteries are being charged.

For more information about charging status, refer to *6.3.2 Product Overview, page 95*.

### 6.3.5 Troubleshooting

- If the POWER LED is off:
  - Check that the charging cable is correctly connected.
  - If the LED still does not light, the battery charger may be faulty. Contact your provider.
- If the CHARGING LED is off:
  - Check that the charging cable is correctly connected.

### 6.3.6 Technical Specifications

Item	• Battery charger
Model	• OH-24V8000MA
Output current (DC)	• 8A
Floating voltage (DC)	• 24 V Rated
Input current (AC)	• max 3.5 A
Input voltage (AC)	• 100 – 240 V, 50 – 60 Hz
Operating temperature	• 0 °C - 40 °C
Battery application	• 26 Ah - 68 Ah (C5)
Dimensions L x W x H [mm]	• 184.1 x 93.6 x 52.5
Weight	• 1172 g

## 7 Transport

### 7.1 Transport — General Information

#### WARNING!

**Risk of death or serious injury to the power wheelchair user and potentially any other nearby occupant of the vehicle, if a power wheelchair is secured using a tie-down system available from a third party supplier and the unladen weight of the power wheelchair exceeds the maximum weight for which the tie-down system is certified**

- Make sure the weight of the power wheelchair does not exceed the weight for which the tie-down system is certified. Consult the tie-down manufacturer's documentation.
- If you are unsure how much your power wheelchair weighs, then you must have it weighed using calibrated scales.

#### WARNING!

#### Risk of injury or damage

If the power wheelchair is fitted with a tray or other auxiliary equipment this could break free during transfer to a vehicle and cause damage or injury to users in the event of a collision.

- When possible, other auxiliary power wheelchair equipment should be either secured to the power wheelchair or removed from the power wheelchair and secured in the vehicle during travel.

#### NOTICE!

- The vehicle should have the floor strength to take the combined weight of the occupant, the power wheelchair and accessories / options.

### 7.2 Transferring Power Wheelchair to Vehicle

#### WARNING!

**Power wheelchair is at risk of tipping over if transferred to a vehicle while user is still seated in power wheelchair**

- Transfer power wheelchair without user whenever possible.
- If power wheelchair with user must be transferred to vehicle using a ramp, ensure that ramp does not exceed rated slope.
- If power wheelchair must be transferred to vehicle using a ramp that does exceed rated slope, a winch must then be used. An attendant can then safely monitor and assist transfer process.
- Alternatively, a platform lift may be used.
- Ensure that total weight of power wheelchair including user does not exceed maximum permitted total weight for ramp or platform lift.
- Power wheelchair should always be transferred to vehicle with backrest in upright position and tilt in upright position (refer to *5.4 Driving up and down Gradients, page 90*).

#### WARNING!

**Risk of injury and damage to power wheelchair and vehicle**

Risk of tipping over or uncontrolled movements of power wheelchair if transferred to vehicle using a ramp that exceeds rated slope.

- Transfer power wheelchair to vehicle without user.
- An attendant must assist transfer process.
- Ensure that all carer fully understand manual of ramp and winch.
- Ensure that winch is suitable for your power wheelchair.
- Use only suitable tie-down points. Do not use removable or movable components of power wheelchair as tie-down points.

#### WARNING!

**Risk of injury and damage to power wheelchair**

If power wheelchair must be transferred to vehicle via a lift, when remote is turned on, there is a risk that device may act erratically and fall off lift.

- Before transferring power wheelchair via lift, turn off product and disconnect either bus cable from remote or batteries from system.

1. Drive or push your power wheelchair into transport vehicle using suitable ramp.
2. Anchor the power wheelchair to the transport vehicle, refer to *7.3 Use Power Wheelchair as Vehicle Seat, page 96* and secure the user in the power wheelchair, refer to *7.3.2 Securing User in Power Wheelchair, page 97*.

### 7.3 Use Power Wheelchair as Vehicle Seat

Not every power wheelchair does automatically have the permission to be used as a vehicle seat. The following labels explain whether the power wheelchair can be used as a vehicle seat or not.

If the power wheelchair may NOT be used as a vehicle seat, this is identified by the following label:

ISO 7176-19

If the power wheelchair can be used as a vehicle seat, the tie-down points are identified by the following label:

ISO 7176-19

In order to use a power wheelchair as a vehicle seat, it must be fitted with tie-down points to enable anchoring in the motor vehicle. These accessories / options may be included in the standard scope of power wheelchair order and delivery in some countries (UK for example), but may also be obtained from Invacare as an option in other countries.

**The following information is only relevant if your power wheelchair can be used as a vehicle seat:**

#### WARNING!

**Risk of serious injury**

The power wheelchair has been designed and tested to conform to the requirements of ISO 7176-19 for use only as forward-facing seat in a motor vehicle.

The power wheelchair has been dynamically tested in a forward-facing orientation with the ATD (anthropomorphic test device, "crash-test dummy") restrained by a three-point belt restraint.

If any of the instructions is not followed, serious injuries or damage can occur in an event of a collision:

- Alterations or substitutions shall not be made to the power wheelchair securement points or to structural and frame parts or components since this can affect the crashworthiness of the power wheelchair, and it can also change the performance of the power wheelchair in normal use. If it is considered necessary to make these kinds of alterations, Invacare shall be consulted.
- Only use spill-proof sealed batteries approved by Invacare.
- It is imperative that the power wheelchair is inspected by an authorised provider for determination whether the power wheelchair is suitable for reuse after any type of vehicle collision.

The power wheelchair may be used as a vehicle seat in connection with an anchoring system that has been checked and approved in accordance with ISO 10542. The transporting vehicle must be professionally converted to anchor the power wheelchair. Contact your vehicle's manufacturer for more information.

If possible, the user should always leave the power wheelchair to use a vehicle seat and the vehicle-manufacturer-installed restraint system. The unoccupied power wheelchair should be stored in a cargo area or secured in the vehicle during travel.

A power wheelchair permitted as vehicle seat has undergone a crash test in accordance to ISO 7176-19 for use in road vehicles and meets the requirements for forward facing transport and head on collisions. The "crash-test dummy" was secured using pelvic and upper body safety belts. Both types of safety belt should be used in order to minimize the risk of injuries to head or upper body.

Invacare tests with tie-down systems, that meet the requirements of ISO 10542-1 and the kerb weight of the power wheelchair. For information concerning the kerb weight, refer to 10.1 *Technical Specifications*, page 103.

### 7.3.1 Anchoring Power Wheelchair Used as Vehicle Seat


The power wheelchair is fitted with tie-down points. Snap hooks or belt loops can be used for fixation.

#### CAUTION!

**There is a risk of injury if the power wheelchair is not properly secured during use as a vehicle seat**

- Always use a tie-down system suitable for the combined weight of occupant and power wheelchair.
- If possible, the user should always leave the power wheelchair to use a vehicle seat and the safety belts provided with the vehicle.
- The power wheelchair should always be anchored facing in the transport vehicle's intended direction of travel.
- The power wheelchair must always be secured in accordance with the power wheelchair and anchoring system manufacturers' operating manual.
- Always remove and secure any accessory parts fixed to the power wheelchair such as chin controls or tables.
- If your power wheelchair is equipped with an angle adjustable backrest, then it must be placed in an upright position.

The power wheelchair in the graphic faces backwards due to better visibility of tie-down points.

The arrows indicate the direction of fixation to the vehicle.

1. Secure forward-facing power wheelchair with tie-down system belts at following locations:
  - a. Four tie-down points on power wheelchair base (two at front (A) and two at rear (B)).
2. Secure power wheelchair by tensioning belts in accordance with tie-down system manufacturer's user manual.

### 7.3.2 Securing User in Power Wheelchair


#### CAUTION!

**Risk of injury during use of the power wheelchair as a vehicle seat if a headrest is wrongly adjusted or not installed**

This can cause the neck to be hyperextended during collisions.

- A headrest must be installed. The headrest optionally supplied for this power wheelchair by Invacare is perfectly suitable for use during transport.
- The headrest must be adjusted to the user's ear height.

The "crash test dummy" was secured using pelvic and upper body safety belts. Both types of safety belt should be used in order to minimize the risk of injuries to head or upper body.


#### CAUTION!

**Risk of injury if the user is not properly secured within the power wheelchair**


- Safety restraint devices must only be used when the wheelchair user's weight is 23 kg or more.
- Even if the power wheelchair is fitted with a posture belt or any other power wheelchair integrated belt system, this is no substitute for a proper safety belt which complies with ISO 10542 in the transport vehicle. Always use the safety belt installed in the transport vehicle.
- Safety belts must be pulled as tightly as possible without causing the user discomfort.
- Safety belts must not be positioned while twisted.
- Ensure that the third seat belt anchorage point is not fixed directly to the vehicle floor, but to one of the vehicle uprights.
- Both pelvic and upper torso restraint belts must be used to restrain the occupant to reduce the possibility of head and chest impacts with the vehicle components. They shall be used together only as designed for.

- Any wheelchair anchored occupant restraint i.e. 3-point belt, harness or postural supports (lap straps, lap belts) should not be used or relied on for occupant restraint in a moving vehicle. Always use a vehicle anchored and certified occupant restraint system instead.
- Care should be taken when applying the occupant restraint to position the seatbelt buckle so that the release button will not be contacted by power wheelchair components during transport and during a crash.
- Safety belts must be in contact with the user's body. They must not be held at a distance from the user's body using parts of the power wheelchair such as armrests or wheels.

The pelvic belt should be positioned in the area between the user's pelvis and thighs so that it is unobstructed and not too loose. The ideal angle of the pelvic belt to the horizontal is between 45° and 75°. The maximum permissible angle is between 30° and 75°. The angle should never be less than 30°!

The safety belt installed in the transporting vehicle should be applied as shown in the illustration above.

- 1) centre line of the body                      2) centre of the sternum

## 7.4 Transporting Power Wheelchair Without Occupant

### CAUTION! Risk of injury

- If you are unable to fasten your power wheelchair securely in a transport vehicle, Invacare recommends that you do not transport it.

Your power wheelchair may be transported without restrictions, whether by road, rail or by air. Individual transport companies have, however, guidelines which can possibly restrict or forbid certain transport procedures. Please ask the transport company regarding each individual case.

- Before transporting your power wheelchair, make sure the motors are engaged and that the remote is switched off.
- Invacare strongly recommends that you additionally disconnect the batteries, refer to 6.2.1 *Connecting / Disconnecting Batteries*, page 92.
- Invacare strongly recommends securing the power wheelchair to the floor of the transporting vehicle.

## 8 Maintenance

### 8.1 Maintenance Introduction

### CAUTION! Risk of injury or damage

- If maintenance work or servicing is done, while the power wheelchair is in use, the occupant can get injured or the power wheelchair can get damaged.
- The occupant should not sit in the power wheelchair during maintenance work or servicing.

The term “Maintenance” means any task performed to ensure that a medical device is in good working order and ready for use as intended. Maintenance encompasses different areas, such as everyday care and cleaning, inspection checks, repair tasks and refurbishment.

It is recommended, to have your power wheelchair checked once a year by an authorised Invacare provider to maintain its driving safety and roadworthiness.

### 8.2 Inspection Checks

The following tables list inspection checks that should be performed by the user and their intervals. If the power wheelchair fails to pass one of the inspection checks, refer to the chapter indicated or contact your authorised Invacare provider. A more comprehensive list of inspection checks and instructions for maintenance work can be found in the service manual for this device, which is available through Invacare. The service manual is intended to be used by trained and authorised service technicians, and describes tasks which are not intended to be performed by the user.

### 8.2.1 Before Each Use of Power Wheelchair

Item	Inspection Check	If Not Passed
Screwed connections	Check all connections, such as backrests and wheels, for tight fit.	<ul style="list-style-type: none"> <li>Contact your provider.</li> </ul>
Signal horn	Check for correct function.	<ul style="list-style-type: none"> <li>Contact your provider.</li> </ul>
Lighting system	Check that all lights, such as turn indicators, head lamps and tail lights, are functioning correctly.	<ul style="list-style-type: none"> <li>Contact your provider.</li> </ul>
Battery box locking system	Check to ensure that the battery box locking system is functioning correctly. Locking pins must be completely engaged in the holes provided for them (refer to 7.4 <i>Transporting Power Wheelchair Without Occupant, page 98</i> ).	<ul style="list-style-type: none"> <li>Contact your provider.</li> </ul>
Batteries	Make sure the batteries are charged. Refer to 3.4.2 <i>Battery Gauge, page 79</i> and 6.2.3 <i>Charging Batteries, page 92</i> for a description of the battery charge indicator.	<ul style="list-style-type: none"> <li>Charge the batteries (refer to 6.2.3 <i>Charging Batteries, page 92</i>).</li> </ul>

### 8.2.2 Weekly

Item	Inspection Check	If Not passed
Armrests / side parts	Check that armrests are firmly attached in their holders and do not wobble.	<ul style="list-style-type: none"> <li>Tighten the screw or clamping lever that holds armrests.</li> <li>Contact your provider.</li> </ul>
Tires (pneumatic)	Check that the tires are undamaged.	<ul style="list-style-type: none"> <li>Contact your provider.</li> </ul>
	Check that the tires are inflated to the correct pressure.	<ul style="list-style-type: none"> <li>Inflate the tire to the correct pressure (refer to 8.4 <i>Wheels and Tyres, page 101</i> and 10.1 <i>Technical Specifications, page 103</i>).</li> </ul>
Tires (puncture-proof)	Check that the tires are undamaged.	<ul style="list-style-type: none"> <li>Contact your provider.</li> </ul>
Anti-tippers	Check that anti-tippers are firmly attached and do not wobble. Check that the spring clips of the anti-tippers are in good order and secure the antitippers correctly.	<ul style="list-style-type: none"> <li>Contact your provider.</li> </ul>

### 8.2.3 Monthly

Item	Inspection Check	If Not Passed
All upholstered parts	Check for damage and wear.	<ul style="list-style-type: none"> <li>Contact your provider.</li> </ul>
Removable legrests	Check whether legrests can be fixed securely and whether loosening mechanism is properly operable.	<ul style="list-style-type: none"> <li>Contact your provider.</li> </ul>
	Check that all adjustment options function properly.	<ul style="list-style-type: none"> <li>Contact your provider.</li> </ul>
Castors	Check that castors rotate and swivel freely.	<ul style="list-style-type: none"> <li>Contact your provider.</li> </ul>
Drive wheels	Check that the tyre pattern is 3 mm minimum. Check that drive wheels rotate without wobbling. To do this, it is easiest to have someone standing behind the power wheelchair and observe the drive wheels as you drive away from the person.	<ul style="list-style-type: none"> <li>Contact your provider.</li> </ul>
Electronics and connectors	Check all cables for damage and all connecting plugs for snug fit.	<ul style="list-style-type: none"> <li>Contact your provider.</li> </ul>

## 8.3 Troubleshooting Remote

### 8.3.1 Fault Diagnosis

If the electronic system shows a fault, use the following fault-finding guide to locate the fault.

Ensure that the drive electronics system is powered up before starting any diagnosis.

#### If the Status Display is OFF:

- Check whether the drive electronics system is powered up.
- Check whether all cables are correctly connected.
- Ensure that the batteries are not discharged.

#### If a Fault Number is Displayed in the Status Display:

- Proceed to the next section.

### 8.3.2 Fault Codes and Diagnosis Codes

If there is a fault with the system when it is powered up, the status indicator flashes red. The number of flashes indicates the type of fault.

The table below describes the fault indication, and a few possible actions that can be taken to rectify the problem. The actions listed are not in any particular order and are suggestions only. The intention is that one of the suggestions may help you clear the problem. If in doubt, contact your provider.

Flash Code	Fault Description	Possible Action
1	Remote fault	<ul style="list-style-type: none"> <li>• Check cables and connectors.</li> <li>• Check other remotes, if installed.</li> <li>• Contact your provider.</li> </ul>
2	Network or configuration fault	<ul style="list-style-type: none"> <li>• Restart power wheelchair.</li> <li>• Check cables and connectors.</li> <li>• Recharge the batteries.</li> <li>• Check charger.</li> <li>• Contact your provider.</li> </ul>
3	Motor 1 <sup>1</sup> fault	<ul style="list-style-type: none"> <li>• Check cables and connectors.</li> <li>• Contact your provider.</li> </ul>
4	Motor 2 <sup>1</sup> fault	<ul style="list-style-type: none"> <li>• Check cables and connectors.</li> <li>• Contact your provider.</li> </ul>
1	For configuration of motors, refer to 10.1 Technical Specifications, page 103.	
5	Motor 1 <sup>1</sup> magnetic brake fault	<ul style="list-style-type: none"> <li>• Check cables and connectors.</li> <li>• Check left magnetic brake is engaged.</li> <li>• Contact your provider.</li> <li>• Refer to 5.7 Pushing the Power Wheelchair in Freewheel Mode, page 91.</li> </ul>
6	Motor 2 <sup>1</sup> magnetic brake fault	<ul style="list-style-type: none"> <li>• Check cables and connectors.</li> <li>• Check right magnetic brake is engaged.</li> <li>• Contact your provider.</li> <li>• Refer to 5.7 Pushing the Power Wheelchair in Freewheel Mode, page 91.</li> </ul>
1	For configuration of motors, refer to 10.1 Technical Specifications, page 103.	
7	Module fault (other than remote module)	<ul style="list-style-type: none"> <li>• Check cables and connectors.</li> <li>• Check modules.</li> <li>• Recharge batteries.</li> <li>• Contact your provider.</li> <li>• If the power wheelchair was stalled, reverse away or remove obstacle.</li> </ul>

### 8.3.3 OON (“Out Of Neutral”)


OON (“Out Of Neutral”) is a safety feature that prevents accidental driving or seating movements, when:

- the system is powering up,
- after a function change or
- when the system comes out of an inhibit or drive lock-out.

#### Drive OON Warning

The joystick must be in the centre position:


- when a system is powering up,
- on a function change or
- when transitioning from a drive lock-out or inhibit state.

Otherwise a drive OON warning is displayed. During a drive OON warning, the battery gauge LEDs and the drive wheel indicator (if fitted) flash continually (all on, followed by all off) to alert the user. In this state the wheelchair does not drive. If the joystick is returned to the centre position, the warning clears and the wheelchair drives normally.


### 8.3.4 Drive Inhibit Indication

The drive inhibit mode ensures that the wheelchair does not drive when connected to the charger.

Drive inhibit mode is indicated by the battery gauge with a right-to-left chase sequence. The chase sequence continues until the fault condition has been cleared.

### 8.3.5 Cut-off Voltage

When the battery voltage decreases below the battery cut-off voltage:

- the status indicator flashes red (Flash code 2, refer to 8.3.2 *Fault Codes and Diagnosis Codes*, page 100),
- the red LED on the battery gauge flashes,
- the horn sounds once every ten seconds.

## 8.4 Wheels and Tyres

### Dealing With Wheel Damages

In case of having a damaged wheel, contact your provider. Because of safety reasons do not have the wheel repaired by yourself or by not authorised persons.

### Dealing With Pneumatic Tyres

#### NOTICE!

Never drive with too low tyre pressure, this could result in damage to tyre.

If tyre pressure is exceeded rim could be damaged.

- Inflate tyres to recommended pressure.

Use tyre gauge to check pressure.

Check weekly that the tyres are inflated to the correct pressure, refer to 8.1 *Maintenance Introduction*, page 98.

For recommended tyre pressure see inscription on tyre/rim or contact Invacare. Compare table below for conversion.

psi														
22	23	25	26	28	29	30	32	33	35	36	38	39	41	44
bar														
1,5	1,6	1,7	1,8	1,9	2,0	2,1	2,2	2,3	2,4	2,5	2,6	2,7	2,8	2,9

## 8.5 Short-term Storage

In case a serious fault is detected, a number of safety mechanisms are built into your power wheelchair and will protect it. The power module prevents your power wheelchair from driving.

When the power wheelchair is in such a condition and while waiting for repair:

1. Switch off power.
2. Disconnect the batteries.  
Depending on the power wheelchair model, you can either remove the battery packs or disconnect the batteries from the power module. Refer to the corresponding chapter about disconnecting the batteries in the service manual, which is available through Invacare.
3. Contact your provider.

## 8.6 Long-term Storage

In case your power wheelchair is not used for a longer period of time, you need to prepare it for storage to ensure a longer life for your power wheelchair and batteries.

### Storing Power Wheelchair and Batteries

- We recommend to store the power wheelchair at a temperature of 15 °C, avoid hot and cold extremes when storing to ensure a long service life of the product and batteries.
- The components are tested and approved for greater temperature ranges as detailed below:
  - Allowable temperature range to store the power wheelchair is -40° up to 65 °C.
  - Allowable temperature range to store batteries is -25° up to 65 °C.
- Even not being used, batteries discharge themselves. Best practice is to disconnect the battery supply from the power module if storing the power wheelchair longer than two weeks, refer to 6.2.1 *Connecting / Disconnecting Batteries*, page 92.
- Batteries should always be fully charged before storing.
- If storing the power wheelchair longer than four weeks, check the batteries once a month and recharge as needed (before gauge reads half full) to avoid damage.
- Store in a dry, well-ventilated environment protected from outer influences.
- Slightly overinflate pneumatic tyres.
- Position the power wheelchair on flooring that is not discoloured by contact with tyre rubber.

### Preparing Power Wheelchair for Use

- Re-connect the battery supply to the power module.
- The batteries must be charged before use.
- Have the power wheelchair checked by an authorised Invacare provider.

## 8.7 Cleaning and Disinfection

### 8.7.1 General Safety Information

#### CAUTION!

##### Risk of contamination

- Take precautions for yourself and use appropriate protective equipment.

#### CAUTION!

##### Risk of electric shock and product damage

- Switch off the device and disconnect from mains, if applicable.
- When cleaning electronic components consider their protection class regarding water ingress.
- Make sure that no water splashes to the plug or the wall outlet.
- Do not touch the power socket with wet hands.

#### NOTICE!

Wrong fluids or methods can harm or damage the product.

- All cleaning agents and disinfectants used must be effective, compatible with one another and must protect the materials they are used to clean.
- Never use corrosive fluids (alkalines, acid etc.) or abrasive cleaning agents. We recommend an ordinary household cleaning agent such as dishwashing liquid, if not specified otherwise in the cleaning instructions.
- Never use a solvent (cellulose thinner, acetone etc.) that changes the structure of the plastic or dissolves the attached labels.
- Always make sure that the product is completely dried before taking into use again.

For cleaning and disinfection in clinical or long-term care environments, follow your in-house procedures.

### 8.7.2 Cleaning Intervals

#### NOTICE!

Regular cleaning and disinfection enhance smooth operation, increases the service life and prevents contamination.

Clean and disinfect the product:

- regularly while in use,
- before and after any service procedure,
- when it has been in contact with any body fluids,
- before using it for a new user.

### 8.7.3 Cleaning

#### NOTICE!

- The product does not tolerate cleaning in automatic washing plants, with high-pressure cleaning equipment or steam.

#### NOTICE!

Dirt, sand and seawater can damage the bearings and steel parts can rust if the surface is damaged.

- Only expose the wheelchair to sand and seawater for short periods and clean it after every trip to the beach.
- If the wheelchair is dirty, wipe off the dirt as soon as possible with a damp cloth and dry it carefully.

1. Remove any installed optional equipment (only optional equipment which does not require tools).
2. Wipe down the individual parts using a cloth or soft brush, ordinary household cleaning agents (pH = 6 - 8) and warm water.
3. Rinse the parts with warm water.

4. Thoroughly dry the parts with a dry cloth.

Car polish and soft wax can be used on painted metal surfaces to remove abrasions and restore gloss.

#### Cleaning Upholstery

For cleaning upholstery refer to the instructions on the labels of the seat, cushion and backrest cover.

If possible, always overlap hook and loop strips (the self-gripping parts) when washing, to minimize lint and thread build-up on hook strips and prevent damage to upholstery fabric by these.

### 8.7.4 Disinfection Instructions

Method: Follow the application notes for the used disinfectant and wipe-disinfect all accessible surfaces.

Disinfectant: Ordinary household disinfectant.

Drying: Allow the product to air-dry.

## 9 After Use

### 9.1 Reconditioning

This product is suitable for reuse. To recondition the product for a new user, carry out the following actions:

- Inspection according to service plan, refer to the service manual, which is available through Invacare.
- Cleaning and disinfection, refer to *8.7 Cleaning and Disinfection, page 102*.
- Adaptation to the new user, refer to *4 Setup, page 83*.

Make sure that the user manual is handed over with the product.

If any damage or malfunction is detected, do not reuse the product.

### 9.2 Disposal

#### WARNING!

##### Environmental hazard

Device contains batteries.

This product may contain substances that could be harmful to the environment if disposed of in places (landfills) that are not appropriate according to legislation.

- DO NOT dispose of batteries in normal household waste.
- DO NOT throw batteries into a fire.
- Batteries MUST be taken to a proper disposal site. The return is required by law and free of charge.
- Only dispose of discharged batteries.
- Cover terminals of batteries prior to disposal.
- For information about the correct handling of damaged batteries, refer to *6.2.11 Handling Damaged Batteries Correctly, page 94*.

Be environmentally responsible and recycle this product through your recycling facility at its end of life.

Disassemble the product and its components, so the different materials can be separated and recycled individually.

The disposal and recycling of used products and packaging must comply with the laws and regulations for waste handling in each country. Contact your local waste management company for information.

## 10 Technical Data

### 10.1 Technical Specifications

The technical information provided hereafter applies to a standard configuration or represents maximum achievable values. These can change if accessories are added. The precise changes to these values are detailed in the sections for the respective accessories.

Note that there may be values in this list, which are not relevant to your product, since this list applies to all available models (on the date of printing). If not otherwise stated, each value in this list refers to all models of the product.

The models and configurations available in your country can be found in the country-specific sales documents.

Note that in some cases the measured values may vary up to  $\pm 10$  mm.

Permissible Operating, Storage and Humidity Conditions	
Temperature range for operation according to ISO 7176-9:	<ul style="list-style-type: none"> <li>-25 °C ... +50 °C</li> </ul>
Recommended storage temperature:	<ul style="list-style-type: none"> <li>15 °C</li> </ul>
Temperature range for storage according to ISO 7176-9:	<ul style="list-style-type: none"> <li>-20 °C – +60 °C with batteries</li> <li>-40 °C – +65 °C without batteries</li> </ul>
Operating humidity range according to ISO 7176-9:	<ul style="list-style-type: none"> <li>0 ... 90% RH</li> </ul>

Electrical System	
Batteries <sup>1</sup>	<ul style="list-style-type: none"> <li>51.0Ah (C20) / 43.3 Ah (C5) Cycle Endurance (IEC60251-1:2005) min 300 Cycles</li> </ul>
<p>1 Usable battery capacity depending on the discharge time. C20: Discharge over a period of 20 hours. C5: Discharge over a period of 5 hours.</p>	
Main Fuse	<ul style="list-style-type: none"> <li>63 A</li> </ul>
Degree of Protection	<ul style="list-style-type: none"> <li>IPX4<sup>2</sup></li> </ul>
<p>2 IPX4 classification means that the electrical system is protected against spray water.</p>	

Electrical Specifications of DLX-REM060/150/210/215				
Parameter	Min.	Nominal	Max.	Units
Operating voltage (V <sub>batt</sub> )	<ul style="list-style-type: none"> <li>17</li> </ul>	<ul style="list-style-type: none"> <li>24</li> </ul>	<ul style="list-style-type: none"> <li>34</li> </ul>	<ul style="list-style-type: none"> <li>V</li> </ul>
Idle current	<ul style="list-style-type: none"> <li>-</li> </ul>	<ul style="list-style-type: none"> <li>56</li> </ul>	<ul style="list-style-type: none"> <li>-</li> </ul>	<ul style="list-style-type: none"> <li>mA at 24V</li> </ul>
Quiescent current (power off)	<ul style="list-style-type: none"> <li>-</li> </ul>	<ul style="list-style-type: none"> <li>-</li> </ul>	<ul style="list-style-type: none"> <li>0.23</li> </ul>	<ul style="list-style-type: none"> <li>mA at 24V</li> </ul>

Charging Device	
Output Current	<ul style="list-style-type: none"> <li>8 A</li> </ul>
Output Current	<ul style="list-style-type: none"> <li>24 V nominal</li> </ul>

Drive Wheel Tyres	
Tyre Type	<ul style="list-style-type: none"> <li>14 inch puncture-proof, pneumatic</li> </ul>
Tyre Pressure	The recommended maximum tyre pressure in bar or kpa is marked on the side wall of the tyre or the rim. If more than one value is listed, the lower one in the corresponding units applies. (Tolerance = -0.3 bar, 1 bar = 100 kpa)

Castor Tyres	
Tyre Type	<ul style="list-style-type: none"> <li>8 inch solid</li> <li>9 inch solid</li> </ul>

Driving Characteristics		
	Seat Width 430 mm	Seat Width 480 mm
Speed		<ul style="list-style-type: none"> <li>• 6 km/h</li> <li>• 10 km/h</li> </ul>
Max. Stopping Distance		
Normal Operation		<ul style="list-style-type: none"> <li>• 870 mm (6 km/h)</li> <li>• 1830 mm (10 km/h)</li> </ul>
Emergency Operation		<ul style="list-style-type: none"> <li>• 850 mm (6 km/h)</li> <li>• 1820 mm (10 km/h)</li> </ul>
Driving Range in Accordance with ISO 7176-4 <sup>3</sup>		
3	Note: The drive range of a power wheelchair is strongly influenced by external factors, such as the speed setting of the wheelchair, the charging state of the batteries, surrounding temperature, local topography, road surface characteristics, tyre pressure, weight of user, drive style and use of batteries for lighting, servos etc. The specified values are theoretical maximum achievable values measured according to ISO 7176-4.	
Continuous Driving Distance Range		<ul style="list-style-type: none"> <li>• 29.2 km (6 km/h)</li> <li>• 25.6 km (10 km/h)</li> </ul>
Manoeuvring Distance Range		<ul style="list-style-type: none"> <li>• 10.4 km (6 km/h)</li> <li>• 8.3 km (10 km/h)</li> </ul>
Max. Climbable Obstacle Height		<ul style="list-style-type: none"> <li>• Forward: 50 mm</li> <li>• Reverse: 50 mm</li> </ul>
Rated Slope <sup>4</sup>		<ul style="list-style-type: none"> <li>• 6° (10.5 %)</li> </ul> <p>according to manufacturer's specifications with 136 kg payload, 10° seat angle, 20° backrest angle</p>
4	Static stability downhill, uphill, and sideways according to ISO 7176-1 = 9° (15.8 %) Dynamic stability according to ISO 7176-2 = 6° (10.5 %)	
Max. Slope with Engaged Parking Brakes		<ul style="list-style-type: none"> <li>• 19.4° (uphill)</li> <li>• 9.7° (downhill)</li> </ul>
Lateral Dynamic Stability:		
Min. Diameter for Turning in Circles at max. Speed		<ul style="list-style-type: none"> <li>• 1785 mm</li> </ul>
Stable While Turning Suddenly		<ul style="list-style-type: none"> <li>• Yes</li> </ul>
Turning Diameter	<ul style="list-style-type: none"> <li>• 1608 mm</li> </ul>	<ul style="list-style-type: none"> <li>• 1619 mm</li> </ul>
Pivot Width	<ul style="list-style-type: none"> <li>• 1195 mm</li> </ul>	<ul style="list-style-type: none"> <li>• 1220 mm</li> </ul>
Reversing Width	<ul style="list-style-type: none"> <li>• 1191 mm</li> </ul>	<ul style="list-style-type: none"> <li>• 1210 mm</li> </ul>
Required Width of Angled Corridor	<ul style="list-style-type: none"> <li>• Forward: 823 mm</li> <li>• Backward: 843 mm</li> </ul>	<ul style="list-style-type: none"> <li>• Forward: 871 mm</li> <li>• Backward: 880 mm</li> </ul>
Required Doorway Entry Depth		<ul style="list-style-type: none"> <li>• 1092 mm</li> </ul>
Required Corridor Width for Side Opening	<ul style="list-style-type: none"> <li>• 825 mm</li> </ul>	<ul style="list-style-type: none"> <li>• 1005 mm</li> </ul>
Operating Forces DLX-REM060/150/210/215		
Joystick		<ul style="list-style-type: none"> <li>• 1.9 N</li> </ul>
Power Button		<ul style="list-style-type: none"> <li>• 2.5 N</li> </ul>
Speed Dial		<ul style="list-style-type: none"> <li>• 1.2 N</li> </ul>
Horn Button		<ul style="list-style-type: none"> <li>• 2.5 N</li> </ul>

Dimensions of Power Wheelchair According to ISO 7176-5		
	Seat Width 430 mm	Seat Width 480 mm
Seat-to-floor Height <sup>5</sup>	<ul style="list-style-type: none"> <li>• 451 mm</li> <li>• 471 mm</li> </ul>	
5 Measured without seat cushion.		
Max. Total Height	<ul style="list-style-type: none"> <li>• 1350 mm</li> </ul>	
Max. Total Width (depending on Seat Width and Base Width)	<ul style="list-style-type: none"> <li>• 669 mm</li> </ul>	<ul style="list-style-type: none"> <li>• 787 mm</li> </ul>
Total Length ( with Standard Legrests)	<ul style="list-style-type: none"> <li>• 1035 mm</li> </ul>	
Stowage Length	<ul style="list-style-type: none"> <li>• 428 mm</li> </ul>	
Stowage Width	<ul style="list-style-type: none"> <li>• 543 mm</li> </ul>	<ul style="list-style-type: none"> <li>• 598 mm</li> </ul>
Stowage Height	<ul style="list-style-type: none"> <li>• 1048 mm</li> </ul>	
Ground Clearance	<ul style="list-style-type: none"> <li>• 70 mm</li> </ul>	
Seat Width	<ul style="list-style-type: none"> <li>• 430 mm</li> </ul>	<ul style="list-style-type: none"> <li>• 480 mm</li> </ul>
Seat Depth	<ul style="list-style-type: none"> <li>• 410 - 510 mm</li> </ul>	
Seat Cushion Thickness	<ul style="list-style-type: none"> <li>• 60 mm</li> </ul>	
Backrest Angle	<ul style="list-style-type: none"> <li>• 90° - 120°</li> </ul>	
Backrest Height <sup>5</sup>	<ul style="list-style-type: none"> <li>• 510 mm</li> </ul>	
5 Measured without seat cushion.		
Armrest Height	<ul style="list-style-type: none"> <li>• 240 - 315 mm</li> </ul>	
Armrest Depth <sup>6</sup>	<ul style="list-style-type: none"> <li>• 305 - 395 mm</li> </ul>	
6 Distance between backrest reference plane and most forward part of armrest assembly.		
Max. Armrest Weight (heaviest part)	<ul style="list-style-type: none"> <li>• 1.2 kg</li> </ul>	
Max. Headrest Weight	<ul style="list-style-type: none"> <li>• 1.0 kg</li> </ul>	
Seat Angle	<ul style="list-style-type: none"> <li>• 0° - 20°</li> </ul>	
Footrest and Legrests		
Standard 80° Legrest	<ul style="list-style-type: none"> <li>• Length: 245 - 370 mm</li> </ul>	
	<ul style="list-style-type: none"> <li>• Angle: +69° - +74.2°</li> </ul>	
	<ul style="list-style-type: none"> <li>• Weight: 0.8 kg</li> </ul>	

**Kerb Weight<sup>7</sup>**

7 The actual kerb weight depends on the fittings your power wheelchair has been supplied with. Every Invacare power wheelchair is weighed when leaving the works. Refer to the nameplate for the kerb weight (including batteries) measured.

- 113.6 kg

**Component Weights**

- |           |   |
|-----------|---|
| Batteries | <ul style="list-style-type: none"> <li>• 16.1 kg</li> </ul> |
|-----------|---|

**Payload**

- |                  |  |
|------------------|--|
| Max. User Weight | <ul style="list-style-type: none"> <li>• 136 kg</li> </ul> |
|------------------|--|

## 11 Service

### 11.1 Inspections Performed

It is confirmed by stamp and signature that all jobs listed in the inspection schedule of the service and repair instructions have been properly performed. The list of the inspection jobs to be performed can be found in the service manual which is available through Invacare.

<b>Delivery Inspection</b>	<b>1st Annual Inspection</b>
Stamp of authorised provider / Date / Signature	Stamp of authorised provider / Date / Signature
<b>2nd Annual Inspection</b>	<b>3rd Annual Inspection</b>
Stamp of authorised provider / Date / Signature	Stamp of authorised provider / Date / Signature
<b>4th Annual Inspection</b>	<b>5th Annual Inspection</b>
Stamp of authorised provider / Date / Signature	Stamp of authorised provider / Date / Signature

**EU Export:**

Invacare Poirier SAS  
Route de St Roch  
F-37230 Fondettes  
Phone: (33) (0) 2 47 62 69 80  
serviceclient\_export@invacare.com  
www.invacare.eu.com

Invacare GmbH  
Am Achener Hof 8  
D-88316 Isny  
Germany

Invacare UK Operations Limited  
Unit 4, Pencoed Technology Park, Pencoed  
Bridgend CF35 5AQ  
UK

1682496-A

2026-02-02

**Yes, you can.**

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1. 移乗面®と電動車椅子座席©の間の隙間を、移乗を行うために必要な最小距離®まで縮小します。これは係員が行う必要があるかもしれません。