


Invacare | **Action⁵**




Did You Know...? The *Invacare Action⁵* has not only one but two patents pending, making it the most stable, adjustable and compact chair in the *Invacare Action* range.



Unique folding mechanism

The **Action⁵** has a patent pending, self-locking, H-shaped folding mechanism providing an greatly improved level of stability against other folding chairs in the marketplace, leading to a smoother driving experience for the client. This level of support makes the wheelchair almost effortless to propel and ensures a secure seating position.



This innovative, simple folding mechanism makes the **Action⁵** compact when folded, allowing for easy handling, storage and transportation.



Only 295mm wide!

Unique adjustable backrest

The **Action⁵** also includes a patent pending backrest system that has been specially designed to allow the user to remain seated whilst backrest adjustments are made to the chair.




With a single Allen Key, the step-less adjustments can be made easily and smoothly, allowing the user to remain in the comfort of their own chair



To make it efficient, we have also added a scale to the backrest attachment to allow for a more accurate angle to be set

This high level of detail in the **Action⁵** allows the wheelchair to combine innovative design with both sophistication and style. The requirements of the user is at the heart of everything we do.

The perfect fit
for your lifestyle



Invacare International Sàrl
Route de Cité-Ouest 2 - 1196 Gland
Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
hgeurope@invacare.com
www.invacare.eu.com



Yes, you can.®