



Product Service

# CERTIFICATE

No. Z1 14 12 52355 019

**Holder of Certificate:** **Küschall AG**  
 Benkenstrasse 260  
 4108 Witterswil  
 SWITZERLAND



**Certification Mark:**



**Product:** **Standard wheelchairs**  
**Manual wheelchairs**

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition the certification holder must not transfer the certificate to third parties. See also notes overleaf.

**Test report no.:** 713045990

**Valid until:** 2019-12-02



**Date,** 2014-12-09

( Torsten Zimmer )

Page 1 of 2



Product Service

**CERTIFICATE**  
**No. Z1 14 12 52355 019**

**Model(s):** **Küschall K-Series attract**

**Parameters:**

total length incl. footrests:	760 mm (90° Frame), 830 mm (75° Frame)
width:	530 - 650 mm
height:	715 - 960 mm
seat width:	360, 380, 400, 420, 440, 460 and 480 mm
seat depth:	375 - 475 mm
seat height:	480 - 510 mm
backrest height:	315 - 480 mm
max. user weight:	120 kg
brake system:	knee lever brake
back wheels:	
- size	24 and 25 x 1"
- tyres	air filled
- adjustable	no
- fixation	quick release
front wheels:	
- size	Ø100, Ø125 x 32 mm
- tyres	puncture proof
- adjustable	no
- fixation	rigid
frame material:	aluminum
armrest:	not applicable
legrest:	height adjustable
anti tip device:	available

**Tested according to:** EN 12183:2014

**Production Facility(ies):** 52355