

**rea**<sup>®</sup>

# Azalea<sup>®</sup> Minor



Designed for small adults  
and teenagers

**INVACARE**<sup>®</sup>  
Yes, you can.<sup>®</sup>



The **Rea Azalea Minor** is designed for users that require a seat with a smaller depth and width. Adapted from the **Rea Azalea**, the **Azalea Minor** boasts all the advantages of a reliable, tilt-able wheelchair with its special weight-shifting mechanism. This model is ideal for all users that require a solid, stable wheelchair. Available in five frame colours and with various upholsteries. A tension adjustable backrest is optional.



### ◀ Tension adjustable backrest

Adjustable in height and width.



### ◀ Table tray

It is transparent and the rail offers gripping opportunities for posture solutions.



### ◀ Shorter Chair

The chassis is 30 mm shorter, the seat depth is 50 mm shorter and the backrest 100 mm lower than on standard models.

## Rea Azalea compatibility continues

**Rea Azalea Minor** designed for users that need a smaller-sized seat, such as petite adults or teenagers. Key features of **Azalea Minor** is the 'Tilt in Space' weight-shifting mechanism, to ensure a stable tilt and makes maneuvering easier. The width and depth of the **Azalea Minor** can be adjusted to suit individual user needs.

## Tension adjustable backrest

The tension adjustable backrest offers individual shaping with the choice of backrest cushions or lateral cover.

For other Rea Azalea accessories, please consult the Rea Azalea brochure or visit [www.invacarearea.com](http://www.invacarearea.com) to see videos of them

## Technical data ▼


For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

|                     | Seat width ▼                 | Seat depth ▼   | Seat height ▼                 | Backrest height ▼ | Armrest height ▼   | Length of footrests ▼ | Seat angle ▼   | Backrest angle ▼                         |
|---------------------|------------------------------|----------------|-------------------------------|-------------------|--------------------|-----------------------|----------------|--|
| <b>Azalea Minor</b> | 340 - 440 mm                 | 380 - 450 mm   | 400/450 mm                    | 550 - 650 mm      | 240 - 360 mm       | 330 - 627 mm          | +1° to 25°     | 0° - 30°                                 |
|                     |                              |                |                               |                   |                    |                       |                |  |
|                     | Driving unit width ▼         | Total height ▼ | Total length inc. Footrests ▼ | Total weight ▼    | Max. User weight ▼ | Transport weight ▼    | Crash tested ▼ | Maximum slope for use of parking brake ▼ |
| <b>Azalea Minor</b> | SW + 220 mm<br>(SW + 230 mm) | 900 - 1250 mm  | 900 - 970 mm                  | 32 kg             | Max 75 kg          | 14,5 kg <sup>1</sup>  | ISO 7176-19    | 7°                                       |

## Frame colour ▼



## Upholstery ▼



1. Transport weight increases with 4,2 Kg when Power tilt & Recline option is used

Tel: +41 61 487 70 70  
hqeurope@invacare.com  
www.invacare.eu.com

Invacare International GmbH  
Benkenstrasse 260  
4108 Witterswil  
Switzerland

Rea, Rea design and DSS (Dual Stability System) design are registered trademarks of Invacare International. The Rea Azalea is designed with DSS (Dual Stability system). It provides the user with enhanced stability and safety when being tilted or the backrest is reclined.



Yes, you can.®