

Dragon[®] Vertic



Access all areas – with the push of
a button





The *Invacare Dragon Vertic* is an agile power chair that not only provides mobility to an active customer, but will bring you face to face with other people. Its comfortable seating system easily moves you into a standing position by a simple joystick command, providing increased access and excellent mobility.

Easy transfer

The knee fixation can be swivelled left or right, or be taken away entirely for easier transfer.

The transfer is also facilitated by swivelling-up armrests and foot plates.



Mobility when upright

The *Invacare Dragon Vertic* greatly enhances activity in your daily life with the vertical function. With mini front castors for extra stability you can even drive whilst in the upright position giving you increased access and mobility. The process of going upright is always from the seating position (driving) into vertical. The speed for moving up and down can be controlled via the joystick.

Comfortable positioning

Seat depth, legrest length, backrest tension and height can be adapted to the needs of the customer, and the angle of the backrest can be electrically adjusted for optimum comfort when driving and when upright.

Features and options



The armrest can be turned upward for an easier transfer.

The headrest is individually adjustable in width, height and horizontal position.

The sliding backrest minimises the shearing forces and is height adjustable.

Features and options



The hip support is width adjustable and removable for transfer.



The adjustable back rest angle optimises comfort in seating and upright position.



The foot plates can be swivelled up for easy transfer.



Puncture proof rear tyres available.



Mini castors allow driving in the upright position.



Tray

Convenient tray for work and for trunk support, can be swivelled down.



Knee support

The knee support is adjustable in depth, width and height it can be swivelled left and right, or can be taken off entirely



Shark controls

Easy controlling with Shark. Variable speed for vertical function and for backrest angle adjustment.



On the road version

Front and rear lights provide safety when driving on public roads. The light is activated via the Shark controller.



Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Seat width	Seat depth	Seat height	Backrest height	Armrest height	Length of footrests	Total height
Dragon® Vertic	385 - 460 mm 450 mm incl. arm pads	380 - 460 mm	500 / 550 mm	600 / 650 mm	260 - 350 mm	310 - 450 mm	1670 / 1720 mm

	Backrest angle	Driving unit width	Total height	Total length inc. Footrests	Total length without footrests	Total weight	Max. User weight
Dragon® Vertic	6° - 22°	610 mm ¹	1160 / 1210 mm	1050 mm ²	850 mm	131 kg ³	100 kg

	Battery capacity	Turning diameter	Max. Climbable obstacle height	Max. Safe slope	Motor capacity	Driving range	Speed
Dragon® Vertic	2 x 12V 60 Ah	1660 mm	40 mm	17.6%	2 x 180 W	up to 25 km ⁴	6 km/h

1. Total width is depending on seat adjustment
2. Incl. legrest and anti-tippers
3. Depending on configuration
4. Acc. to ISO 7176-4:2008

Frame colour

Blue metallic



Scan here to access the website.

Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Benkenstrasse 260
 4108 Witterswil
 Switzerland



Yes, you can.®

