

**Dragon<sup>®</sup>**



**Compact and Versatile**



The *Invacare Dragon* is a solid and cost efficient power wheelchair that provides excellent indoor and outdoor mobility in the suburban environment. It is easy to drive, and the seat can be adjusted to the user's physical requirements.

### Mobility and comfort

The robust frame makes the *Invacare Dragon* a reliable wheelchair in daily use, whereby its small dimensions provide good manoeuvrability. The spring suspension of the motors absorbs shocks for increased comfort, and the large batteries provide a longer range. Driving with the **Shark** controller is powerful and intuitive.

### Adjustability and easy service

With the smart fixation seat height and tilt angle can be adjusted for comfortable seating, without taking the seat apart. The electronic modules and the batteries are easily accessible, making service simple. Powered seat tilt and powered backrest recliner are optionally available.

### Features and options



#### Tilt and seat height

With this new feature seat height and tilt can be easily adapted to the user needs.

#### Contour seat

Contoured upholstered seat for improved comfort.

#### Shark control

The version of the **Shark** remote offers the possibility to control lights with indicators and up to two electric adjustments.



### Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Seat width	Seat depth	Seat height	Backrest height	Armrest height	Length of footrests	Seat angle	Backrest angle
<b>Standard Seat</b>	440 - 490 mm <sup>1</sup>	410, 460, 510 mm	425 - 500 mm	480 / 540 mm	290 - 360 mm	300 - 430 mm	up to 22° <sup>2</sup>	-10° - 30° <sup>3</sup>
<b>Contour Seat</b>	440 - 490 mm <sup>1</sup>	410, 460, 510 mm	425 - 500 mm	560 mm	290 - 360 mm	300 - 430 mm	up to 22° <sup>2</sup>	-10° - 30° <sup>3</sup>
<b>Laguna Seat</b>	440 - 490 mm <sup>1</sup>	410, 460, 510 mm	425 - 500 mm	560 mm	290 - 360 mm	300 - 430 mm	up to 22° <sup>2</sup>	-10° - 30° <sup>3</sup>

	Driving unit width	Total height	Total length inc. Footrests	Total length without footrests	Total weight	Max. User weight	Battery capacity	Turning diameter
<b>Standard Seat</b>	610 mm <sup>4</sup>	980 / 1040 mm	1030 mm <sup>5</sup>	850 mm	100 kg <sup>8</sup>	125 kg	2 x 12 V 60 Ah	1660 mm
<b>Contour Seat</b>	610 mm <sup>4</sup>	1060 mm	1030 mm <sup>5</sup>	850 mm	105 kg <sup>8</sup>	125 kg		
<b>Laguna Seat</b>	610 mm <sup>4</sup>	1060 mm	1030 mm <sup>5</sup>	850 mm	105 kg <sup>8</sup>	125 kg		

	Max. Climbable obstacle height	Max. Safe slope	Motor capacity	Driving range	Speed
	60 mm <sup>6</sup>	15% (8.5°)	2 x 180 W	up to 30 km <sup>7</sup>	6 km/h

1. Seat cushion 430 mm
2. Varies with seat height
3. In 6 steps, 5° - 30° with gas spring, option
4. Depending on seat adjustment
5. Incl. leg rest and antitippers
6. 100 mm with kerb climber
7. Acc. to ISO 7176-4:2008
8. Total weight is depending on configuration




### Frame colour

Blue Metallic



QR Code



Tel: +41 61 487 70 70  
[hqeuropa@invacare.com](mailto:hqeuropa@invacare.com)  
[www.invacare.eu.com](http://www.invacare.eu.com)

Invacare International GmbH  
 Benkenstrasse 260  
 4108 Witterswil  
 Switzerland



Yes, you can.®

