



*Invacare*<sup>®</sup> | **Comet**<sup>®</sup> *HD*

Well-proportioned,  
comfortable, supporting

Trusted Solutions,  
Passionate People



The *Invacare Comet<sup>HD</sup>* is part of the *Invacare Comet* family, renowned for its safe, fast and enjoyable drive. This rugged outdoor vehicle is ideal for users weighing between 150 kg to 220 kg. Packed with safety features, this powerful scooter will take you to your destination in the most direct and secure way. The large contoured extra padded seat together with four individual adjustable suspension units ensures an exceptional driving experience. The *Comet<sup>HD</sup>* has been designed to allow for additional weight capacity.

## Durable design and power in one

Even for users with a weight of up to 220 kg, safety is never a concern with the *Comet<sup>HD</sup>*. All necessary precautions for guaranteeing the stability and safety of the scooter have been taken into consideration during the design process.

## Features and options



### Unique suspension system

It has proven itself in the professional mountain biking industry and allows for simple independent user adjustment.



### Reinforced seat column

It ensures stability and safety while seated.



### Reinforced armrest

It allows for the additional stress of getting on and off the *Comet<sup>HD</sup>*.

## That comfortable feeling

The XL contoured seat has been developed in conjunction with therapists to ensure optimum seating comfort. The seat base and height- and width-adjustable armrests have been reinforced to allow for additional stresses and can easily be used as supports when getting on and off the scooter. Seating comfort can be increased by using the pneumatically adjustable lumbar support.

## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local *Invacare* website.

Backrest width	Seat depth	Seat to floor height	Backrest height	Driving unit width	Total length	Empty weight	The heaviest part
<i>Invacare Comet<sup>HD</sup></i> 660 mm	530 mm	440 - 515 mm	610 mm	660 mm	1475 mm	148 kg	96 kg
Max. User weight	Battery capacity	Turning diameter	Max. Climbable obstacle height	Ground clearance	Max. Safe slope	Max. Driving distance	Max. Driving speed
<i>Invacare Comet<sup>HD</sup></i> 220 kg	2 x 75 Ah	3260 mm	100 mm	135 mm	10°	50 <sup>1</sup> km	10 km/h
Motor capacity			Colour of shrouds		Upholstery		
	Weight capacity 220 kg	1. 10 km/h 2. Peak performance	Slate Grey		Black Vinyl		
<i>Invacare Comet<sup>HD</sup></i> 650 W / 1800 W <sup>2</sup>							

Tel: +41 (0)22 354 60 10  
 Fax: +41 (0)22 354 60 11  
 hqeuropa@invacare.com  
 www.invacare.eu.com

Invacare International Sàrl  
 Route de Cité-Ouest 2 - 1196 Gland  
 Switzerland

**INVACARE**  
 Yes, you can.®

