



INVACARE® TOP END®  
**T-5 7000 SERIES**  
TENNIS WHEELCHAIR



The **Invacare® Top End® T-5 7000 Series Tennis Wheelchair** is faster, lighter, stronger and better than ever. In 2011 it hit the market with new material, an enhanced headtube design and an axle securing system, reducing overall weight by two pounds and improving performance. Each Top End T-5 7000 Series Tennis Wheelchair is custom built to an individual athlete's specifications and is overseen by Paul Schulte, the design engineer for Top End who is a 3 time Paralympian having been a member of the USA National Men's Wheelchair Basketball Team.

## Features

Custom-built fixed frame in 7005 aluminium

Adjustable height fixed swivel anti-tip

58mm fifth wheel

Fixed back with fixed or adjustable height

Adjustable tension seat and back upholstery

Fixed axle position with fixed wheelbase

Fixed camber: 15°, 18° or 20°

High performance Sun Components® wheels:

24", 25", 26" or 700C

Four bed offensive wing

Anodized pushrims

High pressure clincher tires

Precision fork caster system

72mm casters

Adjustable height tubular footrest


Powder coat paint

Rear frame protector

Axle securing system

250 lb. weight capacity

## Model #: T57



1. Custom front frame angle and tubular footrest
2. Adjustable height fixed swivel anti-tip
3. 7005 aluminium heat treating process
4. Precision fork caster system
5. Optional welded fixed tube sideguards
6. Adjustable center-of-gravity option



## Options

Spinergy® Spox or SLX wheels

Sun Components Fusion wheels

Adjustable center-of-gravity

Welded seat rigidizer bars with solid seat

Telescoping removable anti-tip bar

Click straps

Bolt-on fixed side guards

Fixed side guards with welded tube

Fixed aluminium or carbon fiber side guards

Titanium axles

Titanium pushrims

Coated pushrims

Footrest covers

Frame protectors

Custom strapping/toe clips

Alignment gauge

Adjustable angle footplate/toe clips

Recessed platform footplate with toe

stops/footrest strap

Kenda® tires

Close mount/cut pushrim tab

**Invacare Corporation**  
www.invacare.com

## USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

## Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

## Top End

www.topendwheelchair.com  
Sports and Recreation  
Products  
4501 63rd Circle North  
Pinellas Park, Florida  
33781-5914  
(727) 522-8677  
(800) 532-8677  
Fax: (727) 522-1007

©2016

Invacare Corporation.

All rights reserved.

Trademarks are identified

by the symbols ™ and ®.

All trademarks are owned

by or licensed to Invacare

Corporation unless

otherwise noted. Spinergy

is a registered trademark

of Spinergy Inc. Sun and

Equalizer are a registered

trademark of Sun

Components. Kenda is a

registered trademark of

Kenda Rubber Industrial

Co., Ltd. Specifications are

subject to change without

notification.

Form No. 12-508

120930



Yes, you can.®