


Porto



Comfort, variety, functionality
and design



Yes, you can.

The Invacare **Porto** lift-out chair comes with a renewal design, offering different functionalities and fabric choices - something meet every home need and aesthetic perspectives.

Comfort is essential, that is why the **Porto** chair offers several seating positions from upright to reclined, available with a manual version or with one motor.

Features and Options



Relax position

Different relax positions will allow the same comfort in any recliner angle chosen.



Normal position

Different fabrics will allow a perfect fit in all environments.



The **Porto** chair provides a high level of comfort thanks to the continuing of contours between the leg rest, the seat and the lumbar support which maximizes good positioning. **Porto** has a contemporary design and can be easily moved from one place to another, thanks to the small back castors which are simply dismounted by taking out a screw.

- One motor
- Easy-to-use remote control
- Side pocket
- Easy to use and mount

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Seat depth	Seat height	Seat width	Full back seat tilt	Back seat height	Chair width	Depth in upright position	Depth in relaxed position
Invacare Porto	540 mm	480 mm	580 mm	120°	1080 mm	860 mm	850 mm	1580 mm
	Floor to armrests height	Max. user weight	Total product weight					
Invacare Porto	650 mm	125 kg	57 kg					
				Invacare Porto One motor Relaxed position Reclined pos. Lift-out pos. Manual Relaxed position Relaxed pos. Lift-out pos.				

Foam density: Seat 30 Kg/m³, Back 30 kg/m³, Armrest 25 kg/m³

Colours



Spares parts available

Motor Powerlift "PORTO"	1577637
Electric transformer + Cable "PORTO"	1577638
Remote Control W/ Cable (2 Buttons)	1577639
Battery + Cable Porto	1577640

Battery with cable is a spare that allows to use the motor version disconnected from the electricity with same features and is very useful when there is a need to move the chair for a space that is far away of power.

Tel: +44 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Neuhofweg 51
 4147 Aesch BL
 Switzerland



© 2023 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. EU - 07/23

Yes, you can.®