

# Accent



In tune with your needs



With safety and functionality at its core, the **Accent** bed from Invacare offers exceptional value for money. Providing real comfort for the client along with an appealing design, this bed is ideally suited to any environment.

The **Accent** bed is an economical choice that's perfectly in tune with your needs. It meets the highest safety standards, is fully electrically operated, and the design of the mattress platform ensures maximum comfort. With stylish wooden bed ends in natural beech, it provides a homely feel and is perfectly suited to all home and long-term care settings.

### Reliable and safe

A highly affordable choice, the **Accent** bed is reliable and safe. It has been designed using only high-quality components to ensure functionality and durability. Like all Invacare beds, the bed and the side rails meet the strength demands of IEC 60601-2-52 safety standards for beds.

### Reduced shear and friction

The optimised dimensions of the different sections of the mattress platform are designed in accordance with anthropometric statistics and pressure mapping recommendations. Because of this, general shear and friction are reduced and the client can enjoy increased comfort when sitting. Lockable tilt is available as an option to further aid pressure relief and provide a comfortable sitting position.

## Features and Options



### Bella and Aria side rails

Bella, full-length wooden side rail in beech or Aria, full-length side rail in steel are easy to mount without the need for tools.



### Lifting pole

A lifting pole that's light and easy to adjust.













### Hand controls

Hand controls that come with or without locking functions.


## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Width outside Width inside	Length outside Length inside	Height adjustment	Mattress support dimensions	Backrest angle	Thigh angle	Tilt option	Total weight
Invacare <b>Accent</b>	 1020 mm 900 mm	 2240 mm 2080 mm	 400 - 800 mm	 800-230- 300-670 mm	 0 - 70°	 0 - 24°	 11°	 64 kg

	Max. user weight	Heaviest part of the product
Invacare <b>Accent</b>	 145 kg	 17 kg

## Frame colour and wood work



QR  
Code

Tel: +41 61 487 70 70  
hqeurope@invacare.com  
www.invacare.eu.com

Invacare International GmbH  
Benkenstrasse 260  
4108 Witterswil  
Switzerland

