

Taller, wider, bigger.. Invacare SB 755 in Wide Size

The popular **SB 755** from Invacare is now available in wider variants of 1050 mm and 1200 mm, providing even more comfort for taller or bigger clients (weights of 200 kg). These wider models provide extra space and comfort for larger clients and plenty of room to turn. The **SB 755 Wide** has been engineered with an improved auto-regression system to offer total comfort. The bed also offers an optimal seating position and is designed to reduce shear and friction.

The **SB 755 Wide** is attractive and modern; the streamlined design is easy to clean and will blend harmoniously into any homecare or institutional environment.



Trusted Solutions, Passionate People



Yes, you can.®



Features and Options

Comfort for users and carers

The **SB 755 Wide** comes with several key features designed to ensure an excellent sitting position and to reduce shear and friction during movement. As Invacare always strives to cater for all patient needs, we are delighted to now offer an integrated length extension of 100 mm at the head and foot of the bed to meet the requirements of taller users.


Mount it, move it... love it!

Mounting and dismantling the bed has been made particularly easy to facilitate transport and storage. The bed can be stored on a trolley and rolled directly from van to bedroom.



Auto regression function

Client comfort when activating backrest, minimising shear and friction.



Optimised dimensions of mattress platform according to anthropometrical measurements, for increased comfort.



Extendable mattress platform

Additional 100 mm at head and foot. Bed easily extends with 50, 100, 150 or 200 mm to suit user needs.

Technical data

SB 755 Wide 1050 mm	Width outside*: 1150 mm Width inside: 1050 mm	Length outside*: 2130 mm Length inside: 2080 mm Mattress: 2000 mm	Height adjustable: 350 - 820 mm Total height: 1600 - 260 - 680 mm 2070 mm	Mattress support dimensions: 800 - 260 - 260 - 680 mm	Under bed clearance: 170 mm	0° - 73° Backrest angle	0° - 28° Thigh angle
SB 755 Wide 1200 mm	Width outside*: 1300 mm Width inside: 1200 mm	Length outside*: 2130 mm Length inside: 2080 mm Mattress: 2000 mm	Height adjustable: 350 - 820 mm Total height: 1600 - 260 - 680 mm 2070 mm	Mattress support dimensions: 800 - 260 - 260 - 680 mm	Under bed clearance: 170 mm	0° - 73° Backrest angle	0° - 28° Thigh angle
SB 755 Wide	0° - 16° Minimum angle below horizontal	0° - 18° Full legrest lift	25° (Angle below frame) 14° (Angle above frame)	SWL 235 kg Max. user weight 200 kg	Washable in a washing machine		



Bariatric use

Complies with IEC 60601-2-52

NB/ The bed must not be used by patients under 12 years of age, or by patients with body size equivalent to an average 12 year old or smaller.

Invacare International Sàrl
Route de Cité-Ouest 2 - 1196 Gland
Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
hqeurope@invacare.com
www.invacare.eu.com

© 2015 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Frame colour and wood work



Aluminium grey



Beech



Cherry



Light oak



Yes, you can.®

SB755 Wide - EU - 08/2015

