



INVACARE® TOP END®
PRO™
BASKETBALL WHEELCHAIR




The **Invacare® Top End® Pro™ Basketball Wheelchair** delivers custom-built performance with adjustments that take only minutes. Perfect for competitors still dialing in their game.

Features

Custom aluminum frame
Offensive 4 bend wing
Adjustable front and rear seat heights
Anti-tip/fifth wheel: Adjustable height
Fixed back
Adjustable back height
Adjustable tension back upholstery
Adjustable Center-of-Gravity
Rear seat height adjustment range is 6" and the frame is available in a tall or short variety to accommodate every aspiring athlete

Fixed camber: 15°, 18° or 20°
High performance Sun® wheels
Anodized pushrims
High pressure clincher tires
Smart fork caster system
72mm casters
Adjustable platform footrest
Powder coat paint
250 lb. weight capacity
Quick release lever allows the player to change his or her rear seat height in just minutes

Model #: SPBB



1. Front seat height and center-of-gravity adjustments
2. Rear seat height adjustment
3. Quick release lever
4. Single swivel fixed fifth wheel
5. Precision fork caster system
6. Adjustable footrest system



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End

www.topendwheelchair.com
Sports and Recreation Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(727) 522-8677
(800) 532-8677
Fax: (727) 522-1007

©2014

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Spinerger is a registered trademark of Spinerger Inc. Sun and Equalizer are registered trademarks of Sun Components. Kenda is a registered trademark of Kenda Rubber Industrial Co., Ltd. Specifications are subject to change without notification.
Form No. 04-188
Rev. 12/14
141156

Options

Clamp-on fixed aluminum side guards
Spinerger® SLX wheels
Sun® Equalizer™ wheels
12 tab pushrims
Kenda® tires
O-Ring positioning belt
Close mount/cut pushrim tab

Click straps
Spoke guards
Frame protectors
Alignment gauge
Double swivel anti-tip bar
Cushion: soft or hard
Folding back

