

AQUATEC®

Ocean

Series



Shower commode chairs



Aquatec Ocean



- Stable, stainless steel frame
- Individually adjustable seat height: no tools required
- Machine-washable adjustable backrest with quick-release buckles, can be adjusted to the user
- Space-saving push bar to enable optimal positioning over the toilet
- Ergonomic, swing-away armrests can be adjusted to two different height levels
- Height-adjustable, swing-away and detachable footrests for simple transfer
- Fold-away, textured footplates with heel loops provide a safe position

The **Ocean** is available as follows:

- Standard version with castors
- Self-propel model with 24" wheels, quick-release-axles and four castors
- XL version (180 kg load capacity and more space between the armrests) with castors
- All models can be retrofitted with 24" wheels or to the XL version

Features and options

Ocean 24" wheels

Self-propel model with 24" wheels, quick-release axles and four castors.



Ocean XL

More space between the armrests with a load capacity of up to 180 kg.



Individually adjustable seat height: no tools required.



Easy disassembly for mobile use, e.g. on travels.

Aquatec **Ocean Vip**



- Stable, stainless steel frame
- Seat tilt can be individually adjusted (tilt angle 0° - 35°)
- Stainless steel gas spring
- Extremely stable even at maximum tilt angle
- Seat height adjustable in 3 steps, no tools required
- Very comfortable, ergonomically contoured seat
- Machine-washable adjustable backrest with quick-release buckles – can be easily adjusted to the user

The **Ocean Vip** is available as follows:

- Standard version with castors
- Standard version with castors and toilet pan holder
- Standard version with castors, toilet pan holder and soft seat
- **XL** version (more space between the armrests) with castors
- The **Ocean Vip** can be retrofitted with XL armrests

Features and options

Seat tilt

Tilt angle 0° - 35°.



Headrest

Height and depth-adjustable headrests included.



Directional lock

Easy to operate brakes on 3 castors, 1 castor with directional lock to help push the chair down straight floors.



Adjustable seat height

Adjustable in 3 steps, no tools required.

Aquatec ***Ocean Dual-Vip***



- Seat tilt can be individually adjusted (tilt angle 0° - 35°)
- Two additional stainless steel gas springs to lower the backrest (recline angle 10° - 50°)
- Second operating lever with safety lock to lower the backrest
- Particularly stable even at maximum tilt angle, anti-tippers included
- Machine-washable adjustable backrest with quick-release buckles – can be adjusted to the user

The ***Ocean Dual Vip*** is available as follows:

- Standard version with castors
- The ***Ocean Vip*** can be retrofitted with reclining backrests for the ***Ocean Dual Vip*** and be fitted with anti-tippers
- The ***Ocean Dual Vip*** can be retrofitted with XL armrests

Features and options



Seat surface 0°, backrest 25°



Full seat tilt 35°, back rest 50°



Anti-tippers included

Aquatec **Ocean E-Vip**



- Seat height can be steplessly adjusted between 500 and 1050 mm electrically
- Seat tilt can be electrically tilted forwards and backwards (tilt angle -5° to 35°)
- Particularly stable even at maximum tilt angle and maximum seat height
- Clearly laid out hand control for simple adjustment of the seat height and seat tilt (4 buttons only)
- Sealed electronic components: ideal for use in the shower
- Can be cleaned in automatic cleaning stations at temperatures up to 85 °C!
- Detachable battery pack can be recharged easily

The **Ocean E-Vip** is available as follows:

- Standard version with castors
- The **Ocean E-Vip** can be retrofitted with XL armrests

Features and options



Seat surface 0°, backrest 25°



Detachable battery pack



Clearly laid out hand control

Accessories *Aquatec Ocean*



Soft seat

With key-shaped cut out.
16341 for *Ocean* / *Ocean XL* / *Ocean E-Vip*
1470075 for *Ocean Vip* / *Ocean Vip XL* / *Ocean Dual Vip*



Soft seat insert

For soft seat with key-shaped cut out.
16343 for all models of the *Ocean* series



Soft seat

With small, oval opening.
16342 for all models of the *Ocean* series



Splash guard

Simply snap-on, for soft seat with key-shaped cut out.
16511 for all models of the *Ocean* series



Universal soft seat

Closed seat surface.
1526032 for all models of the *Ocean* series



Ergonomic special soft seat

Extra soft, textured surface
1531259 soft for *Ocean*, *Ocean XL*, *Ocean Vip*, *Ocean Vip XL*, *Ocean Dual Vip*



Variable soft seat

Hygiene cut out can be placed in all 4 directions. Easy and tool-less assembly.
1544469 for *Ocean* / *Ocean Vip* / *Ocean Dual Vip*



Safety bar

For additional safety. The surface with soft padding. Easy to assemble.
10167-10 for all models of the *Ocean* series. Length: 540 mm
1470078 XL safety bar for all XL versions of the *Ocean* series. Length: 620 mm



Safety belts

For safe positioning. Individually adjustable, high load capacity and easy to assemble. Machine-washable. For all models of the *Ocean* series.
1470081 Chest belt, grey
1470082 Pelvic belt, grey



Armrest locking device

To lock the armrest.
1470077 for all models of the *Ocean* series



XL back

Increases the distance between the armrests by 80 mm. Retrofit kit for *Ocean*.
1534979 until SN14551780
1558847 from SN14551781



Backrest

Solid backrest with ergonomic shape, easy to clean and dry.
1575299 for *Ocean* / *Ocean Vip* / *Ocean Dual Vip* and *Ocean E-Vip*

Accessories *Aquatec Ocean*



Soft backrest cushion

For additional stability and comfort - also for small users and young adults

1558769 for *Ocean / Ocean XL / Ocean Vip / Ocean Vip XL / Ocean Dual Vip* and *Ocean E-Vip*



Wider XL armrests

Increases the distance between the armrests by 40 mm on each side.

Retrofit kit for *Ocean Vip, Ocean Dual Vip* and *Ocean E-Vip*
1470086 Wider XL armrests (pair)



24" wheels (set)

Incl. brakes and removable axles.

16379 for *Ocean*



Calf-and amputee support

Padded, both height and angle adjustable.

1556894 left for *Ocean / Ocean XL / Ocean Vip / Ocean Vip XL / Ocean Dual Vip* and *Ocean E-Vip*
1556893 right for *Ocean / Ocean XL / Ocean Vip / Ocean Vip XL / Ocean Dual Vip* and *Ocean E-Vip*



Calf support with foot support

Padded, both height and angle adjustable with additional foot support.

1546705 left for *Ocean / Ocean XL / Ocean Vip / Ocean Vip XL / Ocean Dual Vip* and *Ocean E-Vip*
1546439 right for *Ocean / Ocean XL / Ocean Vip / Ocean Vip XL / Ocean Dual Vip* and *Ocean E-Vip*



Lateral

For safe positioning. Width and height adjustable - also for small users and young adults.

1535077 for *Ocean / Ocean XL / Ocean Vip / Ocean Vip XL / Ocean Dual Vip* and *Ocean E-Vip*



Ergonomic head rest holder

Safe and comfort with various adjustment positions (pad not included).

1531613 for *Ocean Vip / Ocean Dual Vip* and *Ocean E-Vip*



Head/Neckrest (neckrest shown)

1531465 Headrest for *Ocean Vip, Ocean Vip XL, Ocean Dual Vip* and *Ocean E-Vip*
1531608 Neckrest for *Ocean Vip, Ocean Vip XL, Ocean Dual Vip* and *Ocean E-Vip*



Front anti tipper

For additional safety

1535936 for *Ocean Vip / Ocean Vip XL / Ocean Dual Vip*



Toilet pan

With handle and cover.
10230 for all models of *Ocean* series. Length: 460 mm
1541289 lid for toilet pan

Toilet pan holder (not shown)

1527985 for *Ocean / Ocean XL*
1470076 for *Ocean Vip / Ocean Vip XL / Ocean Dual Vip*
1471233 for *Ocean E-Vip*


Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.


Order numbers	Aquatec Ocean	Aquatec Ocean 24" wheels	Aquatec Ocean XL	Aquatec Ocean Vip with tilt	Aquatec Ocean Vip XL with tilt	Aquatec Ocean Dual Vip with tilt & recline	Aquatec Ocean E-Vip
Shower commode chair with castors	1534327		1534328	1470000	1470788	1470829	1471208
Shower commode chair with 24" wheels	1534329						
Shower commode chair with castors and toilet pan holder				1470712			1522373
Shower commode chair with castors, toilet pan holder and soft seat	1470713						
Dimensions							
Seat width	480 mm (A)	480 mm (A)	480 mm (A)	500 mm (A)	500 mm (A)	500 mm (A)	480 mm (A)
Seat depth	450 mm (B)	450 mm (B)	450 mm (B)	450 mm (B)	450 mm (B)	450 mm (B)	450 mm (B)
450 mm (B)	475 - 600 mm (C)	475 - 600 mm (C)	475 - 600 mm (C)				
Seat height	-	-	-	510 - 610 mm (C)	510 - 610 mm (C)	510 - 610 mm (C)	500 - 1050 mm (C)
Total width	560 mm (D)	675 mm (D)	640 mm (D)	595 mm (D)	625 mm (D)	625 mm (D)	730 mm (D)
Total depth	900 mm (E)	1060 mm (E)	900 mm (E)	1060 mm (E)	1060 mm (E)	1060 mm (E)	1083 mm (E)
Total height without adjustable headrest	946 - 1094 mm (F)	946 - 1094 mm (F)	946 - 1094 mm (F)	1085 - 1185 mm (F)	1085 - 1185 mm (F)	1085 - 1185 mm (F)	1070 - 1620 mm (F)
Depth without footrests	700 mm (G)	855 mm (G)	700 mm (G)	845 mm (G)	845 mm (G)	845 mm (G)	1010 mm (G)
Footplate depth	155 mm (H)	155 mm (H)	155 mm (H)	160 mm (H)	160 mm (H)	160 mm (H)	155 mm (H)
Height to armrest	680 - 850 mm (I)	680 - 850 mm (I)	680 - 850 mm (I)	720 - 860 mm (I)	720 - 860 mm (I)	720 - 860 mm (I)	730 - 1320 mm (I)
Height, footrests to seat	255 - 495 mm (J)	255 - 495 mm (J)	255 - 495 mm (J)	345 - 465 mm (J)	345 - 465 mm (J)	345 - 465 mm (J)	345 - 465 mm (J)
Height, ground to lower seat edge:							
with holder for toilet pan	390 - 520 mm (K)	390 - 520 mm (K)	390 - 520 mm (K)	420 - 520 mm (K)	420 - 520 mm (K)	420 - 520 mm (K)	420 - 970 mm (K)
without holder for toilet pan	410 - 540 mm (K)	410 - 540 mm (K)	410 - 540 mm (K)	440 - 540 mm (K)	440 - 540 mm (K)	440 - 540 mm (K)	465 - 1015 mm (K)
Height, seat to armrest	210 / 250 mm (L)	210 / 250 mm (L)	210 / 250 mm (L)	215 / 245 mm (L)	215 / 245 mm (L)	215 / 245 mm (L)	230 - 270 mm (L)
Width between the armrests	460 mm (M)	460 mm (M)	540 mm (M)	475 mm (M)	555 mm (M)	475 mm (M)	470 mm (M)
Width, cut out in seat plate	205 mm	205 mm	205 mm	205 mm	205 mm	205 mm	205 mm
Depth, cut out in seat plate	320 mm	320 mm	320 mm	320 mm	320 mm	320 mm	320 mm
Tilt angle seat	-	-	-	0° - 35°	0° - 35°	0° - 35°	-
Tilt angle back	-	-	-	-	-	10° - 50°	-
Graduation, height adjustment	25 mm	25 mm	25 mm	50 mm	50 mm	50 mm	stepless
Ø Castors	5" (127 mm)	5" (127 mm)	5" (127 mm)	5" (127 mm)	5" (127 mm)	5" (127 mm)	5" (127 mm)
Ø Drive wheels	-	24" (610 mm)	-				
Width between the 24" wheels	-	540 mm (N)	-	-	-	-	-
More info							
Weight	14.5 kg	19.5 kg	14.5 kg	23.5 kg	23.5 kg	28.3 kg	39 kg
Load capacity	150 kg	150 kg	180 kg	150 kg	150 kg	150 kg	150 kg
Colour	White	White	White	White	White	White	White

Heel loops are included in the scope of delivery. The toilet pan holder and toilet pan are not.


Ocean with castors




Ocean with 24" wheels




Ocean XL



Ocean Vip / Vip XL / Dual Vip



Ocean E-Vip



Invacare® H263 / H273 Lima



- Stable, corrosion-free aluminium frame
- Easy to operate folding mechanism enables space-saving storage and easy transport
- Variable, padded seat with hygiene cut out can be used in four orientations and thus enables individual personal hygiene
- Swing-away armrests and height-adjustable, swivelling and detachable footrests for easy transfer
- Double stop lock on all 4 castors for reliable stopping and fixing
- Soft back of quick-drying material offers a stable seating position
- Very good positioning over the toilet

The **Lima** is available as follows:

- Standard version with castors
- Self-propel model with 24" wheels
- In both models, the toilet pan holder, toilet pan, padded seat including insert and pelvic belt are included in the scope of delivery

Features and options



H273 / Lima with 24" wheels



Foldable

Variable seat with
hygiene cut out




Pelvic belt included

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Order numbers	Invacare® H263 Lima	Invacare H273 Lima 24" wheels
Shower commode chair with castors	1525732	
Shower commode chair with 22" wheels		1525733
Dimensions		
Seat width (A)	460 mm	460 mm
Seat depth (B)	460 mm	460 mm
Seat height (C)	580 mm	580 mm
Total width (D)	570 mm	680 mm
Total depth (E)	1050 mm	1120 mm
Total height (F)	1000 mm	1000 mm
Total width folded	225 mm	335 mm
Total depth folded	1050 mm	1120 mm
Total height folded	1000 mm	1000 mm
Depth without footrests	760 mm	760 mm
Footplate depth (H)	165 mm	165 mm
Height to armrests	780 mm	780 mm
Height, footrests to seat	325 - 460 mm	325 - 460 mm
Height, ground to lower seat edge: with holder for toilet pan	340 mm	340 mm
without holder for toilet pan	535 mm	535 mm
Height, seat to armrest	210 mm	210 mm
Width between the armrests	475 mm	475 mm
Width, cut out in seat plate	190 mm	190 mm
Depth, cut out in seat plate	360 mm	360 mm
Ø Castors	5" (127 mm)	5" (127 mm)
Ø Drive wheels	-	24" (610 mm)
Width between the 24" wheels	-	550 mm
More info		
Weight	14.0 kg	18.5 kg
Load capacity	135 kg	135 kg
Colour	Dark blue - white	Dark blue - white



Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Benkenstrasse 260
4108 Witterswil
Switzerland

QR
Code



© 2018 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. Ocean Shower chair - EU - 10/2018