

Actio/2 P550/2



Neat and stylish



Invacare **Actio/2** is a neat and stylish rollator. The slim design facilitates accessing narrow spaces. The padded and decorated seat offers great comfort. Enjoy indoor safety with new **Actio/2** suited for an active life-style.

Extra comfort

Enjoy the comfort of the padded and decorated seat. **Actio/2** can also be equipped with an accessory cane holder and soft basket.

Slim design

The frame is narrow to allow you to access all doorways and narrow passes.

Features and Options



Easy folding



Comfortable seat



Cane holder

Makes the cane easily accessible without getting in the way.











Soft basket

Allows to easily transport items. Max. load 2.5kg.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please refer to your local Invacare website.

Width between outside wheels	Depth	Width between wheels	Overall height	Seat height	Width between handles	Wheel diameter	Load capacity	
								
Invacare Actio/2	530 mm	630 mm	380 mm	740 - 1000 mm	550 mm	410 mm	110 mm	5 kg

Max. user weight

Invacare Actio/2
130 kg

Colour

RAL 9006



Art. nr.

Singel pack with seat **1513721**
 Package of 4 with seat **1521861**

QR Code

Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Benkenstrasse 260
 4108 Witterswil
 Switzerland

