

AQUATEC®

Pure Bidet



Modern, luxurious and hygienic



Yes, you can.®

The *Aquatec Pure Bidet* is ideal for those who experience difficulty with personal hygiene, or for those simply seeking added luxury at home. The modern looking bidet fits easily and discreetly to existing toilets and provides a gentle yet effective cleansing solution using a warm water spray. Key features include:

- High functionality
- Inbuilt technology
- Fits to existing toilet

High functionality

- Intuitive hand control with large buttons
- Adjustable water temperature with aerated spray and variable jet settings
- Personal spray head positioning for both posterior and feminine settings
- Adjustable warm air dryer
- Temperature controlled heated seat
- LED night light
- Soft close seat and lid



Inbuilt technology

- Continuous water heating for instant warm water that will not run out during cleaning
- Automatic spray head cleansing before and after use to maintain cleanliness of the system
- Equipped with a thermal protector to stop overheating
- Automatic shut down for energy saving which is activated via its seat sensor
- Self-diagnosis feature to automatically detect malfunction
- Earth Leakage Breaker (ELB) for added safety, giving peace of mind



Fits to existing toilet

- Uses standard pan fixing points
- Quick and easy to assemble and detach
- DIN EN 1717 conformity to ensure safe connection to public water pipes



Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Total width	389 mm	Type of protection, entire device	IP X4
Total length	515 mm	Battery type remote control	2 x AA (LR6) / 1.5 V
Total depth	156 mm	Water supply pressure	79 - 735 kPa
Product weight	3.9 kg	Water supply temperature	3 - 35 °C
Max. user weight	125 kg	Water flow rate CLEANSING	0.7 l/min
Length of power cord	1.4 m	Water flow rate BIDET	0.85 l/min
Power supply	220 - 240 V AC, 50/60 Hz	Water pressure adjustment	3 levels
Max. power consumption	1725 - 2053 W		

Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Neuhofweg 51
 4147 Aesch BL
 Switzerland

