


# Invacare® Comet® Heavy Duty

Comfortable, supportive and well proportioned



### Key Features

- 220kg weight capacity
- Large wheels and powerful motor
- Bright lights incl. indicators and brake light
- Auto speed reduction during turning
- 4 wheel adjustable Heavy Duty suspension
- Built-in adjustable lumbar support for additional comfort

### Colours




Scooter Name	Speed Range	Max Range	Max User Weight	Scooter Weight	Length	Batteries	Colour Options
Colibri	8km/h	16km	136kg	44kg (inc battery) Heaviest part: 15,6kg	101cm	2 x 12 Ah or 2 x 18Ah	Sapphire Blue, Jasper Red, Jade Green, Sunstone Orange, Marble White
Leo (3 wheels)	8km/h	36km	136kg	78kg (inc battery) 83kg (inc battery)	120cm	2 x 36Ah	Moonstone Silver, Onyx Blue, Ruby Red
Leo (4 wheels)					122cm		
Orion (3 wheels)	10 & 12km/h 6, 10 & 12km/h	42km	136kg	103kg (inc battery) 110kg (inc battery)	127cm	2 x 50Ah	Moonstone Silver Onyx Blue
Orion (4 wheels)					129cm		
Comet (3 wheels)	15km/h 15km/h	55km	160kg	133kg (inc battery) 136kg (inc battery)	145cm	2 x 75Ah	Moonstone Silver & Onyx Blue Moonstone Silver & Onyx Blue
Comet (4 wheels)					147.5cm		
Comet Alpine	10km/h	50km	220kg	136kg (inc battery)	147.5cm	2 x 75Ah	Slate Grey
Comet HD	10 km/h				147.5cm		

Whatever your final choice, enjoy the freedom and confidence of being able to get around comfortably and securely.



Invacare Authorised Dealer



CE Distributor  
Invacare International Sàrl  
Route de Cité-Ouest 2  
1196 Gland - Switzerland  
Tel: +41 (0)22 354 60 10  
Fax: +41 (0)22 354 60 11  
hqeurope@invacare.com  
www.invacare.eu.com

© 2015 Invacare International Sàrl  
All rights reserved.  
All information quoted is believed to be correct at time of print.  
Invacare reserves the right to alter product specifications without prior consultation.  
Scooter Range - EU - 01/2015

## Invacare® Scooter Range



Safe and easy to use



Yes, you can.®

## Invacare® Colibri®

Effortless, colourful and light




reddot award 2014 winner



### Key Features

- LiteLock™ for simple and fast dismantling
- Comes apart in less than a minute
- Choice of 5 vibrant colours
- Large wheels and 18Ah battery pack available
- On/off board battery charging
- Fits into most car boots

### Colours



## Invacare® Leo®

Safe, solid and secure



### Key Features

- Compact and maneuverable
- Ergonomic steering with tiller adjustment
- Front and rear lighting includes brake light
- Disassembles into 5 parts
- 2-step disengaging lever to prevent accidental free wheeling
- Integrated splash guards to protect the electronics and motor

### Colours



## Invacare® Orion®

Stable, trustworthy and durable



### Key Features

- Large battery capacity for longer drives
- Bright lights incl. indicators and brake light
- Auto speed reduction during turning
- Adjustable suspension
- Comfortable adjustable seating
- Ergonomic intuitive steering incl. hand brake

### Colours



## Invacare® Comet® / Comet® Alpine

Rugged, robust and reliable



### Key Features

- Fast and responsive
- Large wheels and powerful motor
- Bright lights incl. indicators and brake light
- Auto speed reduction during turning
- 4 wheel adjustable suspension
- Ergonomic intuitive steering incl. hand brake

### Colours



(Comet)

(Comet Alpine)