




Legacy



The original



You want to have a reliable, good quality rollator. A rollator with more possibilities than a basic one but without the most advanced adjustments and accessories.

Dolomite **Legacy** is the original that offers you precisely this. That's why it has also been the inspiration for lots of similar rollator over the years.



Walking area

Large walking area inside the rollator frame is a signum for Dolomite. The flip-up seat gives premium walking space. Makes it easier to reach and access things. The rear forks gives a smooth inside, being friendly to toes and shoes. Structured seat with anti-slip function, rounded rear edge and generous sitting area, means that the user can sit securely and comfortably for a long time.

Safe when folding

Folding handle positioned high-up so that the user avoids having to bend down and loose balance.



Features and options



Basket

Spacious basket with handle. Max load capacity 9 kg.
Art no. 12530-15



One-hand brake

Make it possible to activate both brake with one hand.
Art no. 12505-R (right)
12505-L (left)



Brake booster

Provides extra braking force for weak hands.
Art no. 12651



Kerb-climber

"Pedal" that simplifies passing kerbs or similar. Mounted on the right side.
Art no. 12560-21

Features and options



Cane holder

Makes the cane easily accessible.
Art no. 12606



Soft grips

A softer alternative to the anatomical handle grips.
Art no. 12677 (pair)



Wheel guard

Protects clothing from being dirty from wheel splash.
Art no. 12544 (pair)



Tray

For transportation of a cup of coffee or similar. Can be moved to a table.
Art no. 70311



Hanging bag

Provides extra room when the basket is quickly filled. Can be combined with oxygen cylinder holder.
Art no. 12523



Slow-down brakes

A safety feature preventing the walker from rolling away. Perfect for walking practice and one handed solutions.
Art no. 13041-21



Oxygen cylinder holder

The oxygen cylinder is centrally positioned on the rollator and thereby provides good manoeuvrability. Fits most cylinders.
Art no. 12682-25



Backrest

A comfortable and safe support to lean against. Mounted on the handlebars and therefore follows the height adjustment.
Art no. 12510



Crossbar

Make it possible to handle the rollator with one hand.
Art no. 1531573 Right
Art no. 1536813 Left

Two heights (push handles and seat)



LEGACY 520

Art no. 12052-37-24 DA



LEGACY 600

Art no. 12050-37-24 DA

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please refer to your local Invacare website.

	Width between outside wheels	Depth	Width between wheels	Overall height	Seat height	Width between handles	Wheel diameter
LEGACY 520	600 mm	650 mm	455 mm	760 - 850 mm	540 mm	420 mm	200 mm
LEGACY 600	600 mm	720 mm	455 mm	835 - 920 mm	620 mm	420 mm	200 mm

	Load capacity	Max. user weight
LEGACY 520	7.2 kg	150 kg
LEGACY 600	7.6 kg	150 kg

Frame colour

RAL 5022



Please notice colors may vary slightly from those displayed above.

QR
Code

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Benkenstrasse 260
4108 Witterswil
Switzerland



Yes, you can.®



© 2018 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. Legacy - EU - 10/2018