

AQUATEC®
Ocean



The modular shower
commode to fit all needs



The *Aquatec Ocean / Ocean XL* shower commodes provide highest versatility for institutional and home care. With its height adjustment and the four different varieties it can be easily adapted to different users' body measurements and weights, giving everybody a comfortable and ergonomic sitting position.

Adaptability

The large number of possibilities to adjust single components of the *Ocean / Ocean XL* like the footrests, the backrest or even the height of the frame results in the ability to adapt the chair to many different users' sizes. With the great number of accessories, e.g. the 24"-wheels or the different soft seats, the *Ocean / Ocean XL* gains even more flexibility to exactly fit everybody's needs.

Corrosion free

Every part of the *Aquatec Ocean / Ocean XL* is either made of stainless steel or plastics providing that the chair stays perfectly corrosion free. The tough frame is powder coated which makes the *Ocean / Ocean XL* easy to clean and underlines its high-quality appearance.

Features and Options



Ocean XL

Provides more space between the armrests and a higher load capacity of 150kg.



Knock-down system

The *Ocean / Ocean XL* can be assembled and disassembled in minutes in order to use it for traveling.



Height adjustable - no tools necessary

No tools are needed to adjust the height and adapt the *Ocean / Ocean XL* to the user's needs.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Seat height	Seat width	Total width	Total depth	Total height	Width between armrests	Max. user weight	Total weight
<i>Ocean</i>	477 - 600 mm	480 mm	560 mm	900 mm	964 - 1094 mm	460 mm	150 kg (180 kg for XL-version)	14.5 kg
<i>Ocean 24"</i>	477 - 600 mm	480 mm	675 mm	1000 mm	964 - 1094 mm	460 mm	150 kg	19.5 kg

Colour

White



QR Code

Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Benkenstrasse 260
 4108 Witterswil
 Switzerland

