

# EC-Track™



Discrete, flexible and safe



*Invacare* **EC Track** system completes the versatility of the *Invacare Robin*. To meet the varied transfer needs, we offer a wide range of tracking options for different kind of environments and types of installation. The **EC Track** has a discrete modern design with a low height and long span, ensuring that the complexity of the various transfer needs are met. **EC Track** is constructed from light aluminium material and is easy to install.

Features and options



*Invacare* **EC Track**

Is available in 3 different profiles.



**Curves**

Are available in 30, 45, 60 and 90 degrees as standard. Only S-Profiles.



The *Invacare Robin* trolley



*Invacare Robin*

Connects to the track, without the use of tools.



**Standard bracket**

Available for all track options.



**The Quick bracket**

Single point installation (unavailable in the UK).



**Wall bracket**

An alternative to ceiling mounting.



**Wall support**

To support the wall-mounted track




**Traverse trolley**



**Transit coupling**


To transfer from one room to another, using only one **Robin** motor.

## Features and options




### Room to room transfer

With two *Invacare Robin* motors and a room-to-room trolley, it is possible to transfer the user from room to room.





### Wall support

The tracks can be installed in places where the wall does not have sufficient carrying capacity.




### Transit coupling

To transfer the user from one room to another, using only one *Robin* motor.




### Wall bracket

To install the track along the wall.



### Traverse systems

Permits an unlimited number of lifting places.



### Wall bracket

To install the track in a free span. The tracks are mounted at the ends.



### EC - Gantry

Tracks on „legs“- no need for installing tracks in wall or ceiling. EC - Gantry is very easy to install and dismount and once mounted it is very easy to raise and lower the gantry without using tools.

#### The EC - Gantry is available in two versions:

1. With castors - height adjustable from 195 - 230 cm
2. Without castors - height adjustable from 216 - 251 cm



### EC - Gantry Traverse

Room covering tracks on „legs“- no need for installing tracks in wall or ceiling. EC - Gantry Traverse is very easy to install and dismount and once mounted it is very easy to raise and lower the gantry without using tools. EC- Gantry is height adjustable from 218 to 253 cm.

## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Max. user weight




Invacare **EC-Track**

200 kg

## Frame colours

White

Aluminium



Please notice colors may vary slightly from those displayed above.

Tel: +41 61 487 70 70  
hqeurope@invacare.com  
www.invacare.eu.com

Invacare International GmbH  
Benkenstrasse 260  
4108 Witterswil  
Switzerland

QR  
Code

