




Yes, you can.®

On Track with EC Track!

A discrete track system that enables users to move with freedom around their home.



Move securely from one room to another
The two rooms are connected,
and transferring from one room to another is simply achieved.



In between hung traverse
ECTrack has a discreet
and stylish finish.





Side hung
Side hung tracks can be mounted
very close to the wall.

EC TRACK, EU.10/13 - ©Invacare® International Sàrl 2013. All rights reserved. Patents pending - ASPHÀTIE 334 860 285 R.C.S. TOURS



Invacare International Sàrl
Route de Cité-Ouest 2 - 1196 Gland - Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
hqeurope@invacare.com - www.invacare.eu.com

Yes, you can.®



EC Track brings two rooms together

In order to connect two or more rooms for mobility with the same motor unit, a transit coupling can be mounted.



EC Track system with traverse track under hung
 This configuration allows for more space between the track and the ceiling, which is an advantage if lighting fixtures or other elements are mounted on the ceiling.



EC Track system with end hung wall mounted parallel tracks
 In situations where it is not possible to mount the parallel tracks onto the ceiling, the track or tracks can be directly mounted to the wall.



EC Track operates with electronic ease

The switch for the transit coupling is electronically operated. For optimal security the transit coupling will be locked before transferring from one room to the other.