

AQUATEC®

Sorrento



Bathtub Seat



The *Aquatec Sorrento Bathtub Seat* is an outstanding modular and modern product of the *Aquatec Sorrento* range. This range includes great function with great style and design.

The non-slip soft touch surfaces as well as the soft backrest provide a secure and very comfortable seat for the user. The height adjustable front legs and the adjustable anti-slip attachment legs assure a secure positioning of the bathtub seat and will cover all standard bathtubs.



Sorrento – looks great and perfect in every bathroom

Non-slipping Soft Touch Surface

The *Aquatec Sorrento* is the first Bathroom collection with the innovative TPE Soft Touch Surface on the seat and the armrests. Not only the great look but also the comfortable feeling will persuade the user.

Easy adjustability

Tool-less, fast and easy adjustment – these are only some of the convincing benefits of the *Aquatec Sorrento Bathtub Seat*. The two height adjustable front legs as well as the anti-slip attachment legs are unique in this segment, assure a safe placement in the bathtub and provide great legroom for the user.



Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Sorrento Bathub
Seat	
Bathtub Seat	1532715
Measurements	
Seat width (mm)	400
Seat depth (mm)	400
Total width (mm)	710 - 760
Total depth (mm)	560
Total height (mm)	670 - 725
Width between armrests (mm)	450

Measurements	
Armrest length (mm)	280
Armrest height (mm)	195
Backrest width (mm)	500
Backrest height (mm)	380
Backrest cover height (mm)	280
Width hygiene recess (mm)	140
Depth hygiene recess (mm)	170
High load capacity	130 kg
Weight	4,5 kg
Colour	

Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Benkenstrasse 260
 4108 Witterswil
 Switzerland

