



Invacare[®] | **Sunto**



Invacare **Sunto** offers an excellent support for those users who need to reactivate the ambulatory function and when users require an extra support with a good walking position.

Technical Data

- Painted steel tube Ø 25 mm
- Double cross with plungers
- Push handle
- PU axillary supports adjustable in height, width and depth
- Four pivoting Ø 80 mm wheels
- Also available with brakes on the rear wheels. (Invacare **Sunto-B**)

Foldable / Adjustable

Simple and balanced folding with the help of the release button. Easy to adjust.

Accessories

Comfortable padded seat.
Item No. 1525063
Forearm support.
Item No. 1525064




Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please refer to your local Invacare website.

	Invacare Sunto
Width Inside Base legs	600 mm
Overall depth	Open 450 mm Folded 120 mm
Overall width	Open 680 mm Folded 120 mm
Height adjustment	1180 mm - 1450 mm
Product weight	10.5 kg
Max user weight	100 kg

Colour

Dark Blue



Art. nr.

1524974 Invacare **Sunto**
1524975 Invacare **Sunto-B**

QR
Code

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Benkenstrasse 260
4108 Witterswil
Switzerland



Yes, you can.®

