

AQUATEC®

# Grab handles



Grab handles, support folding handles





Invacare **H140 - H144**

- Grab handles made from high-quality plastic, particularly impact and scratch-resistant
- Available in five different lengths
- Ergonomically shaped for a comfortable grip



Invacare **H335**



Invacare **H345**

- Space-saving: mounted on the wall and can be folded when not in use
- The support folding handle is locked in the upper position
- Also available with height-adjustable supporting leg for additional stability
- Supporting leg automatically folds in when the handle is folded upwards



Invacare **H160**

- Grab handle for the bathtub
- Provides additional safety when getting into and out of the bath
- Easy assembly using quick release screws no tools required
- Safely secured on the edge of the bath using a clamping mechanism made of steel and a pressure plates with anti-slip surface



Invacare **H170 Stileo**



Invacare **H171 Stileo**



Invacare **H172 Stileo**


- Suction grab handle for the bath to support natural movement
- Simple and easy assembly using vacuum lever
- Ergonomical shape and structured surface ensures a secure grip, available in two lengths
- For a short term adoption or as temporary supply

## Technical data


For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Order numbers	Invacare H140	Invacare H141	Invacare H142	Invacare H143	Invacare H144	Invacare H335	Invacare H345	Invacare H160	Invacare H170	Invacare H171	Invacare H172
Grab handle (pack of 1)					H144						
Grab handle (pack of 3)	1525875	1525876	1525877	H143							
Support folding handle						H335					
Rest folding handle with supporting leg							H345				
Bath tub handle (pack of 4)								1522374			
Suction grab handle (pack of 3)									1606999	1526476	1526477
Suction grab handle (pack of 1)									1606979	1546225	1546226
<b>Dimensions</b>											
Length (A)	250 mm	400 mm	500 mm	600 mm	800 mm	-	-	-	290 mm	400 mm	500 mm
Width (B)	35 mm	35 mm	35 mm	35 mm	35 mm	40 mm (A)	40 mm (A)	-	35 mm	35 mm	35 mm
Depth	-	-	-	-	-	620 mm (B)	620 mm (B)	-	-	-	-
Height (C)	105 mm	105 mm	105 mm	105 mm	105 mm	135 mm	135 mm	-	72mm	75 mm	75 mm
Total width (D)	65 mm	65 mm	65 mm	65 mm	65 mm	135 mm	135 mm	215 mm (A)	95 mm (E)	100 mm (E)	100 mm (E)
Total length	-	-	-	-	-	-	-	-	300 mm	450 mm (D)	550 mm (D)
Total depth	-	-	-	-	-	700 mm (E)	700 mm (E)	170 mm (B)	-	-	-
Total height	-	-	-	-	-	220 mm (F)	220 mm (F)	550 mm (C)	-	-	-
Inside length (E)	125 mm	275 mm	375 mm	475 mm	675 mm	-	-	-	105 mm	165 mm (F)	270 mm (F)
Height, supporting leg (G)	-	-	-	-	-	-	635 - 645 mm	-	-	-	-
Depth, folded	-	-	-	-	-	170 mm	170 mm	-	-	-	-
Height, folded	-	-	-	-	-	820 mm	840 - 850 mm	-	-	-	-
Height from the bathtub edge	-	-	-	-	-	-	-	405 mm (D)	-	-	-
Adjustment range, inner dimension (E)	-	-	-	-	-	-	-	75 - 150 mm	-	-	-
Ø Suction cup	-	-	-	-	-	-	-	-	80mm	93 mm	93 mm
<b>More info</b>											
Weight	0.3 kg	0.5 kg	0.6 kg	0.7 kg	0.9 kg	3.0 kg	3,5 kg	2.3 kg	0.35 kg	0.5 kg	0.6 kg
Load capacity	120 kg	120 kg	120 kg	120 kg	120 kg	120 kg	120 kg	120 kg	30 kg	30 kg	30 kg
Colour	White	White	White	White	White	White	White	White	Blue-white	Blue-white	Blue-white


Invacare H140 - H144




Invacare H335/H345



Invacare H160



Invacare H170 / H171 / H172 Stileo





Tel: +41 61 487 70 70  
hqueurope@invacare.com  
www.invacare.eu.com

Invacare International GmbH  
Benkenstrasse 260  
4108 Witterswil  
Switzerland



© 2018 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. Grabhandles - EU - 10/2018