

AQUATEC®

Grab handles



Grab handles, support folding handles



Yes, you can.®



Invacare **H320/1**

- Grab handle made from high-quality plastic, particularly impact and scratch-resistant
- Oval and angled shape for better grip and support



Invacare **H330/1**

- Space-saving: mounted on the wall and can be folded when not in use.
- Toilet paper holder integrated



Invacare **H140 - H142**

- Grab handles made from high-quality plastic, particularly impact and scratch-resistant
- Available in three different lengths
- Ergonomically shaped for a comfortable grip



Invacare **H335**

Invacare **H345**

- Space-saving: mounted on the wall and can be folded when not in use
- The support folding handle is locked in the upper position
- Also available with height-adjustable supporting leg for additional stability
- Supporting leg automatically folds in when the handle is folded upwards
- In different lengths and with different options available. Not all products, models or options may be available in every country. For detailed information contact your local dealer.



Invacare **H170 Stileo**



Invacare **H171 Stileo**



Invacare **H172 Stileo**

- Suction grab handle for the bath to support natural movement
- Simple and easy assembly using vacuum lever
- Ergonomical shape and structured surface ensures a secure grip, available in three lengths
- For a short term adoption or as temporary supply

Technical data


For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

| Order numbers | Invacare H320/1 | Invacare H330/1 | Invacare H140 | Invacare H141 | Invacare H142 | Invacare H335 | Invacare H345 | Invacare H170 | Invacare H171 | Invacare H172 |
|---|-------------------|-----------------|---------------|---------------|---------------|---------------|---------------|---------------|---------------|---------------|
| Grab handle (pack of 3) | H320/1 | | 1525875 | 1525876 | 1525877 | | | | | |
| Support folding handle | | H330/1 | | | | H335 | | | | |
| Rest folding handle with supporting leg | | | | | | | H345 | | | |
| Suction grab handle (pack of 3) | | | | | | | | 1606999 | 1526476 | 1526477 |
| Dimensions | | | | | | | | | | |
| Length (A) | 330 (angled 135°) | - | 250 mm | 400 mm | 500 mm | - | - | 290 mm | 400 mm | 500 mm |
| Width (B) | 40 mm | - | 35 mm | 35 mm | 35 mm | 40 mm (A) | 40 mm (A) | 35 mm | 35 mm | 35 mm |
| Depth | - | - | - | - | - | 620 mm (B) | 620 mm (B) | - | - | - |
| Height (C) | 76 mm | - | 105 mm | 105 mm | 105 mm | 135 mm | 135 mm | 72mm | 75 mm | 75 mm |
| Total width (D) | - | 100 mm | 65 mm | 65 mm | 65 mm | 135 mm | 135 mm | 95 mm (E) | 100 mm (E) | 100 mm (E) |
| Total length | - | 670 mm | - | - | - | - | - | 300 mm | 450 mm (D) | 550 mm (D) |
| Total depth | - | - | - | - | - | 700 mm (E) | 700 mm (E) | - | - | - |
| Total height | - | 205 mm | - | - | - | 220 mm (F) | 220 mm (F) | - | - | - |
| Inside length (E) | 245 mm | - | 125 mm | 275 mm | 375 mm | - | - | 105 mm | 165 mm (F) | 270 mm (F) |
| Height, supporting leg (G) | - | - | - | - | - | - | 635 - 645 mm | - | - | - |
| Depth, folded | - | - | - | - | - | 170 mm | 170 mm | - | - | - |
| Height, folded | - | - | - | - | - | 820 mm | 840 - 850 mm | - | - | - |
| Ø Suction cup | - | - | - | - | - | - | - | 80mm | 93 mm | 93 mm |
| More info | | | | | | | | | | |
| Weight | 0,18 kg | 2,2 kg | 0.3 kg | 0.5 kg | 0.6 kg | 3.0 kg | 3,5 kg | 0.35 kg | 0.5 kg | 0.6 kg |
| Load capacity | 115 kg | 130 kg | 120 kg | 120 kg | 120 kg | 125 kg | 125 kg | 30 kg | 30 kg | 30 kg |
| Colour | White | White | White | White | White | White | White | Blue-white | Blue-white | Blue-white |


Invacare H320/1




Invacare H330/1




Invacare H140 - H142



Invacare H335/H345



Invacare H170 / H171 / H172 Stileo



Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland

