

# Propad Leg Trough



Shaped leg support to facilitate  
comfort and recovery



The *Invacare Propad Leg Trough* has been designed in consultation with clinicians to provide maximum support for the whole leg.

### Offers comfort and support to facilitate recovery

Provides optimum support and protection for the whole leg, reducing pressure and facilitating recovery.

### Relieves pressure on vulnerable heels

Elevates the heels and allows them to be suspended over the edge of the trough, providing pressure relief for the vulnerable heel area.

## Features and Options



### Water resistant cover

Two-way stretch, vapour permeable and anti-static polyurethane (PU) water resistant cover.



### Available in three sizes

Wide range of patients can be accommodated.



### Cover includes a zip, and is removable

Facilitates laundering of the cover.



## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Length	Width	Height	Total product weight	Max. cleaning temperature
<i>Propad Leg Trough</i>	510, 610, 660 mm	110 mm	115 mm	0.8 kg	80°

Manufactured to European Standards:  
 Interior Foam: BS 5852: 1982 Crib 5, BSEN597-1,2  
 Cover: BS7175: 1989 Crib 5

## Colour

Blue



QR Code

Tel: +41 61 487 70 70  
[hqeuropa@invacare.com](mailto:hqeuropa@invacare.com)  
[www.invacare.eu.com](http://www.invacare.eu.com)

Invacare International GmbH  
 Benkenstrasse 260  
 4108 Witterswil  
 Switzerland

