

AQUATEC®

Sorrento



Shower Collection



Yes, you can.

The *Aquatec Sorrento Shower Chair* and Stools are part of the outstanding modular and modern products of the *Aquatec Sorrento* range. The non-slip soft touch surfaces as well as the soft backrest provide a secure and very comfortable seat for the user. The height adjustable, folding legs as well as the fast and easy assembly cover various users' needs.

The all-corrosion free, modular shower chair range is the best choice for users who need assistance in the bathroom for their personal hygiene.

Non-slipping Soft Touch Surface

The *Aquatec Sorrento* is the first shower chair with the innovative TPE Soft Touch Surface on the seat and the armrests. Not only the great look but also the comfortable feeling will persuade the user.

Height adjustable and folding legs

The *Aquatec Sorrento range* can be tool less adjusted to individual body sizes. With its folding mechanism it is ideal to take on vacation or for easy storage. Backrest and armrests are fast and easy to take off.



Sorrento – looks great in every bathroom



Sorrento as Shower Chair



Sorrento as Shower Stool with armrests



Sorrento as Shower Stool



Sorrento as Shower chair without recess

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Sorrento Shower Chair	Sorrento Shower Stool w/armrests	Sorrento Shower Stool	Sorrento Shower Chair w/o recess
Shower Chair	1532633			1532730
Shower Stool w/armrests		1532660		
Shower Stool			1532661	
Measurements				
Seat width (mm)	400	400	400	400
Seat depth (mm)	400	400	400	400
Seat height (mm)	460 - 610	460 - 610	460 - 610	460 - 610
Total width (mm)	520	520	520	520
Total depth (mm)	430	430	430	430
Total height (mm)	850 - 990	850 - 990	850 - 990	850 - 990
Range seating height adjustment (mm)	150	150	150	150
Width between armrests (mm)	450	450		450
Armrest length (mm)	280	280		280
Armrest height (mm)	195	195		195
Backrest width (mm)	500			500
Backrest height (mm)	380			380
Backrest cover height (mm)	250			250
Width hygiene recess (mm)	140	140	140	
Depth hygiene recess (mm)	170	170	170	
High load capacity	135 kg	135 kg	135 kg	135 kg
Weight	5,5 kg	4 kg	3,5 kg	5,5 kg
Colour				



Tel: +44 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland



© 2023 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. EU – 07/23

Yes, you can.®