

AQUATEC®

Bali



The sturdy bath board



Yes, you can.®

The *Aquatec **Bali*** serves for a safe seat on all standard bathtubs. It allows the user to safely and quickly wash with the help of the adjustable shower head whenever it is preferred to a bath.




Sturdy, stable and secure

Load capacity of 175 kg!

Safe and comfortable use

Rounded edges facilitate secure transfer. The non-slip structured surface provides safe sitting.

Features and Options



Adaptable

Secure fixation on different sizes of bathtubs with adjustable stoppers.



Stable handgrip

Offers additional safety.



Compartment for soap

Convenient and practical.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Length	Width	Width	Width	Total weight	Max. User weight
<i>Aquatec Bali</i>	695 mm	300 mm	470 - 610 mm	220 mm	1,8 kg	175 kg

Frame colours

White-Blue



Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Neuhofweg 51
 4147 Aesch BL
 Switzerland



© 2023 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. EU - 07/23

Yes, you can.®