



INVACARE® TOP END®

XLT HANDCYCLE



Yes, you can.®



If recreation is your goal, the Invacare® Top End® XLT Handcycle is great fun and features low maintenance 7 speed hubs with reverse hands-on braking. Just get out there and go! All Invacare Top End recreation handcycles are equipped with a safety flag, full composite chain guard, cushion and seat/chest restraint, as standard, for extra value, comfort and safety.

The **Invacare® Top End® XLT Handcycle** is designed with the recreational handcycling enthusiast in mind. Equipped with seven speeds, the XLT handcycle will easily handle most terrain, whether it's riding around town or taking a spin in the country. The lightweight, oval tube aluminum XLT frame with chrome-moly steel fork features simple geometry for maximum stiffness and easy cornering. The seat and back are totally adjustable for a comfortable fit. The narrow rear-end enhances aerodynamics and the 15° of camber increases stability. If you want a great way to exercise, cross-train or just have fun, choose the XLT handcycle.

Model #: XLT



1. Twist Shifter
2. Quad Twist shifting option
3. QuadGrip handpedals with 3" width extension option



Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End
www.topendwheelchair.com
Sports and Recreation Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(727) 522-8677
(800) 532-8677
Fax: (727) 522-1007

©2014
Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols™ and®.
All trademarks are owned
by or licensed to Invacare
Corporation unless
otherwise noted. Shimano
is a registered trademark of
Shimano Inc.
Form No. 12-510
Rev. 12/14
141163

Features

Choice of 15" wide or 17" wide seat widths
6061 Aluminum frame/chrome-moly fork
Adjustable seat, back angle and crank height
Full chain guard, chain tension idler and safety flag
Adjustable footrest system with safety straps
Seat and back cushions with chest restraint
Shimano® 7 speed hub with reversing drum brake
175 mm crankset with one chain ring
Quad Twist shifter option helps those with hand impairments

Frictionless ovalized ergonomic hand pedals
Parking brake
26" High Performance wheels
High pressure clincher tires
Quick release axles
15° camber
Road crown compensator
Reflectors
250 lb. weight capacity
Many options and aftermarket accessories are available

