

AQUATEC®
Beluga



Independent and safe bathing
in one's own bathtub



Yes, you can.®

The *Aquatec Beluga* is designed to be used in domestic home environments. The bathlift gives people who feel insecure in their movements the possibility to take a bath safely and independently. While keeping the privacy the personal hygiene can be realized in a most relaxing atmosphere.



Maximum bathing comfort

At push of a button of the ergonomic floating hand control the user gets smoothly and safely in and out of the bathtub. The reclining backrest provides a comfortable and relaxing bathing using the full capacity of the bathtub.

Easy handling

Due to its easily detachable components and the self-releasing suction cups the *Aquatec Beluga* can be put in and taken out of the bathtub in no time. Smooth surfaces and the detachable machine-washable covermats (up to 60°) facilitates the cleaning of the maintenance-free bathlift.

Features and Options



Aquatec Beluga without covermat Transfer aid

Due to the closed, non-slip surface the bathlift can also be used without cover.



Various optional transfer aids facilitates the transfer.












Accessories

Extensive accessories for an individual adaptation.


Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Back tilt angle	Seat height	Seat width	Total width	Seat depth	Total weight	Battery weight	Max. User weight
<i>Aquatec Beluga</i>	 38°	 Max. 420 mm Min. 60 mm	 380 mm	 700 mm	 460 mm	 12,9 kg	 1,6 kg	 135 kg



Frame colours



Tel: +44 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland



Yes, you can.®

