



DOLOMITE

Futura



The universal rollator



Dolomite Futura has been developed to be able to satisfy the needs of the vast majority of people who use rollators for operational adaptability. It is available in four different seat heights and it has one of the largest ranges of accessories on the market. With **Futura** the opportunities for personal customization are almost unlimited.



Infinite adjustable push handles

Easily adjusted handles that can be positioned at exactly the right height to compensate different arm lengths, and at exactly the right angle to follow the body's natural motion pattern.

Walking area

Large gate area inside the rollator frame is a signum for Dolomite. The flip-up seat gives premium walking space. And makes it easier to reach and access things. In front of you. The rear forks gives a smooth inside, being friendly to toes and shoes.

Structured seat with anti-slip function, rounded rear edge and generous sitting area, means that the user can sit securely and comfortably for a long time.



User friendly brakes

Long brake handle and brake stopper - provides good leverage, simple to brake even with weak hand strength.

Short distance to park-position - reduces risk of losing balance. Five year warranty on brakes with reinforced cable housing, stainless steel brake stopper.

Features and options



Tray

For transportation of a glass, coffee, mug or similar. Can be moved to a table.
Art no. 70311



Basket

Different baskets are available for different height model.
Art no. 12530-15: 600, 520
Art no. 12582-15: 680
Max load capacity 9 kg
Art no. 12538-15: 450
Max load capacity 7 kg



Cane holder

Makes the cane easily accessible.
Art no. 12606



Backrest

A comfortable and safe support to lean against. Mounted on the handlebars and therefore follows the height adjustment.
Art no. 12600

Features and options



Kerb climber

“Pedal” that simplifies passing kerbs or similar. Mounted on the right side. Art no. 12560-21



Bag

Provides extra room when the basket is quickly filled. Can be combined with oxygen cylinder holder. Art no. 12523



One-hand brake

Makes it possible to activate both brakes with one hand. Art no. 12507-R (right). Art no. 12507-L (left)



Brake booster

Provides extra braking force for weak hands. Art no. 12651



Soft grip

A softer alternative to the anatomical grips. Art no. 12593 (pair)



Slow-down brake

A safety feature preventing the walker from rolling away. Perfect for gait training and one handed solutions. Art no 13041-21 (pair).



Adapt-handle set

Multi adjustable handle to get optimized position. State size of rollator when ordering. Art no. 12613 (pair)

HP-handle set

For adaptation for one-handed/stroke patients. State size of rollator when ordering. “Dolomite **HP**”. Art no. 12611-R; HP right, one-hand brake left. Art no. 12611-L; HP left, one-hand brake right



Oxygen cylinder holder

The oxygen cylinder is centrally positioned on the rollator and thereby provides good maneuverability. Art no. 12682-25 Fits most cylinders.



Wheel guard

Protects clothing from being dirty from wheel splash. Art no. 12544 (pair)



Crossbar

Make it possible to handle the rollator with one hand Art no. 1559174 Right Art no. 1559175 Left



Screw instead of knob

Replace height adjustment knob with the screw when adjustment is not desired. Art no. 60140

Four heights (push handles and seat)



FUTURA 450

Art no. 12074-43-21



FUTURA 520

Art no. 12072-43-21



FUTURA 600

Art no. 12070-4321



FUTURA 680

Art no. 12076-43-21

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please refer to your local Invacare website.

| | Width between outside wheels | Depth | Width between wheels | Overall height | Seat height | Width between handles | Load capacity |
|-------------------|------------------------------|--------|----------------------|----------------|-------------|-----------------------|---------------|
| Futura 450 | 600 mm | 590 mm | 460 mm | 640 - 790 mm | 475 mm | 420 mm | 7.6 kg |
| Futura 520 | 610 mm | 660 mm | 470 mm | 710 - 850 mm | 545 mm | 420 mm | 7.8 kg |
| Futura 600 | 610 mm | 700 mm | 470 mm | 790 - 930 mm | 620 mm | 420 mm | 8 kg |
| Futura 680 | 640 mm | 780 mm | 500 mm | 860 - 1015 mm | 705 mm | 420 mm | 8.23 kg |

| | Wheel diameter | Max. user weight |
|----------------------------------|----------------|------------------|
| Futura 450, 520, 600, 680 | 200 mm | 150 kg |



Frame colours

RAL 7021



Please notice colors may vary slightly from those displayed above.

QR Code

Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Benkenstrasse 260
 4108 Witterswil
 Switzerland



Yes, you can.®

