

Crutches



Standard crutches



Invacare® **Advantik** Invacare® **Advantik Tall**



Invacare® **Advantik Soft**



Invacare® **P2832 - P2831** Invacare® **Big 150**

Advantik is a range of stable light weight crutches with a weight capacity of 130 kg. The click-together function is convenient when quickly needing to get one hand free. Reflector is standard. **Advantik Soft** has an ergo-nomic shape and a soft surface to give good comfort.

Invacare® Childrens crutch has double adjustment in stylish colours to meet the children's needs.

A sturdy crutch with weight capacity for users up to 150 kg. The 1080 mm crutch can be cut to fit any user. Adjustable arm support.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please refer to your local Invacare website.

Product	Code	Floor to handle	Handle to cuff	Weight apacity	Weight/unit
Invacare® Advantik P2800 Grey	P2800	750-980 mm	240 mm	130 kg	0.618 kg
Invacare® Advantik P2802 Yellow	P2802	750-980 mm	240 mm	130 kg	0.618 kg
Invacare® Advantik P2803 Red	P2803	750-980 mm	240 mm	130 kg	0.618 kg
Invacare® Advantik P2804 Bordeaux	P2804	750-980 mm	240 mm	130 kg	0.618 kg
Invacare® Advantik P2805 Blue	P2805	750-980 mm	240 mm	130 kg	0.618 kg
Invacare® Advantik Soft	1529865	750-980 mm	240 mm	130 kg	0.635 kg
Invacare® Advantik Tall	1529863	850-1080 mm	240 mm	130 kg	0.635 kg
Invacare® Advantik children	P2831	495-715 mm	175-225 mm	100 kg	0.4 kg
Invacare® Advantik children	P2832	495-715 mm	175-225 mm	100 kg	0.4 kg
Invacare® Basik	P227S	720-950 mm	240 mm	100 kg	0.515 kg
Invacare® Basik (Red/Green/Blue)	P228 R/G/B	700-1030 mm	230 mm	100 kg	0.61 kg
Invacare® Big 150	1529868	1080 mm Cut to decide length	240-265 mm	150 kg	0.75 kg

QR
Code

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Benkenstrasse 260
4108 Witterswil
Switzerland



Yes, you can.®

