


Softform[®]

Premier Visco

Pressure-redistributing viscoelastic mattress

An exceptionally comfortable support surface
that conforms to the person's body shape.



Stability and pressure redistribution

For personalised pressure injury management



Yes, you can.[®]

Effective pressure injury prevention



Pressure Redistribution

- › High quality mattress facilitates pressure reduction by minimising peak pressures.
- › Deep viscoelastic insert pad is high-density and temperature-sensitive to the body, allowing it to conform to a person's shape.
- › Modular construction enables parts to be replaced easily and increases the longevity.



Shear & Friction Reduction

- › Two-way stretch, water-resistant, polyurethane cover is durable and helps to reduce shear and friction.



Comfort

- › High quality mattress delivering an exceptionally comfortable support surface.
- › Viscoelastic foam conforms around the body for custom support and comfort.
- › Four-sided concealed welded zip reduces potential fluid ingress and allows easy detachment for inspection and laundering.

Technical data ▼

For more information about this product, including the user manual, please visit our website: www.invacare.eu.com

Sizes available ▼

		850 mm	880 mm	900 mm	1000 mm	1050 mm
150 mm height	1970 mm	–	•	–	–	•
	2000 mm	•	•	•	•	–
	2080 mm	–	–	•	–	–



Max. User Weight ▼

247.6 kg



Max. Cleaning Temperature ▼

95 °C

Manufactured to BS EN 597-1 Cig, BS EN 597-2 Match, Foam Interior: BS 5852:1990 (Ignition Source 5); Cover BS 7175:1989 (Ignition Source 5).

Essential nursing care is pivotal in pressure injury management. This mattress may positively contribute to the outcome of a pressure injury prevention plan that includes education, clinical judgment and action planning based on vulnerability, which are fundamental factors in the prevention of pressure injuries. A range of assessment scales can be used as a formal method of assessing risk from pressure injury development, and should be used in conjunction with an informal assessment (informed nursing judgement) as part of a comprehensive and individualised plan.

Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland

