



Soprano



The big wheeler in tune



The *Dolomite Soprano* is a lightweight outdoor walker with large wheels, providing excellent skills on rough uneven ground. Easy to deal with as it can be folded and stands in the folded position. **Soprano** has heavy-duty functions, in a lightweight suit.



Takes you anywhere

Easy to push thanks to large, light and solid wheels. The excellent driving skills on uneven surfaces and the sturdy frame taking 160 kg gives mobility for everyone and everywhere.

Stands in folded position

Makes the walker easy to park on a small surface. Ideal when visiting crowded places, the theatre or for storing at home.

Forward facing

Handles and a flip up seat. Together with a wide frame, this gives the user a large gait area and stable rollator for outdoor use.

Features and options



Tray

For transportation of a cup of coffee or similar. Can be moved to a table.
Art no. 70311



Cane holder

Makes the cane easily accessible.
Art no. 12606



Backrest

A comfortable and safe support to lean against. Mounted on the handles and therefore follows the height adjustment.
Art no. 12628

Features and options



Kerb-climber

“Pedal” that simplifies passing kerbs or similar.
Art no. 12563



Wheel guard

Protects clothing from being dirty from wheel splash.
Art no. 12544 (pair)



One-hand brake

Makes it possible to activate both brakes with one hand. Can only be ordered as complete product
Art no. 12030-44-21L (left)
12030-44-21R (right)



Drip-feed holder

Holds drip feed on place. Height adjustable. No tools needed.
Art no. 12547



Adapt handle set

Multi adjustable to get optimized hand position. Can only be ordered as complete product.
Art no. 12030-44-21 ADAPT



Screw instead of knob

Replace height adjustment knob with the screw when adjustment is not desired.
Art no. 13655 (pair)



Hanging bag

Provides extra room when the basket is quickly filled. Can be combined with oxygen cylinder holder.
Art no. 12523



HP-handle set

Adaptation for one-handed/ stroke patients. Can only be ordered as complete product.
Art no.
12030-44-21 HPL (left)
12030-44-21 HPR (right)



Slow-down brake

A safety feature preventing the walker from rolling away. Perfect for gait training and one handed solutions.
Art no. 1452396-L (left)
1452396-R (right)



Oxygen cylinder holder

The oxygen cylinder is centrally positioned on the rollator and thereby provides good maneuverability. Fits most cylinders.
Art no. 12683



Basket

Spacious basket with handles. Max load capacity 12 kg.
Art no. 12541

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please refer to your local Invacare website.

 <p>Width between outside wheels</p>	 <p>Depth</p>	 <p>Width between wheels</p>	 <p>Overall height</p>	 <p>Seat height</p>	 <p>Width between handles</p>	 <p>Wheel diameter</p>
Dolomite Soprano 655 mm	750 mm	545 mm	790 - 990 mm	620 mm	440 mm	260 mm
 <p>Product weight</p>	 <p>Max. user weight</p>	 <p>Width folded</p>				
Dolomite Soprano 10,8 kg	160 kg	330 mm				



Frame colour

RAL 1035



Please notice colors may vary slightly from those displayed above.

Art. nr.

Soprano 12030-44-21

QR
Code

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Benkenstrasse 260
4108 Witterswil
Switzerland



Yes, you can.®

