

Flo-tech Lumbar Pad



Improving posture,
comfort and support



The *Invacare Flo-tech Lumbar Pad* is designed to improve posture, comfort and support; the **Lumbar Pad** is simple and easy to use. It is fitted with adjustable straps and is available in a black fabric cover.

Features and Options



Adjustable Straps

Provides the user with improved comfort, adjustability and positioning.





Low Memory Visco Foam

The **Lumbar Pad** is made up of two foams; Low density Visco foam on one side and Blue foam on the other. This combination of foams provide the user with added support on their lower back, as the Blue foam increases the density of the product whilst the Low density Visco foam absorbs the pressure from the user.

Technical data


For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Size	Max. cleaning temperature
		
Flo-tech Lumbar Pad	420 x 200 x 50 mm	80°


The Lumbar Pad is manufactured to European Standards: BS5852: 1990 Crib 5

Colour

Black fabric



QR Code



Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Benkenstrasse 260
4108 Witterswil
Switzerland

