

9630 Styxo



The *Invacare 9630 Styxo* combines two functions in one product. It can either be used as a toilet frame over the toilet or as a toilet commode with usage of the toilet pan. The toilet pan with lid and the splash guard are included in the delivery.




Product description:

- Easy to handle seating height adjustment in 7 steps from 41 cm to 56 cm
- Complete with elongated toilet seat and lid, pan with lid and splash guard
- Tool-less assembly and easy to clean
- High load capacity: 160 kg

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Order numbers	9630 Styxo
Toilet Commode	9630E
Dimensions	
Seat width	360 mm
Seat depth	420 mm
Seat height	410-560 mm
Total width	620 mm
Total depth	490 mm
Total height	600-750 mm
More info	
Load capacity	160 kg
Weight	5,5 kg
Colour	Grey



Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Benkenstrasse 260
 4108 Witterswil
 Switzerland



Yes, you can.®