

iPortal2™



Media and mobility that
keeps you connected



iPortal2 is our most advanced wheelchair media-connectivity system yet. This pioneering solution allows wheelchair users the ability to connect with personal computers and portable computing devices from their power wheelchair joystick or special input device.

The new **iPortal2** module includes:

- A real-time display for a wheelchair's settings and battery status.
- Bluetooth pairing to PC and iOS devices.
- Simple switching between devices.



Dashboard App

The **iPortal2** Dashboard Application displays wheelchair information in real-time using Bluetooth technology to connect with the wheelchair's operating system. It displays the speed range and seat adjustment indicators as well as a large battery status gauge. When connected to the iPhone 3GS (or later versions) users can also access a compass heading.

Accessibility

iPortal2's ground-breaking functionality makes a significant difference to users' lives - accessibility. It's a pioneering solution that gives wheelchair users the ability to control an iDevice via their wheelchair joystick or special input device. Apple's Voiceover and Assistive Touch enables wheelchair users to control the Apple device. Surfing the net, making calls, playing music is now easy and accessible without the need to touch the Apple device.

Features and options



Mouse Mover

With the Mouse Mover feature enabled with a simple Bluetooth pairing, the wheelchair's joystick can be used as a computer mouse pointer, allowing the user to navigate and interact with their PC or laptop.

Compatibility

The **iPortal2** uses standard Bluetooth technology, whether built-in or via a dongle (Bluetooth 2.1 or higher) and is compatible with the iOS (iPod, iPhone, iPad), Mac, Windows and Linux operating systems as well as tablet computers. It is available to Dynamic Controls' Shark, DX & DX2 power wheelchairs controllers.



QR Code

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Benkenstrasse 260
 4108 Witterswil
 Switzerland

