

# H340 Izzo



The *Invacare H340 Izzo* can be used over the toilet as toilet frame or toilet seat raiser as well as a commode with the toilet pan with lid – depending on your individual requirements. The easy to operate flip down armrests and the height adjustable seat support the transfer.



### Product Description:

Easy to handle seating height adjustment in 5 steps from 45 cm to 55 cm

Complete with elongated toilet seat with lid and pan with lid


Tool-less assembly

High load capacity: 113 kg

### Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

| Order numbers  | H340 Izzo          |
|----------------|--------------------|
| Toilet Commode | H340               |
| Dimensions     |                    |
| Seat width     | 360 mm             |
| Seat depth     | 430 mm             |
| Seat height    | 450-550 mm         |
| Total width    | 470 mm             |
| Total depth    | 530 mm             |
| Total height   | 820-920 mm         |
| More info      |                    |
| Load capacity  | 113 kg             |
| Weight         | 8,8 kg             |
| Colour         | White / light grey |



Tel: +41 61 487 70 70  
[hqeuropa@invacare.com](mailto:hqeuropa@invacare.com)  
[www.invacare.eu.com](http://www.invacare.eu.com)

Invacare International GmbH  
 Benkenstrasse 260  
 4108 Witterswil  
 Switzerland



Yes, you can.®